

# Advances in GC Column Selectivity for Environmental, Food Sample and Pharma Applications.


Gustavo Serrano, Ph.D.  
Global Product Manager – GC Columns  
Chemistries and Supplies Division  
Agilent CrossLab Group

# Outline


- **Volatile Organic Compounds (VOCs)**  
Pharmaceutical and Environmental Applications
- **Polycyclic Aromatic Hydrocarbon (PAHs)**  
Food and Environmental Applications
- **Polychlorinated Biphenyls (PCBs)**  
Environmental Application
- **Dioxins**  
Food and Environmental Applications

# VOCs & Volatiles: Pharma and Environmental Applications


## Agilent JW 624 GC Columns



# VOCs & Volatiles Analysis



# J&W Portfolio --- Cyanopropylphenyl Liquid Phases




6% cyanopropylphenyl 94% dimethylpolysiloxane

## 624-type phases

- DB-624 UI
- DB-624
- VF-624ms
- CP-Select 624 CB
- DB-Select 624 UI for <467>
- CP-Select for Hexanes

## 1301-type phases


- DB-1301
- VF-1301ms



14% cyanopropylphenyl 86% dimethylpolysiloxane

## 1701-type phases

- DB-1701
- DB-1701P
- VF-1701ms
- CP-Sil 19 CB



50% cyanopropylphenyl 50% dimethylpolysiloxane

## 225-type phases

- DB-225
- DB-225ms
- CP-Sil 43 CB

# The 624 is the Most Popular Column Choice for Volatiles Analysis

Typical GC phases

DB-5ms, 30m x 0.25 mm i.d. x 0.25  $\mu\text{m}$

DB-624, 30m x 0.25 mm i.d. x 1.4  $\mu\text{m}$

$$R_s = 1/4 \cdot \frac{\alpha - 1}{\alpha} \cdot \frac{k}{k+1} \cdot \sqrt{N_{th}}$$

$N_{th}$  Plate count

$k'$  Capacity factor

$\alpha$  Selectivity factor

Many VOCs & solvents have low BPs with therefore little retention/ $k$ . The thick-film creates higher retention/ $k$  and thus improves separation.

VOCs and solvents contain various functional groups. Cyano/phenyl/methyl mix of 624 provides different levels of interaction mechanisms: Dispersion, pi-pi stacking, dipole,  $\text{H}_2$  bonding. This provides excellent selectivity with improved separations across many analytes.



# 624 Applications & Workflows

## Waters



VOCs / Purgeables  
Nitrosamines

## Soils



VOCs

## Alcoholic Beverages



Alcohols & aromas

## Air



Volatiles & Hydrocarbons  
Halogenates / CFCs

## Industrial



Solvents & Volatiles

## Pharmaceuticals



Residual solvents

# Residual Solvents

- Most are headspace applications
  - Hand sanitizer could also be done by headspace
- Can use FID or MSD
- There are Intuvo and 8890 applications
- Variety of columns (select and regular 624 UI) and dimensions
  - DB-624 UI
  - DB-Select 624 UI for 467





# USP <467> Application Compares Select 624 to Innowax

- HS/GC/dual FID
  - Split using CFT
- Class I and Class II (A & B) Residual Solvents

## Class 2A solvents




Figure 4. Chromatogram of the USP residual solvents class 2A standard solution resolved on a J&W DB-Select 624 UI for 467 GC column.

Application Note  
Residual Solvent



Analysis of USP Method <467>  
Residual Solvents on the Agilent 8890  
GC System

### Author

Lukas Wieder, Jie Pan, and  
Rebecca Veeneman  
Agilent Technologies, Inc.  
2850 Centerville Road  
Wilmington DE 19808

### Abstract

This Application Note highlights the use of the Agilent 8890 GC and the Agilent J&W DB-Select 624 UI for 467 and Agilent J&W HP-INNOWax columns in the detection and confirmation of <467> residual solvents. The system meets all specifications required in USP Method <467>, and demonstrates excellent repeatability across several injections.

5994-0442EN

### 5994-0442EN Instrument Conditions

Column	J&W DB-Select 624 UI for 467, 30 m x 0.32 mm, 1.8 $\mu$ m (p/n 123-0334UI)
Carrier Column 1	Helium, constant flow, 2.0 mL/min
Oven	40 $^{\circ}$ C (5.0 min), ramp 18 $^{\circ}$ C/min to 180 $^{\circ}$ C (3.0 min)
Inlet	Split mode, 140 $^{\circ}$ C, split ratio 5:1

# Class I, II, & III RSA

- GC/FID/MS
- Headspace for 624 UI
- Liquid injection on Wax




Figure 2. GC/MS-SCAN and FID Chromatogram of 52 compounds of standard solution (Level 7) on an Agilent DB-624 column. The effluent split ratio is MSD: FID= 1:1.

No.	Name
1	Methanol
2	Pentane
3	Ethanol
4	Ethyl ether
5	1,1-Dichloroethene
6	Acetone
7	Isopropanol
8	Ethyl formate
9	Acetonitrile
10	Methyl acetate
11	Methylene chloride
12	tert-Butylmethyl ether
13	trans-1,2-Dichloroethene
14	Hexane
15	1-Propanol
16	Nitromethane
17	cis-1,2-Dichloroethene
18	2-Butanone
19	Ethyl acetate
20	2-Butanol
21	Tetrahydrofuran
22	Chloroform
23	1,1,1-Trichloroethane
24	Cyclohexane
25	Carbon tetrachloride
26	2-Methyl-1-propanol
27	1,2-Dimethoxyethane
28	Benzene
29	1,2-Dichloroethane
30	Isopropyl acetate
31	Heptane
32	1-Butanol
33	Trichloroethylene
34	Methylcyclohexane
35	1,4-Dioxane
36	Propyl acetate
37	4-Methyl-2-pentanone
38	Isoamyl alcohol
39	Pyridine
40	Toluene
41	Isobutyl acetate
42	1-Pentanol
43	2-Hexanone
44	Butyl acetate
45	Tetrahydrothiophene
46	Chlorobenzene
47	Ethylbenzene
48	m,p-Xylene
49	o-Xylene
50	Isopropylbenzene
51	Anisole
52	1,2,3,4-Tetrahydronaphthalene

Application Note  
Pharma & Biopharma



## Analysis of USP <467> Residual Solvents of Class 1, Class 2, and Class 3 using the Agilent 8890 GC/ FID /5977B MSD System

### Authors

Youjuan Zhang, Shun Na  
Agilent Technologies  
(Shanghai) Co. Ltd. Shanghai  
200131 P. R. China

### Abstract

An Agilent 8890 Gas Chromatograph (GC) configured with a flame ionization detector (FID) and a mass spectrometric detector (MSD) was used for the analysis of USP <467> residual solvents of Class 1, Class 2, and Class 3. There were 52 compounds with low boiling points introduced by an Agilent 7697A headspace sampler onto a DB-624 column, while 10 compounds with relatively high boiling points were introduced by the automatic liquid sampler onto a DB-WAX column. A purged two-way capillary flow technology (OFT) device was used to split the sample 1:1 onto FID and MSD.


This Application Note demonstrates excellent peak shape, resolution and great repeatability, which shows this FID and MSD dual-channel system is a powerful tool for qualitative and quantitative analysis of residual solvents.

# 5994-1488EN

### 5994-1488EN Instrument Conditions

Column	J&W DB-624 UI for 467, 60 m x 0.25 mm, 1.4 μm (p/n 122-1364UI)
Column 1	Helium, constant flow, 1.0 mL/min
Oven	40 °C (10 min), ramp 5 °C/min to 80 °C (0 min), ramp 12 °C/min to 200 °C (10 min)
Inlet	Split mode, 250 °C, split ratio 10:1

# Hydrochlorothiazide



C<sub>7</sub>H<sub>8</sub>ClN<sub>3</sub>O<sub>4</sub>S<sub>2</sub> MW 297.74

## HS:


Oven Temperature: 80 °C  
Loop Temperature: 100 °C  
Transfer Line Temperature: 120 °C  
Vial Equilibration (min): 30  
Sample Loop : 1mL

Standards: 0.075mg/mL MeoH, 0.125mg/mL  
EtoH, 0.02mg/ML isopropanol(ISTD)  
Sample :25mg/mL Hydrochlorothiazide

## GC:

Column: **DB- select 624 UI**  
30 m x 0.530 mm x 3 µm, P/N125-0334UI  
Carrier : He, constant flow mode, 4.5mL/min  
Inlet: 200 °C, split ratio 5:1  
OVEN: 40 °C for 8 min  
45 °C/min to 200 °C  
200 °C for 1 min  
Detector: FID 250 °C

# Hydrochlorothiazide




NO.	Compound	Retention time (min)		Area (mAU*s)		Theoretical Plate (N)	Resolution (Rs)	USP Tf
		RT	RSD (%) (n=6)	Area	RSD (%) (n=6)			
1	Methanol	2.43	0.05	165.75	1.60	15848	/	1.14
2	Ethanol	3.21	0.05	446.8	0.43	19565	9.47	1.05
3	Isopropanol (IS)	3.98	0.04	112.85	0.96	16902	7.28	1.02

# Ethylene oxide sterilization residuals in PPE

## DB-Select 624 UI, 30m x 0.53mm x 3um (P/N 125-0334UI)

- Engineered to optimize pharmaceutical residual solvents analysis per USP Method <467>
- Ultra inertness and low bleed
- Recommended columns for ethylene oxide detection



### Headspace:

Oven: 80°C

Loop: 90°C

Transfer Line: 100°C

Equilibration time: 30min

### GC:

Inlet: 200°C, split ratio: 3:1

Inlet liner: Ultra Inert, splitless, straight, 2mm id (p/n 5190-6168)

Column: DB-Select 624UI, 30mx0.53mmx3um (p/n125-0334UI)




Carrier gas: Nitrogen, 5mL/min

Oven: 40°C(1min), 15°C/min to 130°C(1min)

Detector: FID, 250°C

# Capabilities of 8860 and HP-5ms UI/DB-624 UI for Environmental Applications

- Volatile halogenated and aromatic hydrocarbons
- Flexibility of 8860
  - Allows three detectors on one system
  - FID, ECD, FPD



## Analysis of Drinking Water with the Agilent 8860 GC and 7697A Headspace Sampler

Application Note  
Environmental



### Authors

Ian Eisele  
Agilent Technologies, Inc.  
Wilmington, DE, USA  
  
Youjuan Zhang  
Agilent Technologies  
(Shanghai) Co. Ltd.  
Shanghai, P. R. China

### Abstract

The flexibility of the Agilent 8860 GC is demonstrated by an analysis requiring multiple types of detection and sample introduction. Volatile halogenated and aromatic hydrocarbons in water are sampled by headspace for detection by FID and ECD, and organophosphorus pesticides are detected by FPD. A single 8860 GC can be configured to run both analyses.

5994-1239EN

### 5994-1239EN Instrument Conditions

Column	J&W DB-624 UI for 467, 30 m x 0.25 mm, 1.4 μm (p/n 122-1334UI)
Column 1	Helium, constant flow, 1.0 mL/min
Oven	40 °C (2.0 min), ramp 6 °C/min to 120 °C (4.0 min)
Inlet	Split mode, 250 °C, split ratio 50:1

# VOC and the stability of the 8890 with the DB-624 UI GC Column

- The Chinese Ministry of Environmental Protection
  - HJ642-2013
- Demonstrates performance using a difficult matrix
  - Soil




Figure 1. Overlaid TIC SIM traces for six replicates of 20 µg/L samples.

Application Note  
Environmental



## Determination of Volatile Organic Compounds in Soil and Sediments

Using an Agilent 7697A Headspace Sampler, 8890 GC, and 5977B GC/MSD combined platform

### Author

Zhang Jie  
Agilent Technologies Co. Ltd  
Shanghai

### Abstract

The accurate detection of volatile organic compounds (VOCs) in contaminated soil and sediment is of particular importance. The Chinese Ministry of Environmental Protection has produced a headspace GC/MS method, HJ642-2013, for the analysis of VOCs in soil and sediments. This Application Note has followed China method HJ642-2013 and demonstrated the excellent instrument performance of Agilent 7697A headspace, 8890 GC, and 5977B MSD combined platform for the targeted analysis.

5994-0662EN

### 5994-0662EN Instrument Conditions

Column	J&W DB-624 UI for 467, 60 m x 0.25 mm, 1.4 µm (p/n 122-1364UI)
Column 1	Helium, constant flow, 1.2 mL/min
Oven	40 °C (2.0 min), ramp 8 °C/min to 90 °C (4.0 min), ramp 6 °C/min to 200°C (10 min)
Inlet	Split mode, 250 °C, split ratio 10:1

# Recent 624 Phase Applications

## **5994-1926EN**

Novel Residual Solvents Analysis of Cannabinoid Products with the Agilent Headspace-GC/MS System

## **5994-0442EN**

Analysis of USP Method <467> Residual Solvents on the Agilent 8890 GC System

## **5994-1488EN**

“Analysis of USP <467> Residual Solvents of Class 1, Class 2, and Class 3 using the Agilent 8890 GC/ FID /5977B MSD System”

## **5994-1239EN**

“Analysis of Drinking Water with the Agilent 8860 GC and 7697A Headspace Sampler”

## **5994-0662EN**

“Determination of Volatile Organic Compounds in Soil and Sediments:  
Using an Agilent 7697A Headspace Sampler, 8890 GC, and 5977B GC/MSD combined platform”


## **5994-1923EN**

“Analysis of Ethanol and Isopropanol in Alcohol based hand sanitizers by direct injection GC/FID”




# Food and Environmental Applications: PAHs


## GC Column Updates & Solutions




Benzo[1,2,3-cd]fluoranthene




Dibenz[a,e]pyrene




Benzo[a]pyrene




Benz[a]anthracene



Dibenz[a,h]pyrene

















Benzo[b]fluoranthene



Benzo[ghi]perylene


# Agilent GC-PAH Solutions

Sample Containment	Sample preparation	Standards	Inlet Supplies	Separation	System Supplies	Application Notes
						
						
	<p>EMR-Lipids</p>	<p>EPA 610 EPA 8270</p>		<p>DB-5ms UI DB-EUPAH Select PAH</p>		

# Consumables Selection Toolkit for PAHs in Food Workflow


5994-2016EN

Includes Checkout Standards, Column, Liners, Unions, Vials, Caps, filament, gas filters etc. needed to get started with PAH testing with suggested Agilent methods and compatible with PAH Workflow Kits

  
From Insight to Outcome

## Polycyclic Aromatic Hydrocarbons Analysis in Food

Using triple quadrupole GC/MS/MS: Consumables workflow ordering guide



Item	Part Number
<b>View <a href="#">MyList</a> of sample preparation supplies</b>	
Captiva EMR-Lipid cartridge, 3 mL, 300 mg	5190-1003
Ceramic homogenizers	5982-9312
Agilent positive pressure manifold processor (PPM-48)	5191-4101
Captiva EMR-Lipid cartridge, 6 mL, 600 mg	5190-1004
Bond Elut Jr PSA, 500 mg	12162042B
<b>View <a href="#">MyList</a> of standards</b>	
EU PAH(15+1) standard kit; 250 µg/mL	5190-0487
EPA PAH standard (calibration standard); 500 µg/mL	8500-6035
Please go to <a href="http://www.agilent.com/chem/standards">www.agilent.com/chem/standards</a> for additional volume and concentration standard options	
<b>View <a href="#">MyList</a> of Inlet Supplies</b>	
Inlet septa, Advanced green, nonstick, 11 mm, 50/pk	5183-4759
Inlet septa, Advanced green, nonstick, 11 mm, 100/pk	5183-4759-100
Ultra Inert Splitless, single taper, glass wool	5190-2293
Ultra Inert splitless single-taper liner with glass frit	5190-5112
Ultra Inert Gold seal, with washer, 1/pk	5190-6144
Ultra Inert Gold seal, with washer, 10/pk	5190-6145
Self-Tightening column nut, collared, inlet	G3440-81011
Self-Tightening column nut, collared, MSD	G3440-81013
Replacement collar for self tightening nut	G3440-81012
15% Graphite / 85% Vespel Ferrules, 0.4 mm i.d., 10/pk	5181-3323
5 µL ALS syringe, fixed needle, 23-26s/42/cone	5181-1273
5 µL ALS syringe, fixed needle, 23-26s/42/cone 6/pk	5181-8810
10 µL ALS syringe, fixed needle, 23-26s/42/cone	5181-1267
10 µL ALS syringe, fixed needle, 23-26s/42/cone 6/pk	5181-3360
20x magnifier loop	430-1020
<b>View <a href="#">MyList</a> of GC Columns- 7890, 8890 and 8860</b>	
Agilent J&W DB-EUPAH, 20 m x 0.18 mm, 0.14 µm	121-9627
Agilent J&W DB-5ms 20 m x 0.18, 0.18 µm	121-5522UI
Agilent J&W Select PAH, 30 m x 0.25 mm, 0.15 µm	CP7462



# Consumables Selection Toolkit for PAHs in Environmental Workflow


## 5994-2060EN

Includes Checkout Standards, Column, Liners, Unions, Vials, Caps, filament, gas filters etc. needed to get started with PAH testing with suggested Agilent methods and compatible with PAH Workflow Kits

**Agilent CrossLab**  
From Insight to Outcome

## Polycyclic Aromatic Hydrocarbons Analysis in Environmental Samples

Using single quadrupole GC/MS and triple quadrupole GC/MS/MS:  
Consumables workflow ordering guide



View <a href="#">Mylist</a> of GC column	
Agilent J&W DB-5ms UI 20 m x 0.18, 0.18 µm	121-5522UI
Agilent J&W Select PAH, 30 m x 0.25 mm, 0.15 µm	CP7462
Agilent J&W Select PAH, 15 m x 0.15 mm, 0.10 µm	CP7461
Agilent DB-UI8270D, 30 m x 0.25 mm, 0.25 µm	122-9732
Agilent DB-UI8270D, 20 m x 0.18 mm, 0.36 µm	121-9723
Agilent J&W DB-EUPAH, 20 m x 0.18 mm, 0.14 µm	121-9627
View <a href="#">Mylist</a> of Intuvo GC columns	
Agilent J&W DB-5ms UI 20 m x 0.18, 0.18 µm	121-5522UI-INT
Agilent J&W Select PAH, Intuvo GC column, 30 m x 0.25 mm, 0.15 µm	CP7462-INT
Agilent J&W Select PAH, 15 m x 0.15 mm, 0.10 µm	CP7461-INT
Agilent DB-UI8270D INT, 30 m x 0.25 mm, 0.25 µm	122-9732-INT
Agilent DB-UI8270D INT, 20 m x 0.18 mm, 0.36 µm	121-9723-INT
Agilent J&W DB-EUPAH, Intuvo GC column 20 m x 0.18 mm, 0.14 µm	121-9627-INT
View <a href="#">Mylist</a> of Intuvo Supplies	
Guard Chip, Intuvo Split/Splitless	G4587-60565
Intuvo inlet chip	G4581-60031
Flow Chip, Intuvo, D2-MS	G4581-60033
Flow Chip, Intuvo, swaged HES MS tall	G4590-60109
Inlet/MSD (Intuvo) Polyimide gasket	5190-9072
View <a href="#">Mylist</a> of MS Supplies	
EI filament (for 7000A/B/C/D, 5977B Inert Plus, 5977A extractor, inert or stainless steel and 5975 systems)	G7005-60061
HES Filament for 7010 Triple Quadrupole GC/MS	G7002-60001
Drawout plate, 9 mm, inert source	G3440-20022
Drawout plate, 9 mm, extractor source	G3870-20449
View <a href="#">Mylist</a> of Gas filters	
Gas Clean carrier gas Kit for 7890	CP17988
Gas Clean Carrier Gas Kit for 8890 and 8860	CP179880
Gas Clean carrier gas purifier replacement cartridge	CP17973
Gas Clean Filter Kit for Intuvo	CP17995
View <a href="#">Mylist</a> of Vials and caps	
2 mL screw top amber, write-on spot, deactivated, certified, 100 pc	5183-2072
Screw caps, blue, certified, PTFE/silicone/PTFE septa	5182-0723
100 µL vial insert, glass with polymer feet	5181-8872



# Priority PAHs in Food & Environment

	16 EPA	EU PAH (15+1)	EU PAH8	EU PAH4	AOAC	TEF
Napthalene	x				x	0.001
Acenaphthylene	x					0.001
Acenaphthene	x					0.001
Fluorene	x				x	0.001
Phenanthrene	x				x	0.001
Anthracene	x				x	0.01
Fluoranthene	x				x	0.001
Pyrene	x				x	0.001
Benzo[a]anthracene	x	x	x		x	0.1
Chrysene	x	x	x	x	x	0.01
Benzo[b]fluoranthene	x	x	x	x	x	0.1
Benzo[k]fluoranthene	x	x	x		x	0.1
Benzo[j]fluoranthene		x				0.1
Benzo[a]pyrene	x	x	x	x	x	1
Indeno[123cd]pyrene	x	x	x		x	0.1
Dibenzo[ah]anthracene	x	x	x	x	x	1
Benzo[ghi]perylene	x	x	x		x	0.01
Dibenzo[al]pyrene		x				10
Dibenzo[ai]pyrene		x				10
Dibenzo[ah]pyrene		x				10
Dibenzo[ae]pyrene		x				1
Cyclopenta[cd]pyrene		x				0.1
Methylchrysene		x				
Benzo[c]fluorene		x				


} PAHs in high conc.  
Low toxicity

} PAHs in low conc.  
High toxicity


Gradual shift in focus from occurrence to toxicity & health risk → Changing analytical requirements

TEF = Toxicity Equivalence Factor


# 15+1 EU Priority PAHs & PAH4




Benz[a]anthracene  
MW: 228




Cyclopenta[cd]pyrene  
MW: 226



Chrysene  
MW 228




Indeno[1,2,3-cd]pyrene  
MW 276




Dibenzo[a,h]anthracene  
MW 278


## Isomeric PAHs



Benzo[b]fluoranthene  
MW: 252




Benzo[k]fluoranthene  
MW: 252




Benzo[j]fluoranthene  
MW: 252


## Late Eluters/Heavy PAHs




Dibenzo[a,l]pyrene  
MW 302




Dibenzo[a,e]pyrene  
MW 302




Benzo[c]fluorene  
MW: 216




Benzo[e]pyrene  
MW 252




Benzo[a]pyrene  
MW 252



Benzo[ghi]perylene  
MW 276



Dibenzo[a,i]pyrene  
MW 302



Dibenzo[a,h]pyrene  
MW 302

# Analytical Challenges & Solutions


- || Food related matrices:
  - Fats, oils, proteins → QuEChers + EMR Lipids
  - Inlet and column contaminating → Liner selection & preventive maintenance SOP
  - Increased background interferences → Sample prep & column backflush
- || PAH isobars
  - Many PAHs have identical or quite similar Mw: no differentiation by Mass Spec
  - Chromatographic separation becomes essential**
- || GC column selectivity
  - Some PAHs with similar Mw are challenging to separate chromatographically
  - Dedicated, specially engineered GC liquid phases are required**

# Column Selection for PAHs -- Environmental and Food

## Environmental Analysis

**DB-5ms UI**  
DB-5ms

16 EPA PAHs  
5ms is robust  
option



DB-5ms UI  
in case polars  
are included


## Food Analysis

**DB-EUPAH**  
DB-17ms

16 EPA PAHs  
+  
Benzo b, k, j

EUPAH  
for  
More separation  
detail

20 x 018 x 014  
EUPAH optimized  
for  
Separation/  
Speed/productivity



**Select PAH**

16 EPA PAHs  
+  
Benzo b, k, j

4 EUPAHs in Food  
Benzo b, k, j  
&  
Chrysene/TriP

15m x 0.15mm  
for  
Speed/productivity




# Agilent Specific-application GC column for PAH analysis

## DB-EUPAH

- Mid-polar column that resolves benzo (b,j,k) fluoroanthenes
- Resolution of 24 combined regulated PAHs can be achieved under 28 min

App Notes:  
5990-5872EN  
5990-6155EN



- |                        |                     |                          |                            |
|------------------------|---------------------|--------------------------|----------------------------|
| 1. Naphthalene-d8      | 8. Fluorene         | 15. Benzo(a)anthracene   | 22. Benzo(a)pyrene         |
| 2. Naphthalene         | 9. Phenanthrene-d10 | 16. Chrysene-d12         | 23. Perylene-d12           |
| 3. 2-methylnaphthalene | 10. Phenanthrene    | 17. Chrysene             | 24. Indeno(1,2,3-cd)pyrene |
| 4. 1-methylnaphthalene | 11. Anthracene      | 18. Benzo(b)fluoranthene | 25. Dibenzo(a,h)anthracen  |
| 5. Acenaphthylene      | 12. Fluoranthene    | 19. Benzo(k)fluoranthene | 26. Benzo(g,h,i)perylene   |
| 6. Acenaphthene-d10    | 13. Pyrene          | 20. Benzo(j)fluoranthene |                            |
| 7. Acenaphthene        | 14. p-Terphenyl-d14 | 21. Benzo(e)pyrene       |                            |


# Agilent Specific-application GC column for PAH analysis

## Select PAH


- Mid-polar column that optimize chrysene/triphenylene resolution
- Resolution of 54 PAHs under 40 min

App Notes:  
SI-02232  
SI-02263

## Full resolution of benzo [b] fluoroanthene



## Optimal resolution of triphenylene/chrysene




# What happens if we increase the column aspect ratio?


## Custom Select PAH

- 40m x 0.18mm i.d. x 0.07  $\mu\text{m}$  df
- Full resolution of chrysene/triphenylene and benzo [b,j,k] fluoroanthenes
- Keep in mind limitations of thin-film phases (column lifetime, loading capacity)

### Best resolution of triphenylene/chrysene



### Full resolution of benzo [b,j,k] fluoroanthene



\*Courtesy of Carlos Bueno – Application Scientists, Agilent Spain

# Factors to consider in optimizing EU 15+1 and US-EPA PAH analyses

- Choose Agilent J&W DB-EUPAH when resolution of benzo (b,j,k) fluoroanthene isomers required.
- Select an Agilent J&W Select-PAH when Triphenylene is presented as an interference.

Offers the best resolution value between Triphenylene and Chrysene


- Choose an Agilent J&W DB-5ms Ultra Inert when benzo (b,j,k) fluoroanthene isomers can be reported as a sum of the isomers.

DB-5ms Ultra Inert resolves 23 of 24 regulated PAHs in 27% faster cycle time than DB-EUPAH

- Consider use of retention gaps and inlet backflushing to reduce cycle time and maintenance
- Achieve faster analysis times with no loss of resolution using 0.18mm i.d. GC columns
- Optimize injection volume, temperature, purge time, and solvent focusing for best results
- Minimize inlet and system dwell time with higher linear velocities

# A typical customer setup for PAHs in Food samples

## GC/MS/MS Configuration with Backflush, Jetclean and RT Locking



### Liner

- 4-mm single-tapered liner with glass wool

### Column

- Select PAH
- 30 m x 250 µm x 0.25 µm
- Column 1 Flow at 1.2 mL/min
- Column 2 Flow at 1.5 mL/min

### Oven program

- Initial: 80 °C (0.5 min)
- Ramp at 120 °C/min to 120 °C
- Ramp at 40 °C/min to 180 °C
- Ramp at 3 °C/min to 280 °C
- Ramp at 120 °C/min to 325 (18 min)


### Modifications

- Backflush (BF)
- JetClean
- Retention Time Locking (RTL)


# BackFlush

To maintain column lifetime

How it works? At the end of every run...



Reduce cycle time and eliminates carryover



# Key PAH App Notes incl. Sample Prep

## **PAHs in Salmon**

SI-02424: <https://www.agilent.com/cs/library/applications/SI-02424.pdf>

Sample prep: 7.8 x 300 mm PLgel 100A  
Column: J&W Select PAH, 15m x 0.15mm x 0.10µm

## **PAHs in Peanuts and chocolate with Select PAH column**

5991-2299EN: <https://www.agilent.com/cs/library/applications/5991-2299EN.pdf>

Sample prep: Bond-Elut Si  
Column: J&W Select PAH, 15m x 0.15mm x 0.10µm

## **PAH Analysis in Salmon with Enhanced Matrix Removal**

5991-6088EN: <https://www.agilent.com/cs/library/applications/5991-6088EN.pdf>

Sample prep: Bond-Elut EMR Lipids & Bond-Elut Polish  
Column: DB-5ms UI, 20m x 0.18mm x 0.18µm


## **Determination of 19 Polycyclic Aromatic Hydrocarbon Compounds in Salmon and Beef**

5994-0553EN: <https://gcms.cz/labrulez-bucket-strapi-h3hsga3/paper/application-pahs-salmon-beef-5994-0553en-agilent.pdf>

Sample prep: 2 step SoLE using 20:80 EtOAc/ACN and EMR Lipids  
Column: DB-EUPAH, 30m x 0.25mm 0.25µm


# Environmental applications: Polychlorinated Biphenyls (PCBs)

Column Updates & Solutions





# Polychlorinated biphenyls



- A polychlorinated biphenyl (PCB) is an organic chlorine compound with the formula  $C_{12}H_{10-x}Cl_x$ . They were widely used as dielectric and coolant fluids in electrical apparatus, carbonless copy paper and in heat transfer fluids
- Around 209 PCB's
- List of reported PCB congeners are determined by different organizations. Examples of PCB congener lists include:
  - 7 Indicator Congeners (marker PCB's) in Europe: 28, 52, 101, 118, 138, 153, 180, indicate presence/absence of PCB's. Related to PCB occurrence in Aroclors 1242, 1248
  - 12 Most Toxic Congeners (WHO), 77, 81, 126, 169, 105, 114, 118, 123, 156, 157, 167, 189
  - 18 NOAA Status & Trend Congeners
  - 36 CDC PCB congener list (found in biological sample)
  - 159 Congeners Quantified in Aroclor mixes

# Polychlorinated Biphenyls (PCBs) analysis by GC

## Aroclors

- 'Old' way to measure PCBs. Chromatographic 'fingerprint' indicative of one of Aroclors mixture. No concentration/identification of specific PCB
- 9 common PCB Aroclor mixtures (1221, 1232, 1242, 1016, 1248, 1254, 1260, 1262, 1268), each with distinctive GC pattern.
- **Method:** GC/ECD, SW-846, EPA 8082A, GC column: **DB-5ms, 30m x 0.25mm x 1 um**, DB- 608, 30m x 0.25mm x 1 um

## Homologs


- Grouping PCB congeners by number of chlorine atoms, this can vary from one to ten. All PCBs with same number of chlorine atoms belong to same homolog group. Testing for PCB homologs provides more reliable results than the Aroclor method
- Laboratory results for PCB homologs will list the amount of PCB present in the sample by the number of chlorine atoms.
- **Method:** GC/MS, EPA Methods 8270D, 625, and 8081 A/B, GC Column: **DB-UI 8270D 30m x 0.25mm x 0.25 um**, DB-5ms UI, 30m x 0.25mm x 0.25 um

## Congeners

- More congener detail can be needed, generally a subset of full 209 PCB congeners.
- NOAA PCB congener method cites 20 congeners , SW-846 8082 method cites 19 congeners. WHO list cites 12 congeners.
- **Method:** GC/MS EPA Method 1668C, GC column: **DB-XLB 30m x 0.18mm x 0.18 um meets method requirements.** Column must resolve congeners 34/23, and 187/182. Congener 156/157 must co-elute within 2 seconds at peak maximum

# Analysis of Aroclors 1016-1268 (without 1221) with DB-XLB

GC System: Agilent 7890  
Column: **DB-XLB**, 30m x 0.18 mm i.d. x 0.18  $\mu$ m df (121-1232)  
Carrier: Helium at 37 cm/sec, measured at 150 °C  
Oven: 100 °C (1 min), 1.2 °C/min to 265 °C  
Injector: Hot on-column, 250 °C, 1  $\mu$ L in isoctane, 12.5 ppm  
Detector: MSD, SIM, 340 °C transfer line



IUPAC	Retention time
PCB 23	55.54 min
PCB 34	55.33 min
PCB 156	105.3 min
PCB 157	105.72 min


Other GC columns:  
DB-5ms UI, DB-8270D UI  
DB-35ms UI  
CP-Sil 8 for PCB  
CP-Sil 8/C18 CB for PCB

Source: <https://www.agilent.com/cs/library/chromatograms/361.pdf>

# Analytical Challenges & Solutions for PCBs

- For PCB's, -5 type columns can be used to quantify PCB homologs, **but polarity is too low to resolve key PCB congeners**.

## Key PCB separations using GC columns with different selectivity




Select PAH can partially separate 28/32 and 138/163 PCB pairs, making it a good option for multi-residue analysis

# How about DB-XLB?

- As mentioned earlier, DB-XLB offers the best resolution for PCB congeners
- $R_s = 1.0$  for PCB 31/28, and baseline resolution for other key PCBs

## EU Monitor PCBs on DB-XLB



**Column:** DB-XLB  
30 m x 0.25 mm I.D., 0.5  $\mu$ m  
1 m x 0.53 mm I.D. Retention Gap

**J&W P/N:** 122-1236

**Carrier:** Helium at 34.2 cm/sec, measured at 150°C


**Oven:** 100°C for 1 min  
100-320°C at 5.6°/min

**Injector:** Hot On-column, 250°C, Split Flow  
100 mL/min

**Detector:** 2  $\mu$ L dilute Aroclor mixture  
MSD, SIM of 221.9, 255.9, 291.9,  
325.8, 359.8, 395.8, 429.7, 463.7  
300°C transfer line


# Agilent offers the only one-column solution for OCPs, PAHs and PCBs congeners: Select PAH

## GC/MS analysis of 16 PAHs, 17 PCBs, 24 OCPs with Select PAH GC Column



App Note: SI-02438


# Environmental applications: Dioxins and Furans by GC/MS/MS and a DB-5ms UI



# PCDD and PCDF Nomenclature and Isomers

## Dioxins

Polychlorinated dibenzodioxin (PCDD)




2,3,7,8 – TCDD

(2,3,7,8-Tetrachlorinated dibenzodioxin)

## Furans

Polychlorinated dibenzofurans (PCDF)



2,3,4,7,8 – PCDF

(2,3,4,7,8-Pentachlorinated dibenzofurans)

Chlorine Atoms	PCDD Isomers	PCDF Isomers
1	2	4
2	10	16
3	14	28
4	22	38
5	14	28
6	10	16
7	2	4
8	1	1
<b>Total</b>	<b>75</b>	<b>135</b>

75 total PCDD isomers

135 total PCDF isomers

**17 toxic isomer**




## EPA has accepted GC/TQ as an equivalent Technology to analyze dioxins in Environmental Samples

**This method meets all requirements for measurements of 2,3,7,8 substituted tetra-Octa CDDs and CDFs in wastewater... The performance of this method is substantially similar to methods listed in 40 CFR Part 136 for measurement of PCDDs/PCDFs in wastewater... EPA will recommend this method be included in list of approved methods at 40 CFR Part 136**

-EPA Office of Water

# System Configuration Used in the Application 5994-3029EN (16130 Method)


Most labs use CTC PAL for sample introduction (1 mL vial)




# Dioxins/Furans – Chromatogram

Excellent separation of the difficult  
*hexa-dioxin/furan isomers*

GC Column – DB5MS UI (60m x 0.25mm x 0.25 μm)



# TCDD Isomer Specificity



Method 1613B calls for calculation of the percent valley between the GC peaks that elute most closely to the 2,3,7,8-TCDD and TCDF isomers.

The height of the valley between the isomers most closely eluting to the 2,3,7,8-TCDD labeled “x” **does not exceed 25% of the 2,3,7,8-TCDD peak height “y”**.

This parameter can be set as an outlier in the MassHunter Quantitative Analysis method. If the valley exceeds 25%, the analytical conditions need to be adjusted or the analysis repeated using a different GC column.

Method Table									
Time Segment: <All>				Compound: 2,3,7,8-TCDD		Reset Table View			
Quantifier									
Name	TS	Transition	Scan	Type	Resolution Calculation Type	Resolution Limit			
2,3,7,8-TCDD	1	319.9 -> 256.9...	MRM	Target	Valley Height/Peak Heig...	25.0			

Sample							2,3,7,8-TCDD Results							
Name	Data File	Type	Level	Vial	Acq. Date-Time	Acq. Method File	RT	Resp.	MI	Calc. Conc.	Accuracy	S/N	Resolution F.	Resolution R.
DX041D-CAL_01-73	DX9Z0444.D	Cal	CS3	7	8/22/2019 2:31 AM	TQEI_DB5_DX_11	26.351	221151		9.0183	90.2	2339.21	20.4	7.8

# Method detection Limits (MDLs) with 7890/7010 GC/TQ Far Surpassed Method 1613B MDLs

The GC/TQ MDL results for the aqueous (1 L), solid (10 g), and tissue (10 g) samples:

Compound	Aqueous	Solid	Tissue
	MDL and (MRL) in pg/L	MDL and (MRL) in pg/g	MDL and (MRL) in pg/g
2,3,7,8-TCDD	1.1 (10) *	0.029 (1)	0.057 (0.5)
1,2,3,7,8-PeCDD	1.39 (50)	0.037 (5)	0.051 (2.5)
1,2,3,4,7,8-HxCDD	1.05 (50)	0.042 (5)	0.061 (2.5)
1,2,3,6,7,8-HxCDD	1.08 (50)	0.045 (5)	0.033 (2.5)
1,2,3,7,8,9-HxCDD	1.78 (50)	0.064 (5)	0.067 (2.5)
1,2,3,4,6,7,8-HpCDD	1.19 (50)	0.070 (5)	0.032 (2.5)
OCDD	9.4 (100)	0.311 (10)	0.085 (5)
2,3,7,8-TCDF	0.56 (10)	0.60 (1)	0.056 (0.5)
1,2,3,7,8-PeCDF	1.0 (50)	0.037 (5)	0.046 (2.5)
2,3,4,7,8-PeCDF	1.25 (50)	0.039 (5)	0.033 (2.5)
1,2,3,4,7,8-HxCDF	0.89 (50)	0.032 (5)	0.029 (2.5)
1,2,3,6,7,8-HxCDF	1.11 (50)	0.031 (5)	0.046 (2.5)
1,2,3,7,8,9-HxCDF	1.22 (50)	0.048 (5)	0.084 (2.5)
2,3,4,6,7,8-HxCDF	1.26 (50)	0.026 (5)	0.034 (2.5)
1,2,3,4,6,7,8-HpCDF	0.92 (50)	0.255 (5)	0.064 (2.5)
1,2,3,4,7,8,9-HpCDF	1.35 (50)	0.028 (5)	0.043 (2.5)
OCDF	2.81 (100)	0.365 (10)	0.113 (5)

**MDL ~ 1.1 ppt for 2,3,7,8-TCDD**

\*The current MDL for 2,3,7,8-TCDD with the GC/HRMS 1613B is 4.4 pg/L. The MDL achieved with GC/TQ is 4 times better.

# Dioxin Analysis Method for Food & Feed



## Validation of a Confirmatory GC/MS/MS Method for Dioxins and Dioxin-like PCBS to Meet the Requirements of EU Regulation 709/2014

### Application Note

Food

#### Author

Joerg Riener  
Agilent Technologies GmbH  
Waldbronn  
Germany

#### Abstract

Using the Agilent 7890B GC and the Agilent 7000C Series Triple Quadrupole GC/MS System, a GC/MS/MS method has been developed and fully validated to meet the requirements of EU Regulation 709/2014 for the monitoring of polychlorinated dibenzo-*p*-dioxins (PCDDs), polychlorinated dibenzofurans (PCDFs), and polychlorinated biphenyls (PCBs) in animal feedstuffs. It provides similar performance to GC/HRMS, the analytical platform required by previous EU regulations, in spite of the difference in mass analyzer technologies.



## Determination of Polychlorinated Dibenzo-*p*-dioxins (PCDD) and Polychlorinated Dibenzofurans (PCDF) in Foodstuffs and Animal Feed using the Agilent 7000 Triple Quadrupole GC/MS System

### Application Note

Food Safety

#### Authors

Prof. Dr. Peter Fürst  
Dr. Thorsten Bernsmann  
Dominik Baumeister  
Chemical and Veterinary Analytical  
Institute  
Münsterland-Emscher-Lippe  
(Chemisches und Veterinär-  
untersuchungsamt, CVUA-MEL)  
Joseph-König-Str. 40  
48147 Münster  
Germany

#### Abstract

A method has been developed on the Agilent 7000 GC Triple Quadrupole GC/MS system for the analysis of polychlorinated dibenzo-*p*-dioxins (PCDD) and polychlorinated dibenzofurans (PCDF) in foodstuffs and animal feed. The method was shown to give linear response over the required concentration range, good repeatability of response and quantitation down to low pg TEQ/g levels.



# Consumables Selection Toolkit for Dioxins Workflow

## 5994-3649EN

Includes Checkout Standards, Column, Liners, Unions, Vials, Caps, filament, gas filters etc. needed to get started with Dioxin testing with suggested Agilent methods and compatible with Dioxin Workflow Kit




From Insight to Outcome

## Analysis of Dioxins in Food and Feed using GC/MS

Consumable Workflow Ordering Guide





PCDDs



MyList of Chemical Standards	
Dioxin Analyzer Checkout Standard Kit	G3440-85039
Chlorinated Dibenzop-dioxin Mixture, 10 µg/mL, Toluene, 1 mL	RPE-065M-1
Chlorinated Dibenzofuran Mixture, 10 µg/mL, Toluene, 1 mL	RPE-045M-1
2,3,7,8-Tetrachlorodibenzo-p-dioxin, 50 µg/mL, Toluene, 1 mL	RPE-029S-1
Octachlorodibenzo-p-dioxin, 50 µg/mL, Toluene, 1 mL	RPE-017S-1
2,3,7,8-Tetrachlorodibenzofuran, 50 µg/mL, Toluene, 1 mL	RPE-037S-1
1,2,3,7,8-Pentachlorodibenzofuran, 50 µg/mL, Toluene, 1 mL	RPE-042S-1
1,2,3,4,7,8-Hexachlorodibenzofuran, 50 µg/mL, Toluene, 1 mL	RPE-043S-1
Octachlorodibenzofuran, 50 µg/mL, Toluene, 1 mL	RPE-019S-1
MyList of GC Column and retention gap	
DB-5ms Ultra Inert, 60 m x 0.25 mm, 0.25 µm	I22-5562UI
Ultimate Plus Deactivated fused silica tubing, 5 m, 0.25 mm <sup>o</sup>	CP802505
MyList of GC Supplies	
Purged Ultimate Union**	G3186B
Liner, Ultra Inert, splitless, single taper, glass wool (4 mm id)	5190-2293
Inlet Septa Non-Stick Adv Green, 11 mm, 50/pk	5183-4759
Ultra Inert Gold seal, with washer, 1/pk	5190-6144
Ultra Inert Gold seal, with washer, 10/pk	5190-6145
Self-Tightening column nut, collared, inlet	G3440-81011
Self-Tightening column nut, collared, MSD	G3440-81013
ALS syringe, Blue Line, 10 µL, fixed needle, 23/42/cone, PTFE-tip plunger	G4513-80220
15%Graphite/85% Vespel Ferrules, 0.4 mm id, 10/pk	5181-3323
QuickPick preventative maintenance kit, for split vent trap, includes 1 cartridge**	5181-6495
Ceramic scribe column cutter, 4/pk	5181-8836
Column installation tool	G1099-20030
MyList of Vials and Caps	
Vial, screw top, amber, silanized, write-on spot, certified, 2 mL, 100/pk	5183-2072
ALS Vial Cap, screw, green, PTFE/red silicone septa, 100/pk	5182-0718
Vial insert, 250 µL, glass conical with polymer feet, deactivated, 100/pk Insert size: 5.6 x 30 mm	5181-8872
MyList of MS supplies	
EI Filament (for 7000A/B/C/D, 5977B Inert Plus, 5977A Extractor, Inert or Stainless steel and 5975 systems)	G7005-60061

This guide provides recommendations for Agilent products used in this analysis, so you can find what you're looking for quickly. Click the MyList\* links in the header to add items to your "Favorite Products" list at the Agilent online store. Then, enter the quantities for the products you need. Your list will remain under "Favorite Products" for your use with future orders.



# Conclusions

- **Agilent have the most complete portfolio of GC columns for general purpose and specific applications**
- Agilent is committed to continue innovation that provides better lab productivity and economic value to our customers
- Agilent continues to develop new technical notes to solve the toughest challenges in gas chromatography

**We are here to help!**



# Thank you!