



Agilent 8890 Features to Improve Lab Productivity

Jamie Wetuski
Product Support
22 November 2022



Agilent
CrossLab
From Insight to Outcome






Communicating With Your GC



Touch screen



Agilent software



Browser



Touch Screen Home

The screenshot displays the Agilent CrossLab Touch Screen Home interface. At the top, there is a navigation bar with icons for Home, Method, Sequence, Diagnostics, Maintenance, Logs, Settings, and Help. Below this, the main area shows two method configurations. Each configuration consists of three components: a pump (S/SL) at 250.00 °C, a column (Col 1 or Col 3) at 70.00 °C, and a detector (FID) at 300.00 °C. The flow rates are 5.876 mL/min. The status bar at the bottom indicates 'STATUS: WAITING FOR PREP RUN' and 'Est. Remaining 2.00'. There are also icons for a sequence, a play button, and a stop button.

Component	Value
Pump (S/SL)	23.029 psi, 250.00 °C
Column (Col 1)	5.876 mL/min, 70.00 °C
Detector (FID)	2.3 pA, 300.00 °C
Pump (S/SL)	23.029 psi, 250.01 °C
Column (Col 3)	5.876 mL/min, 70.00 °C
Detector (FID)	9.6 pA, 300.00 °C

STATUS: WAITING FOR PREP RUN

Sequence Method Sample Est. Remaining **2.00**



Touch Screen Method

The screenshot displays the Agilent CrossLab touch screen interface. At the top, there is a navigation bar with icons for Home, Method, Sequence, Diagnostics, Maintenance, Logs, Settings, and Help. Below this, the main area shows two parallel method paths. Each path starts with a pump icon (S/SL) at 23.030 psi and 250.00 °C. The first path goes through a column (Col 1) at 5.876 mL/min and 70.01 °C, and ends with a detector (FID) at 2.4 pA and 300.00 °C. The second path goes through a column (Col 3) at 5.876 mL/min and 70.01 °C, and ends with a detector (FID) at 9.7 pA and 300.00 °C. A green status bar at the bottom indicates 'STATUS: WAITING FOR PREP RUN'. Below the status bar, there is a table with columns for Sequence, Method, Sample, and Est. Remaining (2.00). At the bottom right, there are icons for a pump, a play button, and a stop button.

Sequence	Method	Sample	Est. Remaining
			2.00

Touch Screen Method, Continued

The screenshot displays the 'Method' configuration screen in the Agilent CrossLab software. The top navigation bar includes 'Method', 'Sequence', 'Diagnostics', 'Maintenance', 'Logs', a settings gear icon, and a help icon. A left-hand menu lists various system components: ALS/Tray, Valves, Inlets, Columns (highlighted), Oven, Detectors, Analog Out, and Events. The main area is titled 'Column 1' and shows the 'Control Mode' set to 'Constant Flow'. Below this, a table displays 'Setpoint' and 'Actual' values for 'Flow' and 'Pressure', along with 'Avg. Velocity' and 'Holdup Time'.

	Setpoint	Actual
<input checked="" type="checkbox"/> Flow	5.876 mL/min	5.876 mL/min
<input checked="" type="checkbox"/> Pressure	21.461 psi	21.461 psi
Avg. Velocity	72.234 cm/sec	
Holdup Time	0.692 min	



Touch Screen Sequence

The screenshot displays the Agilent CrossLab touch screen interface. At the top, there is a navigation bar with icons for Home, Method, Sequence, Diagnostics, Maintenance, Logs (with a green circle containing the number 1), Settings, and Help. Below the navigation bar is a sidebar menu with the following items: ALS/Tray, Valves, Inlets, Columns, Oven, Detectors, Analog Out, and Events (which is highlighted in blue). The main content area is titled "Events" and contains a "Run Events Table". The table has four columns: Time, Event, Position, and Setpoint. Below the table are two buttons: "Add Event" and "Remove Event".

Time	Event	Position	Setpoint
------	-------	----------	----------

Touch Screen Diagnostics

The screenshot displays the Agilent CrossLab touch screen interface. At the top, a navigation bar includes icons for Home, Method, Sequence, Diagnostics (with a '2' notification), Maintenance, Logs (with a '1' notification), Settings, and Help. Below the navigation bar, the 'Sequence Name' is set to 'Front_BackFID'. To the right of the dropdown are four control buttons: 'Run' (green), 'Pause' (red), 'Resume' (red), and 'Abort' (red). Below these controls is a table with the following data:

Line	Sample Name	Vial Number	No. of Injections	Method
1	Front FID	1	1	Front_SSL_FID
2	Back FID	1	1	Back_SSL_FID



Gas Identification Test

The screenshot displays a software interface for diagnostic tests. At the top left, there is a back arrow and the text "< Diagnostic Tests". The main title is "Inlets Diagnostic Tests" in bold, with a help icon (a question mark in a circle) to its right. Below the title, there are two tabs: "Front Inlet - S/SL" (which is selected and underlined) and "Back Inlet - S/SL". Under the "Front Inlet - S/SL" tab, a list of diagnostic tests is shown:

- Gas Identification Test
- Gas Supply Pressure Check
- Leak and Restriction Test
- Pressure Decay Test
- Septum Purge Test
- Split Vent Restriction Test



Gas Supply Pressure Check

< Diagnostic Tests **Inlets Diagnostic Tests** ⓘ

Front Inlet - S/SL	Back Inlet - S/SL
Gas Identification Test	
Gas Supply Pressure Check	
Leak and Restriction Test	
Pressure Decay Test	
Septum Purge Test	
Split Vent Restriction Test	



Leak and Restriction Test

< Diagnostic Tests **Inlets Diagnostic Tests** ?

Front Inlet - S/SL	Back Inlet - S/SL
Gas Identification Test	
Gas Supply Pressure Check	
Leak and Restriction Test	
Pressure Decay Test	
Septum Purge Test	
Split Vent Restriction Test	




Pressure Decay Test

< Diagnostic Tests **Inlets Diagnostic Tests** ?

Front Inlet - S/SL	Back Inlet - S/SL
Gas Identification Test	
Gas Supply Pressure Check	
Leak and Restriction Test	
Pressure Decay Test	
Septum Purge Test	
Split Vent Restriction Test	



Septum Purge Test

< Diagnostic Tests **Inlets Diagnostic Tests** 

Front Inlet - S/SL	Back Inlet - S/SL
Gas Identification Test	
Gas Supply Pressure Check	
Leak and Restriction Test	
Pressure Decay Test	
Septum Purge Test	
Split Vent Restriction Test	

Split Vent Restriction Test


The screenshot shows a software interface for diagnostic tests. At the top left, there is a back arrow and the text "< Diagnostic Tests". The main title is "Inlets Diagnostic Tests" with a help icon (question mark in a circle) to its right. Below the title, there are two tabs: "Front Inlet - S/SL" and "Back Inlet - S/SL". The "Back Inlet - S/SL" tab is currently selected. Underneath the tabs, a list of diagnostic tests is displayed:

- Gas Identification Test
- Gas Supply Pressure Check
- Leak and Restriction Test
- Pressure Decay Test
- Septum Purge Test
- Split Vent Restriction Test

Split Vent Restriction Test

Back Inlet : Split Vent Restriction Test

Close Cancel ?



State	In Progress
Time Remaining	0 seconds
Inlet Pressure Target	12.661 psi
Inlet Pressure Actual	12.664 psi

Verify Basic Control

50% Complete



Touch Screen Diagnostics

The screenshot displays a touch screen interface for diagnostic tests. At the top left, there is a back arrow and the text "< Diagnostic Tests". The main title is "Inlets Diagnostic Tests" with a help icon (question mark in a circle) to its right. Below the title, there are two tabs: "Front Inlet - S/SL" and "Back Inlet - S/SL". The "Back Inlet - S/SL" tab is currently selected. A list of diagnostic tests is shown below the tabs:

- Gas Identification Test
- Gas Supply Pressure Check
- Leak and Restriction Test
- Pressure Decay Test
- Septum Purge Test
- Split Vent Restriction Test




Troubleshooting

< Diagnostics		Troubleshooting Report		Start Troubleshooting	?
Issue Date	Issue Name	Issue resolution accepted	Diagnosis		
Wed 19 Oct 2022 09:47:05 AM	No Peaks	Cancelled	-----		
Wed 19 Oct 2022 09:44:44 AM		Cancelled	-----	^	
Wed 19 Oct 2022 09:10:55 AM		Cancelled	-----		
Wed 19 Oct 2022 09:09:34 AM		Cancelled	-----		
Wed 19 Oct 2022 09:06:42 AM		Cancelled	-----		
Wed 19 Oct 2022 09:06:18 AM		Cancelled	-----		v
Mon 29 Aug 2022 03:05:49 PM		Cancelled	-----		

Diagnostics Gas and Power Usage

The screenshot displays the Agilent CrossLab software interface. At the top, a navigation bar includes icons for Home, Method, Sequence, Diagnostics (highlighted in blue), Maintenance, Logs, Settings, and Help. Below this, a grid of diagnostic options is shown: Warnings and Errors, Diagnostic Tests, System Health Report, Detector Evaluation, Blank Evaluation, Peak Evaluation, Troubleshooting, and Gas and Power Usage. A green status bar at the bottom indicates 'STATUS: WAITING FOR PREP RUN' with an upward arrow and a user icon. Below the status bar, a control bar shows 'Sequence', 'Method', 'Sample', 'Est. Remaining 6.50', and icons for a gas flow controller, play, and stop.

Touch Screen Maintenance



Touch Screen Inlet Maintenance

The screenshot displays the 'Inlets Maintenance' screen. At the top, there is a navigation bar with '< Overview', 'Inlets Maintenance', 'Perform Maintenance', and a help icon. Below this, there are two tabs: 'Front Inlet - SS' and 'Back Inlet - SS'. The 'Front Inlet - SS' tab is active, showing a table with columns for 'Part', 'Status', and 'Reset All'. The table lists four items: 'Gold seal age' (2 mth 4 days), 'Gold seal injections' (6 injections), 'Liner age' (11 min 26 sec), and 'Liner injections' (0 injections). A yellow status bar at the bottom indicates 'STATUS: NOT READY'. Below the status bar, there is a row of controls including 'Sequence', 'Method', 'Sample', 'Est. Remaining 6.50', and three icons: a flask with 'P', a play button, and a square stop button.

Part	Status	Reset All
✓ Gold seal age	2 mth 4 days	
✓ Gold seal injections	6 injections	^
✓ Liner age	11 min 26 sec	
✓ Liner injections	0 injections	v

STATUS: NOT READY

Sequence Method Sample Est. Remaining **6.50** [Flask P] [Play] [Stop]

Maintenance Counter

The screenshot displays the 'Inlets Maintenance' screen in the Agilent CrossLab software. At the top, there is a navigation bar with '< Overview', 'Inlets Maintenance', 'Perform Maintenance', and a help icon. Below this, two tabs are visible: 'Front Inlet - SS' (selected) and 'Back Inlet - SS'. The main content area is a table with columns for 'Part', 'Status', and 'Reset All'. The table lists four maintenance items, all with green checkmarks indicating they are completed or within limits. A yellow status bar at the bottom of the table area reads 'STATUS: NOT READY'. Below the status bar, there is a row of controls including 'Sequence', 'Method', 'Sample', 'Est. Remaining 6.50', and three icons: a flask with 'P', a play button, and a square stop button.

Part	Status	Reset All
✓ Gold seal age	2 mth 4 days	
✓ Gold seal injections	6 injections	^
✓ Liner age	28 min 56 sec	
✓ Liner injections	0 injections	v

STATUS: NOT READY

Sequence Method Sample Est. Remaining **6.50**


Maintenance Counter Column

< Overview **Columns Maintenance** Perform Maintenance ?

Column 1		Column 2
Part	Status	
✓ Injections	23 injections	^ v
✓ Length	30 m	
✓ Max applied temperature	299 °C	
✓ On time	2 mth 3 days	

STATUS: WAITING FOR PREP RUN ^

Sequence Method Sample Est. Remaining **6.50** [P] [▶] [■]



Maintenance Column Key



Maintenance Column Key

The screenshot displays the Agilent CrossLab software interface. At the top, a navigation bar includes icons for Home, Method, Diagnostics (with a green circle containing the number 1), Maintenance, Logs, Settings, and a Help icon (?). Below the navigation bar are three buttons: "View System Health Report", "Start Troubleshooting", and "Diagnostic Tests". A prominent alert banner features a yellow warning triangle icon, the text "Instrument: New Column Smart Key", and a timestamp "2021-11-04 16:26:42.". The main area of the interface is mostly empty, with a red circle highlighting a small black dot in the center. At the bottom, a yellow status bar displays "STATUS: NOT READY" on the left, an upward-pointing arrow in the center, and a user profile icon on the right. Below the status bar is a control panel with columns for "Sequence", "Method", "Sample", and "Est. Remaining" (showing "0.53"). To the right of the control panel are three circular icons: a blue circle with a white 'P' (Pause), a blue circle with a white play button (Run), and a blue circle with a white square (Stop).

Maintenance Detector

The screenshot displays the 'Maintenance' section of the Agilent CrossLab interface. The navigation bar includes 'Method', 'Sequence', 'Diagnostics', 'Maintenance' (active), 'Logs' (with a '1' notification), a settings gear, and a help icon. Two buttons, 'View Logs' and 'Install Column', are located in the top right of the main area. The central part of the screen features six white cards, each representing a component: ALS, Valves, Inlets, Columns, Detectors, and Instrument. Each card has a green checkmark in the top right corner, indicating that all components are in good status. At the bottom, a green status bar shows 'STATUS: WAITING FOR PREP RUN' with an upward arrow, and a control bar with 'Sequence', 'Method', 'Sample', 'Est. Remaining 6.50', and icons for a pipette, play, and stop.

Detector Diagnostic

The screenshot displays the 'Maintenance' section of the Agilent CrossLab interface. The top navigation bar includes 'Method', 'Sequence', 'Diagnostics', 'Maintenance' (selected), 'Logs 1', a settings gear, and a help icon. Below the navigation bar are 'View Logs' and 'Install Column' buttons. The main area contains six diagnostic cards: ALS, Valves, Inlets, Columns, Detectors, and Instrument, each with a green checkmark indicating a successful status. A green status bar at the bottom shows 'STATUS: WAITING FOR PREP RUN' and 'Est. Remaining 6.50'. The bottom-most bar contains labels for 'Sequence', 'Method', 'Sample', and 'Est. Remaining 6.50', along with icons for a pipette, play, and stop.


Touch Screen Logs

The screenshot displays the 'Detectors Maintenance' screen. At the top, there is a navigation bar with '< Overview', 'Detectors Maintenance', 'Perform Maintenance', and a help icon. Below this, there are three tabs: 'Front Detector - FID' (selected), 'Back Detector - FID', and 'Aux Detector 2 - TCD'. The main content area is a table with the following data:


Part	Status	Reset All
✓ Collector assembly injections	5 injections	
✓ Ignitor ignitions	21 ignitions	
✓ Jet injections	5 injections	

At the bottom of the screen, there is a yellow status bar that reads 'STATUS: NOT READY' with an upward arrow icon. Below this is a control bar with labels 'Sequence', 'Method', 'Sample', and 'Est. Remaining 6.50'. To the right of the 'Est. Remaining' label are three icons: a 'P' in a circle, a play button, and a square stop button.


Touch Screen Settings



Touch Screen Settings, Continued




Touch Screen Help



Browser Interface

Chemical Analysis, Life Sciences, x Agilent 8890 x +

← → ↻ 🏠 ⚠ Not secure | 130.30.250.249/#login 🔗 ☆ 🛠 🗄 Ⓜ Paused ⋮


Agilent
8890
Gas Chromatograph

Enter PIN

Enter

Browser Interface, Continued

The screenshot displays the Agilent CrossLab browser interface. The browser window title is "Agilent 8890" and the address bar shows "130.30.250.249/#method". The interface has a dark blue header with navigation tabs: Home, Method, Sequence (selected), Diagnostics, Maintenance, Logs, Settings, and Help. Below the header is a toolbar with icons for back, forward, refresh, home, and various file management actions. The main content area is divided into two sections: "Injection" and "Washes and Pumps".

Injection Section:

- Syringe:** Syringe Size (µL) is 10; Injection Volume (µL)* is 1.
- Dwell Time:** Pre-Injection (min)* is 0; Post-Injection (min)* is 0.
- Multiple Injections:** Number of Injections* is 1; Injection Delay (sec)* is 0.
- Sample Depth:** Enable checkbox is unchecked; Setpoint (mm)* is 0.

Washes and Pumps Section:

- Solvent A:** Pre-Injection* is 0; Post-Injection* is 0; Volume (µL) is empty.
- Solvent B:** Pre-Injection* is 0; Post-Injection* is 0; Volume (µL) is empty.
- Sample:** Washes* is 0; Volume (µL) is Max (8); Pumps* is empty.

A left-hand sidebar contains a tree view of system components: ALS/Tray (Front, Back, Other Injectors; Valves), Inlets (Front, Back Inlet - S/SL), Columns, Oven, and Detectors (Front, Back, Aux Detectors; Analog Out).

At the bottom, a green status bar shows "STATUS: WAITING FOR PREP RUN" and a user profile icon. Below this is a table with the following data:

Sequence	Method	Sample Name	Est. Remaining
			6.5

Browser Interface, Continued

The screenshot shows a web browser window with the following details:

- Browser tabs: "Chemical Analysis, Life Sciences, ..." and "Agilent 8890".
- Address bar: "Not secure | 130.30.250.249/#settings/main/0".
- Navigation icons: Back, Forward, Refresh, Home, and a "Paused" button.
- Menu bar: Home, Method, Sequences, Diagnostics, Maintenance, Logs, Settings (highlighted), and Help.
- Main content area: A grid of six dark blue buttons labeled "Service Mode", "Calibration", "About", "System Settings", "Tools", and "Scheduler".
- Footer: A green bar with "STATUS: WAITING FOR PREP RUN", an up arrow, and a user icon.
- Bottom status bar: A light green bar with a table header.

Sequence	Method	Sample Name	Est. Remaining
			6.5

Firmware

Firmware on
www.agilent.com

Download the
latest firmware
update tool today!

The screenshot shows the Agilent website homepage. At the top, there is a navigation bar with the Agilent logo and the tagline "Trusted Answers". The navigation menu includes: PRODUCTS, APPLICATIONS & INDUSTRIES, TRAINING & EVENTS, SERVICES, SUPPORT & RESOURCES, and ORDER CENTER. There is also a search bar and a shopping cart icon. The main banner features a large image of a molecular structure with the text "Cancer Genomics: The long path to truth" and a "READ MORE >" link. Below the banner is a horizontal menu with categories: Inoculating Labs, IVDR, The Visionary Lab, Cancer Genomics (highlighted), Vacuum Tech, Net Zero Commitment, and Diversity Matters. The "HIGHLIGHTS" section contains three articles: "Introducing Cary 3500 UV-Vis with Agilent OpenLab", "High-Resolution Ion Mobility MS", and "Customer insights from AGQ Labs". The "NEWS" section lists three recent news items: "Agilent Announces the Innovative Seahorse XF Pro Analyzer", "Agilent to Webcast First-Quarter Fiscal Year 2022 Financial Results Presentation", and "Agilent Expands Manufacturing Capacity in Shanghai".



Agilent CrossLab

From Insight to Outcome



Contact Agilent Chemistries and Supplies Technical Support



1-800-227-9770 Option 3, Option 3:

Option 1 for GC and GC/MS columns and supplies

Option 2 for LC and LC/MS columns and supplies

Option 3 for sample preparation, filtration, and QuEChERS

Option 4 for spectroscopy supplies

Option 5 for chemical standards

Available in the USA and Canada 8–5, all time zones



gc-column-support@agilent.com

lc-column-support@agilent.com

spp-support@agilent.com

spectro-supplies-support@agilent.com

chem-standards-support@agilent.com