


# Extraction of Anabolic Steroids from Horse Urine Using ISOLUTE® SLE+ Prior to LC-MS/MS Analysis



**Figure 1.** Structures of Ethylestrenol and Stanozolol

## Introduction


Ethylestrenol and stanozolol are anabolic steroids which can be used to increase muscle mass and enhance performance. These drugs have been linked to instances of doping in race horses. This application note describes a Supported Liquid Extraction (SLE) protocol for the extraction of ethylestrenol and stanozolol from horse urine prior to LC-MS/MS analysis.

The method described in this application note achieves high reproducible analyte recoveries from both gelding and filly urine. Protocols for 48-well, 96-well plate and column formats are described.

ISOLUTE® SLE+ Supported Liquid Extraction plates and columns offer an efficient alternative to traditional liquid-liquid extraction (LLE) for bioanalytical sample preparation, providing high analyte recoveries, no emulsion formation, and significantly reduced sample preparation.

## Analytes

### 4-Ethylestrenol and Stanozolol.



**Figure 1.** Supported liquid extraction mechanism.

## Sample Preparation Procedure

**Format: ISOLUTE® SLE+ 1 mL Supported Liquid Extraction plate, 48-well, part number 820-1000-Q01.** Protocols for 96-well plates and column formats are also included, see **Table 1**.

### Sample Pre-treatment

Take appropriate volume of urine and add same volume of H<sub>2</sub>O. Mix.

### Sample Loading

ISOLUTE SLE+ 48-well plate: Load the pre-treated sample (800 µL) to each well of the 48-well plate followed by a pulse of vacuum or positive pressure to initiate flow. Leave to absorb for 5 minutes.

### Elution

Ensure a suitable collection vessel is in place.

Apply 1 mL of MTBE and allow to flow under gravity.

Apply a second 1 mL of MTBE and allow to flow under gravity.

Apply a third 1 mL of MTBE and allow to flow under gravity until the solvent reaches the top frit. Pull through the remaining solvent with vacuum or positive pressure for 10–20 seconds.

**Note:** DCM is a suitable alternative elution solvent

### Post Elution

Evaporate to dryness at 40 °C in a stream of air or nitrogen.

### Reconstitution\*

Reconstitute using 500 µL of 20/80 H<sub>2</sub>O/ACN with 0.1% Formic acid. Mix gently.

\*Recovery and reproducibility for ethylestrenol was found to be affected by non-specific binding and/or losses on evaporation and reconstitution. We recommend that collection vessels, evaporation vessels and reconstitution solvents are investigated during method development to minimize losses/variability.

**Table 1.** Loading and elution details for alternative ISOLUTE SLE+ formats

Format	Load Volume	Elution Protocol	Reconstitution Solvent
ISOLUTE SLE+ 400 µL capacity plate	300 µL of pre-treated urine	2 x 600 µL MTBE or DCM	200 µL 20/80 H <sub>2</sub> O/MeOH 0.1% Formic acid.
ISOLUTE SLE+ 400 µL capacity column	300 µL of pre-treated urine	2 x 600 µL MTBE or DCM	200 µL 20/80 H <sub>2</sub> O/MeOH 0.1% Formic acid.
ISOLUTE SLE+ 1 mL capacity column	800 µL of pre-treated urine	1 x 3 mL MTBE or DCM	200 µL 20/80 H <sub>2</sub> O/MeOH 0.1% Formic acid.

## UPLC Conditions

### Instrument

Waters Acquity UPLC (Waters Assoc., Milford, MA, USA)

### Column

ACE Excel 2 C18 column (50 x 2.1 mm id)

### Mobile Phase

A: 2 mM ammonium acetate 0.1% Formic Acid (aq)

B: 2 mM ammonium acetate 0.1% Formic Acid (MeOH)

### Flow Rate

0.4 mL/min

**Table 2.** UPLC Gradient Conditions.

Time	%A	%B	Curve
0	25	75	
0.25	25	75	6
1.50	10	90	6
1.51	0	100	6
3.50	0	100	6
3.51	25	75	1
4.00	25	75	1

**Curve 6:** Linear Gradient

### Injection Volume

10 µL (partial loop with overfill)

### Sample Temperature

20 °C

### Column Temperature

40 °C

## Mass Spectrometry Conditions

### Instrument

Premier XE triple quadrupole mass spectrometer equipped with an electrospray interface for mass analysis.

### Desolvation Temperature

450 °C

### Ion Source Temperature

120 °C

Positive ions acquired in the multiple reaction monitoring (MRM) mode:

**Table 3.** Mass Spectrometry Conditions.

Compound	MRM Transition	Cone Voltage (V)	Collision Energy (eV)
Stanozolol	329.0 > 94.8	50	40
Stanozolol	329.0 > 106.8	50	40
Ethylestrenol	271.1 > 120.8	35	15
Ethylestrenol	274.1 > 146.8	35	15


## Results

### Analyte Recovery and Reproducibility

Horse urine (gelding and filly) was spiked (n=7 for each matrix) with both analytes at a sample concentration of 40 ng/mL. Both dichloromethane and MTBE were evaluated as elution solvents. Typical recovery and reproducibility for each format are shown in **Table 4** below.

**Table 4.** Typical analyte recovery and reproducibility using various ISOLUTE® SLE+ formats

Elution Solvent	Format	Stanozolol (40 ng/ mL)		Ethylestrenol (40 ng/ mL)	
		Recovery	% RSD	Recovery	% RSD
DCM	ISOLUTE SLE+ 400 µL capacity plate	>85%	<3%	>84%	<7%
DCM	ISOLUTE SLE+ 400 µL capacity column	>87%	<7%	>82%	<12%
MTBE	ISOLUTE SLE+ 1 mL capacity plate (48-well)	>90%	<9%	>88%	<11%
DCM	ISOLUTE SLE+ 1 mL capacity column	>89%	<11%	>85%	<16%



**Figure 2.** Recovery profile for anabolic steroids with MTBE elution, 20:80 H<sub>2</sub>O:ACN 0.1% Formic acid reconstitution solution using 48-well plate format.




Figure 3. LC/MS chromatography

## Ordering Information

Part Number	Description	Quantity
<b>820-0400-P01</b>	ISOLUTE® SLE+ 400 Supported Liquid Extraction Plate	1
<b>820-1000-Q01</b>	ISOLUTE® SLE+ 1 mL Supported Liquid Extraction Plate, 48-well	1
<b>820-0055-B</b>	ISOLUTE® SLE+ 400 µL Sample Volume Columns	50
<b>820-0140-C</b>	ISOLUTE® SLE+ 1 Sample Volume Columns	30
<b>121-9600</b>	Biotage® VacMater™-96 Sample Processing Manifold	1
<b>PPM-96</b>	Biotage® PRESSURE+ 96 Positive Pressure Manifold	1
<b>SD-9600-DHS-EU</b>	Biotage® SPE Dry Sample Concentrator System 220/240 V	1
<b>SD-9600-DHS-NA</b>	Biotage® SPE Dry Sample Concentrator System 100/120 V	1
<b>C103263</b>	TurboVap® 96, Evaporator 100/120V	1
<b>C103264</b>	TurboVap® 96, Evaporator 220/240V	1

## Additional Notes

### Buffer Preparation

- 2 mM ammonium acetate aq: Weigh 0.15416 g and dissolve in H<sub>2</sub>O. Dilute and make up to 1 L in H<sub>2</sub>O then add 0.1 mL Formic acid.
- 2 mM ammonium acetate in methanol: Weigh 0.15416 g and dissolve in methanol. Dilute and make up to 1 L in methanol then add 0.1 mL Formic acid.

### EUROPE

Main Office: +46 18 565900  
 Toll Free: +800 18 565710  
 Fax: +46 18 591922  
 Order Tel: +46 18 565710  
 Order Fax: +46 18 565705  
 order@biotage.com  
 Support Tel: +46 18 56 59 11  
 Support Fax: + 46 18 56 57 11  
 eu-1-pointsupport@biotage.com

### NORTH & LATIN AMERICA

Main Office: +1 704 654 4900  
 Toll Free: +1 800 446 4752  
 Fax: +1 704 654 4917  
 Order Tel: +1 704 654 4900  
 Order Fax: +1 434 296 8217  
 ordermailbox@biotage.com  
 Support Tel: +1 800 446 4752  
 Outside US: +1 704 654 4900  
 us-1-pointsupport@biotage.com

### JAPAN

Tel: +81 3 5627 3123  
 Fax: +81 3 5627 3121  
 jp\_order@biotage.com  
 jp-1-pointsupport@biotage.com

### CHINA

Tel: +86 21 2898 6655  
 Fax: +86 21 2898 6153  
 cn\_order@biotage.com  
 cn-1-pointsupport@biotage.com

To locate a distributor,  
 please visit our website  
 www.biotage.com

### Part Number: AN843

© 2015 Biotage. All rights reserved. No material may be reproduced or published without the written permission of Biotage. Information in this document is subject to change without notice and does not represent any commitment from Biotage. E&OE. A list of all trademarks owned by Biotage AB is available at [www.biotage.com/legal](http://www.biotage.com/legal). Other product and company names mentioned herein may be trademarks or registered trademarks and/or service marks of their respective owners, and are used only for explanation and to the owners' benefit, without intent to infringe.