

Application News

Observation of Shock Waves Using the Schlieren Method with the HPV-X3

Yuki Nishikawa

User Benefits

- ◆ HPV-X3 offers three times the resolution of conventional models, enabling high-speed imaging at much higher resolutions.
- ◆ With a maximum frame rate of 20 Mfps, the HPV-X3 is well-suited for observing high-speed phenomena such as shock waves.

Introduction

Shock waves are characterized by abrupt changes in pressure and density and are generated when an object travels faster than the speed of sound or during explosive phenomena. They are important research subjects in fields such as aerospace and fluid dynamics. In particular, accurately characterizing the structure and temporal evolution of shock waves generated around projectiles is a crucial issue directly related to vehicle design, behavior prediction, and safety evaluation. Optical techniques capable of detecting minute changes in refractive index are effective for visualizing such phenomena. Among these, the Schlieren method is widely used to visualize density variations in fluids.¹⁾

In this paper, the HPV-X3 high-speed video camera (Fig. 1), which provides three times the resolution of conventional models, was used to observe shock waves generated by a pellet launched from a gas gun using the Schlieren method.



Fig. 1 HyperVision™ HPV-X3 High-Speed Video Camera




Fig. 2 Projectile

Table 1 Imaging Devices

High-Speed Video Camera:	HPV-X3
Lens:	180 mm macro lens ×2 teleconverter
Light:	Laser illumination

Table 2 Imaging Conditions

	Recording Speed	Objects placed along the projectile path
Imaging (1)	1,150 ns	Two rubber bands
Imaging (2)	1,150 ns	Obstacle (with a circular hole)

Imaging System

The pellet shown at the center of Fig. 2 was placed in the firing chamber of the gas gun, sandwiched between two plastic holders. The holders were used to position the pellet in the firing chamber; after firing, the holders fell away, allowing only the pellet to be captured by the camera.

The imaging system is detailed in Table 1, the imaging conditions are shown in Table 2, and the experimental setup is shown in Fig. 3. Light from the laser source was reflected by Concave mirror (1) and collimated between Concave mirrors (1) and (2). The HPV-X3 captured images of shock waves generated along the projectile path through a mirror and a knife-edge.

The experiment was conducted twice: in the first trial, two rubber bands were placed in the projectile path; in the second trial, an obstacle with ten circular holes (Fig. 4) was placed in the path. This allowed us to observe phenomena such as fragmentation and shock wave reflection.




Fig. 3 Experimental Setup



Fig. 4 View of the Pellet Path

■ Imaging Results

The results of the Imaging (1) are shown in Fig. 5. In image (1), the flying pellet and a shock wave in front of it can be observed. The shock wave ahead is then reflected by the rubber band in image (2), splitting into a transmitted shock wave and a reflected shock wave. In image (3), the pellet collides with the rubber band, and shock waves are generated both ahead of and behind the collision point. At the second rubber band, no clear reflection of the shock wave was observed, possibly because the pellet's speed had decreased.

Next, Fig. 6 shows the imaging results when the obstacle shown in Fig. 4 was placed in the path. In Image (1), the pellet is positioned immediately before the obstacle. In Image (2), shock waves generation can be observed. In Image (3), the pellet collides with the obstacle, fragments into small pieces, and multiple fragments become visible. It was also confirmed that shock waves were generated from each fragment.




Fig. 5 Observation Images of the Pellet Passing through the Rubber Band (Time Interval between Images: Approximately 15 μ s)




Fig. 6 Observation Images of the Pellet Impacting and Passing through the Obstacle (Time Interval between Images: Approximately 23 μ s)

■ Conclusion

Using the HPV-X3 high-speed video camera, shock waves generated by pellets fired from a gas gun were observed using the Schlieren method. This approach enabled close observation of phenomena, including shock waves generated as the pellet passed through or was reflected by rubber bands placed in the projectile path, as well as pellet fragmentation by impact with an obstacle and the subsequent generation of shock waves from each fragment.

The HPV-X3 high-speed video camera has three times the resolution of the previous, HPV-X2, model and is well-suited for detailed observation of high-speed phenomena such as shock waves.

HyperVision and HPV are trademarks of Shimadzu Corporation or its affiliated companies in Japan and/or other countries.

<References>

- 1) K. Ohtani, "High Speed Visualization Measurement of Underwater Shock Wave and Cavitation Bubble Generation," 17-22, 2023.



Shimadzu Corporation

www.shimadzu.com/an/

For Research Use Only. Not for use in diagnostic procedures.

This publication may contain references to products that are not available in your country. Please contact us to check the availability of these products in your country.

The content of this publication shall not be reproduced, altered or sold for any commercial purpose without the written approval of Shimadzu. See <https://www.shimadzu.com/about/trademarks/index.html> for details.

Third party trademarks and trade names may be used in this publication to refer to either the entities or their products/services, whether or not they are used with trademark symbol "TM" or "®".

Shimadzu disclaims any proprietary interest in trademarks and trade names other than its own.

The information contained herein is provided to you "as is" without warranty of any kind including without limitation warranties as to its accuracy or completeness. Shimadzu does not assume any responsibility or liability for any damage, whether direct or indirect, relating to the use of this publication. This publication is based upon the information available to Shimadzu on or before the date of publication, and subject to change without notice.

01-00983-EN

First Edition: Feb. 2026

[> Please fill out the survey](#)

Related Products

Some products may be updated to newer models.



[> HyperVision HPV-X3](#)
High-Speed Video Camera

Related Solutions

[> Engineering Materials](#)

[> Rubber](#)

[> Price Inquiry](#)

[> Product Inquiry](#)

[> Technical Service / Support Inquiry](#)

[> Other Inquiry](#)