

Application News

HyperVision™ HPV™-X3 High-Speed Video Camera
Autograph™ AGX™-V2 Series Precision Universal Testing Machine

3D-DIC Analysis in Compression Test of Glass Tube

Fumiaki Yano, Yuki Nishikawa

User Benefits

- ◆ The HPV-X3 high-speed video camera has resolution three times higher than the conventional device, resulting in improved analysis performance in DIC analyses.
- ◆ The HPV-X3 high-speed video camera allows recording at a maximum framerate of 20 Mfps, two times faster than the conventional device.
- ◆ Three-dimensional displacement measurement is possible by synchronous recording using two HPV-X3 high-speed video cameras and 3D-DIC analysis.

■ Introduction

Glass is used in various industries, such as electronic devices and automobiles. Depending on the intended application, mechanical properties are measured by material testing, but to achieve further improvements in performance, it is necessary to observe the origin of fracture and the condition of crack propagation. However, glass displays brittle fracture behavior, which is difficult to observe visually. In our previous reports, we observed the fracture behavior of glass specimens in ring-on-ring tests^{1,2}, and those observation results indicated that a high-speed video camera with a recording speed (framerate) of several Mfps or higher is required to observe the fracture behavior of glass. In this experiment, the fracture behavior of a glass tube during a compression test was observed with an HPV-X3 high-speed video camera. In addition, synchronous videorecording was carried out using two HPV-X3 cameras, and the strain distribution during fracture was visualized by a 3D-DIC (three-dimensional digital correlation) analysis.

■ Specimen Information and Measurement System

The glass tube used as the specimen was compressed in the axial direction using an AGX-V2 precision universal testing machine. Table 1 and Table 2 show the equipment configuration and the specimen dimensions, respectively. In this article, fracture observation in the compression test and visualization of the strain distribution by 3D-DIC analysis were carried out. Fig. 1 shows the condition of the test for the 3D-DIC analysis. When an observation target such as a glass tube has a curved shape, a 3-dimensional strain analysis is necessary. Therefore, synchronous recording from two directions was conducted using two HPV-X3 high-speed video cameras, with a random pattern drawn on the specimen surface, as shown in Fig. 2. Table 3 shows the recording conditions.




Fig. 1 Condition of Test (for 3D-DIC Analysis)




Fig. 2 Test Jig and Specimen (for 3D-DIC Analysis)

Table 1 Equipment Configuration

High-speed video camera	: HPV-X3, 2 units
Lens	: 100 mm macro lens, 2 units
Illumination	: SLG-600G
Precision universal testing machine	: AGX-100kNV2
Load cell	: 100 kN
Test jig	: Compression plate (φ50)
DIC analysis software	: VIC-3D

Table 2 Specimen Dimensions

Outer diameter	: 20 mm
Inner diameter	: 18 mm
Height	: 40 mm

Table 3 Recording Conditions

Recording target	: Recording 1 Fracture observation Recording 2 3D-DIC analysis
Recording speed	: 5 Mfps
Recording trigger	: Acceleration sensor

■ Recording Results

Fig. 3 shows the recording results of fracture observation in the compression test of the glass tube. In image (1) in Fig. 3, an axial crack has already initiated, and in image (2), a new crack different from the one in image (1) has also occurred. In images (3) to (5), axial propagation of the crack in image (2) was observed. Subsequently, in images (6) to (8), a similar crack propagated around the circumference of the glass tube, and in images (9) to (12), the condition when the tube shatters into pieces was observed.

Fig. 4 shows the results of the 3D-DIC analysis. From image (1) in Fig. 4, the origin of fracture is on left rear side, and the crack propagates around the circumference in images (2) to (9). Compressive strain concentrates in the unfractured area in front of the crack, as shown in purple. Since the crack propagates around the circumference, the purple area gradually moves to the right. On the other hand, because compressive strain cannot be maintained after crack propagation, strain decreases, as shown in red, and with crack propagation, the decrease in strain also spreads to the upper right.




Fig. 3 Condition of Fracture Observation (Time Interval Between Images: 2 μ s)




Fig. 4 Condition of Fracture Observation (Time Interval Between Images: 400 ns)

■ Conclusion

The fracture of a glass tube during a compression test was observed using two HPV-X3 high-speed video cameras, and a 3D-DIC analysis was carried out. In fracture observation, the crack that became the origin of fracture could be identified, and the condition of crack propagation was clearly captured. In the 3D-DIC analysis, it was possible to observe the concentration of strain in front of the crack and the decrease in strain after crack propagation. Since glass fracture is brittle and is extremely rapid, a high-speed video camera with a recording speed of 5 Mfps or higher, such as the HPV-X3, is an effective instrument for fracture observation.

<Related Application News Articles>

1. Fracture Observation of Glass in Ring-on-Ring Bending Test, [Application News No. V30](#)
2. Fracture Observation in Glass Ring-on-Ring Bending Tests Using HPV-X3, [Application News No. 01-00881](#)

The authors wish to thank Energy Cycle Systems Research, K.Sato Lab., Institute of Advanced Materials Strength, Graduate School of Engineering, Tohoku University for his generous cooperation in the recording.

HyperVision, HPV, Autograph, and AGX are trademarks of Shimadzu Corporation or its affiliated companies in Japan and/or other countries.



SHIMADZU

Shimadzu Corporation

www.shimadzu.com/an/

For Research Use Only. Not for use in diagnostic procedures.

This publication may contain references to products that are not available in your country. Please contact us to check the availability of these products in your country.

The content of this publication shall not be reproduced, altered or sold for any commercial purpose without the written approval of Shimadzu. See <https://www.shimadzu.com/about/trademarks/index.html> for details.

Third party trademarks and trade names may be used in this publication to refer to either the entities or their products/services, whether or not they are used with trademark symbol "TM" or "®".

Shimadzu disclaims any proprietary interest in trademarks and trade names other than its own.

The information contained herein is provided to you "as is" without warranty of any kind including without limitation warranties as to its accuracy or completeness. Shimadzu does not assume any responsibility or liability for any damage, whether direct or indirect, relating to the use of this publication. This publication is based upon the information available to Shimadzu on or before the date of publication, and subject to change without notice.

01-01052-EN

First Edition: Jan. 2026

➤ Please fill out the survey

Related Products

Some products may be updated to newer models.



➤ **HyperVision HPV-X3**
High-Speed Video Camera



➤ **AUTOGRAPH**
AGX™-V2 Series
Precision Universal Testing Machines

Related Solutions

➤ Engineering
Materials

➤ Glass, Ceramics

➤ Price Inquiry

➤ Product Inquiry

➤ Technical Service /
Support Inquiry

➤ Other Inquiry