

Quantitation of Pesticide Residues in Milk Using the Agilent 6470 Triple Quadrupole LC/MS

Analysis of pesticide residues in milk with reporting of lipophilic pesticides based on milk fat percentage



Figure 1. Agilent 1290 Infinity II LC coupled to an Agilent 6470 triple quadrupole LC/MS.

Authors

Prasanth Joseph,
Saikat Banerjee, and
Samir Vyas
Agilent Technologies, Inc.

Abstract

Pesticides are used to support the growth of food crops, control vector-borne diseases, or to control animal parasites. Such pesticides may transfer to livestock animals primarily through feed and fodder. After ingestion, some pesticides are metabolized and deposited throughout the fat and muscle tissue of the animal's body. Milk from these animals contains a considerable amount of fat. Many regulatory agencies specify the maximum residue limits (MRLs) for pesticides based on fat (%) in milk and milk-based products. To detect the MRLs required in milk and milk-based products, the instrument must deliver good sensitivity for analytes in matrix. This application note provides a method for the analysis of pesticide residue based on fat percentage in milk samples. Additionally, this application note describes a detailed methodology for the quantitation of pesticide residues carbofuran, carbendazim, benomyl, edifenphos, and 2,4-D in milk matrix using LC/MS/MS. Quantitation of pesticides in milk is based on QuEChERS extraction, followed by instrument analysis without adopting extensive cleanup procedures.

Introduction

Pesticides used in agriculture may enter an animal's body through feed and drinking water. Some of these pesticides are fat soluble and deposit into the animal's fat and muscle tissue. Milk contains a considerable amount of fat; therefore, regulatory agencies mention a maximum residue limit (MRL) for lipophilic pesticides, based on milk fat – for safe consumption.

QuEChERS methodology has been a widely accepted method to extract pesticide residues from various food commodities such as fruits and vegetables. For milk, a QuEChERS-based extraction of pesticide residues was found to be a straightforward sample preparation to adopt. For detection and quantitation, a highly selective multiple reaction monitoring (MRM)-based LC/MS/MS method was developed using the Agilent 6470 triple quadrupole LC/MS (LC/TQ). The sensitivity of the 6470 LC/TQ can provide analysis of milk samples at the required limits of detection, and report the results based on milk fat percentage.

Experimental

Chemicals and reagents

LC/MS pesticide standards were purchased from Sigma-Aldrich. These standards were used for method development and analysis of milk samples. Agilent Bond Elut QuEChERS extraction kit (part number 5982-5650) was used for sample preparation. LC/MS-grade solvents such as acetonitrile and water were purchased from Honeywell (Charlotte, NC, USA). Acetic acid of MS grade was purchased from Fluka (now of Honeywell).

Instrument configuration

- Agilent 1290 Infinity II high-speed pump (G7120A)
- Agilent 1290 Infinity II multisampler (G7167B)
- Agilent 1290 Infinity II multicolumn thermostat (G7116B)
- Agilent 6470 triple quadrupole LC/MS (G6470B)

Data acquisition and data analysis

All samples were acquired using the Agilent MassHunter Data Acquisition software version 10.1. Chromatograms were reviewed through MassHunter qualitative analysis software version 10.0. Quantitation of each batch was carried out using MassHunter quantitative analysis software version 10.1.

Sample preparation




Figure 2. Flowchart for sample preparation.

Table 4. MRM parameters.

Compound ID	Precursor Ion (m/z)	Product Ion (m/z)	Frag.	CE	CAV	Ionization
Edifenphos	311	111	88	24	4	ESI Positive
	311	109	88	40	4	ESI Positive
Carbofuran	222	165.1	84	12	4	ESI Positive
	222	123	84	24	4	ESI Positive
Carbendazim	192	160	78	20	4	ESI Positive
	192	132	78	36	4	ESI Positive
2,4-D	221	163	56	12	4	ESI Negative
	219	161	56	12	4	ESI Negative

Table 1. Chromatography conditions.

Parameter	Value
Mobile Phase A	0.1% acetic acid in water
Mobile Phase B	Acetonitrile
Flow Rate	0.3 mL/min
Injection Volume	5 µL
Column Temperature	40 °C
Sample Diluent	Acetonitrile/water (60/40)
Needle Wash	MeOH/acetonitrile/water (25/50/25)
Column	Agilent XDB C18, 3.0 × 150 mm, 3.5 µm (p/n 963954-302)

Table 2. Gradient.

Time (min)	%A	%B
0	95	5
1.5	95	5
4	50	50
9	0	80
12	0	80
Post Run	2 minutes	

Table 3. MS source parameters.

Ionization Source	AJS ESI
Ionization Mode	ESI Positive
Gas Temperature	200 °C
Gas Flow	12 L/min
Nebulizer	35 psi
Sheath Gas	390 °C
Sheath Gas Flow	11.5 L/min
Capillary Voltage	2,000 V
Nozzle Voltage	2,000 V

Table 5. Dilution chart for prespike calibration curve.

Working Standard Concentration	Volume Taken (μL)	Volume of Milk (mL)	Obtained Concentration (ng/mL)
100 ppb	40	10	0.4
100 ppb	100	10	1
100 ppb	200	10	2
1 ppm	50	10	5
1 ppm	100	10	10
1 ppm	200	10	20
10 ppm	50	10	50
10 ppm	75	10	75
10 ppm	100	10	100
10 ppm	200	10	200

Results and discussion

Individual stock solutions were made and mixed to appropriate volumes of stock solutions to make working standards. Working standards prepared were of concentrations 1 ppm, 10 ppm, and 100 ppb. Individual volumes of working standards were spiked in 10 mL of milk samples, then followed the sample preparation flowchart to generate matrix-based calibration (prespike) points.

Matrix-based calibration curves (Figure 4) were found to be linear for carbofuran and edifenphos from 0.4 ng/mL to 100 ng/mL with both linear regression and 1/x. A linear regression line was plotted for carbendazim and 2,4-D with 1/x² weighing. Regression coefficient values were above 0.995 for all four compounds. Benomyl was found to be very unstable and converted to carbendazim, therefore, no calibration curve or measurement was made. In this case, benomyl was reported as carbendazim.




Figure 3. Extracted ion chromatogram (EIC) at 10 ng/mL prespiked in milk matrix.




Figure 4. Matrix-based calibration curves from 0.4 ng/mL to 100 ng/mL.

Spiking was made per the maximum residue limits mentioned in BIS (Bureau of Indian Standards) regulation documents. Recovery percentage was calculated based on this experiment.

Spiking was made for edifenphos using the maximum residue limit of 10 ppb fat

- This analysis considered milk fat percentage to be 4%.
- Fat level in milk: 4 g in 100 mL
- Sample volume: 10 mL
- Fat (g) in 10 mL: 0.4 g
- To get a spike level of 10 ppb (fat basis), a spike concentration of 4 ng in 10 mL milk sample is needed.

- The calculation is as follows:
 $4 \text{ ng}/10 \text{ mL} = 4 \text{ ng}/0.4 \text{ g fat} = 10 \text{ ng/g fat} = 10 \text{ ppb fat}$.
- To spike 4 ng in 10 mL, 40 μL of 100 ppb in 10 mL milk must be spiked.
- A spike sample at the 0.4 ppb level showed a recovery of more than 84%.

Batch Table												
Sample: MB 100 PPB		Sample Type: <All>		Compound: Ediphenphos		ISTD:						
Sample					Ediphenphos		Ediphenphos Results				Qualifi...	
ID	Name	Data File	Type	Level	Exp. Conc.	RT	Resp.	MI	Final Conc.	Accuracy	Ratio	
	MB BLANK	MB BLANK-r003...	Sample			10.992	1896	<input type="checkbox"/>	0.0000		14.2	
	MB 0.4 PPB	MB 0.4 PPB.d	Cal	1	0.4000	10.981	35672	<input type="checkbox"/>	0.4416	110.4	62.1	
	MB 1 PPB	MB 1.0 PPB.d	Cal	2	1.0000	10.969	71705	<input type="checkbox"/>	1.0355	103.5	63.4	
	MB 5 PPB	MB 5.0 PPB.d	Cal	3	5.0000	10.957	319565	<input type="checkbox"/>	5.1209	102.4	74.2	
	MB 10 PPB	MB 10.0 PPB.d	Cal	4	10.0000	10.957	635246	<input type="checkbox"/>	10.3241	103.2	75.1	
	MB 20 PPB	MB 20.0 PPB.d	Cal	5	20.0000	10.957	1201605	<input type="checkbox"/>	19.6592	98.3	77.6	
	MB 50 PPB	MB 50.0 PPB.d	Cal	6	50.0000	10.957	2998108	<input type="checkbox"/>	49.2702	98.5	77.9	
	MB 75 PPB	MB 75.0 PPB.d	Cal	7	75.0000	10.969	4489234	<input type="checkbox"/>	73.8479	98.5	78.6	
	MB 100 PPB	MB 100.0 PPB.d	Cal	8	100.0000	10.969	6183108	<input type="checkbox"/>	101.7673	101.8	78.2	
	MB SPIKE	MB SPIKE.d	QC	9	0.4000	10.969	29451	<input type="checkbox"/>	0.3390	84.8	63.3	

Figure 5. Edifenphos results using the maximum residue limit of 10 ppb fat.

Spiking was made for carbofuran at the maximum residue limit of 50 ppb fat

- To get a spike level of 50 ppb (fat basis), a spike concentration of 20 ng in 10 mL milk sample is needed.
- The calculation is as follows:
 $20 \text{ ng}/10 \text{ mL} = 20 \text{ ng}/0.4 \text{ g fat} = 50 \text{ ng/g fat} = 50 \text{ ppb fat}$.
- To spike 20 ng in 10 mL, 200 μL of 100 ppb in 10 mL milk must be spiked.
- A spike sample at the 20 ppb level showed a recovery of more than 92%.

Spiking was made for carbendazim at the maximum residue limit of 100 ppb fat

- To get a spike level of 100 ppb (fat basis), a spike concentration of 40 ng in 10 mL milk sample is needed.
- The calculation is as follows:
 $40 \text{ ng}/10 \text{ mL} = 40 \text{ ng}/0.4 \text{ g fat} = 100 \text{ ng/g fat} = 100 \text{ ppb fat}$.
- To spike 40 ng in 10 mL, 40 μL of 1 ppm in 10 mL milk must be spiked.
- A spike sample at the 40 ppb level showed a recovery of more than 95%.

Spiking was made for 2,4-D at the maximum residue limit of 50 ppb

- To get a spike level of 50 ppb, a spike concentration of 50 ng in 1 mL milk sample is needed.
- The calculation is as follows:
 $500 \text{ ng}/10 \text{ mL} = 50 \text{ ppb}$.
- To spike 500 ng in 10 mL, 50 μL of 10 ppm in 10 mL milk must be spiked.
- A spike sample at the 50 ppb level showed a recovery of more than 92%.

Batch Table												
Sample: MB 100 PPB		Sample Type: <All>		Compound: Carbofuran		ISTD:						
Sample					Carbofur...		Carbofuran Results					Qualifi...
?	▼	Name	Data File	Type	Level	Exp. Conc.	RT	Resp.	MI	Final Conc.	Accuracy	Ratio
		MB BLANK	MB BLANK-r003...	Sample			7.806	229	<input type="checkbox"/>	0.1357		
		MB 0.4 PPB	MB 0.4 PPB.d	Cal	1	0.4000	7.818	77814	<input type="checkbox"/>	0.4615	115.4	75.2
		MB 1 PPB	MB 1.0 PPB.d	Cal	2	1.0000	7.818	205061	<input type="checkbox"/>	0.9957	99.6	77.8
		MB 5 PPB	MB 5.0 PPB.d	Cal	3	5.0000	7.818	1088528	<input type="checkbox"/>	4.7051	94.1	76.8
		MB 10 PPB	MB 10.0 PPB.d	Cal	4	10.0000	7.818	2213712	<input type="checkbox"/>	9.4294	94.3	78.5
		MB 20 PPB	MB 20.0 PPB.d	Cal	5	20.0000	7.818	4522421	<input type="checkbox"/>	19.1230	95.6	78.7
		MB 50 PPB	MB 50.0 PPB.d	Cal	6	50.0000	7.818	11654697	<input type="checkbox"/>	49.0691	98.1	79.5
		MB 75 PPB	MB 75.0 PPB.d	Cal	7	75.0000	7.818	18041329	<input type="checkbox"/>	75.8846	101.2	78.2
		MB 100 PPB	MB 100.0 PPB.d	Cal	8	100.0000	7.818	24197302	<input type="checkbox"/>	101.7316	101.7	80.3
		MB SPIKE	MB SPIKE.d	QC	9	20.0000	7.818	4387832	<input type="checkbox"/>	18.5579	92.8	80.0

Figure 6. Carbofuran results using the maximum residue limit of 50 ppb fat.

Batch Table												
Sample: MB 100 PPB		Sample Type: <All>		Compound: Carbendazim		ISTD:						
Sample					Carbendaz...		Carbendazim Results					Qualifi...
?	▼	Name	Data File	Type	Level	Exp. Conc.	RT	Resp.	MI	Final Conc.	Accuracy	Ratio
		MB BLANK	MB BLANK-r003...	Sample			5.414	4606	<input type="checkbox"/>	0.0831		24.9
		MB 0.4 PPB	MB 0.4 PPB.d	Cal	1	0.4000	5.414	27874	<input type="checkbox"/>	0.4116	102.9	23.1
		MB 1 PPB	MB 1.0 PPB.d	Cal	2	1.0000	5.414	65455	<input type="checkbox"/>	0.9422	94.2	22.1
		MB 5 PPB	MB 5.0 PPB.d	Cal	3	5.0000	5.414	332331	<input type="checkbox"/>	4.7099	94.2	19.7
		MB 10 PPB	MB 10.0 PPB.d	Cal	4	10.0000	5.414	680043	<input type="checkbox"/>	9.6189	96.2	19.8
		MB 20 PPB	MB 20.0 PPB.d	Cal	5	20.0000	5.414	1379767	<input type="checkbox"/>	19.4975	97.5	20.3
		MB 50 PPB	MB 50.0 PPB.d	Cal	6	50.0000	5.402	3629817	<input type="checkbox"/>	51.2634	102.5	20.8
		MB 75 PPB	MB 75.0 PPB.d	Cal	7	75.0000	5.414	5540561	<input type="checkbox"/>	78.2391	104.3	21.1
		MB 100 PPB	MB 100.0 PPB.d	Cal	8	100.0000	5.414	7659812	<input type="checkbox"/>	108.1584	108.2	21.0
		MB SPIKE	MB SPIKE.d	QC	9	40.0000	5.414	2718107	<input type="checkbox"/>	38.3920	96.0	20.8

Figure 7. Carbendazim results using the maximum residue limit of 100 ppb fat.

Batch Table												
Sample: MB 100 PPB		Sample Type: <All>		Compound: 2,4-D		ISTD:						
Sample					2,4-D Met...		2,4-D Results					Qualifi...
?	▼	Name	Data File	Type	Level	Exp. Conc.	RT	Resp.	MI	Final Conc.	Accuracy	Ratio
		MB BLANK	MB BLANK-r003...	Sample			8.140	54	<input type="checkbox"/>	0.0811		
		MB 0.4 PPB	MB 0.4 PPB.d	Cal	1	0.4000	8.187	9467	<input type="checkbox"/>	0.4133	103.3	65.0
		MB 1 PPB	MB 1.0 PPB.d	Cal	2	1.0000	8.187	24200	<input type="checkbox"/>	0.9331	93.3	60.8
		MB 5 PPB	MB 5.0 PPB.d	Cal	3	5.0000	8.187	130273	<input type="checkbox"/>	4.6760	93.5	62.4
		MB 10 PPB	MB 10.0 PPB.d	Cal	4	10.0000	8.175	269695	<input type="checkbox"/>	9.5955	96.0	61.1
		MB 20 PPB	MB 20.0 PPB.d	Cal	5	20.0000	8.175	549477	<input type="checkbox"/>	19.4677	97.3	62.1
		MB 50 PPB	MB 50.0 PPB.d	Cal	6	50.0000	8.175	1447667	<input type="checkbox"/>	51.1607	102.3	62.5
		MB 75 PPB	MB 75.0 PPB.d	Cal	7	75.0000	8.175	2216920	<input type="checkbox"/>	78.3041	104.4	62.7
		MB 100 PPB	MB 100.0 PPB.d	Cal	8	100.0000	8.175	3110061	<input type="checkbox"/>	109.8188	109.8	62.2
		MB SPIKE	MB SPIKE.d	QC	9	50.0000	8.175	1316350	<input type="checkbox"/>	46.5271	93.1	63.1

Figure 8. 2,4-D results using the maximum residue limit of 50 ppb.

A matrix-based standard, prepared by prespiking the pesticide standards at the 10 ppb level in milk matrix, was injected with six replicates. Consistency of the chromatogram with respect to the response and retention time demonstrates the repeatability of the instrument result.

Conclusion

A highly sensitive and robust polarity-switching method based on Multiple Reaction Monitoring was developed to quantify pesticide residues in milk matrix.

Method performance was evaluated for sensitivity, specificity, linearity, reproducibility, and recovery. This note describes the reporting of pesticides based on the percentage of milk fat. This method can be adopted for routine pesticide quantitation in milk samples.

References

1. Comprehensive LC/MS/MS Workflow of Pesticide Residues in Food Using the Agilent 6470 Triple Quadrupole LC/MS System. *Agilent Technologies application note*, publication number 5994-2370EN, **2020**.
2. Bedi, J. S. *et al.* Pesticide residues in milk and their relationship with pesticide contamination of feedstuffs supplied to dairy cattle in Punjab (India). *Journal of Animal and Feed Sciences* **2018**, *27*, 18–25.




Figure 9. Overlaid chromatograms of six replicate injections at 10 ppb concentration spiked in milk matrix.




Figure 10. Repeatability in area response for six replicate injections at 10 ppb in milk matrix.

www.agilent.com/chem

DE44260.6023148148

This information is subject to change without notice.

© Agilent Technologies, Inc. 2021
Printed in the USA, April 7, 2021
5994-3225EN