

# Fast and Reliable Analysis of Soil and Sediments using ICP-MS with an Innovative Cell

High-throughput and low detection limits with  
the Agilent 9500 ICP-QQQ and a Dual-Cell System



## Author

Tetsuo Kubota  
Agilent Technologies, Inc.

## Introduction

Ever since ICP-MS instruments were first commercialized, their high sensitivity and multi-elemental capabilities have made them a suitable technique for quantifying trace elements in soil and sediment, particularly in support of environmental monitoring and regulatory compliance. To help laboratories transition to ICP-MS for these important environmental applications, the US Environmental Protection Agency (EPA) first introduced the performance-based test Method 6020 in 2007, and revised it as Method 6020B in 2014.<sup>1</sup> Method 6020B remains widely used worldwide for the quantitative measurement of a broad range of elements, including arsenic (As), cadmium (Cd), mercury (Hg), selenium (Se), and lead (Pb), in solid and aqueous wastes by ICP-MS.

Traditionally, single-quadrupole ICP-MS (SQ ICP-MS) instruments equipped with collision/reaction cells (CRCs) have been used to meet the requirements of Method 6020B for sample analysis. Laboratories have typically used a combination of acquisition strategies to control interferences, including no gas mode for

interference-free analytes and low mass analytes, helium (He) collision mode for analytes affected by polyatomic ion interferences, and half-mass mode to correct doubly charged ion ( $M^{2+}$ ) interferences. This multi-mode approach has enabled laboratories to meet Method 6020B performance and quality control (QC) criteria over many years.

Now, for the first time, the Agilent 9500 Triple Quadrupole ICP-MS featuring the Dual-Cell System (DCS) cell introduces two innovative cell gas modes, Advanced Helium Mode (AHM) and Air cell mode.<sup>2,3</sup> These modes simplify interference removal approaches, enhance analytical performance, and accelerate analysis times.

- AHM
  - Significantly enhances sensitivity for low-mass analytes while providing robust interference removal across the entire mass range.
  - Eliminates the need for multiple gas modes, replacing no gas, conventional He, and high energy (HE) He mode with a single gas mode.
- Air cell mode
  - Uses oxygen ( $O_2$ ) from ambient air as a reactive cell gas to improve detection limits (DLs) for some analytes.
  - Effectively eliminates interferences, including  $M^{2+}$  that are challenging to resolve in He collision mode.
  - Eliminates the need for  $O_2$  gas cylinders, reducing operating costs and simplifying workflows.

In this study, a range of environmentally certified reference materials (CRMs) was selected to represent sample types routinely analyzed in high-throughput environmental testing laboratories. As these laboratories prioritize productivity across all aspects of their operations to remain competitive, the 9500 ICP-QQQ configuration and analytical method were optimized to maximize sample throughput. The performance of the method was evaluated according to criteria set out in Method 6020B.

## Experimental

### Instrumentation

All measurements were performed using an Agilent 9500 ICP-QQQ that includes the following standard components: Ultra High Matrix Introduction (UHMI), nickel (Ni) interface cones, and a unique DCS cell. The design features and functionality of the 9500 and DCS are explained elsewhere.<sup>2,3</sup>

The 9500 was also equipped with the optional Agilent Advanced Valve System (AVS MS) discrete sampler for improved productivity. The AVS MS, which uses a high-speed piston pump, integrates seamlessly with the Agilent SPS 4 and SPS 6 autosamplers for sample delivery, enabling the high sample throughput required by most environmental testing laboratories. The UHMI guarantees high matrix tolerance of the 9500 using aerosol dilution. The AVS MS further improves robustness for large sample batches of high matrix samples by minimizing the time that the instrument is exposed to the sample during each measurement.

Method parameters (recommended list of elements, preferred isotopes, integration times, internal standards, etc.) were loaded from the "EPA 6020B" preset method in the Agilent OpenLab ICP-MS software, allowing the analyst to quickly create a new batch method. All the plasma parameters listed in Table 1 were automatically loaded by selecting UHMI-4.

**Table 1.** Agilent 9500 ICP-QQQ operating parameters.

|                                 | Advanced He Mode (AHM) | Air Cell Mode |
|---------------------------------|------------------------|---------------|
| Plasma Mode                     | UHMI-4                 |               |
| RF Power (W)                    | 1600                   |               |
| Sampling Depth (mm)             | 12                     |               |
| Nebulizer Gas (L/min)           | 0.94                   |               |
| Dilution Gas (L/min)            | 0.16                   |               |
| Extract 1 (V)                   | -11.2                  |               |
| Extract 2 (V)                   | -300                   |               |
| Omega Bias                      | -200                   |               |
| Omega Lens                      | -5.0                   |               |
| Cell Gas Flow Rate (mL/min)     | 14                     | 0.4           |
| Front Energy Discrimination (V) | 10                     | -1            |

**Table 2.** Agilent AVS MS operating parameters.

|                                 | Time (s) | Speed (%) |
|---------------------------------|----------|-----------|
| Sample Load                     | 11       | 50        |
| Stabilization                   | 18       | 5         |
| Probe Rinse                     | 15       | 5         |
| Rinse 1                         | 6        | 80        |
| Rinse 2                         | 35       | 5         |
| Optional Loop Probe Wash        | 10       | 50        |
| Optional Loop Wash              | 1        | 50        |
| Lower Sample Consumption Option | 3        | -         |
| Sample-to-Sample Analysis Time  | 122      | -         |

## Chemicals and reagents

A solution containing high-purity (EL grade) 1% nitric acid (HNO<sub>3</sub>) and 0.5% hydrochloric acid (HCl) from Kanto Chemical Co., Inc., Japan, was used as the blank, diluent, and rinse solutions. For most elements, calibration standards, spikes, and continuous calibration verification (CCV) solutions were prepared from the Agilent multi-element Environmental Calibration Standard, and initial calibration verification (ICV) solutions were prepared from the Agilent mixed multi-element ICV standard. Calibration standards for Li, Al, Mn, Zn, Sr, and Hg were prepared using 1000 ppm single-element stocks (Kanto Chemical Co., Inc.). A mixed internal standard (ISTD) comprising Sc, Ge, Rh, In, Tb, and Lu was prepared from 1000 ppm single-element stocks (Kanto Chemical Co., Inc.).


## Standard and sample preparation

All standards and samples were prepared in an acid matrix of 1% HNO<sub>3</sub> and 0.5% HCl. HCl is routinely added to samples for ICP-MS analysis to ensure the stability of elements such as antimony (Sb) and mercury (Hg). The 9500 ICP-QQQ in AHM collision-cell mode removes any Cl-based interferences that would be formed from the added chloride matrix.

Six-point calibrations, including the calibration blank, were prepared at the concentration ranges given in Figure 1. The term “mineral elements” refers to Na, Mg, Al, K, Ca, Mn, Fe, and Zn.

For in-run QC, the ICV and CCV standards were used.

The CRMs included Trace Elements in Water (NIST 1643f), River Sediment A, River Sediment B, Soil A, and Soil B (High Purity Standards). The CRMs were all supplied as predigested solutions. NIST 1643f was introduced with no dilution, and the other CRMs were diluted 10x using the 1% HNO<sub>3</sub> and 0.5% HCl acid matrix. The ISTD solution containing 1 ppm Sc, Ge, Rh, In, Tb, and Lu was automatically added online via the dedicated seventh port of the AVS MS valve. The flow rate of the ISTD solution was about 1/15 of the sample flow rate, achieved using pump tubing with a smaller internal diameter.



**Figure 1.** Analytical sequence of calibrants, QC solutions, and samples. The sample block was analyzed repeatedly by the Agilent 9500 ICP-QQQ with automatic insertion of the periodic QC block after every 10 samples.

## Results and discussion

### Limits of Detection (LODs)

All 28 analytes were measured by the 9500 ICP-QQQ in AHM, while selected elements (K, Ca, V, As, and Se) were also measured in Air cell mode using the acquisition parameters listed in Table 1. Limits of Detection (LODs) were calculated from three times the standard deviation of 10 measurements of the reagent blank solution (Table 3).

Operating the DCS of the 9500 ICP-QQQ in AHM achieved low-level LODs for low-mass elements such as Li and Be, which are traditionally analyzed in no gas mode. Additionally, AHM delivered an LOD for Se approximately one order of magnitude lower than that of the conventional He mode, even with a shorter integration time (reduced from 3 to 1 s).

Operating the 9500 ICP-QQQ DCS in Air cell mode reduced the K background, improving the LOD to 0.029 µg/L compared with 0.97 µg/L in AHM. For elements such as V, As, and Se, which exhibit high reactivity with O<sub>2</sub>, Air cell mode utilizes oxygen-based mass-shift reactions (for example, As<sup>+</sup> → AsO<sup>+</sup>), resulting in enhanced sensitivity and reduced background compared to conventional He mode. Ca was measured in Air cell mode to demonstrate the M<sup>2+</sup> removal capability of the DCS (removal of <sup>88</sup>Sr<sup>2+</sup> interference on <sup>44</sup>Ca<sup>+</sup>).

LODs for the trace analytes were mostly in the low ng/L (ppt) range, confirming the high sensitivity of the 9500 ICP-QQQ for trace analysis of complex environmental samples.

**Table 3.** Integration times and LODs of analytes in the blank.

| Analyte     | DCS Cell Mode | Integration Time (s) | LOD (µg/L) |
|-------------|---------------|----------------------|------------|
| 7 Li        | AHM           | 0.3                  | 0.0013     |
| 9 Be        | AHM           | 0.3                  | 0.0006     |
| 23 Na       | AHM           | 0.1                  | 0.096      |
| 24 Mg       | AHM           | 0.1                  | 0.020      |
| 27 Al       | AHM           | 0.3                  | 0.024      |
| 39 K        | AHM           | 0.1                  | 0.970      |
| 39 -> 39 K  | Air           | 0.1                  | 0.029      |
| 44 Ca       | AHM           | 0.1                  | 0.260      |
| 44 -> 60 Ca | Air           | 0.1                  | 0.410      |
| 51 V        | AHM           | 0.3                  | 0.0014     |
| 51 -> 67 V  | Air           | 0.3                  | 0.0004     |
| 52 Cr       | AHM           | 0.3                  | 0.0035     |
| 55 Mn       | AHM           | 0.1                  | 0.0051     |
| 56 Fe       | AHM           | 0.1                  | 0.011      |
| 59 Co       | AHM           | 0.3                  | 0.0003     |
| 60 Ni       | AHM           | 0.3                  | 0.0025     |
| 63 Cu       | AHM           | 0.3                  | 0.0040     |
| 66 Zn       | AHM           | 0.1                  | 0.0044     |
| 75 As       | AHM           | 1                    | 0.0006     |
| 75 -> 91 As | Air           | 1                    | 0.0004     |
| 78 Se       | AHM           | 1                    | 0.0057     |
| 78 -> 94 Se | Air           | 1                    | 0.0068     |
| 88 Sr       | AHM           | 0.3                  | 0.0028     |
| 95 Mo       | AHM           | 0.3                  | 0.0009     |
| 107 Ag      | AHM           | 0.3                  | 0.0002     |
| 111 Cd      | AHM           | 0.3                  | 0.0003     |
| 121 Sb      | AHM           | 0.3                  | 0.0049     |
| 137 Ba      | AHM           | 0.3                  | 0.0004     |
| 201 Hg      | AHM           | 1                    | 0.0001     |
| 205 Tl      | AHM           | 0.3                  | 0.00004    |
| Pb*         | AHM           | 0.3                  | 0.0005     |
| 232 Th      | AHM           | 0.3                  | 0.00004    |
| 238 U       | AHM           | 0.3                  | 0.00002    |

\*Reported as the sum of the lead isotopes 206+207+208.

### ICV

The recommended control limits for the ICV in Method 6020B are ±10% of the actual value. All analytes were measured within these control limits (Table 4). The ICV solution contained 1000 µg/L Sr, which caused a Sr<sup>2+</sup> interference on Ca, resulting in a high recovery in AHM. However, the Sr<sup>2+</sup> interference was effectively resolved using a mass-shift reaction with O<sub>2</sub> in Air cell mode, resulting in a recovery of 93.9%.

### Long-term stability


To demonstrate the robustness of the 9500 ICP-QQQ, 180 samples and 38 QC solutions were analyzed for approximately eight hours.

Figure 2 shows the percent recovery of the raw ISTD responses for the entire eight-hour sequence normalized to the calibration blanks. All samples passed the >70% test, with minimal drift. The recovery data confirms that the robust, matrix tolerant plasma of the 9500 ensured that minimal matrix deposition occurred over eight hours. The ISTD plot also shows that the 9500 provides excellent control of signal suppression, with ISTD signals being consistent for elements covering a range of masses and ionization potentials. The results demonstrate the long-term robustness and high matrix tolerance of the 9500 ICP-QQQ with UHMI.


**Table 4.** ICV recoveries.

| Analyte     | DCS Cell Mode | ICV                  |                    |              |
|-------------|---------------|----------------------|--------------------|--------------|
|             |               | Measured Conc (µg/L) | Actual Conc (µg/L) | Recovery (%) |
| 7 Li        | AHM           | 9.84                 | 10                 | 98           |
| 9 Be        | AHM           | 102                  | 100                | 102          |
| 23 Na       | AHM           | 10,000               | 10,000             | 100          |
| 24 Mg       | AHM           | 9,900                | 10,000             | 99           |
| 27 Al       | AHM           | 9,520                | 10,000             | 95           |
| 39 K        | AHM           | 10,200               | 10,000             | 102          |
| 39 -> 39 K  | Air           | 10,100               | 10,000             | 101          |
| 44 Ca       | AHM           | 12,600               | 10,000             | 126          |
| 44 -> 60 Ca | Air           | 9,390                | 10,000             | 94           |
| 51 V        | AHM           | 96.7                 | 100                | 97           |
| 51 -> 67 V  | Air           | 98.7                 | 100                | 99           |
| 52 Cr       | AHM           | 97.8                 | 100                | 98           |
| 55 Mn       | AHM           | 9,570                | 10,000             | 96           |
| 56 Fe       | AHM           | 9,940                | 10,000             | 99           |
| 59 Co       | AHM           | 99.5                 | 100                | 100          |
| 60 Ni       | AHM           | 99.6                 | 100                | 100          |
| 63 Cu       | AHM           | 103                  | 100                | 103          |
| 66 Zn       | AHM           | 9,830                | 10,000             | 98           |
| 75 As       | AHM           | 98.0                 | 100                | 98           |
| 75 -> 91 As | Air           | 98.0                 | 100                | 98           |
| 78 Se       | AHM           | 98.7                 | 100                | 99           |
| 78 -> 94 Se | Air           | 99.2                 | 100                | 99           |
| 88 Sr       | AHM           | 988                  | 1,000              | 99           |
| 95 Mo       | AHM           | 93.7                 | 100                | 94           |
| 107 Ag      | AHM           | 99.4                 | 100                | 99           |
| 111 Cd      | AHM           | 99.7                 | 100                | 100          |
| 121 Sb      | AHM           | 96.4                 | 100                | 96           |
| 137 Ba      | AHM           | 99.1                 | 100                | 99           |
| 201 Hg      | AHM           | 0.524                | 0.5                | 105          |
| 205 Tl      | AHM           | 106                  | 100                | 106          |
| *Pb         | AHM           | 104                  | 100                | 104          |
| 232 Th      | AHM           | 99.4                 | 100                | 99           |
| 238 U       | AHM           | 101                  | 100                | 101          |

\*Reported as the sum of the lead isotopes 206+207+208.



**Figure 2.** ISTD stability during the analysis of 180 samples over approximately eight hours using the Agilent 9500 ICP-QQ. ISTD recoveries normalized to the calibration blank for all samples. Due to limited space, not all sample names are shown.



**Figure 3.** CCV stability of 18 CCVs measured using the Agilent 9500 ICP-MS.

## Accurate analysis of CRMs

The five CRMs were each analyzed 18 times during the sequence, and the mean concentration and recovery were calculated for each analyte (Table 5). Most of the certified elements, which included major elements at 100s of mg/L (ppm) and trace elements at levels as low as 10 µg/L (ppb), were recovered within ±10%. Not all analytes have certified

values in every reference material, therefore, any blank cells in Table 5 indicate the absence of a certified value. Regarding Co in River Sediment A, the spike recovery test showed a good recovery (94.6% for recovery of a 20 µg/L spike), suggesting that the CRM may be contaminated.

**Table 5.** Recovery of certified values for CRMs. Concentration units µg/L.

| n=18        |               | NIST 1643f        |                       |              | HPS River Sediment A |                       |              |
|-------------|---------------|-------------------|-----------------------|--------------|----------------------|-----------------------|--------------|
| Analyte     | DCS Cell Mode | Mean Conc. (µg/L) | Expected Conc. (µg/L) | Recovery (%) | Mean Conc. (µg/L)    | Expected Conc. (µg/L) | Recovery (%) |
| 7 Li        | AHM           | 15.4 ± 0.178      | 16.6                  | 93           | 0.267 ± 0.172        | -                     | -            |
| 9 Be        | AHM           | 12.9 ± 0.159      | 14                    | 94           | 0.346 ± 0.235        | -                     | -            |
| 23 Na       | AHM           | 18,900 ± 198      | 18,830                | 100          | 51,500 ± 753         | 50,000                | 103          |
| 24 Mg       | AHM           | 7,320 ± 70.9      | 7,454                 | 98           | 70,600 ± 881         | 70,000                | 101          |
| 27 Al       | AHM           | 121 ± 16.5        | 134                   | 90           | 257,000 ± 2,700      | 250,000               | 103          |
| 39 K        | AHM           | 1,840 ± 12.9      | 1,933                 | 95           | 143,000 ± 1,100      | 150,000               | 96           |
| 39 -> 39 K  | Air           | 2,000 ± 25.0      | 1,933                 | 104          | 151,000 ± 2,770      | 150,000               | 101          |
| 44 Ca       | AHM           | 29,200 ± 219      | 29,430                | 99           | 282,000 ± 3,030      | 300,000               | 94           |
| 44 -> 60 Ca | Air           | 29,100 ± 473      | 29,430                | 99           | 311,000 ± 6,680      | 300,000               | 104          |
| 51 V        | AHM           | 33.7 ± 0.274      | 36                    | 93           | 232 ± 2.04           | 250                   | 93           |
| 51 -> 67 V  | Air           | 33.7 ± 0.435      | 36                    | 93           | 226 ± 3.86           | 250                   | 90           |
| 52 Cr       | AHM           | 17.4 ± 1.41       | 19                    | 94           | 286,000 ± 2,440      | 300,000               | 95           |
| 55 Mn       | AHM           | 34.3 ± 1.70       | 37                    | 92           | 7,440 ± 72.7         | 8,000                 | 93           |
| 56 Fe       | AHM           | 80.3 ± 12.6       | 93                    | 86           | 1,180,000 ± 10,900   | 1,200,000             | 98           |
| 59 Co       | AHM           | 24.0 ± 0.259      | 25                    | 95           | 118 ± 1.61           | 100                   | 118          |
| 60 Ni       | AHM           | 55.5 ± 0.667      | 60                    | 93           | 497 ± 7.97           | 500                   | 99           |
| 63 Cu       | AHM           | 20.7 ± 0.336      | 22                    | 95           | 964 ± 14.9           | 1,000                 | 96           |
| 66 Zn       | AHM           | 70.3 ± 1.22       | 74                    | 94           | 14,100 ± 135         | 15,000                | 94           |
| 75 As       | AHM           | 52.8 ± 0.417      | 57                    | 92           | 536 ± 6.80           | 600                   | 89           |
| 75 -> 91 As | Air           | 53.5 ± 0.625      | 57                    | 93           | 543 ± 6.28           | 600                   | 91           |
| 78 Se       | AHM           | 10.8 ± 0.146      | 12                    | 92           | 18.3 ± 0.568         | 20                    | 92           |
| 78 -> 94 Se | Air           | 10.9 ± 0.182      | 12                    | 93           | 18.3 ± 1.37          | 20                    | 92           |
| 88 Sr       | AHM           | 304 ± 1.94        | 314                   | 97           | 2.77 ± 1.65          | -                     | -            |
| 95 Mo       | AHM           | 104 ± 0.532       | 115                   | 90           | 1.87 ± 0.707         | -                     | -            |
| 107 Ag      | AHM           | 0.913 ± 0.062     | 1.0                   | 94           | 0.668 ± 0.448        | -                     | -            |
| 111 Cd      | AHM           | 5.42 ± 0.081      | 6.0                   | 92           | 95.0 ± 1.37          | 100                   | 91           |
| 121 Sb      | AHM           | 50.0 ± 0.475      | 55                    | 90           | 448 ± 4.17           | 0                     | -            |
| 137 Ba      | AHM           | 477 ± 3.14        | 518                   | 92           | 457 ± 4.27           | 500                   | 91           |
| 201 Hg      | AHM           | 0.004 ± 0.001     | -                     | -            | 0.0371 ± 0.0094      | -                     | -            |
| 205 Tl      | AHM           | 6.71 ± 0.079      | -                     | -            | 9.51 ± 0.0720        | 10                    | 95           |
| *Pb         | AHM           | 18.6 ± 0.327      | 18                    | 101          | 6,950 ± 39.5         | 7,000                 | 99           |
| 232 Th      | AHM           | 0.006 ± 0.006     | -                     | -            | 18.6 ± 0.113         | 20                    | 93           |
| 238 U       | AHM           | 0.010 ± 0.008     | -                     | -            | 9.50 ± 0.0972        | 10                    | 95           |

\*Reported as the sum of the lead isotopes 206+207+208.

Table 5 continued. Recovery of certified values for CRMs. Concentration units µg/L.

| n=18        | DCS Cell Mode | HPS River Sediment B |                       |              | HPS Soil A        |                       |              | HPS Soil B        |                       |              |
|-------------|---------------|----------------------|-----------------------|--------------|-------------------|-----------------------|--------------|-------------------|-----------------------|--------------|
|             |               | Mean Conc. (µg/L)    | Expected Conc. (µg/L) | Recovery (%) | Mean Conc. (µg/L) | Expected Conc. (µg/L) | Recovery (%) | Mean Conc. (µg/L) | Expected Conc. (µg/L) | Recovery (%) |
| 7 Li        | AHM           | 1.11 ± 0.196         | -                     | -            | 0.608 ± 0.135     | -                     | -            | 0.617 ± 0.146     | -                     | -            |
| 9 Be        | AHM           | 0.270 ± 0.172        | -                     | -            | 0.304 ± 0.188     | -                     | -            | 0.256 ± 0.161     | -                     | -            |
| 23 Na       | AHM           | 50,600 ± 510         | 5,000                 | 101          | 70,500 ± 658      | 70,000                | 101          | 100,000 ± 855     | 100,000               | 100          |
| 24 Mg       | AHM           | 121,000 ± 737        | 12,000                | 101          | 70,900 ± 508      | 70,000                | 101          | 81,200 ± 506      | 80,000                | 101          |
| 27 Al       | AHM           | 633,000 ± 3,610      | 60,000                | 105          | 526,000 ± 3,360   | 500,000               | 105          | 731,000 ± 4,720   | 700,000               | 104          |
| 39 K        | AHM           | 189,000 ± 1,170      | 20,000                | 94           | 186,000 ± 1,670   | 200,000               | 93           | 195,000 ± 2,010   | 210,000               | 93           |
| 39 -> 39 K  | Air           | 199,000 ± 3,360      | 20,000                | 100          | 202,000 ± 1,810   | 200,000               | 101          | 207,000 ± 2,370   | 210,000               | 99           |
| 44 Ca       | AHM           | 286,000 ± 2,730      | 30,000                | 95           | 336,000 ± 1,750   | 350,000               | 96           | 115,000 ± 471     | 125,000               | 92           |
| 44 -> 60 Ca | Air           | 298,000 ± 5,070      | 30,000                | 99           | 353,000 ± 8,550   | 350,000               | 101          | 128,000 ± 3,110   | 125,000               | 103          |
| 51 V        | AHM           | 997 ± 6.85           | 100                   | 100          | 91.9 ± 0.710      | 100                   | 92           | 726 ± 5.37        | 800                   | 91           |
| 51 -> 67 V  | Air           | 966 ± 15.3           | 100                   | 100          | 91.6 ± 1.06       | 100                   | 92           | 720 ± 11.2        | 800                   | 90           |
| 52 Cr       | AHM           | 13,800 ± 170         | 1,500                 | 92           | 1.04 ± 0.814      | -                     | -            | 364 ± 6.78        | 400                   | 91           |
| 55 Mn       | AHM           | 5,470 ± 55.0         | 600                   | 91           | 103 ± 42.3        | 100                   | 103          | 90,100 ± 661      | 100,000               | 90           |
| 56 Fe       | AHM           | 373,000 ± 2,960      | 40,000                | 93           | 187,000 ± 1,540   | 200,000               | 93           | 326,000 ± 2,990   | 350,000               | 93           |
| 59 Co       | AHM           | 141 ± 2.22           | 15                    | 94           | 2.54 ± 0.149      | -                     | -            | 95.5 ± 1.35       | 100                   | 95           |
| 60 Ni       | AHM           | 480 ± 7.65           | 50                    | 96           | 273 ± 3.17        | 300                   | 91           | 180 ± 2.98        | 200                   | 90           |
| 63 Cu       | AHM           | 925 ± 15.4           | 100                   | 92           | 284 ± 4.34        | 300                   | 94           | 2,840 ± 31.8      | 3,000                 | 94           |
| 66 Zn       | AHM           | 4,530 ± 77.2         | 500                   | 91           | 962 ± 38.3        | 1,000                 | 96           | 63,800 ± 662      | 70,000                | 91           |
| 75 As       | AHM           | 196 ± 2.56           | 20                    | 98           | 182 ± 2.28        | 200                   | 91           | 5,300 ± 75.2      | 6,000                 | 88           |
| 75 -> 91 As | Air           | 194 ± 1.64           | 20                    | 98           | 184 ± 1.42        | 200                   | 92           | 5,420 ± 57.7      | 6,000                 | 90           |
| 78 Se       | AHM           | 10.1 ± 0.418         | 1.0                   | 101          | 9.18 ± 0.355      | 10                    | 92           | 1.49 ± 0.166      | -                     | -            |
| 78 -> 94 Se | Air           | 9.50 ± 0.583         | 1.0                   | 95           | 9.06 ± 0.557      | 10                    | 91           | 1.50 ± 0.300      | -                     | -            |
| 88 Sr       | AHM           | 3.30 ± 0.117         | -                     | -            | 3.37 ± 0.0832     | -                     | -            | 1.35 ± 0.0555     | -                     | -            |
| 95 Mo       | AHM           | 4.10 ± 0.158         | -                     | -            | 0.379 ± 0.100     | -                     | -            | 1.10 ± 0.105      | -                     | -            |
| 107 Ag      | AHM           | 0.300 ± 0.115        | -                     | -            | 0.551 ± 0.225     | -                     | -            | 0.266 ± 0.156     | -                     | -            |
| 111 Cd      | AHM           | 28.4 ± 0.323         | 3.0                   | 91           | 3.17 ± 0.0855     | 3                     | 96           | 191 ± 2.61        | 200                   | 92           |
| 121 Sb      | AHM           | 36.7 ± 0.377         | 4.0                   | 92           | 29.5 ± 0.209      | 30                    | 98           | 362 ± 3.02        | 400                   | 91           |
| 137 Ba      | AHM           | 3,690 ± 27.8         | 400                   | 92           | 4,630 ± 37.2      | 5,000                 | 93           | 6,440 ± 52.1      | 7,000                 | 92           |
| 201 Hg      | AHM           | 0.033 ± 0.008        | -                     | -            | 0.0280 ± 0.00870  | -                     | -            | 0.0470 ± 0.0148   | -                     | -            |
| 205 Tl      | AHM           | 9.50 ± 0.103         | 1.0                   | 95           | 0.0640 ± 0.0204   | -                     | -            | 0.491 ± 0.0432    | -                     | -            |
| *Pb         | AHM           | 1,930 ± 14.9         | 200                   | 97           | 389 ± 3.65        | 400                   | 97           | 58,600 ± 269      | 60,000                | 98           |
| 232 Th      | AHM           | 92.8 ± 0.636         | 10                    | 93           | 94.0 ± 0.666      | 100                   | 94           | 95.3 ± 0.543      | 100                   | 95           |
| 238 U       | AHM           | 27.7 ± 0.231         | 3.0                   | 92           | 9.43 ± 0.124      | 10                    | 94           | 235 ± 3.33        | 250                   | 94           |

\*Reported as the sum of the lead isotopes 206+207+208.

## Matrix spike recoveries

Matrix Spike (MS), Matrix Spike Duplicate (MSD), and relative percent difference (RPD) results for 10x diluted Soil A are given in Table 6. Soil A contains Na, Mg, Al, K, Ca, Fe, and Ba at a higher concentration than the 100 µg/L spike amount, so no spike recoveries are given for these elements. However, the CRM recovery results for Soil A in Table 5 show that all these major elements can be analyzed with good accuracy. Recoveries for the other elements were generally within ±10% of the spike level, and the RPD between MS and MSD was less than 3.5%, demonstrating the method's accuracy.

**Table 6.** Matrix Spike (MS)/Matrix Spike Duplicate (MSD) spike concentration recoveries and relative percent difference (RPD) for 10x diluted Soil A CRM.

| Analyte     | DCS Mode | Soil A Measured Conc (µg/L) | Soil A MS Measured Conc (µg/L) | Soil A MSD Measured Conc (µg/L) | Soil A MSD Recovery (%) | RPD MS/MSD (%) |
|-------------|----------|-----------------------------|--------------------------------|---------------------------------|-------------------------|----------------|
| 7 Li        | AHM      | 0.061                       | 9.650                          | 9.710                           | 96                      | 0.6            |
| 9 Be        | AHM      | 0.030                       | 18.900                         | 19.000                          | 95                      | 0.6            |
| 51 V        | AHM      | 9.19                        | 28.8                           | 28.7                            | 98                      | 0.1            |
| 51 -> 67 V  | Air      | 9.16                        | 29.1                           | 29.2                            | 100                     | 0.6            |
| 52 Cr       | AHM      | 0.104                       | 19.3                           | 19.3                            | 96                      | 0.3            |
| 55 Mn       | AHM      | 13.3                        | 9270                           | 9440                            | 94                      | 1.8            |
| 59 Co       | AHM      | 0.254                       | 19.4                           | 19.3                            | 95                      | 0.9            |
| 60 Ni       | AHM      | 27.3                        | 47.8                           | 47.5                            | 101                     | 0.6            |
| 63 Cu       | AHM      | 28.4                        | 49.7                           | 49.6                            | 106                     | 0.3            |
| 66 Zn       | AHM      | 96.2                        | 9590                           | 9560                            | 95                      | 0.3            |
| 75 As       | AHM      | 18.2                        | 38.6                           | 38.3                            | 100                     | 0.9            |
| 75 -> 91 As | Air      | 18.4                        | 38.7                           | 38.5                            | 100                     | 0.4            |
| 78 Se       | AHM      | 0.918                       | 20.0                           | 19.8                            | 94                      | 1.1            |
| 78 -> 94 Se | Air      | 0.906                       | 20.1                           | 19.5                            | 93                      | 3.1            |
| 88 Sr       | AHM      | 0.337                       | 10.1                           | 9.92                            | 96                      | 2.2            |
| 95 Mo       | AHM      | 0.038                       | 18.7                           | 18.6                            | 93                      | 0.8            |
| 107 Ag      | AHM      | 0.055                       | 19.1                           | 19.0                            | 95                      | 0.3            |
| 111 Cd      | AHM      | 0.289                       | 19.5                           | 19.5                            | 96                      | 0.1            |
| 121 Sb      | AHM      | 2.95                        | 22.3                           | 22.4                            | 97                      | 0.6            |
| 201 Hg      | AHM      | 0.003                       | 0.483                          | 0.481                           | 96                      | 0.3            |
| 205 Tl      | AHM      | 0.006                       | 19.8                           | 19.8                            | 99                      | 0.1            |
| *Pb         | AHM      | 38.9                        | 60.5                           | 60.7                            | 109                     | 0.4            |
| 232 Th      | AHM      | 9.40                        | 29.4                           | 29.9                            | 103                     | 1.7            |
| 238 U       | AHM      | 0.943                       | 20.3                           | 20.4                            | 97                      | 0.5            |

\*Reported as the sum of the lead isotopes 206+207+208.

## Measurement of Ca, As, and Se in an M<sup>2+</sup> matrix

If present at high enough concentrations, elements with low second ionization energies can form M<sup>2+</sup> ions that can interfere with analytes at half their mass (for example, <sup>88</sup>Sr<sup>2+</sup> on <sup>44</sup>Ca<sup>+</sup>, <sup>150</sup>Nd<sup>2+</sup> and <sup>150</sup>Sm<sup>2+</sup> on <sup>75</sup>As<sup>+</sup>, <sup>156</sup>Gd<sup>2+</sup> and <sup>156</sup>Dy<sup>2+</sup> on <sup>78</sup>Se<sup>+</sup>). However, operating the DCS in Air cell mode minimizes the M<sup>2+</sup> background on the analyte ions via an O<sub>2</sub> mass-shift reaction.<sup>4</sup> Table 7 shows the recovery results for Ca, As, and Se with no added M<sup>2+</sup> interferences and with added interferences. High recoveries for all three analytes were observed when operating the DCS in AHM due to M<sup>2+</sup> interferences, whereas excellent recoveries of 100 ± 10% were achieved in Air cell mode.

**Table 7.** Recoveries of Ca, As, and Se in 10x diluted Soil A with and without M<sup>2+</sup> interferences

| Analyte [DCS Mode] | Soil A Certified Value (µg/L) | Soil A (No Interferent) Recovery (%) | M <sup>2+</sup> Interferent | M <sup>2+</sup> Interferent Conc (µg/L) | Soil A (With Interferent) Recovery (%) |
|--------------------|-------------------------------|--------------------------------------|-----------------------------|---|--|
| 44 Ca [AHM]        | 35,000                        | 94                                   | Sr                          | 10,000                                  | 198                                    |
| 44 -> 60 Ca [Air]  |                               | 98                                   |                             |   | 107                                    |
| 75 As [AHM]        | 20                            | 95                                   | Nd & Sm                     | 400 (each)                              | 132                                    |
| 75 -> 91 As [Air]  |                               | 94                                   |                             |   | 94                                     |
| 78 Se [AHM]        | 1.0                           | 92                                   | Gd & Dy                     | 20 (each)                               | 288                                    |
| 78 -> 94 Se [Air]  |                               | 105                                  |                             |   | 100                                    |

## Conclusion

The study has shown that the Agilent 9500 ICP-QQQ, equipped with the unique Dual-Cell System (DCS) and the optional AVS MS discrete sampling system, is ideal for high-throughput, multi-element analysis of soil and sediment samples.

The use of the EPA 6020B Preset Method and selection of UHMI-4 plasma conditions in the Agilent OpenLab ICP-MS software enabled rapid setup of the instrument's operating parameters, ensuring consistent performance across operators.

The instrument's matrix-tolerant plasma (low CeO/Ce ratio), UHMI, and DCS effectively addressed the challenges of analyzing high-matrix samples, including signal drift, ion suppression, and spectral interferences, thereby delivering an optimized analytical workflow.

The method's productivity was enhanced by operating the DCS in Advanced Helium Mode (AHM) and Air cell mode and by using short integration times, saving 5 to 10 s per sample. AHM replaces conventional no gas, helium (He)

collision, and high energy He (HEHe) in a single mode, while Air cell mode uses oxygen from the laboratory, alleviating the need for gas cylinders. The 9500 maintains exceptionally high sensitivity in AHM and Air cell mode, minimizing background signals.

Overall, the 9500 ICP-QQQ with DCS achieved excellent analytical data:

- Low detection limits in the low ng/L (ppt) range for most analytes.
- Removal of  $M^{2+}$  interferences on Ca, As, and Se generated by Sr and rare earth elements (REEs) in the soil samples using Air cell mode.
- Recoveries of all certified elements in the five environmental CRMs or matrix spike samples within  $100 \pm 10\%$ , confirming the effective control of interferences using AHM and Air cell mode.
- Method robustness and reproducibility over more than eight hours of continuous measurements, as demonstrated by the recoveries of ISTDs and CCVs within the limits specified in EPA Method 6020.

The 9500 ICP-QQQ with DCS and an integrated AVS MS enables high-throughput, routine analysis of complex environmental matrices. The EPA Method 6020-compliant method provided accurate data by effectively handling polyatomic ion and  $M^{2+}$  interferences using AHM and Air cell mode.

## References

1. U.S. EPA. 2014. "Method 6020B (SW-846): Inductively Coupled Plasma-Mass Spectrometry," Revision 2. Washington, DC, <https://www.epa.gov/esam/epa-method-6020b-sw-846-inductively-coupled-plasma-mass-spectrometry> (accessed April 2026)
2. Dual-Cell System (DCS) and Advanced Helium Mode (AHM), Agilent publication, [5994-8985EN](#)
3. Air Cell Mode of the Agilent 9500 ICP-QQQ with Dual-Cell System, Agilent publication, [5994-8987EN](#)
4. Sugiyama, N. Solving Doubly Charged Ion Interferences using an Agilent 8900 ICP-QQQ, Agilent publication, [5994-1155EN](#)

## Products used in this application

### Agilent products

| Product Type               | Description  | Part Number                 |
|----------------------------|--|-----------------------------|
| Sample Introduction System | Quartz sample introduction system for 9500 ICP-MS  | <a href="#">M5150-67107</a> |
|                            | MicroMist nebulizer for 9500 ICP-MS  | <a href="#">M5150-67024</a> |
| Interface                  | ICP-MS sampler cone for 9500 ICP-MS, Ni tip with Cu base   | <a href="#">M5150-67000</a> |
|                            | Nickel skimmer cone for 9500 ICP-MS with u-lens  | <a href="#">M5150-67005</a> |
|                            | Extraction-Omega lens assembly, u-lens, stainless steel base   | <a href="#">M5150-67022</a> |
| Tubing Kits                | Easy-fit peristaltic-pump tubing, PVC, white/white, 1.02 mm id, for sample                                     | <a href="#">5005-0020</a>   |
|                            | Easy-fit peristaltic-pump tubing, PVC, blue/orange, 0.25 mm id, for internal standard                          | <a href="#">5005-0021</a>   |
|                            | Easy-fit peristaltic-pump tubing, beige thermoplastic, yellow/blue, 1.52 mm id, for drain                      | <a href="#">5005-0022</a>   |
|                            | Sample loop for ADS 2/AVS MS, 1.50 mL 1.00 mm id   | <a href="#">5005-0425</a>   |
|                            | AVS MS preconfigured tubing kit for 9500 ICP-MS  | <a href="#">M5171-67001</a> |
| Bottle Kits                | Diluent/carrier 6 L bottle kit, includes a 6 L can, GL45 StaySafe cap, fittings, and venting valve             | <a href="#">5005-0435</a>   |
|                            | Diluent 2 L PFA bottle kit for ICP-MS, includes 2 L PFA bottle, GL45 StaySafe cap, fittings, and venting valve | <a href="#">5005-0436</a>   |
|                            | Waste container kit, includes a 10 L waste can, S60 StaySafe cap, fittings, and acid vapor filter              | <a href="#">5005-0437</a>   |
| Chemical Standards         | Initial calibration verification standard  | <a href="#">5183-4682</a>   |
|                            | Environmental calibration standard   | <a href="#">5183-4688</a>   |

[www.agilent.com/chem/9500icpqqq](http://www.agilent.com/chem/9500icpqqq)

DE-013879

This information is subject to change without notice.

© Agilent Technologies, Inc. 2026  
Published in the USA, June 01, 2026  
5994-9128EN