

Analysis of Permanent Gases More Challenging Than You Might Think

Mark Sinnott
May 17, 2023
Application Engineer

Analysis of Permanent Gases

- Reduce the list to only mission critical solutes
- Sample introduction – syringes and valves
- General discussion on PLOT columns
 - The molesieve column is at the heart of permanent gas separations
- Techniques when $\text{CO}_2 + \text{C}_2, \text{C}_3 +$ is also needed
 - Column isolation
 - Parallel columns
 - Cryogenic separations
 - Unique selectivity packed columns
- Techniques for low level detection of hydrogen
- Micro GC
- Ammonia
- MSD?



Reduce the List First!

- Reduce your list to only mission critical solutes
 - This is a critical step and can greatly simplify things
 - Eliminating even one solute can help (i.e. do I really need CO₂?)



Gas Tight Syringes



Syringes for Gas Injection

Point Style

Description

Application

<http://pconlab.net/Hamilton-Syringes.html>

2



10–12° sharp, beveled, curved non-coring

Gas chromatography, septum piercing

3



Blunt, electro-polished

High performance liquid chromatography (HPLC) injection, thin layer chromatography (TLC), general liquid handling, controlled animal injections

3T



Blunt, electro-polished, coated with PTFE 19 mm from the tip

Thin layer chromatography (TLC) applications

4



Sharp 10–12° beveled needle

Life science/animal injections

5



Conical with side port for penetration without coring

Headspace, applications prone to needle clogging, causes minimal septum damage


AS



Conical, non-coring designed to withstand multiple injections

Autosampler injection, pre-pierced septa

Gas Sampling Valve



Gas sampling valve

PLOT Columns

Column Type	Phase Type	Chromatographic Process	Stationary Phases
WCOT	Liquid	Gas – Liquid partitioning	Polysiloxanes PEG
PLOT	Solid	Gas – Solid adsorption	Porous polymers Al ₂ O ₃ , zeolites

WCOT = Wall Coated Open Tubular


PLOT = Porous Layer Open Tubular

Why Use a PLOT Column?


- Highly retentive
 - Can resolve gases at non-cryo temperatures
- Unique selectivity
- Permits higher initial oven temperatures than WCOT



PLOT Columns



WCOT Ethylene Analysis



ETHYLENE
1. Methane
2. Ethylene
3. Ethane
4. Propylene

Column: DB-Petro 100
100 m x 0.25 mm I.D. , 0.5 μm
Carrier: Helium at 28.5 cm/s
Oven: -40°C for 20 min
-40°C - 120°C at 5°/min
Injector: Split 1:16, 100 μL, 200°C
Detector: FID, 250°C

PLOT Ethylene Analysis

Ethylene

Column: GS-Alumina

50 m x 0.53 mm I.D.

J&W P/N: 115-3552

Carrier: Helium at 11 mL/min, measured at 35°C

Oven: 35°C for 2 min

35-190°C at 6°/min

190°C for 3 min


Injector: Split 1:30, 200°C

0.2 mL of trace hydrocarbons in ethylene

Detector: FID, 200°C

Nitrogen makeup gas at 20 mL/min

1. Methane
2. Ethane
3. Ethylene
4. Propane
5. Cyclopropane
6. Propylene
7. Isobutane
8. Acetylene
9. *n*-Butane
10. Propadiene
11. *trans*-2-Butene
12. 1-Butene
13. Isobutylene
14. *cis*-2-Butene
15. Isopentane
16. *n*-Pentane
17. Propyne
18. 1,3-Butadiene
19. 1-Pentene
20. *n*-Hexane



C585

JW Column PLOT column Portfolio- DB, HP, CP

HP-PLOT Al₂O₃

GS-Alumina

GS-Alumina KCl

HP-PLOT Al₂O₃ S

HP-PLOT Al₂O₃ M

CP- Al₂O₃/KCl

CP- Al₂O₃/Na₂SO₄

CP-SilcaPLOT

CP-CarboBOND

CP-CarboPLOT P7

GS-CarbonPLOT

CP-PoraPLOT Amines

CP-PoraPLOT S

CP-PoraBOND Q

CP-PoraPLOT Q

CP-PoraPLOT Q-HT

HP-PLOT Q

CP-PoraBOND U

HP-PLOT U

GS-Q

CP-PoraPLOT U

HP-PLOT molesieve

CP-molesieve 5A


GS-GasPro

ShinCarbon ST

Select Permanent Gas Column

But I Also Need to Analyze for CO₂...It's Getting Complicated

- Recall that CO₂ cannot be done on a molesieve column (absorbed)
 - Alternative options:
 - Two unique injections on two columns
 - Molesieve and PLOT-Q
 - Be aware how retentive the molesieve is!
 - One injection on one column at Cryo temps (-80 °C)
 - GasPro
 - One injection on one column – ShinCarbon ST (packed column)
 - One injection with two columns and a valve to “direct traffic”
 - Column isolation – molesieve and PLOT-Q
 - One injection with two columns in parallel
 - Select permanent gas column – molesieve and PLOT-Q
- ...did I mention this would get complicated?



Molesieve Column

- Separation based mostly on molecular size and shape
- Excellent at what it does but very limited
 - H_2 , O_2 , N_2 , CH_4 , CO
 - Noble gases
 - NO_2
- Limitations...cannot use for...
 - CO_2
 - Water
 - C_2 HCs and larger


Molesieve Separations

Conditions

Technique : GC-wide-bore
Column : Agilent CP-Molsieve 5Å, 0.53 mm x 50 m fused silica
PLOT CP-Molsieve 5Å (df = 50 µm)
(Part no. CP7539)
Temperature : 30 °C
Carrier Gas : H₂, 55 kPa (0.55 bar, 8 psi), 25.6 cm/s
Injector : Splitter, 200 mL/min
T = 30 °C
Detector : TCD
T = 180 °C
Concentration Range : 5 - 10%

Peak identification

1. helium
2. neon
3. argon
4. oxygen
5. nitrogen
6. methane




GC: 6890
Sample: 250 ml, split (75:1)
Carrier: Helium, 2 µl/min
Column: HP-PLOT/MoleSieve, 30 m x 0.32 mm x 12 µm
(Part No. 19091P-MS4)
Oven: 40°C Isothermal
Detector: TCD

Pub No.: 5964-2129E

Permanent Gases

1. Neon
2. Argon
3. Oxygen
4. Nitrogen
5. Methane
6. Carbon monoxide



Molesieve Separations

Separation of hydrogen and helium


Technique : GC
Column : CP-Molsieve 5A fused silica
50 m x 0.53 mm df = 50 µm, Part nr. CP7539

Temperature : 40°C
Carrier Gas : Nitrogen, 50kPa (7.2 psi)
Injector : Splitter, 40 ml/min
Detector : µ-TCD, 200°C
Sample Size : 40 µl
Concentration range : 1% in nitrogen

Reference : C. Duvekot, Varian R&D laboratories, Middelburg, The Netherlands



Peak identification

- 1 Helium
- 2 Hydrogen



Noble Gases

1. Neon
2. Argon
3. Oxygen
4. Nitrogen
5. Krypton
6. Xenon



GC: 6890
Sample: 250 µl, split (50:1)
Carrier: Helium, 4 ml/min
Column: HP-PLOT/MoleSieve, 30 m x 0.53 mm x 50 µm
(Part No. 19095P-MS0)
Oven: Temperature program listed above
Detector: TCD


Pub No.: 5964-2129E

Tips For Using a Molesieve Column



- Molesieve material is commonly used as a moisture trap
- Samples must be reasonably dry
- Separation will change as column absorbs water over time
 - Resolution loss between CH₄ and CO indicates that the column should be reconditioned (300 °C 8 hours+)

Permanent Gases

1. Neon
2. Argon
3. Oxygen
4. Nitrogen
5. Methane
6. Carbon monoxide




PLOT-Q




- | | |
|---------------------|-------------------|
| 1. CO/air | 7. Propylene |
| 2. Methane | 8. Propane |
| 3. Carbon dioxide | 9. Dimethyl ether |
| 4. Ethylene | 10. Methanol |
| 5. Ethane | 11. Butylene |
| 6. Hydrogen sulfide | 12. Butane |

Permanent Gases in Hydrocarbon Blend





Column: GS-GasPro, 60 m x 0.32 mm I.D.
J&W P/N: 113-4362
Carrier: Helium, 20 psig (constant pressure)
Inlet: 200°C, 1:15 split
Sample: 100 µL, manual inject
Detector: TCD, 250°C, background compensated
Oven: -80°C for 3 min, -80° to 260°C at 25°/min, hold at 260°C



ShinCarbon ST






Valves



One Injection, Two Columns and Valve...Column Isolation




System must be flow balanced
Use flow meter G6691A




- 1. Carbon Dioxide
- 2. Ethylene
- 3. Ethane
- 4. Hydrogen
- 5. Argon
- 6. Oxygen
- 7. Nitrogen
- 8. Methane
- 9. Carbon monoxide
- * Signal of valve switching

Column Isolation




- 1. Carbon Dioxide
- 2. Ethylene
- 3. Ethane
- 4. Hydrogen
- 5. Argon
- 6. Oxygen
- 7. Nitrogen
- 8. Methane
- 9. Carbon monoxide
- * Signal of valve switching

Column Isolation: Setting Up Valve Timing



- | | |
|---------------------|-------------------|
| 1. CO/air | 7. Propylene |
| 2. Methane | 8. Propane |
| 3. Carbon dioxide | 9. Dimethyl ether |
| 4. Ethylene | 10. Methanol |
| 5. Ethane | 11. Butylene |
| 6. Hydrogen sulfide | 12. Butane |



Column Isolation – Flexibility of Elution Order

Carrier gas: Helium

Inlet: Purged packed 55 °C

Sample loop: 0.1 mL

Oven: 50 °C (10 min) with 10 °C/min
to 120 °C (5 min)

TCD: 180 °C

1 H₂: 50%

2 CO₂: 10%


3 O₂: 172 ppm

4 N₂: Balance


5 CH₄: 5%

6 CO: 50 ppm

* Signal of valve switching



One Injection, Two Columns in Parallel – Select Permanent Gas Column




Set-up for Select Permanent Gases/CO2

Select Permanent Gas Column

Technique : GC
Column : Select Permanent Gases/CO2
Part No. CP7429
Temperature : 50°C
Carrier Gas : Helium, 100 kPa
Injector : Split 50 ml/min
Detector : TCD
Sample Size : 20 µl
Concentration range : % level

Peak identification:


1. Oxygen
2. Nitrogen
3. Methane
4. Carbon monoxide
5. Carbon dioxide



Select Permanent Gas Column

Peak Identification

1. Oxygen
2. Nitrogen
3. Methane
4. Oxygen + Nitrogen
5. Methane
6. CO₂
7. Ethane
8. Propane
9. Ethane
10. Isobutane
11. Butane



Technique : GC
Column : Fused Silica, Select Permanent Gases/CO₂
Part No. CP7429
Temperature : 50°C, 4 min, 250°C, 15°C/min
Carrier Gas : Helium, 100 kPa
Injector : Split 50 ml/min
Detector : μ -TCD
Sample Size : 10 μ l natural gas
Concentration range : % level

AN 2029 - GC

So, You Think That Was Complicated, What if I Need to Analyze for Hydrogen Too...

- The trouble with hydrogen
 - Must use TCD**
 - Sensitivity is carrier gas dependent
 - Creates unique issues related sensitivity
 - Recall the negative hydrogen peak from earlier



***Can also use HID/DID, but we will not be covering this here*

Thermal Conductivity Detector (TCD)

- When using He carrier sensitivity for hydrogen is on the order of only ~10%
- No problem, I'll use argon or nitrogen carrier and get down to ppm levels for hydrogen...
- This is true; however, the sensitivity for all the other compounds is now very poor...a real "catch 22"...
- Let's see why this is the case...



Thermal Conductivity Detector – TCD


Name	Thermal Conductivity at 400 K (mW m ⁻¹ K ⁻¹)
Argon	22.4
Hydrogen	230.9
Helium	189.6
Nitrogen	32.8
Carbon monoxide	32.3
Carbon dioxide	25.2
Acetylene	33.3
Ethylene	34.7
Ethane	36.0
Propane	31.0
Butane	28.3
Pentane	24.9
Hexane	23.4

		Thermal conductivity in mW m ⁻¹ K ⁻¹						Ref.
		100 K	200 K	300 K	400 K	500 K	600 K	
	Air	9.5	18.5	26.4	33.5	39.9	46.0	1
Ar	Argon (<i>P</i> = 0)	6.3	12.4	17.7	22.4	26.5	30.3	2, 3*
BF ₃	Boron trifluoride			19.0	24.6			4
ClH	Hydrogen chloride		9.2	14.5	19.5	24.0	28.1	4
F ₂ S	Sulfur hexafluoride (<i>P</i> = 0)			13.0	20.6	27.5	33.8	5
H ₂	Normal hydrogen (<i>P</i> = 0)	68.2	132.8	186.6	230.9	270.9	309.1	6
H ₂ O	Water (<i>P</i> = 0)			18.6	26.1	35.6	46.2	7
D ₂ O	Deuterium oxide (<i>P</i> = 0)			18.2	26.6	36.3	47.6	8
H ₂ S	Hydrogen sulfide			14.6	20.5	26.4	32.4	4
H ₃ N	Ammonia			25.1	37.2	53.1	68.6	9
He	Helium (<i>P</i> = 0)	74.7	118.3	155.7	189.6	221.4	251.6	10
Kr	Krypton (<i>P</i> = 0)		6.5	9.5	12.3	14.8	17.1	11
NO	Nitric oxide		17.8	25.9	33.1	39.6	46.2	4
N ₂	Nitrogen	9.4	18.3	26.0	32.8	39.0	44.8	1
N ₂ O	Nitrous oxide		9.8	17.4	26.0	34.1	41.8	4
Ne	Neon (<i>P</i> = 0)	22.3	37.4	49.4	59.9	69.5	78.5	12
O ₂	Oxygen	9.1	18.2	26.5	34.0	41.0	47.7	1
O ₂ S	Sulfur dioxide			9.6	14.3	20.0	25.6	4
Xe	Xenon (<i>P</i> = 0)		3.7	5.5	7.2	8.8	10.3	3*, 11
CCl ₂ F ₂	Dichlorodifluoromethane			9.9	15.0	20.1	25.2	13
CF ₄	Tetrafluoromethane (<i>P</i> = 0)			16.0	24.1	32.2	39.9	5
CO	Carbon monoxide (<i>P</i> = 0)			25.0	32.3	39.2	45.7	14
CO ₂	Carbon dioxide		9.6	16.8	25.2	33.5	41.6	15
CHCl ₃	Trichloromethane			7.5	11.1	15.1		4
CH ₄	Methane (<i>P</i> = 0)	10.4	21.8	34.4	50.0	68.4	88.6	16
CH ₃ O	Methanol				26.2	38.6	53.0	4
C ₂ Cl ₂ F ₄	1,2-Dichloro-1,1,2,2-tetrafluoroethane			10.3	15.7	21.1		13
C ₂ Cl ₂ F ₃	1,1,2-Trichloro-1,2,2-trifluoroethane			9.0	13.6	18.3		13
C ₂ H ₂	Acetylene			21.4	33.3	45.4	56.8	4
C ₂ H ₄	Ethylene		11.3	20.6	34.7	49.9	68.6	17
C ₂ H ₆	Ethane		10.7	21.2	36.0	53.8	73.3	18
C ₂ H ₅ O	Ethanol			14.4	25.8	38.4	53.2	4
C ₃ H ₆ O	Acetone			11.5	20.2	30.6	42.7	4
C ₃ H ₈	Propane			18.5	31.0	46.4	64.6	19
C ₄ F ₈	Perfluorocyclobutane			12.5	19.5			13
C ₄ H ₁₀	Butane			16.7	28.3	43.0	60.9	20
C ₄ H ₁₀	Isobutane			17.1	28.9	43.2	60.2	21
C ₄ H ₁₀ O	Diethyl ether			15.1	25.0	37.1		4
C ₅ H ₁₂	Pentane				24.9	37.8	52.7	4
C ₆ H ₁₄	Hexane				23.4	35.4	48.7	4

* More accurate data covering a restricted temperature range.

https://ws680.nist.gov/publication/get_pdf.cfm?pub_id=907540

ShinCarbon ST



Workarounds for Low Level Hydrogen Detection


- Two unique injections:
 - One for hydrogen using Argon or N₂ carrier second injection for all other gases using helium carrier
- One injection using column isolation and argon carrier + methanizer
 - TCD coupled to FID
 - Methanizer placed after TCD but before FID to convert CO/CO₂ to methane for enhanced detection of these gases by FID
- Sabatier Reaction (1897):
$$\text{CO}_2 + 4 \text{H}_2 \xrightarrow[\text{Ni}]{400 \text{ }^\circ\text{C}} \text{CH}_4 + 2\text{H}_2\text{O}$$
- Transformer Oil Gas Analyzer (TOGA)

Conversion of CO/CO2 to Methane


- Discovered by French chemist Paul Sabatier in 1897 (Sabatier Reaction)




Transformer Oil Gas Analyzer (TOGA) Chromatograms



Transformer Oil Gas Analyzer (TOGA) Chromatograms



Transformer Oil Gas Analyzer (TOGA) Chromatograms



Micro GC

- Up to four separate channels
- 14 column chemistries available
- Gas only injections
- Single injection simultaneously distributed to each channel
- TCD detection (helium carrier)**
 - 0.5 ppm WCOT
 - 2 ppm PLOT
 - 10 ppm micropacked
 - 10^5 dynamic range
- Field portable option



***Best case scenario, some solutes will not meet these. All TCD detection levels are carrier gas dependent*





Micro GC



Single injection

4

s
i
g
n
a
l
s



Peak identification

1. hydrogen
2. oxygen
3. nitrogen
4. methane
5. carbon monoxide
6. carbon dioxide
7. ethylene
8. ethane
9. acetylene
10. hydrogen sulfide
11. propane
12. propylene
13. iso-butane
14. normal-butane
15. trans-2-butene
16. 1-butene
17. iso-butene
18. cis-2-butene
19. neo-pentane
20. iso-pentane
21. normal-pentane
22. 1,3-butadiene
23. hexane
24. iso-octane

Micro GC



Single injection

4
s
i
g
n
a
l
s

4 channels

molesieve – Helium Carrier
Low sensitivity H₂ (+10%)
High sensitivity all others (~ppm)

molesieve – Argon or N₂ Carrier
High sensitivity H₂ (~ppm)
Low sensitivity all others (~1-5%)


Plot-Q – Helium Carrier
CO₂ + C₂+ HC's (~ppm)

CP-Sil 5 – Helium Carrier
~C₅+ HC's (~ppm)

Ammonia


Technique	GC capillary
Column	Agilent CP-Volamine, 0.32 mm x 60 m fused silica
Temperature	WCOT (p/n CP7448)
Carrier Gas	40 °C (10 min) → 250 °C, 20 °C/min : He, 100 kPa (1 bar, 14 psi)
Injector	Split, 1:15 T = 180 °C
Detector	FID/TCD
Sample Size	T = 250 °C
Solvent Sample	1 µL, liquid Bulk DMA
Courtesy	Dr. F. de Boever, UCB research center Drogenbos, Dr. G. Baele, UCB Gent Belgium

1. Ammonia (NH³) 700 ppm
2. Water 2000 ppm
3. Methylamine (MMA) 585 ppm
4. Methanol 575 ppm
5. Dimethylamine (DMA) bulk
6. Ethylamine 630 ppm
7. Trimethylamine (TMA) 980 ppm



Ammonia is very reactive!
Lowest practical DL ~100 ppm (TCD)
Will depend on sample handling

Ammonia – HP-PLOT U




Volatile Primary Amines

Column: HP-PLOT U
0.53 mm x 30 m I.D., 20 µm.
P/N: 19095P-U04
Oven: 150°C for 1 min
150-190°C at 10°C/min
190°C for 5 min
Carrier: Hydrogen, 7mL/min
Inlet: 180°C, Split ratio 30:1
Detector: TCD 200°C
Sample: 2 µL amines in methanol, 1-10%.

- 1. Air
- 2. Ammonia
- 3. Water
- 4. Methanol
- 5. Isopropylamine
- 6. Triethylamine
- 7. Propylamine
- 8. Butylamine

Ammonia – Select Low Ammonia NCD



NCD=Nitrogen
Chemiluminescence detector

Very short column results in...

Very short column residence
time and therefore...

Less reactivity and an improved
response

Precision: ammonia 50 ppbv.

Permanent gases on MSD – Generally, not suitable



- Molecules are simply too small for suitable sensitivity (poor signal-to-noise)
- MSD can be modified for smaller mass range, but requires extensive hardware changes
- Most PLOT columns common to permanent gases generate particles which should be avoided

Summary

Analytes	Column	Technique
H ₂ /O ₂ /N ₂ /CH ₄ /CO	molesieve	<ul style="list-style-type: none"> • One injection one column • Micro GC
H ₂ /O ₂ /N ₂ /CH ₄ /CO + Ar	molesieve (thick film)	<ul style="list-style-type: none"> • One injection one column • Micro GC
H ₂ /O ₂ /N ₂ /CH ₄ /CO + CO ₂ + C2s	molesieve + PLOT-Q	<ul style="list-style-type: none"> • Two injections two columns • One injection and valve (column isolation) • One injection parallel columns (Select perm gas column) • Micro GC
H ₂ /O ₂ /N ₂ /CH ₄ /CO + CO ₂ + C2s	GasPro	<ul style="list-style-type: none"> • Cryo (–80 °C)
H ₂ /O ₂ /N ₂ /CH ₄ /CO + CO ₂ + C2s	ShinCarbon	<ul style="list-style-type: none"> • Packed only
H ₂ /O ₂ /N ₂ /CH ₄ /CO + CO ₂ + C2s + Low level hydrogen	molesieve + PLOT-Q	<ul style="list-style-type: none"> • Single injection onto molesieve for H₂ detection only • Argon carrier + methanizer (TOGA) • Micro GC

Contact Agilent Chemistries and Supplies Technical Support



1-800-227-9770 Option 3, Option 3:
Option 1 for GC and GC/MS columns and supplies
Option 2 for LC and LC/MS columns and supplies
Option 3 for sample preparation, filtration, and QuEChERS
Option 4 for spectroscopy supplies
Option 5 for chemical standards
Available in the USA and Canada 8–5, all time zones



gc-column-support@agilent.com
lc-column-support@agilent.com
spp-support@agilent.com
spectro-supplies-support@agilent.com
chem-standards-support@agilent.com



Agilent Technologies

Valve System

Date :

S/N :

Order No :

Item :

FLOW SOURCE

OVEN LEFT SIDE

INJECTION PORT

VALVE COMPARTMENT


BLOCK 2

BLOCK 1

DETECTOR

○ = #872

1/702/3x872



2/701/504/872

Options

SP1 options

Remarks

210,220,307,309,504,701,702,752,4x872

7890-0261

