

PLOT Columns

Separation Solutions for Light Hydrocarbons & Gases


Porous Layer Open Tubulars

Johan Kuipers
Training & Development
Sept 2021




Table of Boiling Point Fractions

Carbon #	C ₁	C ₂	C ₃	C ₄	C ₅	C ₆	C ₇	C ₈	C ₉	C ₁₀	C ₁₁	C ₁₂	C ₁₃	C ₁₄	C ₁₅	C ₁₆
Bpt of n-Paraffin @ 760 mm Hg																
Centigrade	-161	-89	-42	-0.5	+36	69	98	126	151	174	196	216	235	253	270	287
Fahrenheit	-259	-127	-44	+31	97	156	209	258	303	345	384	421	421	488	519	548




Wall Coated Open Tubular (WCOT) Columns



$K_c \Rightarrow$ Large


$K_c \Rightarrow$ Small

Ethylene Analysis on non-polar WCOT



- Column Nth = 350000
- Oven start @ -40°C
- Poor separation of C2s
- Poor separation of C3s

Ethylene Analysis on Alumina PLOT



Column: GS-Alumina
50 m x 0.53 mm I.D.

J&W P/N: 115-3552

Carrier: Helium at 11 mL/min, measured at 35°C

Oven: 35°C for 2 min

35-190°C at 6°/min

190°C for 3 min

Injector: Split 1:30, 200°C

0.2 mL of trace hydrocarbons in ethylene

Detector: FID, 200°C

Nitrogen makeup gas at 20 mL/min


1. Methane
2. Ethane
3. Ethylene
4. Propane
5. Cyclopropane
6. Propylene
7. Isobutane
8. Acetylene
9. *n*-Butane
10. Propadiene
11. *trans*-2-Butene
12. 1-Butene
13. Isobutylene
14. *cis*-2-Butene
15. Isopentane
16. *n*-Pentane
17. Propyne
18. 1,3-Butadiene
19. 1-Pentene
20. *n*-Hexane

What is a PLOT column?

- Porous layer (surface) on the inner wall capillary
- Deposition of porous particles on the wall from suspension
- Porous layer serves as stationary phase
- Gas-Solid Chromatography - PLOT
Differences in analyte distribution between carrier gas & surface of the adsorbent (stationary phase)
- Gas-Liquid Chromatography - WCOT
Differences solubility in liquid phase (stationary phase)

Advantages of GSC:

- Improved selectivity (α)
- Retention (k) of highly volatile solutes



Shape / Size





Surface

Zeolites

Porous Polymers

Al₂O₃

**Bonded Carbon
Molecular Sieves**



Bonded Silica



*Molecular sieve by Bearfruitidea from the Noun Project

Creative Commons: Daniel Kulinski

Porous Layer Open Tubular (PLOT) Columns


-- Advantages --

- Retention for high vapor pressure solutes
 - No cryogenics needed
 - Useful for multi-column & valve switching & heart cutting techniques
- Variety adsorbents
 - Porous polymers, alumina, molecular sieve, silica
 - Selectivity for isomeric compounds (HC isomers)
 - Gases & solvents separations
- Efficiency of a capillary
 - Preferred over packed
 - Packed columns have advantages in specific applications



PLOT -- Some Limitations & Considerations --

- Lower sample capacity than thick film (1 – 5 μ m) liquid phase columns
0.53mm PLOT offers highest loadability
- More limited maximum temperatures than siloxanes – limiting analyte scope
- Highly active -- some compounds do not elute
or require extensive column conditioning (CO₂ / H₂O)
- Molecular rearrangements above MAOT – selectivity might change
Do not exceed MAOT!
- Particles may dislodge under mechanical / flow stress – but solutions do exist



Porous Polymers

HP-PLOT Q & U
PoraBOND Q & U
PoraPLOT Q, U & S
GS-Q
PoraPLOT Q-HT

Zeolites

HP-Molesieve
CP-Molsieve 5A

Oxygenates

CP-Lowox
GS-OxyPLOT

Alumina

Al₂O₃/KCl, HP-PLOT Al₂O₃/KCl
Al₂O₃/Na₂SO₄, HP-PLOT Al₂O₃ S
GS-Alumina
Select MAPD

Porous Silica

GS-GasPro
SilicaPLOT

Graphitised Carbon


CarboPLOT P7
CarboBOND
GS-CarbonPLOT

Compounds & PLOT Column Selection

<ul style="list-style-type: none">○ Permanent gases<ul style="list-style-type: none">- Noble gases, O₂, N₂, H₂, CO, CO₂, CH₄	Molsieve	Select Permanent Gases
<ul style="list-style-type: none">○ Light hydrocarbons<ul style="list-style-type: none">- C1 – C8, saturated, unsaturated,- Aromatics, natural gas, C2, C3, C4 streams	Alumina Silica	Porous carbon Porous polymers
<ul style="list-style-type: none">○ Volatile sulfur compounds<ul style="list-style-type: none">- H₂S, COS, mercaptanes, sulfides, disulfides	Select Low Sulfur Silica	Porous polymers
<ul style="list-style-type: none">○ Oxygenates<ul style="list-style-type: none">- Alcohols, ketones, ethers	Lowox	OxyPLOT
<ul style="list-style-type: none">○ Solvents<ul style="list-style-type: none">- Oxygenates, aromatics, alkanes,- Chlorinated hydrocarbons	Porous polymers	
<ul style="list-style-type: none">○ Chlorinated and fluorinated hydrocarbons	Silica	Porous polymers

Porous Polymer PLOT Columns

- Analyses of polar and non-polar volatile compounds
Oxygenates, gases, halogenated compounds,
hydrocarbons C1 –C6, ketones, solvents
- Most versatile of adsorbent materials
- Suitable for aqueous injections
(in split & direct mode through inlet)
- Elution of water as a sharp peak + quantification of water
(on the polar U type)
- Recommended for column switching systems
- Styrene/divinyl benzene polymer = Non-polar “Q type”
- Styrene-glycol methacrylate co-polymer = Polar “U type”
- Divinylbenzene/vinylpyridine polymer = Mid-Polar “S type”




HP-PLOT Q, GS-Q & PoraPLOT Q

- Separation ethane/ethylene, propane/propylene
 - % levels, not ppm impurities in C2, C3 matrix
- H₂S and COS
 - Refinery gas analysis (RGA)
 - 100 ppm - %, **but not low ppm due to strong surface interactions for H₂S**
- Pre-column for O₂, N₂, CH₄, CO separations on Molsieve column
- Poor C4 isomer separation
 - Alumina preferred
- Not sensitive to moisture -- Stable retention -- No selectivity changes
- Good choice for solvents and CFCs


Refinery fuel gas, porous polymer separation

-- HP- PLOT Q --




Aqueous Injections on Porous Polymer Q-PLOT


Before



After 5 x 5 μ l water, Splitless at 80°C



Retention times
unchanged
for all compounds




Samples containing water
can be analyzed

Isothermal & short cycle
time


Agilent "BOND" Technology & Traditional Particle Traps

Particle PLOT Q




Particle size: 0.1- 2 μm
Number: $> 10^{12}$

PoraBOND Q




Bonded polymer layer
Reduced need particle trap



Solvent Analysis on PoraBOND Q

Conditions


Technique	: GC-wide-bore
Column	: Agilent PoraBOND Q, 0.53 mm x 25 m fused silica PLOT (df = 10 μm) (Part no. CP7354)
Temperature	: 100 °C (2 min) → 300 °C, 5 °C/min
Carrier Gas	: He, 25 kPa (0.25 bar, 3.5 psi)
Injector	: Split, T = 250 °C
Detector	: FID, T = 250 °C
Sample Size	: 0.5 μL
Concentration Range	: 0.1% per compound
Solvent Sample	: DMSO



Peak identification

1. methane
2. methanol
3. ethanol
4. acetonitrile
5. acetone
6. dichloromethane
7. 2-propanol (isopropanol)
8. dimethyl sulfide
9. diethyl ether
10. 1-propanol
11. pentane
12. 2-butanone
13. trichloromethane
14. tetrahydrofuran
15. ethyl acetate
16. 2-methoxyethanol
17. isobutanol
18. butanol
19. hexane
20. benzene
21. trichloroethylene
22. cyclohexane
23. 1,4-dioxane
24. pyridine
25. N,N-dimethylformamide
26. dimethyl sulfoxide

Integrated Particle Trap (PT)





Integrated particle trap at front & back of the column

PLOT column

Assurance & Safety
for valves & detectors
in
for- and backflush modes

Proof of Agilent PLOT-PT Column Performance

- 150°C → 250°C @ 20 °C/min
- Pressure 3x higher than optimal
- 15 cycles
- 10 x Carrier gas off / on



Agilent PLOT-PT

- Zero particle shedding
- No unions and fittings – No blockage
- Multi column valve switching
- Particle traps integrated on both ends – supports back flush & CFT technology
- Similar selectivity, efficiency & peak shape performance to existing Agilent porous polymer PLOT columns
- Minimum method redevelopment required
- Available for most popular porous polymer PLOT columns

PLOT-PT Chromatographic Performance

Column: PoraBond Q PT, 30m×0.25mm,3um

Carrier : Helium, 43cm/s @ 90 °C

Oven: : 90 °C- 140 °C at 10°C/min

140 °C for 6min


140 °C - 200 °C @ 5°C/min

200 °C for 10 min

Injection: Split, 250 °C, split ratio1:160


Detector: MSD, 280 °C Transfer line, full Scan @ m/z 30-350

- | | |
|-----------------------|---------------------------|
| 1. Methyl Alcohol | 14. Trichloromethane |
| 2. Acetaldehyde | 15. 2-Butanone (MEK) |
| 3. Ethanol | 16. Ethyl Acetate |
| 4. Acetonitrile | 17. sec-Butyl alcohol |
| 5. Acetone | 18. MTBE |
| 6. Methylene Chloride | 19. 2-chlorobutane |
| 7. Isopropyl Alcohol | 20. 1-Butanol |
| 8. 2-Propanamine | 21. Benzene |
| 9. Ethyl Formate | 22. 1,1,1-Trichloroethane |
| 10. 1-Propanol | 23. 1-chlorobutane |
| 11. Ethyl ether | 24. Carbon Tetrachloride |
| 12. t-Butyl alcohol | 25. Hexane |
| 13. 1,2-Ethanediol | 26. 1,4-Dioxane |



- | |
|-------------------------------|
| 27. Pyridine |
| 28. Dimethyl Formamide (DMF) |
| 29. Isoamyl Alcohol |
| 30. Dimethyl Sulfoxide (DMSO) |
| 31. Toluene |
| 32. Heptane |
| 33. Paraldehyde |
| 34. Chlorobenzene |
| 35. Ethylbenzene |
| 36. m-Xylene |
| 37. p-Xylene |
| 38. o-xylene |

Selectivity Comparison of Porous Polymer PLOT



GC Conditions	
Columns	Agilent J&W PoraPLOT Q, 25 m × 0.53 mm × 20 μm (p/n CP7554) Agilent J&W PoraPLOT S, 25 m × 0.53 mm × 20 μm (p/n CP7574) Agilent J&W PoraPLOT U, 25 m × 0.53 mm × 20 μm (p/n CP7584) Agilent J&W PoraPLOT Q, 10 m × 0.53 mm × 20 μm (p/n CP7553) Agilent J&W PoraBOND Q, 25m × 0.53 mm × 10 μm (p/n CP7354)
Carrier	Helium, constant flow, 5 mL/min
Oven	50 °C (1.0 min), Ramp 10 °C/min to 200 °C (2 min) *Ramp to 180 °C for the J&W PoraPLOT U
Inlet	SSL Inlet, split mode, 200 °C, split ratio 20:1
Inlet Liner	Ultra Inert, split, low pressure drop, glass wool (p/n 5190-2295)
GC/FID	Agilent 7890B GC equipped with FID
Sampler	Agilent 7693 automatic liquid sampler
FID Conditions	
Temperature	200 °C
Hydrogen	30 mL/min
Air	350 mL/min
Col + Make Up	25 mL/min

Alumina Adsorbent and PLOT Columns

-- Best selectivity for light hydrocarbon separations --

Applications

- General C1 – C6 (C9) hydrocarbons
- Natural Gas
- Ethylene streams, impurities
- Impurities in propylene
- Butylene streams, impurities
complex C4 composition
- Environmental hydrocarbon distributions

Separation & Conditions


- Degree of saturation
Elution: alkane, alkene, alkyne, (dialkenes)
- Selectivity: types of deactivation
KCl, Na₂SO₄ and Proprietary
- Selectivity: column flow & temperature



Moisture carrier gas (& sample)
Retention & selectivity impacted
Gas Clean moisture filters recommended!

Selectivity Difference between KCl and Na₂SO₄

- 1. methane
- 2. ethane
- 3. ethene
- 4. propane
- 5. propene
- 6. cyclopropane
- 7. ethyne
- 8. iso-butane
- 9. propadiene
- 10. n-butane
- 11. t-2-butene
- 12. 1-butene
- 13. iso-butene
- 14. c-2-butene
- 15. 1,3-butadiene
- 16. propyne



KCl for
Impurities in Ethylene
Impurities Propylene
Impurities in 1,3 Butadiene


Na₂SO₄ for C4 feeds

Alumina columns tend to exhibit
overload for highly unsaturated analytes
more quickly

Impurities in 1,3-Butadiene

50 m x 0.32 mm Al₂O₃ / KCl, 5 μm

100 °C → 200 °C, 6 °C/min



Absorption and Response of Hydrocarbons


-- The Case for Al₂O₃ MAPD --



BE AWARE

Alumina columns exhibit partial irreversible absorption of highly unsaturated HC

No 100% level calibration!




Inorganic Gases on GC-GasPro

30m x 0.32mm

Carrier: Helium at 53 cm/sec
Oven: 25°C for 3 min
25-200°C at 10°/min
200°C Hold
Injector: Split 1:50, 200°C, 50 µL
Detector: TCD, 250°C

1. Nitrogen
2. CO₂
3. SF₆
4. COS
5. H₂S
6. Ethylene oxide
7. SO₂



- GS-GasPro & CP-SilicaPLOT
- Light hydrocarbon separation: C1 – C4
- Extended hydrocarbon range compared to other PLOTs
- Inert enough for low ppm light sulfurs, H₂S, COS, mercaptans
- CFC's
- Not sensitive to moisture in carrier gas
- MSD compatible
- MAOT of 200°C

CP-SilicaPLOT -- Hydrocarbons

Column : Agilent CP-SilicaPLOT, 0.53 mm x 30 m, fused silica
PLOT CP-SilicaPLOT (df = 6 µm) (Part no. CP8570)

Temperature : 50 °C (5 min) → 225 °C, 5 °C/min

Carrier Gas : He, 20 kPa (0.2 bar, 2.9 psi)

Injector : Direct
T = 225 °C


Detector : FID
T = 250 °C

Sample Size : 2 µL



Concentration Range : ppm %

Peak identification

1. methane
2. ethane
3. ethylene
4. acetylene
5. propane
6. cyclopropane ●
7. propylene
8. isobutane
9. butane
10. 1-butene
11. propyne (methylacetylene)
12. 1,3-butadiene
13. isobutene
14. cis-2-butene
15. hexane



Permanent Gases on Parallel Column Solution



- **Application:** All permanent gases incl. CO₂, water, ethane, ethylene
- Optimized dimensions for fast analysis
- PoraBOND Q allows high temperature bake-out (300°C) of CO₂ & H₂O


Oxygenates blended in gasoline


- Additives to boost octane content, prevent engine “knocking”
- MTBE, ETBE, Ethanol
- **% level oxygenates**
- **GC analyses on Wax or TCEP polar liquid phase columns**

Oxygenate in intermediates (monomers, naphtha's)

- Lower effectiveness catalysts - lower yield
- Higher costs of catalysts
- More refinery downtime
- **ppm level oxygenates**
- **GC analyses on Lowox, OxyPLOT columns, Agilent exclusives**

Analyzing Oxygenates in Hydrocarbon Matrix


- Low ppm concentration level oxygenates – Low analyte capacity of stationary phase
- FID detection (MSD uncommon)
- Agilent columns with high selectivity hydrocarbons/oxygenates: CP-Lowox / OxyPLOT
- Very moisture sensitive  Use Gas Clean moisture filters in carrier gas




ASTM Trace Oxygenate Analysis Methods

-- Valve Configuration --


**Transfer of Oxygenates
Valve Off**



**Venting Hydrocarbons
Valve On**



Oxygenate Analysis




Select Low Sulfur - PLOT

- Super permeable & porous stationary phase
Proprietary material
- Unique selectivity characteristics for volatile sulfurs/hydrocarbons
No detector quenching for propane/propylene matrix
- Highly inert for volatile sulfurs, H₂S, Methyl mercaptane
- One part number: 60m x 0.32mm, p/n CP8575
- Near 100% recovery for H₂S at 20 – 100 ppb levels
- Zero particle loss
- FPD, PFPD, SCD, AED Compatible

Sulfur components in propylene by PFPD

- 500 ppbv H₂S & COS
- Same response
- No absorption H₂S



Technique: GC-PFPD
Column: Agilent J&W Select Low Sulfur, 60 m × 0.32 mm (p/n CP8575)
Oven: 65 °C isotherm
Carrier gas: Helium, constant flow, 2.0 mL/min
Injector: 200 °C, split 1:20
Detector: PFPD, 200 °C
Sample: Propylene matrix containing ~500 ppb H₂S, COS, and CH₃SH
Injection volume: 1 mL
Injection: Gas sampling valve


Considerations for PLOT Column Analysis

Columns

- Selectivity – Consider elution pattern
- Capacity – Overloading – 0.53mm --
- Inertness (low conc. sulfurs)
- Temperature limits

Contamination

- Efficiency loss; “ghost peaks”; increase in bleed
- Water, CO₂, high molecular weight hydrocarbons
- Carrier gas purifiers highly recommended



- Agilent supplies largest selection of PLOT columns in the market for all gases and volatiles applications with dedicated columns for challenging analyses in the petrochemical industry.
- Fully QC tested to assure column to column reproducibility with the excellent peak shape performance and separation for the best data accuracy.
- Agilent PLOT columns come with the lowest levels of particle shedding for better baseline stability and trouble-free analyses.
- Exclusive to Agilent, are the PLOT-PT columns with integrated particle traps to assure “spike free” detection, mass spectrometer compatibility and improved system performance with complex valve applications.

Agilent CrossLab

From Insight to Outcome

