

Development and Validation of Liquid Chromatography-Time of Flight Method of Bacopaside-I and Metabolites for Pharmacokinetic Study

IMSC 2012 PTh-086

Sontaya Sookying¹, Kornkanok Ingkaninan^{1*},
Dumrongsak Pekthong², Sarawut Oo-puthinan²,
Jie Xing^{3**}, Zhaoqi Zhan³

¹Department of Pharmaceutical Chemistry and
Pharmacognosy, Faculty of Pharmaceutical Sciences,
Naresuan University, Phitsanulok 65000, Thailand

²Department of Pharmacy Practice, Faculty of
Pharmaceutical Sciences, Naresuan University,
Phitsanulok 65000, Thailand

³Customer Support Centre, Shimadzu (Asia Pacific)
Pte Ltd, 79 Science Park Drive, #02-01/08, SINTECH
IV, Singapore Science Park 1, Singapore 118264

*Corresponding author, **Presenter

Development and Validation of Liquid Chromatography-Time of Flight Method of Bacopaside-I and Metabolites for Pharmacokinetic Study

Introduction

Bacopa monnieri (L.) Wettst. (Brahmi) is an Ayurvedic medicine plant used for centuries as a memory enhancer. There have been several studies on biological effects of this plant for a therapeutic potential in treatment and prevention of neurological diseases and improvement of cognitive processes. Bacopaside I (C₄₆H₇₄O₂O₅), a dammarane-type triterpenoid saponin, is one of the main active components of the plant and, therefore, it has been produced and used in biological and therapeutic studies. Various analytical method such as HPTLC, ELISA and

HPLC methods were used to determine bacopasides I and other active components in the plant extracts as well as biological samples. However, so far LC/MS method for quantitative analysis of Bacopaside I in biological samples has not been reported. We report for the first time the development and validation of a high sensitivity LCMS (TOF) method, aiming for pharmacokinetic study of bacopasides I and identification of its metabolites produced in vitro and in vivo.

Experimental

Rat urine and feces were used as the biological matrix for method development and validation of bacopaside I for quantitative analysis. Liquid-liquid extraction method was employed in sample extraction and purification. The LC-MS method used is based on high resolution TOF which is suitable for both quantitative analysis of Bacopaside I in biological samples and identification of metabolites formed. An LCMS-IT-TOF (Shimadzu Corporation, Japan) was used

in this study. A gradient elution method using a Kinetex C18 column (1.7 μm) was developed and optimized. Both target compound Bacopaside I and the internal standard Digitoxin showed very strong mass peaks (*m/z* 977.4387 and 763.4243) in negative ESI mode. A high resolution MS and MS/MS method were established to achieve high sensitivity and reliability of analysis.

Table 1 LCMS-IT-TOF conditions for analysis of Bacopaside-I and metabolites

Column	KinetexC18, 1.7 μm, 50 mm L x 2.1 mm ID	Interface	ESI
Mobile phase	A: Water (0.1% formic acid) B: Acetonitrile (0.1% formic acid)	MS mode	MS (-) and MS ² (-) <i>m/z</i> 200~1200 CID with Ar gas
Elution mode	Gradient, 16 min	Block temp.	200°C
Flow rate	0.3 mL/min	CDL temp.	250°C
Column temp.	45°C	Nebulizing gas	Nitrogen, 1.5 L/min
Injection vol.	10 μL	Drying gas	Nitrogen, 10 L/min

Results and Discussion

This study was carried out in three stages, development of a LC-MS (TOF) method, validation trials and identification of metabolites of bacopaside I from in vitro enzymatic

hydrolysis and faces samples using a MetID Solution program.

1. Method Development: Fig. 1 shows the chromatograms and MS, MS/MS spectra (TOF) of bacopaside I. The internal standard (IS) digitoxin eluted immediately

Development and Validation of Liquid Chromatography-Time of Flight Method of Bacopaside-I and Metabolites for Pharmacokinetic Study

immediately after the target peak. Both target compound and IS were detected in negative ESI mode as $[M-H]^-$ ions. The LCMS-IT-TOF is a high resolution MS with mass error less than 5 ppm. The spectrum results in Fig. 1 are -3.8 ppm and -3.0 ppm for Bacopaside I and digitoxin, respectively. The advantage of high resolution MS is the superior mass selectivity of target peaks. Therefore, it is suitable for quantitation of trace level compounds in complex matrix. The fragmentation pathway of bacopasidel in negative mode ($m/z977.4387 \rightarrow m/z845.3962$) is shown in Fig. 2.




Fig 1. Chromatograms and mass spectra of bacopaside I and internal standard (digitoin, $C_{41}H_{64}O_{13}$) on high resolution LCMS-IT-TOF.




Fig. 2 Fragmentation scheme of bacopaside I (left) to main product ion under CID in negative mode

2. Validation of method: The performance of the quantitation method based on extract ion chromatogram (EIC) of MS spectrum was evaluated systemically using spiked samples prepared from rat urine and feces by L-L extraction. The calibration curve shows good linearity ($r^2=0.9985$) for the testing range from 4.8 ng/mL to 380 ng/mL (Fig. 3). The specificity and matrix effect were studied and it was found that the method specificity was highly reliable due to the accurate mass (<5 ppm) (Fig. 4). Matrix effect was not obvious under the conditions. The LOQ of the method was lower than 4.8 ng/mL. The reproducibility at this level (spiked in urine) was 8.7% (RSD of concentration, $n=6$).

The recoveries of bacopaside I spiked in urine were between 91.1% and 107.1% for three concentration levels (50, 80 and 100 ng/mL) as shown in Table 2. However, the recovery of feces extract samples was found at 40~70%. Further study to improve the recovery of bacopaside I in feces samples is undergoing.




Fig. 3 Calibration curve of bacopaside I (IS method) based on EIC (m/z 977.4393) in negative mode.

Development and Validation of Liquid Chromatography-Time of Flight Method of Bacopaside-I and Metabolites for Pharmacokinetic Study




Table 2 Recovery of bacopaside I spiked in urine samples

Sample Name	Spiked conc. (ng/mL)	Measure Conc. (ng/mL)	Recovery (%)
3/1 n1	50	45.6	91.1
3/2 n2	50	49.8	99.6
3/3 n3	50	53.8	107.5
3/4 n1	80	83.9	104.9
3/5 n2	80	80.5	100.6
3/6 n3	80	81.1	101.4
3/7 n1	100	98.7	98.7
3/8 n2	100	107.1	107.1
3/9 n3	100	105.5	105.5

Fig. 4 EICs of bacopaside I and IS in blank urine (top) and spiked sample (bottom).

3. Identification of metabolite: The LCMS-IT-TOF method established with its full spectrum (m/z 200~1200) and accurate mass measurement (<5 ppm) was applied to *in vitro* hydrolysis and feces samples to detect metabolites formed. The method was modified slightly to include MS positive mode (m/z 200~1200) in order to detect both positive and negative ions formed from metabolites. Three samples used are control, enzymatic hydrolyzed samples (20 mins) and rat feces collected after 12 hrs of oral

administration of bacopaside I. In this preliminary study, we found seven metabolite candidates including $C_{30}H_{48}O_4$ (observed (+) m/z 473.3594, RT3.98 min), $C_{30}H_{46}O_3$ (observed (+) m/z 455.3493, RT3.76 min) and $C_{35}H_{54}O_7$ (observed (+) m/z 587.3891, RT3.56 min). These three candidates are likely corresponding to hydrolysis of bacopaside I on the glycoside with leaving the pentacyclic structure intact (see Fig. 5).




Fig. 5 Structures of metabolite candidates M1, M2 and M3 from bacopaside I

Development and Validation of Liquid Chromatography-Time of Flight Method of Bacopaside-I and Metabolites for Pharmacokinetic Study

Conclusions

For the first time a high resolution LCMS-IT-TOF method for quantitative analysis of bacopaside I in biological samples has been established and validated using rat urine/feces spiked samples for pharmacokinetics study. The method with its full spectrum (m/z 200~1200) and accurate mass measurement ability (<5 ppm) was also used to detect and identify metabolites formed.