

Development of High Performance Liquid Chromatography Tandem Mass Spectrometry Method for Quantitative Analysis of Bacopaside-I in Rat Urine and Feces Samples

### IMSC 2012 PWe-133

Sontaya Sookying<sup>1</sup>, Kornkanok Ingkaninan<sup>1\*</sup>, Dumrongsak Pekthong<sup>2</sup>, Sarawut Oo-puthinan<sup>2</sup>, Jie Xing<sup>3</sup>, Zhaoqi Zhan<sup>3\*\*</sup>

PO-CON1249E

<sup>1</sup>Department of Pharmaceutical Chemistry and Pharmacognosy, Faculty of Pharmaceutical Sciences, Naresuan University, Phitsanulok 65000, Thailand <sup>2</sup>Department of Pharmacy Practice, Faculty of Pharmaceutical Sciences, Naresuan University, Phitsanulok 65000, Thailand

<sup>3</sup>Customer Support Centre, Shimadzu (Asia Pacific) Pte Ltd, 79 Science Park Drive, #02-01/08, SINTECH IV, Singapore Science Park 1, Singapore 118264 \*Corresponding author; \*\*Presenter

# Introduction

Bacopaside I (C46H74O20S) is one of active components isolated from Bacopa monnieri (L.) Wettst. (Brahmi), which has been used as an Ayurvedic medicine plant for centuries <sup>[1]</sup>. It has been reported that Brahmi has potential therapeutic effect in treatment and prevention of neurological diseases and improvement of cognitive processes<sup>[2,3]</sup>. Memory enhancing effect of Brahmi has been established through animal experiments and in healthy volunteers. The damarane triterpenoid saponin, Bacopaside I, is a main active compound in Brahmi. Reversed phase HPLC method has been established and used to determine bacopaside I and other active components in the plant extracts as well as biological samples<sup>[4]</sup>. However, so far LC/MS method for quantitative analysis of Bacopaside I in biological samples has not been reported. We report for the first time the development and validation of a LCMS/ MS method, aiming for pharmacokinetic study of bacopasides I.

## Experimental

A LCMS-8030 (Shimadzu Corporation, Japan) triple quadrupole LC-MS/MS system was used in this method development study. A fast gradient elution separation program using a Kinetex C18 HUPLC column (1.7 um, 50 mmL × 2.1 mmID) was developed and optimized. The MRM transitions employed were 979.4>473.4 for bacopaside I and 609.3>195.0 for reserpine as internal standard (IS) in ESI positive mode. Rat urine and feces were used as biological matrix in the development and validation of a quantitative analysis method for bacopaside I. Liquid-liquid extraction method was employed in sample extraction and purification.

Table 1 LCMS-8030 conditions for quantitative analysis of Bacopaside-I


LC conditions		MS conditions & MRM parameters			
Column	Kinetex C18 column (1.7 um), 50 mm L x 2.1mm ID	Interface	ESI		
Mobile phase	A: Water (0.1% formic acid) B: Acetonitrile (0.1% formic acid)	MRM parameters (positive)	Bacopaside I: 979.4>483.4, CE:-26V, Q1B: -25V; Q3B: -25V Reserpine (IS): 609.3>195.0		
Elution mode	Gradient elution	Die els terrer exeture	CE:-37V 300°C		
	10 minutes	Block temperature	300 C		
Flow rate	0.3 mL/min	CDL temperature	250°C		
Column temp.	40°C	Nebulizing gas flow	Nitrogen, 2.0 L/min		
Injection vol.	15 uL	Drying gas flow	Nitrogen, 15 L/min		

## Results and Discussion

1. Method development: Bacopaside I is ionized by ESI in both positive and negative modes to form  $[M+H]^+$  (*m/z* 979.4) and  $[M-H]^-$  (*m/z* 977.4). We have established an LCMS-IT-TOF method based on negative ESI mode due to its higher intensity in MS mode [5]. It was found that for MRM method, positive ESI method (979.4>483.4) showed better sensitivity than negative mode. The structure and fragmentation of bacopaside I under MRM conditions are shown Fig 1.

#### **SHIMADZU** Excellence in Science Development of High Performance Liquid Chromatography Tandem Mass Spectrometry Method for Quantitative Analysis of Bacopaside-I in Rat Urine and Feces Samples

Fig. 2 shows the chromatograms of a mixed sample of bacopaside I and reserpine. It can be seen that reserpine is an ideal internal standard (IS) in this method. The elution of the two compounds were closed, but not overlapped. Therefore, ion suppression of the IS to bacopaside I was avoided.



2. Performance evaluation: The performance of the MRM quantitation method was evaluated systemically using spiked samples prepared from rat urine and feces by LL extraction method. The quantitation calibration curve of




Fig 1. Structure and Fragmentation pathway of Bacopaside-I on LC/MS/MS in positive mode.

Fig 2. MRM chromatograms of Bacopaside-I and IS on LCMS-8030

bacopaside I is shown in Fig. 2. The linearity  $(r^2)$  was 0.9999 for a range from 8.6 ng/mL to 900 ng/mL.




Fig 3. Calibration curve of Bacopaside-I by internal standard method on LCMS-8030. Conc. of bacopaside I: 8.6, 25.6, 85.5, 342 and 900 ng/mL (IS: 10 ng/mL).

Due to extremely flat baseline, the LOD and LOQ could not be calculated from S/N ratio. The lowest concentration of bacopaside I detectible by this method was about 4.3 ng/mL (spiked in urine) with poorer reproducibility (Fig. 4). The peak area reproducibility at 8.6, 25.7 and 85.5 ng/mL (spiked in urine) was 17.2%, 6.3% and 2.9% (RSD, n=6), respectively. Based on these results, the LOQ of the method is estimated to be about 15 ng/min and LOD is about 5 ng/mL. The recoveries of bacopaside I in urine (spiked) were between 104% and 129% for three concentration levels (45, 72 and 90 ng/mL) (Table 2). However, the recovery




Fig 4. Chromatograms of urine blank and spiked standard samples in urine.

Т

Table 2 Recovery of Bacopaside I in urine samples

No	Sample Name	Conc. Spiked (ng/mL)	Ret. Time (min)	Conc. (ng/mL)	Recovery (%)	Ave Recovery (%)
1	3/1 urine n1	45	3.532	45.42	100.9	
2	3/2 urine n2	45	3.534	45.36	100.8	103.6
3	3/3 urine n3	45	3.531	49.12	109.2	
4	3/4 urine n1	72	3.531	75.48	104.8	
5	3/5 urine n2	72	3.534	87.12	121.0	115.2
6	3/6 urine n3	72	3.535	86.24	119.8	
7	3/7 urine n1	90	3.529	116.52	129.5	
8	3/8 urine n2	90	3.545	108.73	120.8	129.3
9	3/9 urine n3	90	3.534	123.77	137.5	

Table	3	Recovery	of	Bacopaside	I	in	feces	samples	
-------	---	----------	----	------------	---	----	-------	---------	--

No	Sample Name	Conc. Spiked (ng/mL)	Ret. Time (min)	Conc. (ng/mL)	Recovery (%)	Ave Recovery (%)
1	3/1 feces n1	45	3.519	20.37	45.3	
2	3/2 feces n2	45	3.515	15.8	35.1	42.1
3	3/3 feces n3	45	3.522	20.65	45.9	
4	3/4 feces n1	72	3.504	32.71	45.4	
5	3/5 feces n2	72	3.501	32.38	45.0	44.4
6	3/6 feces n3	72	3.506	30.9	42.9	
7	3/7 feces n1	90	3.511	26.58	29.5	
8	3/8 feces n2	90	3.502	46.41	51.6	41.2
9	3/9 feces n3	90	3.503	38.37	42.6	

## Conclusions

For the first time, a high sensitivity LC/MS/MS method for quantitative analysis of bacopaside I in rat urine and feces has been developed and validated. The LOD and LOQ of the method were estimated to be 5 and 15 ng/mL. The recovery of the method for urine samples was at 104~129%, but much lower for feces samples. Further study to improve extraction recovery from feces sample is undergoing.

#### References

- 1. H.K. Singh, B.N., Dhawan, Indian J. Pharmacol. 29 (1977) \$359
- 2. A. Russo, F. Borrelli, Phytomedicine 12 (2005) 305
- 3. D. Vohora, S.N Pal, K.K. Pillai, J. Ethnopharmacol. 71 (2000) 383
- 4. M. Deepak, G.K Sanggli, P.C. Arun, A Amit, Phytochem Anal 16 (2005) 24
- 5. Sontaya Sookying et al., IMSC 2012 (Kyoto), Poster, PTh-086 (2012)

First Edition: September, 2012



Shimadzu Corporation www.shimadzu.com/an/

For Research Use Only. Not for use in diagnostic procedures.

The content of this publication shall not be reproduced, altered or sold for any commercial purpose without the written approval of Shimadzu. The information contained herein is provided to you "as is" without warranty of any kind including without limitation warranties as to its accuracy or completeness. Shimadzu does not assume any responsibility or liability for any damage, whether direct or indirect, relating to the use of this publication. This publication is based upon the information available to Shimadzu on or before the date of publication, and subject to change without notice.

© Shimadzu Corporation, 2012