

APPLICATIONS

Orthogonal Selectivity for Separation of Light and Heavy Petroleum Hydrocarbons by 2D GC using Zebron™ ZB-1HT Inferno™ and ZB-35HT Inferno GC Columns

Badaoui Omais¹, Ramkumar Dhandapani¹, Tim Nelson¹, Chantal Lorentz², and Christophe Geantet²

¹ Phenomenex, Inc., 411 Madrid Ave., Torrance, CA 90501, USA

² IRCELYON (Institute of Research on Catalysis and the Environment of Lyon), ECI2D (Energy, Fuels and Intermediaries for Sustainable Development), Université de Lyon, UMR 5256 France



Badaoui Omais

has completed his engineering studies at Ecole des Mines by a Ph.D. at IFP Energies Nouvelles in collaboration with ESPCI. He is now the sales manager of Phenomenex SAS and passionately supports GC labs with trainings, seminars, and webinars.

Introduction

Complex sample mixtures, such as petroleum hydrocarbon fractions, require a detailed separation of the individual organic compound types. A traditional GC approach will lead to multiple coelutions, however in this present work a 2D GC approach was utilized to explore the multi-dimensional separation of complex light and heavier hydrocarbons. The complex hydrocarbon mixture can be thoroughly resolved due to the combination of the high peak capacity in 2D GC, along with the orthogonal selectivity of the Zebron GC columns that provide a complementary selectivity in two-dimensions. This complete resolution includes a range of compounds, for example linear alkanes (paraffins), linear alkenes (olefins), cyclic alkanes, aromatics (one - tetralin ring derivatives), and aromatics (one, two, three, four, and five - benzene ring derivatives). The combination of Zebron ZB-1HT Inferno with ZB-35HT Inferno provides a high temperature capability to further elute all the fractions and to provide a targeted selectivity that separates individual aliphatic and aromatic analytes from each other.

GC Conditions for Analysis

Column 1: Zebron ZB-1HT Inferno

Phase: 100 % Dimethylpolysiloxane

Dimensions: 15 meter x 0.10 mm x 0.10 µm

Part No.: [7FB-G014-02-C](#) (15 m length column, cut from a 20 m column)

Column 2: Zebron™ ZB-35HT Inferno

Phase: 35 % Phenyl 65 % Dimethylpolysiloxane

Dimensions: 5 meter x 0.25 mm x 0.18 µm

Part No.: [7FG-G025-08-C](#) (5 m length column, cut from a longer one)

Injection: Split 20:1 @ 250 °C, 1 µL

Recommended Liner: ZebronPLUS Straight Z-Liner™

Part No.: [AG2-0A03-05](#)

Instrument: 2D GC - Agilent system 7890A modified with gas phase Microfluidics

Carrier Gas (column 1): Hydrogen @ 1 mL/min (constant flow)

Carrier Gas (column 2): Hydrogen @ 2.4 mL/min (constant flow)

Detector: Flame Ionization (FID) @ 360 °C, 100 Hz acquisition rate

Oven Program for Figure 1

Column 1: 100 °C for 10 min, to 400 °C @ 2.2 °C/min, hold for 20 min



Column 2: 100 °C for 10 min, to 360 °C @ 2.2 °C/min, hold for 20 min

Oven Program for Figure 2

Column 1: 40 °C for 8 min, to 360 °C @ 8 °C/min, hold for 5 min

Column 2: 50 °C for 8 min, to 360 °C @ 2.25 °C/min, hold for 27 min


Figure 1.
Gas Oil Petroleum Fraction Separation by 2D GC




Having trouble reproducing this method? We would love to help!
Visit www.phenomenex.com/LiveChat to get in touch with one of our Technical Specialists

Figure 2.

Heavy Petroleum Fraction Separation by 2D GC

**Figure 3.**

The Zebron HT Inferno Technical Advantage



Results and Discussion

A two-dimensional GC analysis of both a gas oil petroleum fraction and a heavy petroleum fraction were performed using a Zebron™ ZB-1HT Inferno™ coupled to a Zebron ZB-35HT Inferno on an Agilent® 7890A GC system that was modified with a gas phase microfluidics system. The Zebron ZB-35HT Inferno column with its unique combination of a high temperature capability and a polar selectivity provides a detailed orthogonal separation along with the ZB-1HT that provides a true boiling-point based separation. For this study hydrogen was used as the carrier gas for a faster and more efficient analysis. The FID detector provides a great sensitivity and wider dynamic range for hydrocarbons. A high data acquisition rate of 100 points per second was utilized to capture all the narrow peaks eluting out from the second dimension GC column.

The Zebron HT Inferno columns have a special high temperature resistant polyimide coating. In addition, this stationary phase has an Engineered Cross-linking™ (ESC™) bonding technology which provides extreme low bleed even at high temperature (**Figure 3**). **Figure 1** shows the carbon range for the gas oil fraction up to C42 and all the compound classes were successfully separated and eluted using this 2D GC Zebron Inferno HT column set. The use of the Zebron ZB-1HT Inferno on the first dimension and the Zebron ZB-35HT Inferno on the second dimension provides a traditional, yet much improved approach in 2D GC. This selectivity combination enables an orthogonal separation of the analytes along both the axes to provide an improved targeted separation. The Zebron ZB-1HT Inferno separates by boiling point on the first dimension and the Zebron ZB-35HT Inferno then separates orthogonally with pi-pi interactions on the second dimension.

There is a clear separation of the following classes of compounds along both dimensions:

1. Linear alkanes (n- and iso-paraffins)
2. Linear alkenes (olefins)
3. Cyclic alkanes
4. Aromatics (one-ring, tetralin derivatives)
5. Aromatics (one-ring, benzene derivatives)
6. Aromatics (2-rings, naphthalene derivatives)
7. Aromatics (2-rings, di-benzene derivatives)
8. Aromatics (3-rings)
9. Aromatics (4-rings)
10. Aromatics (5-rings)

In similar fashion, detailed separation of heavier fractions was done by using a 2D GC approach as shown in **Figure 2**. The detailed molecular mapping in **Figures 1** and **2** are crucial to optimize processes in petroleum refining.

Conclusion

2D GC is a versatile technique that resolves complex mixtures. The column combination of the Zebron ZB-1HT Inferno and the unique Zebron ZB-35 HT Inferno not only provided a high temperature capability to elute all the fractions, but also provided complementary orthogonal selectivity to separate individual aliphatic and aromatics from each other.

Acknowledgement: This study was performed at IRCELYON GC lab.

We thank IRCELYON for their collaboration on this project. Phenomenex is not affiliated with IRCELYON

APPLICATIONS

Ordering Information

Zebron ZB-1HT Inferno GC Columns

ID(mm)	df(µm)	Temp. Limits °C	Part No.
5-Meter			
0.53	0.10	-60 to 400/430	7AK-G014-02
10-Meter			
0.32	0.25	-60 to 400/430	7CM-G014-11
15-Meter			
0.25	0.10	-60 to 400/430	7EG-G014-02
0.25	0.25	-60 to 400/430	7EG-G014-11
0.32	0.10	-60 to 400/430	7EM-G014-02
0.32	0.25	-60 to 400/430	7EM-G014-11
0.53	0.15	-60 to 400	7EK-G014-05
20-Meter			
0.18	0.18	-60 to 400/430	7FD-G014-08
30-Meter			
0.25	0.10	-60 to 400/430	7HG-G014-02
0.25	0.25	-60 to 400/430	7HG-G014-11
0.32	0.10	-60 to 400/430	7HM-G014-02
0.32	0.25	-60 to 400/430	7HM-G014-11
0.53	0.15	-60 to 400	7HK-G014-05

Note: If you need a 5 in. cage, simply add a (-B) after the part number, e.g., [7HG-G014-11-B](#). Some exceptions may apply. Agilent 6850 and some SRI and process GC systems use only 5 in. cages.

Australia
t: +61 (0)2 9428-6444
[auinfo@phenomenex.com](#)

Austria
t: +43 (0)1-319-1301
[anfrage@phenomenex.com](#)

Belgium
t: +32 (0)2 503 4015 (French)
t: +32 (0)2 511 8666 (Dutch)
[beinfo@phenomenex.com](#)

Canada
t: +1 (800) 543-3681
[info@phenomenex.com](#)

China
t: +86 400-606-8099
[cninfo@phenomenex.com](#)

Denmark
t: +45 4824 8048
[nordicinfo@phenomenex.com](#)

Finland
t: +358 (0)9 4789 0063
[nordicinfo@phenomenex.com](#)

France
t: +33 (0)1 30 09 21 10
[franceinfo@phenomenex.com](#)

Germany
t: +49 (0)6021-58830-0
[anfrage@phenomenex.com](#)

India
t: +91 (0)40-3012 2400
[indiainfo@phenomenex.com](#)

Ireland
t: +353 (0)1 247 5405
[ireinfo@phenomenex.com](#)

Italy
t: +39 051 6327511
[italiainfo@phenomenex.com](#)

Luxembourg
t: +31 (0)30-2418700
[nlinfo@phenomenex.com](#)

[www.phenomenex.com](#)

Phenomenex products are available worldwide. For the distributor in your country, contact Phenomenex USA, International Department at [international@phenomenex.com](#)

Ordering Information

Zebron ZB-35HT Inferno GC Columns

ID(mm)	df(µm)	Temp. Limits °C	Part No.
15-Meter			
0.25	0.10	40 to 400	7EG-G025-02
0.25	0.25	40 to 400	7EG-G025-11
0.32	0.25	40 to 400	7EM-G025-11
20-Meter			
0.18	0.18	40 to 400	7FD-G025-08
30-Meter			
0.25	0.10	40 to 400	7HG-G025-02
0.25	0.25	40 to 400	7HG-G025-11
0.32	0.25	40 to 400	7HM-G025-11

Note: If you need a 5 in. cage, simply add a (-B) after the part number, e.g., [7HG-G025-11-B](#). Some exceptions may apply. Agilent 6850 and some SRI and process GC systems use only 5 in. cages.

Mexico
t: 01-800-844-5226
[tecnicomx@phenomenex.com](#)

The Netherlands
t: +31 (0)30-2418700
[nlinfo@phenomenex.com](#)

New Zealand
t: +64 (0)9-4780951
[nzinfo@phenomenex.com](#)

Norway
t: +47 810 02 005
[nordicinfo@phenomenex.com](#)

Portugal
t: +351 221 450 488
[ptinfo@phenomenex.com](#)

Singapore
t: +65 800-852-3944
[sginfo@phenomenex.com](#)

Spain
t: +34 91-413-8613
[espinfo@phenomenex.com](#)

Sweden
t: +46 (0)8 611 6950
[nordicinfo@phenomenex.com](#)

Switzerland
t: +41 (0)61 692 20 20
[swissinfo@phenomenex.com](#)

Taiwan
t: +886 (0) 0801-49-1246
[twinfo@phenomenex.com](#)

United Kingdom
t: +44 (0)1625-501367
[ukinfo@phenomenex.com](#)

USA
t: +1 (310) 212-0555
[info@phenomenex.com](#)

All other countries 
Corporate Office USA
t: +1 (310) 212-0555
[info@phenomenex.com](#)

Live Chat today!

Our Technical Nerds are
available nearly 24/7.

[www.phenomenex.com/chat](#)

BE-HAPPY™
guarantee

Your happiness is our mission. Take
45 days to try our products. If you
are not happy, we'll make it right.
[www.phenomenex.com/behappy](#)

Terms and Conditions

Subject to Phenomenex Standard Terms & Conditions, which may be viewed at [www.phenomenex.com/TermsAndConditions](#).

Trademarks

Zebtron, Inferno, ESC, and Z-Liner are trademarks of Phenomenex. Agilent is a registered trademark of Agilent Technologies, Inc.

FOR RESEARCH USE ONLY. Not for use in clinical diagnostics procedures.

© 2019 Phenomenex, Inc. All rights reserved.



Having trouble reproducing this method? We would love to help!

Visit [www.phenomenex.com/LiveChat](#) to get in touch with one of our Technical Specialists