

# **Packed Column GC Application Guide**

 **SUPELCO**

ISO 9001  
REGISTERED

## Packed Column GC Application Guide

*This guide to packed column applications provides information on the best approach to designing your application. It includes helpful information about packings, phases, tubing, and different types of deactivations. Nearly 200 different applications are organized according to type of analysis. In most cases, the catalog number for the packing is included with the application chromatogram.*

### Key Words

- packed column
- gas chromatography

### What Makes a Good Column?

Tubing, support, and stationary phase are all that is needed to prepare a column, but to prepare a column that will provide the desired separation may require much more. First, the proper materials must be selected; this can be the most difficult part. Also, these ingredients must be combined in such a manner that optimum performance can be obtained; this can turn out to be more than one bargains for! Any compromise in materials or procedures yields a compromised performance.

### The Rohrschneider and McReynolds Constants

An important feature is the use of the Rohrschneider and McReynolds Constants. Using these systems, most of the stationary phases now in use have been characterized by a series of constants. These constants provide the worker with a wealth of practical information to help solve problems in column selection. For more information and a listing of McReynolds values, request Bulletin 880.

### The Support

The support plays a critical role in several ways in the performance of the column. First, it governs the efficiency of the column (narrowness of peaks). The structure of the support, and the manner in which it is coated, both contribute to the efficiency. Secondly, the wrong choice of support can interact with the sample to cause the chromatographic peaks to tail, i.e., they can be highly asymmetrical and consequently difficult or impossible to measure. Ideally, the support should not interact with the sample but, in practice, this does occur. By careful selection of the support and conditions, one can minimize this problem.

The tailing phenomenon is caused by active sites on the surface of the support. These sites are ones that can form a hydrogen bond; consequently samples that form a strong hydrogen bond tail badly. In practice, compounds such as water, glycols, alcohols, acids, and amines tail severely while carbonyl

compounds such as esters, ketones, and aldehydes tail to a lesser degree. Hydrocarbons that do not hydrogen bond are not bothered by tailing.

To eliminate or reduce the tailing problem one modifies the support surface by: 1) removing the active sites by acid and/or base washing, 2) modifying the surface by silanization, or 3) covering the active sites with the stationary phase having polar functional groups in it.

Acid washing is effective in removing mineral impurities from the support surface as well as miscellaneous extraneous material. Base washing does not impart any special advantage to the support that is not obtained with a well acid washed support. Acid washing by itself is not effective in reducing tailing but it is recommended where a polar phase is used, such as the polyesters and polyglycols.

Silanization, particularly with dimethyldichlorosilane (DMCS) is very effective in reducing tailing. Combined with acid washing and DMCS treatment, the resulting support is recommended for most columns. Silane treatment is a very difficult process to control and should not be attempted unless you wish to make a research project out of it. When silicone stationary phases are used, it is mandatory that an acid-washed and DMCS-treated support be used. The silicone stationary phases, particularly when used in the 1-5% level, are not effective in deactivating the support and require a silane-treated support. SUPELCOPORT is recommended where an acid-washed and DMCS-treated support is required.

The third procedure for deactivation, using a polar stationary phase, requires that the phase contain functional groups such as an ester, an ether, a hydroxyl, and an amine group. These functional groups have strong hydrogen bonding characteristics and tie up the active sites on the support surface. These phases do not require a silanized support, although, when they are used at a level of 5% or less, silanization can be useful. Normally we recommend using an acid-washed support with these phases.

When analyzing acids it is necessary that the stationary phase contain an acid to deactivate the support. See Bulletin 856.

When working with basic compounds such as amines, the stationary phase must contain a base to deactivate the support; otherwise severe tailing will result. KOH is frequently used at a 1-2% level for the purpose. A basic stationary phase such as polyethyleneimines also may be used.

Most of the GC supports in current use are made from diatomaceous earth, also called diatomite. The diatomite is processed in several ways producing two basic types of supports. These are conveniently recognized by their color.

The particle size of supports are generally expressed in terms of screen openings since screens are normally used to prepare them. The particle sizes normally used in GC are as follows:

60/80 mesh	250-177 microns
80/100 mesh	177-149 microns
100/120 mesh	149-125 microns

The designation 60/80 mesh means that the particles have passed through a 60 mesh screen (-60) and will not pass through the 80 mesh screen (+80). It then means that the particles are between 250-177 microns in size. The column efficiency improves with decreasing particle size. At present the 80/100 mesh is the most popular size, but 100/120 mesh is used with increasing frequency when more efficiency columns are desired.

## Column Tubing

Glass, stainless steel, aluminum, and copper are the tubing materials commonly used for columns. While glass is the most inert of the tubing, stainless steel is the most widely used.

Glass is used in situations where the sample might interact with the walls of the tube. It is standard operating procedure to use glass columns when working with pesticides and biochemicals such as steroids and hormones. Glass is more inert than the metals and rarely causes tailing or decomposition of the samples. Glass columns are also used generally in situations where it is desirable that the sample be injected directly into the column. An alternative is to use a glass lined injection port to avoid the contact with metal since metal in the inlet can destroy certain samples or cause tailing. Glass columns also permit you to see how well a column has been packed.

Fused silica lined stainless steel tubing is a suitable substitute for glass tubing. When properly deactivated, it meets or exceeds the deactivation of conventional deactivated glass columns.

Metal columns are used where glass is not required. Stainless steel is generally considered more inert than aluminum or copper. We recommend stainless steel as the metal of choice because it is easier to prepare columns of higher efficiency with stainless steel than with aluminum or copper. The hardness of the material appears to be important in transmitting the vibrational energy when the column is vibrated or tapped. However, if in doubt . . . use glass . . . and silane treat it! (All Supelco glass columns are silane-treated.)

SP-Alloy T-1, our nickel alloy, approaches glass in inertness. Most newer instruments are designed to handle 1/8" OD metal columns. When glass columns are to be used, the instruments are usually equipped to handle 1/4" OD columns. Here the glass column can be made with a heavy glass wall cutting down considerably on the problem of breakage. Where 1/8" OD glass columns are used, the breakage problems are usually severe.

From our own experience we recommend using 2mm ID rather than 4mm ID glass columns because of the inherently higher efficiency. If you shift from 4mm to 2mm columns, remember to reduce the carrier gas flow rate to 1/4, not 1/2. A reasonable flow rate for 4mm columns is 80mL/min, while for 2mm columns a good rate is 20mL/min.

## Upper Temperature Limit

Each stationary phase has an upper temperature limit above which the column should not be operated. Most stationary phases are polymers that consist of materials having a range of molecular weights. As the column temperature is increased, the more volatile portion of the polymer is swept out of the column by the carrier gas. The volatile products could also be formed by thermal degradation of the stationary phase while the column is being used. This is called bleed, and is seen on the recorder as a rise in the baseline or as noise. Above the maximum upper limit the bleed rate is very high and the column will have a relatively short life. In some cases the bleed rate is so high that it will not be possible to move the recorder pen off of full scale. Note that the reported maximum upper limits are approximate values. Each column produced is shipped with a conditioning sheet, and minimum/maximum operating temperature is indicated. For a listing of phases and their maximum/minimum operating temperatures, refer to our general catalog.

## Conditioning the Column

Before a column is used for the first time, it must be conditioned for a period of hours to rid it of the very volatile portion of the stationary phase and also the last traces of solvent which were used in the coating step.

We recommend conditioning the column at the maximum isothermal column temperature *which you expect to reach*. One can start to use a column after it has been conditioned at that temperature after only 3-4 hours, but conditioning overnight is usually better. We recommend that you always condition the column with a flow of carrier gas. Others have recommended "no-flow" conditioning for such materials as UC W98 silicone. This might be necessary for technical grade silicone, but chromatographic grade silicones such as the SP™ and OV® phases, and most other phases, should not be subjected to this drastic treatment.

Many stationary phases are subject to oxidation and should not be exposed to air while hot. At all times when the column is hot, the carrier gas should be allowed to flow through it. Cutting back on the carrier gas overnight can be detrimental to the column. It may also be necessary to insert an oxygen scrubber in the carrier gas line.

## The Lower Temperature Limit

Observing the lower temperature limit is also important to obtain good column performance. Some stationary phases are solids and are ineffective as a stationary phase in this form. They must be operated at a temperature above their melting point.

## Selecting a Column

The simplest method of obtaining a column to handle a particular problem is to duplicate one that is referred to in the applications pages of this guide. Another method is to use the literature as a guide and then to substitute some of the newer, improved stationary phases and solid supports. If you cannot find the solution to your separation problem in the literature, the problem is more complex. When ordering the column, it is essential that you list all of the following:

- a. the stationary phase
- b. the support
- c. the mesh range of the support
- d. the treatment given the support
- e. the percent coating of stationary phase
- f. the tube composition, i.e., glass, stainless, etc.
- g. the dimensions of the tube, i.e., length and diameter, both OD and ID if possible
- h. the configuration of the column or instrument it is to be used with

**Example:** 3% SP-2100 on 100/120 mesh SUPELCOPORT™, 6' x 1/4" OD, 2mm ID, Glass for HP6890.

For a complete listing of our current line of columns, packings, and accessories, please refer to our general catalog.

## Index to Chromatograms

Compound Class	Page
Acetaldehyde .....	6
Acids .....	6
Acids, Methyl Esters (also see Fatty Acid Methyl Esters) .....	11
Alcohols .....	13
Amines .....	19
Amino Acids .....	21
Aroclors .....	22
Ammonia .....	24
Aromatics .....	24
Chlorinated Hydrocarbons .....	28
Creosote .....	31
Drugs .....	33
EPA Methods (600 Series) .....	37
Ethylene Oxide and Residues .....	46
Fatty Acid Methyl Esters .....	47
Fragrances .....	51
Freons .....	53
Formaldehyde .....	54
Hydrocarbons	
C1-C2 Hydrocarbons & Permanent Gases .....	55
C1-C3 Hydrocarbons .....	59
C1-C4 Hydrocarbons .....	60
C1-C5 Hydrocarbons .....	60
C1-C6 Hydrocarbons .....	62
Permanent Gases and Hydrocarbons .....	64
Pesticides .....	68
Phenols .....	69
Pheromones .....	71
Simulated Distillation .....	72
Solvents .....	78
Sugars .....	84
Sulfur and Mercaptans .....	85
Steroids .....	92
Sterols .....	94
Triglycerides .....	95

## Trademarks


Aroclor — Monsanto Co.  
 Bentone — National Lead Co., Baroid Sales Div.  
 Carbopack — Sigma-Aldrich Co.  
 Carbotrap — Sigma-Aldrich Co.  
 Carbowax — Union Carbide Corp.  
 Carboxen\* — Sigma-Aldrich Co.  
 Cellosolve — Union Carbide Corp.  
 Chromosorb — Celite Corp.  
 Dexsil — Dexsil Chemical Corp.  
 Fluorcol — Sigma-Aldrich Co.  
 Freon — E.I. du Pont de Nemours & Co., Inc.  
 Hall — Tracor Instruments, Austin, Inc.  
 HayeSep — Hayes Separations Inc.

ORBO — Sigma-Aldrich Co.  
 Oronite — Uniroyal, Inc.  
 OV — Ohio Valley Specialty Chemical Company  
 Petrocol — Sigma-Aldrich Co.  
 Qualmix — Sigma-Aldrich Co.  
 Scott — Scott Specialty Gases, Inc.  
 SP — Sigma-Aldrich Co.  
 SUPELCOPORT — Sigma-Aldrich Co.  
 Supelpak — Sigma-Aldrich Co.  
 Teflon — E.I. du Pont de Nemours & Co., Inc.  
 TightSpec — Sigma-Aldrich Co.  
 UCON — Union Carbide Corporation  
 Valco — Valco Instruments Co., Inc.  
 \*US Pat. No. 4,839,331.

## Figure 1. Trace Acetaldehyde in Air (2500ppm)

Packing: **Hayesep® D, 100/120 mesh**  
Cat. No.: **10293** (packing)  
Column: 3' x 1/8" ID stainless steel  
Oven: 100°C  
Carrier: helium, 30cc/min  
Det.: P.E. 900 TC, 225 ma, 140°C  
Inj.: Valco® valve, 100mL, 140°C

1. Air
2. Water
3. Acetaldehyde (vol. 2500ppm)




For more information, request Bulletins 790 and 816.  
Chromatogram provided courtesy of Hayes Separations, Inc.

795-0144

## Figure 2. Volatile Fatty Acids (C1-C7)

Packing: **10% SP-1000/1% H<sub>3</sub>PO<sub>4</sub> on 100/120 Chromosorb® W AW**  
Cat. No.: **11841** (packing)  
Column: 6' x 4mm ID glass  
Oven: Fisher Model 2400, 147°C  
Carrier: helium, 86mL/min  
Det.: TC  
Inj.: 14µL

1. Acetic acid
2. Formic acid
3. Propionic acid
4. Isobutyric acid
5. n-Butyric acid
6. Isovaleric acid
7. n-Valeric acid
8. Isocaproic acid
9. n-Caproic acid
10. Heptanoic acid




For more information, request Bulletin 856.

795-0077

**Figure 3. Volatile Acid Standard Mix (C1-C7)**

Packing: 15% SP-1220/1% H<sub>3</sub>PO<sub>4</sub> on 100/120 Chromosorb W AW  
 Cat. No.: 12144 (packing)  
 Column: 6' x 4mm ID glass  
 Oven: Fisher Model 2400, 145°C  
 Carrier: helium, 70mL/min  
 Det.: TC  
 Inj.: 14µL




For more information, request Bulletin 856.

795-0076

**Figure 4. Volatile Fatty Acids in Sewage Effluents (ppm Levels)**

Packing: GP 60/80 Carbowax® C/0.3% Carbowax® 20M/0.1% H<sub>3</sub>PO<sub>4</sub>  
 Cat. No.: 11825-U (packing, 15g/bottle)  
 Column: 30' x 4mm ID glass  
 Oven: 120°C  
 Carrier: nitrogen, 50mL/min  
 Det.: FID, 200°C  
 Inj.: 1µL water containing 50ppm each analyte, 200°C



Useful for analyses of sewage treatment effluents and other wastewater.

For trace analyses, use a glass column.


For more information, request Bulletin 856.

795-0078



## Figure 5. Free Fatty Acids (C2-C10)

Packing: **GP 10% SP-1200/1% H<sub>3</sub>PO<sub>4</sub>** on 80/100 Chromosorb W AW  
Cat. No.: **11965** (packing, 20g/bottle)  
Column: 6' x 1/8" stainless steel  
Oven: 125°C to 175°C at 15°C/min  
Carrier: nitrogen, 20mL/min  
Det.: FID, 200°C  
Inj.: 0.5µL carbon disulfide containing 0.1% each analyte, 150°C




713-0963

## Figure 6. Free Fatty Acids (C14:0 to C18:3)

Packing: **5% DEGS-PS on 100/120 SUPELCOPORT**  
Cat. No.: **11870-U** (packing)  
Column: 3' x 2mm ID glass  
Oven: 200°C  
Carrier: nitrogen, 20mL/min  
Det.: FID, 210°C  
Inj.: 1µL Qualmix™ FA


1. Myristic acid (C14:0)
2. Palmitic acid (C16:0)
3. Palmitoleic acid (C16:1)
4. Stearic acid (C18:0)
5. Oleic acid (C18:1)
6. Linoleic acid (C18:2)
7. Arachidic acid (C20:0)
8. Linolenic acid (C18:3)



795-0079

### Figure 7. Whole Milk Free Acids

Packing: **10% SP-216-PS on 100/120 SUPELCOPORT**  
 Cat. No.: **11879** (packing)  
 Column: 3' x 2mm ID glass  
 Det: FID, 200°C  
 Oven: 130°C to 200°C at 15°C/min  
 Carrier: nitrogen, 20mL/min  
 Inj.: 0.5µL carbon disulfide containing 0.1% each analyte, 150°C




For more information, request Bulletin 856.

713-0962

### Figure 8. Lactic Acid in Corn Silage (C2-C5)

Packing: **80/120 Carbowax B-DA/4% Carbowax 20M**  
 Cat. No.: **11889** (packing)  
 Column: 2m x 2mm ID TightSpec™ glass  
 (stock column available)  
 Oven: 175°C  
 Carrier: nitrogen, 24mL/min  
 Det.: FID, 200°C  
 Inj.: 1µL aqueous mixture, 200°C

	Conc. (ppm)
U Unknown	
1. Acetic acid	370
2. Propionic acid	78
3. Isobutyric acid	282
4. Butyric acid	1125
IS Trimethylacetic acid (int. std.)	50
5. 2-Methylbutyric acid	220
6. Isovaleric acid	284
7. Lactic acid	710
8. Valeric acid	465




For more information, request Bulletin 856.

795-0080

## Figure 9. Lactic Acid in Rumen Fluid

Packing: **80/120 Carbowax B-DA/4% Carbowax 20M**  
 Cat. No.: **11889** (packing)  
 Column: 2m x 2mm ID TightSpec glass  
 (stock column available)  
 Oven: 175°C  
 Carrier: nitrogen, 24mL/min  
 Det.: FID, 200°C  
 Inj.: 1µL aqueous mixture, 200°C

	Conc. (ppm)
U Unknown	
1. Acetic acid	1300
2. Propionic acid	1060
3. Isobutyric acid	43
4. Butyric acid	313
IS Trimethylacetic acid (int. std.)	50
5. 2-Methylbutyric acid	60
6. Isovaleric acid	41
7. Valeric acid	70




For more information, request Bulletin 856.

795-0082

## Figure 10. Acrylic/Propionic Separation in C2-C5 Acids

Packing: **80/120 Carbowax B-DA/4% Carbowax 20M**  
 Cat. No.: **11889** (packing, 15g/bottle)  
 Column: 2m x 2mm ID  
 TightSpec glass  
 (stock column available)  
 Oven: 175°C  
 Carrier: nitrogen, 24mL/min  
 Det.: FID, 200°C  
 Inj.: 1µL water, 10ppm each  
 acid, 100ppm each  
 internal standard, 200°C




1. Methyl ethyl ketone (int. std.)
2. Acetic acid
3. Propionic acid
4. Acrylic acid (5000ppm)
5. Isobutyric acid
6. Butyric acid
7. Trimethylacetic acid (int. std.)
8. 2-Methyl butyric acid
9. Isovaleric acid
10. Valeric acid

For more information, request Bulletin 856.

713-0964

### Figure 11. Nonvolatile Methylated Acids

Packing: 10% SP-1000/1% H<sub>3</sub>PO<sub>4</sub> on  
100/120 Chromosorb W AW  
Cat. No.: 11841 (packing)  
Column: 6' x 4mm ID glass  
Oven: 155°C  
Carrier: helium, 86mL/min




For more information, request Bulletin 856.

795-0081

### Figure 12. N-TFA-O-Methyl Esters of D- & L-Valine Amino Acids

Packing: 5% SP-300 on 100/120 SUPELCOPORT  
Cat. No.: 11835-U (packing)  
Column: 6' x 2mm ID glass  
Oven: 120°C  
Carrier: nitrogen, 20mL/min  
Inj.: 1µL (1µg/µL conc.)




795-0158

### Figure 13. Methyl Esters of Dibasic Acids

Packing: 10% SP-2340 on 100/120 Chromosorb W AW  
 Cat. No.: 11852 (packing)  
 Column: 6' x 1/8" ID stainless steel  
 Oven: 170°C to 220°C at 4°C/min  
 Carrier: nitrogen, 20mL/min  
 Inj.: 0.75µL, 10mg/cc each acid

1. Dimethyl malonate (C3)
2. Dimethyl succinate (C4)
3. Dimethyl glutarate (C5)
4. Dimethyl adipate (C6)
5. Dimethyl pimelate (C7)
6. Dimethyl suberate (C8)
7. Dimethyl azelate (C9)
8. Dimethyl sebacate (C10)
9. Dimethyl terephthalate
10. Dimethyl isophthalate
11. Dimethyl phthalate




713-0976

### Figure 14. Bacterial Acid Standard

Packing: 3% SP-2100 DOH on 100/120 SUPELCOPORT  
 Cat. No.: 12101 (packing)  
 Column: 10' x 2mm ID glass  
 Oven: 150°C to 225°C at 4°C/min  
 Carrier: nitrogen, 20mL/min  
 Det.: FID  
 Inj.: 1µL

- |               |               |
|---------------|---------------|
| 1. 11:0       | 13. 16:0      |
| 2. 2-OH 10:0  | 14. □-17:0    |
| 3. 12:0       | 15. 17:0      |
| 4. 13:0       | 16. 17:0      |
| 5. 2-OH 12:0  | 17. 2-OH 16:0 |
| 6. 3-OH 12:0  | 18. 3-OH 16:0 |
| 7. 14:0       | 19. 18:1      |
| 8. □-15:0     | 20. 18:0      |
| 9. 15:0       | 21. 19:0      |
| 10. 2-OH 14:0 | 22. 19:0      |
| 11. 3-OH 14:0 | 23. 20:0      |
| 12. 16:1      |               |




For more information, request Bulletin 856 and Application Note 50.

795-0159

### Figure 15. Bile Acid Methyl Esters

Packing: **3% SP-2250 on 100/120 SUPELCOPORT**  
 Cat. No.: **11878** (packing, 20g/bottle)  
 Column: **3' x 2mm ID glass**  
 Oven: **275°C**  
 Carrier: **helium, 40mL/min**  
 Det.: **FID**

1. Methyl lithocholate (0.94µg/µL)
2. Methyl deoxycholate (1.25µg/µL)
3. Methyl chenodeoxycholate (1.57µg/µL)
4. Methyl cholate (2.50µg/µL)




For more information, request Bulletin 856.

713-0977

### Figure 16. Alcoholic Beverage/Fermentation Products

Packing: **80/120 Carboxpack B AW/6.6% PEG 20M**  
 Cat. No.: **11814** (packing, 15g/bottle)  
 Column: **2m x 2mm ID glass**  
 Oven: **80°C to 200°C at 4°C/min**  
 Carrier: **nitrogen**  
 Det.: **FID**  
 Inj.: **1µL water:ethanol (50:50), 40-60ppm each analyte**

- |                      |                         |
|----------------------|-------------------------|
| 1. Acetaldehyde      | 20. Butanol             |
| 2. Methanol          | 21. Ethyl isobutyrate   |
| 3. Propanol          | 22. 3-Methylbutan-2-ol  |
| 4. Acetone           | 23. 3-Pentanol          |
| 5. Methyl acetate    | 24. 2-Pentanol          |
| 6. Ethanol           | 25. Isobutyl acetate    |
| 7. Isobutanol        | 26. Ethyl butyrate      |
| 8. Butanol           | 27. Butyl acetate       |
| 9. Isopropanol       | 28. 2-Methylpentan-1-ol |
| 10. Ethyl acetate    | 29. 1-Methylpentan-1-ol |
| 11. Diacetal         | 30. Acetic acid         |
| 12. n-Propanol       | 31. Pentanol            |
| 13. Isopentanol      | 32. Isoamyl acetate     |
| 14. sec-Butanol      | 33. Furfural            |
| 15. Pentanol         | 34. Propionic acid      |
| 16. Ethyl propionate | 35. Hexanol             |
| 17. Propyl acetate   | 36. Isobutyric acid     |
| 18. Isobutanol       | 37. Butyric acid        |
| 19. Acetal           |                         |



Use this packing in glass columns only.

For more information, request Bulletin 790.


Figure provided by Dr. Antonio DiCorcia and colleagues, University of Rome (used with permission of Elsevier Scientific Publishing Co., Amsterdam).

713-0711

## Figure 17. Tequila Headspace

Packing: **HayeSep D, 80/120 mesh**  
 Cat. No.: **10293** (packing)  
 Column: 3' x 1/8" ID stainless steel  
 Oven: 100°C  
 Carrier: helium, 30cc/min  
 Det.: P.E. 900 T.C., 225 ma, 140°C  
 Inj.: Valco valve, 100µL, 140°C

- |                   |           |
|-------------------|-----------|
| 1. Air            | (Att. x1) |
| 2. Carbon dioxide | (Att. x1) |
| 3. Water          | (Att. x1) |
| 4. Methanol       | (Att. x1) |
| 5. Acetaldehyde   | (Att. x1) |
| 6. Ethanol        | (Att. x8) |



For more information, request Bulletin 790.


Chromatogram provided courtesy of Hayes Separations, Inc.

795-0160

## Figure 18. Cherry Brandy

Packing: **80/100 Carbowax B AW/5% Carbowax 20M**  
 Cat. No.: **11812-U** (packing)  
 Column: 2m x 2mm ID glass  
 (stock column available)  
 Oven: 70°C to 170°C at 5°C/min  
 Carrier: helium, 20mL/min  
 Det.: FID  
 Inj.: 0.5µL cherry brandy

- |                        |
|------------------------|
| 1. Acetaldehyde        |
| 2. Methanol            |
| 3. Ethanol             |
| 4. Ethyl acetate       |
| 5. n-Propanol          |
| U. Unknown             |
| 6. Isobutanol          |
| 7. Acetal              |
| 8. Acetic acid         |
| 9. Active amyl alcohol |
| 10. Isoamyl alcohol    |
| U. Unknown             |




For more information, request Bulletin 790.

795-0083

## Figure 19. Bourbon

Packing: **80/120 Carbowax B AW/5% Carbowax 20M**  
 Cat. No.: **11812-U** (packing)  
 Column: 2m x 2mm ID glass  
 (stock column available)  
 Oven: 70°C to 170°C at 5°C/min  
 Carrier: helium, 20mL/min  
 Det.: FID  
 Inj.: 0.5µL bourbon whiskey

1. Acetaldehyde
2. Methanol
3. Ethanol
4. Ethyl acetate
5. n-Propanol
6. Isobutanol
7. Acetic acid
8. Active amyl alcohol
9. Isoamyl alcohol




For more information, request Bulletin 790.

795-0084

## Figure 20. Scotch Whiskey

Packing: **80/120 Carbowax B/6.6% Carbowax 20M**  
 Cat. No.: **11814** (packing)  
 Column: 2m x 2mm ID glass  
 Oven: 80°C to 150°C at 4°C/min  
 Carrier: nitrogen, 20mL/min  
 Det.: FID  
 Inj.: 0.5µL Scotch whiskey

1. Acetaldehyde
2. Methanol
3. Ethanol
4. Ethyl acetate
5. n-Propanol
6. Isobutanol
7. Active amyl alcohol
8. Isoamyl alcohol
9. Acetic acid



For more information, request Bulletin 790.


795-0085



## Figure 21. Denaturants in Alcohol

Packing: **GP 80/100 Carbowax C/0.1% SP-1000**  
Cat. No.: **11820** (packing)  
Column: **6' x 2mm ID glass**  
(stock column available)  
Oven: **80°C**  
Carrier: **nitrogen, 20mL/min**  
Det.: **FID**  
Inj.: **70% ethanol & 2% of each denaturant**

1. Methanol
2. Ethanol
3. Acetone
4. Isopropanol
5. n-Propanol
6. tert-Butanol
7. Chloroform, methyl ester
8. Methyl Cellosolve®
9. sec-Butanol
10. Isobutanol
11. Ethyl acetate
12. n-Butanol
13. Ethyl Cellosolve
14. Benzene
15. Methyl isobutyl ketone




For more information, request Bulletin 790.

795-0161

## Figure 22. Blood Alcohols

Packing: **60/80 Carbowax B/5% Carbowax 20M**  
Cat. No.: **11766** ((packing, 15g/bottle)  
Column: **6' x 2mm ID glass**  
(stock column available)  
Oven: **85°C**  
Carrier: **helium, 20mL/min**  
Det.: **FID**  
Inj.: **1µL water,**  
**0.05-0.10% each analyte**

1. Acetaldehyde
2. Methanol
3. Acetone
4. Ethanol
5. Isopropanol
6. n-Propanol




For more information, request Application Note 72.

713-0622

### Figure 23. Alcohols (C1-C5) and Acetaldehyde

Packing: **80/120 Carbowax B AW/5% Carbowax 20M**  
 Cat. No.: **11812-U** (packing, 15g/bottle)  
 Column: 2m x 2mm ID glass  
 (stock column available)  
 Oven: 70°C to 170°C at 5°C/min  
 Carrier: nitrogen, 20mL/min  
 Det.: FID  
 Inj.: 0.5µL Scotch whiskey

1. Acetaldehyde
2. Methanol
3. Ethanol
4. Ethyl acetate
5. n-Propanol
6. Isobutanol
7. Acetic acid
8. Active amyl alcohol
9. Isoamyl alcohol




Use this packing in glass columns only.  
 For more information, request Bulletin 790.

713-0710

### Figure 24. Isomeric Alcohols (C1-C5)

Packing: **80/100 Carbowax C/0.2% Oronite® NIW**  
 Column: 2m x 2mm ID glass  
 Oven: 125°C  
 Carrier: nitrogen, 20mL/min  
 Det.: FID  
 Inj.: 0.1µL

1. Methanol
2. Ethanol
3. 2-Propanol
4. 1-Propanol
5. 2-Methyl-2-propanol (tert-butyl)
6. 2-Butanol (sec-butyl)
7. 2-Methyl-1-propanol (isobutyl)
8. 1-Butanol
9. 2-Methyl-2-butanol (tert-amyl)
10. 2,2-Dimethyl-1-propanol
11. 3-Methyl-2-butanol
12. 3-Pentanol
13. 2-Pentanol
14. 2-Methyl-1-butanol (active)
15. 3-Methyl-1-butanol
16. 1-Pentanol




For more information, request Bulletin 790.

795-0162

## Figure 25. Isomeric Alcohols (C1-C5)

Packing: **80/100 Carbopack F-SL**  
Column: 2m x 2mm ID glass  
Oven: 120°C  
Carrier: nitrogen, 15mL/min  
Det.: FID  
Inj.: 0.1µL

1. Methanol
2. Ethanol
3. 2-Propanol
4. 1-Propanol
5. 2-Methyl-2-propanol (tert-butyl)
6. 2-Butanol (sec-butyl)
7. 2-Methyl-1-propanol (isobutyl)
8. 1-Butanol
9. 2-Methyl-2-butanol (tert-amyl)
10. 2,2-Dimethyl-1-propanol
11. 3-Methyl-2-butanol
12. 3-Pentanol
13. 2-Pentanol
14. 2-Methyl-1-butanol (active)
15. 3-Methyl-1-butanol
16. 1-Pentanol




For more information, request Bulletin 790.

713-1069

## Figure 26. Alcohols (C10-C18)

Packing: **10% SP-2100 on 80/100 SUPELCOPORT**  
Cat. No.: **12140** (packing)  
Column: 6' x 1/8" ID stainless steel  
(stock column available)  
Oven: 225°C  
Carrier: nitrogen, 20mL/min  
Det.: FID  
Inj.: 0.5µL, 0.1% each  
in chloroform




1. 1-Decanol (C10)
2. 1-Dodecanol (C12)
3. 1-Tetradecanol (C14)
4. 1-Hexadecanol (C16)
5. 1-Octadecanol (C18)

For more information, request Bulletin 790.

795-0163

**Figure 27. Alcohols (C14:0-C18:0, C18:1, & C18:2)**


Packing: 10% SP-2300 on 80/100 SUPELCOPORT  
 Cat. No.: 12091 (packing)  
 Column: 6' x 1/8" ID stainless steel  
 Oven: 200°C  
 Carrier: nitrogen, 20mL/min  
 Det.: FID  
 Inj.: 1µL, 1% total sample in chloroform



795-0164

**Figure 28. Aliphatic Amines in Water (100ppm)**

Packing: Carbowax B/4% Carbowax 20M/0.8% KOH  
 Cat. No.: 11887 (packing)  
 Column: 6' x 2mm ID glass  
 (stock column available)  
 Oven: 90°C to 150°C at 4°C/min  
 Carrier: nitrogen, 20mL/min  
 Det.: FID  
 Inj.: 0.5µL




For more information, request Bulletin 737.

795-0165

## Figure 29. Amines

Packing: **80/100 HaysSep B**  
Cat. No.: **10286** (packing, 75cc/bottle)  
Column: 5' x 1/8" ID stainless steel  
Oven: 140°C to 190°C at 16°C/min  
Carrier: helium, 30cc/min  
Det.: thermal conductivity (175ma, 180°C)  
Inj.: 0.2µL, on-column injection, 150°C




For more information, request Bulletin 737.  
Chromatogram provided courtesy of Hayes Separations, Inc.

713-0717

## Figure 30. Heterocyclic Amines in Water (1000ppm)

Packing: **Carbopack B/4% Carbowax 20M/0.8% KOH**  
Cat. No.: **11887** (packing)  
Column: 6' x 2mm ID glass  
(stock column available)  
Oven: 140°C  
Carrier: nitrogen, 20mL/min  
Det.: FID  
Inj.: 0.4µL



For more information, request Bulletin 737.


795-0166

### Figure 31. Amines

Packing: **60/80 Carbowax B/4% Carbowax 20M/0.8% KOH**  
 Cat. No.: **11887** (packing, 15g/bottle)  
 Column: **2m x 2mm ID glass**  
 (stock column available)  
 Oven: **100°C (2 min) to 220°C**  
 at 16°C/min  
 Carrier: **helium, 30mL/min**  
 Det.: **thermionic (nitrogen-specific)**

#### 11ng Each Amine

1. Butylamine
- IS Triethylamine (int. std.)
2. n-Amylamine
3. n-Hexylamine
4. n-Heptylamine
5. n-Octylamine




For more information, request Bulletin 737.

713-0715

### Figure 32. TMS Amino Acids

Packing: **10% OV-11 on 100/120 SUPELCOPORT**  
 Cat. No.: **11961** (packing)  
 Column: **6' x 2mm ID glass**  
 Oven: **100°C (5 min) to 230°C at 3°C/min**  
 Carrier: **nitrogen, 30mL/min**


1. ALA
2. GLY<sub>2</sub>
3. IMP<sub>2</sub>
4. VAL
5. LEU
6. GLY<sub>3</sub>
7. PRO
8. SER
9. THR
10. HYPRO
11. ASP
12. MET
13. GLU<sub>2</sub>
14. GLU<sub>3</sub>
15. PHE<sub>3</sub>
16. ARG
17. LYS
18. TYR
19. HIS
20. TRP
21. CYS



795-0167

### Figure 33. Aroclors

Packing: **1.5% SP-2250/1.95% SP-2401 on 100/120 SUPELCOPORT**  
Cat. No.: **11947** (packing)  
Column: 2m x 1/4" OD x 4mm ID  
TightSpec glass  
(stock column available)  
Oven: 160°C (Aroclor 1016, 1221, 1232, 1242, 1248)  
200°C (Aroclor 1254, 1260)  
Carrier: nitrogen, 60mL/min  
Det.: ECD  
Inj.: 5ng each Aroclor in 5µL isooctane




For more information, request Bulletin 817 and Application Note 67.

795-0170

### Figure 34. Aroclors

Packing: **3% SP-2100 on 100/120 SUPELCOPORT**  
 Cat. No.: **11738** (packing)  
 Column: **2m x 1/4" OD x 2mm ID TightSpec glass**  
 (stock column available)  
 Oven: **170°C (Aroclor 1016, 1242, 1248, 1254, 1260)**  
**140°C (Aroclor 1221, 1232)**  
 Carrier: **nitrogen, 40mL/min**  
 Det.: **ECD**  
 Inj.: **3ng each Aroclor in 3µL isoctane**




For more information, request Bulletin 817 and Application Note 67.



### Figure 35. Ammonia

Packing: **HayeSep P, 60/80 mesh**  
Cat. No.: **10297** (packing)  
Column: 8' x 1/8" ID stainless steel  
Oven: 80°C  
Carrier: helium, 30cc/min  
Det.: T.C.(175ma), 200°C  
Inj.: 0.1µL of NH<sub>4</sub>OH, on-column, 150°C




For more information, request Bulletin 786.

795-0168

### Figure 36. Aromatics and Aliphatics

Packing: **10% TCEP on 100/120 Chromosorb P AW**  
Cat. No.: **12106-U** (packing, 20g/bottle)  
Column: 8' x 1/8" ID stainless steel  
Oven: 80°C  
Carrier: nitrogen, 20mL/min  
Det.: FID  
Inj.: 0.2µL, approx. equal volumes each analyte




For more information,  
request Bulletin 743.

713-0724

### Figure 37. Aromatics in Gasoline

Packing: **35% BC-150 on 100/120 Chromosorb P AW**  
 Cat. No.: **11840-U** (packing)  
 Column: 10' x 1/8" ID stainless steel  
 Oven: 150°C (5 min) to 200°C at 16°C/min  
 Carrier: nitrogen, 20mL/min  
 Det.: FID  
 Inj.: 0.5µL gasoline




For more information, request Bulletin 743.

795-0169

### Figure 38. Aromatic Impurities in Ethylbenzene

Packing: **80/100 Carbopack C/0.1% SP-1000**  
 Cat. No.: **11820** (packing)  
 Column: 6' x 1/8" ID stainless steel  
 Cat. No.: **12495-U** (general configuration stock column;  
 other stock columns available)  
 Oven: 225°C  
 Carrier: nitrogen, 20mL/min at 42psi




For more information, request Bulletin 743.

711-0097

### Figure 39. Aromatics

Packing: **80/100 Carboxpack C/0.1% SP-1000**  
Cat. No.: **11820** (packing)  
Column: 6' x 1/8" ID stainless steel  
Cat. No.: **12495-U** (general configuration stock column;  
other stock columns available)  
Oven: 200°C  
Carrier: nitrogen, 20mL/min at 42psi




For more information, request Bulletin 743.

711-0096

### Figure 40. Aromatics

Packing: **5.0% SP-1200/1.75% Bentone® 34 on 100/120 SUPELCOPORT**  
Cat. No.: **12134** (packing)  
Column: 6' x 1/8" ID stainless steel  
Cat. No.: **12721** (general configuration stock column;  
other stock columns available)  
Oven: 75°C  
Carrier: nitrogen, 20mL/min  
Inj.: 0.1µL




For more information, request Bulletin 743.

711-0094

**Figure 41. Butylbenzenes**

Packing: **80/100 Carbo-pack C/0.1% SP-1000**  
 Cat. No.: **11820** (packing, 15g/bottle)  
 Column: 6' x 1/8" ID stainless steel  
 Cat. No.: **12495-U** (general configuration stock column,  
 other stock columns available)  
 Oven: 225°C  
 Carrier: nitrogen, 20mL/min  
 Det.: FID  
 Inj.: 1µL chloroform  
 (~0.1% each analyte)




For more information, request Bulletin 743.

713-0721

**Figure 42. Butylbenzenes**

Packing: **60/80 Carbo-pack F-TA**  
 Cat. No.: packing available in packed columns only  
 Column: 2m x 2mm ID glass  
 Oven: 225°C  
 Carrier: nitrogen, 20mL/min  
 Det.: FID  
 Inj.: 0.05µL methyl ethyl ketone  
 (5% each analyte)




For more information, request Bulletin 743.

713-0720

### Figure 43. Tri- and Tetramethylbenzenes

Packing: **10% SP-2250 on 100/120 SUPELCOPORT**  
Cat. No.: **12132** (packing)  
Column: **10' x 1/8" ID stainless steel**  
Oven: **110°C**  
Carrier: **nitrogen, 20mL/min**  
Inj.: **0.6µL, solution 0.1% each component in CHCl<sub>3</sub>**




For more information, request Bulletin 743.

711-0100

### Figure 44. Vinyl Chloride

Packing: **80/100 Carbopack C/0.19% picric acid**  
Cat. No.: **11824** (packing)  
Column: **6' x 1/8" ID stainless steel**  
Cat. No.: **13867** (general configuration stock column only)  
Oven: **50°C**  
Carrier: **nitrogen, 40mL/min**  
Det.: **FID, Varian 3700**  
Inj.: **1µL, 1ppm vinyl chloride in carbon disulfide**

1. Vinyl chloride (1ppm)  
2. Carbon disulfide




For more information, request Bulletin 846.

795-0172

**Figure 45. Vinyl Chloride**

Packing: **60/80 Carbowax C/0.2% Carbowax 1500**  
 Cat. No.: **11826** (packing)  
 Column: **6' x 1/8" ID stainless steel**  
 Cat. No.: **13860-U** (general configuration stock column;  
 other stock columns available)  
 Oven: room temperature  
 Carrier: nitrogen, 20mL/min  
 Det.: FID  
 Inj.: 1mL nitrogen containing 1ppm VC




For more information, request  
 Bulletins 846, and 866.

795-0173

**Figure 46. Methylene Chloride and Vinyl Chloride by OSHA Method 51**

Sample: Carboxen™-564 spiked with vinyl acetate (18ng) and desorbed with methanol:methylene chloride (5:95)  
 Sampling Tube: **ORBO™-92**  
 Cat. No.: **20362**  
 Packing: **80/100 Carbowax C/0.2% Carbowax 1500**  
 Cat. No.: **11827** (packing, 15g/bottle)  
 Column: **6' x 2mm ID glass**  
 Oven: **90°C**  
 Carrier: nitrogen, 20mL/min  
 Det.: FID  
 Inj.: 1µL




For more information, request Bulletins 846 and 861.

713-0694,0695

## Figure 47. Chlorinated Volatile Hydrocarbons by US EPA Method TO-1

Sample: 2ng each component  
 Desorption Tube: **Carbotrap™ 300, 4mm ID**  
 Cat. No.: **20379**  
 Thermal Desorption: 330°C, 4 min  
 Transfer Line Temp.: 220°C  
 Carrier (desorption tube to column): helium, 20mL/min  
 Packing: **80/100 Carbo-pack C/0.1% SP-1000**  
 Cat. No.: **11820** (packing, 15g/bottle)  
 Column: 6' x 1/8" stainless steel  
 Cat. No.: **12495-U** (general configuration stock column; other stock columns available)  
 Oven: 35°C (2 min) to 220°C at 16°C/min  
 Carrier: helium, 20mL/min  
 Det.: FID

1. Methanol (carrier solvent)
2. Chloroform
3. 1,2-Dichloropropane
4. 1,1,1-Trichloroethane
5. Carbon tetrachloride
6. 1,2-Dichloroethane
7. Trichloroethylene
8. Bromoform
9. Tetrachloroethylene
10. Impurity
11. Chlorobenzene




For more information, request Bulletins 846 and 861.

713-0702

## Figure 48. Volatile Compounds by US EPA Methods TO-2 and TO-3

Sample: 1-5ng each component  
 Desorption Tube: **Carbotrap 300, 4mm ID**  
 Cat. No.: **20379**  
 Carrier: helium, 20mL/min  
 Packing: **60/80 Carbo-pack B/1% SP-1000**  
 Cat. No.: **11815** (packing, 15g/bottle)  
 Column: 8' x 1/8" stainless steel  
 Cat. No.: **12543** (general configuration stock column; other stock columns available)  
 Oven: 35°C (2 min) to 220°C at 16°C/min  
 Det.: FID

1. Vinyl chloride
2. Dichloromethane
3. Acrylonitrile
4. 1,1-Dichloroethylene
5. 3-Chloropropene
6. Chloroform
7. 1,2-Dichloroethane
8. 1,1,1-Trichloroethane
9. Carbon tetrachloride
10. Benzene
11. Tetrachloroethylene
12. Toluene
13. Chlorobenzene




For more information, request Bulletins 846 and 861.

713-0709

## Figure 49. Chlorinated Solvents

Packing: **20% SP-2100/0.1% Carbowax 1500 on 100/120 SUPELCOPORT**  
 Cat. No.: **11821** (packing)  
 Column: **10' x 1/8" ID stainless steel**  
 Cat. No.: **12718-U** (general configuration stock column;  
 other stock columns available)  
 Oven: **100°C**  
 Carrier: **nitrogen, 20mL/min**  
 Det.: **FID**




1. Methylene chloride
2. 1,1-Dichloroethane
3. Chloroform
4. 1,2-Dichloroethane
5. 1,1,1-Trichloroethane
6. Carbon tetrachloride
7. Trichloroethylene
8. 1,1,2-Trichloroethane
9. Perchloroethylene
10. 1,1,1,2-Tetrachloroethane
11. 1,1,2,2-Tetrachloroethane

For more information,  
 request Bulletin 824.

795-0174

## Figure 50. Creosote

Packing: **10% SP-2100 on 100/120 SUPELCOPORT**  
 Cat. No.: **11989** (packing)  
 Column: **10' x 1/8" ID stainless steel**  
 Cat. No.: **12717** (general configuration stock column;  
 other stock columns available)  
 Oven: **100°C to 300°C at 6°C/min**  
 Carrier: **nitrogen, 20mL/min**  
 Det.: **FID**  
 Inj.: **0.1µL**



1. Naphthalene
2. 2-Methylnaphthalene
3. 1-Methylnaphthalene
4. Biphenyl
5. 2,6-Dimethylnaphthalene
6. 2,3-Dimethylnaphthalene
7. Acenaphthalene
8. Dibenzofuran
9. Fluorene
10. Methyl fluorenes
11. Phenanthrene & Anthracene
12. Carbazole
13. Fluoranthene
14. Pyrene
15. Chrysene


For more information, request Bulletin 743 and Application Note 108.

713-0749



## Figure 51. Creosote

Packing: **10% SP-2250 on 100/120 SUPELCOPORT**  
Cat. No.: **12132 (packing)**  
Column: **10' x 1/8" ID stainless steel**  
Oven: **100°C to 300°C at 6°C/min**  
Carrier: **nitrogen, 20mL/min**  
Det.: **FID**  
Inj.: **0.1µL**




For more information, request Bulletin 743 and Application Note 108.

711-0099

## Figure 52. Creosote with Liquid Crystal Phase

Packing: **1.5% SP-2100/1% BMBT on 100/120 SUPELCOPORT**  
Cat. No.: **custom**  
Column: **10' x 1/8" ID stainless steel**  
Oven: **85°C to 225°C at 6°C/min**  
Carrier: **nitrogen, 20mL/min**  
Det.: **FID**  
Inj.: **0.1µL creosote**




For more information, request Bulletin 743.

795-0175

**Figure 53. Cannabinoids as TMS Derivatives**


Packing: **3% SP-2250 on 100/120 SUPELCOPORT**  
 Cat. No.: **11875 (packing)**  
 Column: **6' x 2mm ID glass**  
 Oven: **235°C**  
 Carrier: **nitrogen, 20mL/min**  
 Det.: **FID, 275°C**



795-0176

**Figure 54. Nicotine**

Packing: **10% Carbowax 20M/2% KOH on 80/100 Chromosorb W AW**  
 Cat. No.: **11805 (packing)**  
 Column: **6' x 2mm ID**  
 Oven: **200°C**  
 Carrier: **20mL/min**  
 Inj.: **1µL**




For more information, request Bulletin 737.

795-0177

## Figure 55. Anticonvulsant Drugs, Underivatized

Packing: **2% SP-2110/1% SP-2510-DA on 100/120 SUPELCOPORT**  
 Cat. No.: **11768** (packing, 20g/bottle)  
 Column: 3' x 2mm ID glass  
 Oven: 120°C to 250°C at 16°C/min, hold 4 min  
 Carrier: helium, 50mL/min  
 Det.: FID  
 Inj.: 2µL, analyte amounts shown on figure




1. Ethosuximide, 1.0µg
2. p-Methyl, n-propyl succinimide (int. std.), 1.0µg
3. Phenylethylmalonamide (PEMA), 0.2µg
4. p-Methyl PEMA (int. std.), 0.2µg
5. Phenobarbital, 0.5µg
6. p-Methyl phenobarbital (int. std.), 0.5µg
7. Carbamazepine, 0.2µg
8. Promidone, 0.2µg
9. p-Methyl primidone (int. std.), 0.2µg
10. Phenytoin, 0.3µg
11. p-Methyl phenytoin (int. std.), 0.5µg

713-0625

## Figure 56. Anticonvulsant Drugs, Underivatized

Packing: **2% SP-2510-DA on 100/120 SUPELCOPORT**  
 Cat. No.: **11776** (packing, 20g/bottle)  
 Column: 3' x 2mm ID glass  
 Oven: 150°C to 265°C at 16°C/min, hold 4 min  
 Carrier: nitrogen, 50mL/min  
 Det.: FID  
 Inj.: 2µL




1. Ethosuximide
2. Methsuximide
3. Phensuximide
4. Barbitol (int. std.)
5. N-Desmethyl methsuximide
6. Mephenytoin
7. N-Desmethyl phensuximide
8. Phenylethylmalonamide
9. 5-Methyl-5-phenythydantoin (int. std.)
10. Phenobarbital
11. Carbamazepine
12. Primidone
13. Diphenythydantoin
14. 5-p-Methylphenyl-5-phenythydantoin (int. std.)

713-0626

**Figure 57. Anticonvulsant Drugs**


Packing: **10% SP-1000 on 80/100 SUPELCOPORT**  
 Cat. No.: **11872** (packing)  
 Column: **3' x 2mm ID glass**  
 Oven: **190°C**  
 Carrier: **helium, 40mL/min**  
 Det.: **FID, 250°C**  
 Inj.: **0.25µg/mL each drug**



713-0675

**Figure 58. Psychotropic Drugs, Antidepressants, Underivatized**


Packing: **3% SP-2250 on 80/100 SUPELCOPORT**  
 Cat. No.: **11767** (packing, 20g/bottle)  
 Column: **6' x 2mm ID glass**  
 Oven: **250°C**  
 Carrier: **helium, 40mL/min**  
 Det.: **NPD**  
 Inj.: **5µL methanol, 10ng each drug**



713-0676

### Figure 59. Barbiturates


Packing: 3% SP-2250 DA on 100/120 SUPELCOPORT  
Cat. No.: 11981 (packing)  
Column: 6' x 2mm ID glass  
Oven: 230°C  
Carrier: nitrogen, 20mL/min  
Inj.: 1µL, 0.1µg/µL each component



795-0178

### Figure 60. Alkaloids

Packing: 3% SP-2250 DB on 100/120 SUPELCOPORT  
Cat. No.: 11983 (packing, 20g/bottle)  
Column: 3' x 2mm ID glass  
Oven: 230°C  
Carrier: nitrogen, 40mL/min  
Det.: FID  
Inj.: 1µL methanol, analyte quantities shown on figure




713-0664

### Figure 61. US EPA Method 601, Purgeable Halocarbons

Packing: **60/80 Carboxpack B/1% SP-1000**  
 Cat. No.: **11815** (packing)  
 Column: 8' x 1/8" OD stainless steel  
 Cat. No.: **12543** (general configuration stock column; other stock columns available)  
 Oven: 45°C (3 min) to 220°C at 8°C/min, 15 min hold  
 Carrier: nitrogen, 40mL/min  
 Det.: Hall®, 85°C  
 Inj.: 1µL of a synthetic mixture of volatile pollutants in dodecane, 200-500ng/µL each component, 200°C

- |   |  |
|---|--|
| 1. Chloromethane                                    | 17. cis-1,3-Dichloropropene, 1,1,2-Trichloroethane, & Chlorodibromomethane |
| 2. Bromomethane                                     | IS 1-Chloro-2-bromopropane (int. std.)                                     |
| 3. Vinyl chloride & Dichlorofluoromethane           | 18. 2-Chloroethyl vinyl ether  |
| 4. Chloroethane                                     | 19. Bromoform  |
| 5. Methylene chloride                               | 20. 1,1,2,2-Tetrachloroethylene & 1,1,2,2-Tetrachloroethane                |
| 6. Trichlorofluoromethane                           | IS 1,4-Dichlorobutane (int. std.)  |
| 7. 1,1-Dichloroethylene                             | 21. Chlorobenzene  |
| IS Bromochloromethane (int. std.)                   |  |
| 8. 1,1-Dichloroethane                               |  |
| 9. trans-1,2-Dichloroethylene                       |  |
| 10. Chloroform                                      |  |
| 11. 1,2-Dichloroethane                              |  |
| 12. 1,1,1-Trichloroethane                           |  |
| 13. Carbon tetrachloride                            |  |
| 14. Bromodichloromethane                            |  |
| 15. 1,2-Dichloropropane & trans-1,3-Dichloropropene |  |
| 16. Trichloroethylene                               |  |




For more information, request Bulletins 775 and 865.

795-0179

### Figure 62. US EPA Method 602, Purgeable Aromatics

Packing: **5% SP-1200/1.75% Bentone 34 on 100/120 SUPELCOPORT**  
 Cat. No.: **12134** (packing)  
 Column: 6' x 1/8" ID stainless steel  
 Cat. No.: **12721** (general configuration stock column; other stock columns available)  
 Oven: 50°C (2 min) to 90°C at 6°C/min and hold  
 Carrier: helium, 36mL/min  
 Det.: FID  
 Inj.: 1µL of 0.2mg/mL standard of each purgeable aromatic in methanol




For more information, request Bulletins 775 and 865.

795-0180

### Figure 63. US EPA Method 603, Acrolein and Acrylonitrile

Packing: **80/100 Chromosorb 101**  
Cat. No.: **20214** (packing)  
Column: **6' x 1/8" OD stainless steel**  
Cat. No.: **12712** (general configuration stock column; other stock columns available)  
Oven: **100°C (5 min) to 140°C at 10°C/min and hold**  
Carrier: **helium, 45mL/min**  
Det.: **FID**  
Inj.: **2µL synthetic standard in water**




For more information, request Bulletins 775 and 865.

795-0181

### Figure 64. US EPA Method 604, Derivatives of Phenols

Packing: **5% SP-2250 on 80/100 SUPELCOPORT**  
Cat. No.: **11737** (packing)  
Column: **1.8m x 2mm ID glass**  
Oven: **160°C (2 min) to 240°C at 8°C/min**  
Carrier: **nitrogen, 30mL/min**  
Det.: **EC**  
Inj.: **0.1µL with 0.2µg/mL of each component in methanol**




For more information, request Bulletins 775 and 865.

795-0182

**Figure 65. US EPA Method 604, Phenols**

Packing: **1% SP-1240-DA on 100/120 SUPELCOPORT**  
 Cat. No.: **11832** (packing, 20g/bottle)  
 Column: 2m x 2mm ID glass  
 Oven: 70°C (2 min) to 200°C at 8°C/min and hold  
 Carrier: helium, 30mL/min  
 Det.: FID, 250°C  
 Inj.: 0.5µL EPA Phenols Mix (Cat. No. 48859)




1. Methanol
2. o-Chlorophenol (0.25µg)
3. 2-Nitrophenol (0.25µg)
4. Phenol (0.25µg)
5. 2,4-Dimethylphenol (0.25µg)
6. 2,4-Dichlorophenol (0.25µg)
7. 2,4,6-Trichlorophenol (0.75µg)
8. 4-Chloro-m-cresol (1.25µg)
9. 2,4-Dinitrophenol (0.75µg)
10. 4,6-Dinitro-p-cresol (1.25µg)
11. Pentachlorophenol (1.25µg)
12. Nitrophenol (1.25µg)

For more information, request Bulletins 775 and 865.

713-1032

**Figure 66. US EPA Method 606, Phthalates**

Packing: **1.5% SP-2250/1.95% SP-2401 on 100/120 SUPELCOPORT**  
 Cat. No.: **11947** (packing)  
 Column: 2m x 4mm ID glass (stock column available)  
 Oven: 180°C and 220°C  
 Carrier: nitrogen, 60mL/min  
 Det.: ECD, 250°C  
 Inj.: 1µL of phthalate mixture at concentrations indicated



1. Dimethyl phthalate (1.0ppm)
2. Diethyl phthalate (0.2ppm)
3. Di-n-butyl phthalate (1.0ppm)
4. Butyl benzyl phthalate (1.0ppm)
5. bis(2-Ethylhexyl) phthalate (0.4ppm)
6. Di-n-octyl phthalate (1.0ppm)


For more information, request Bulletins 775 and 865.

795-0183



### Figure 67. Phthalates

Packing: **5% SP-2100 on 100/120 SUPELCOPORT**  
Cat. No.: **11782-U** (packing)  
Column: **6' x 1/8" stainless steel**  
Oven: **200°C (10 min) to 250°C at 4°C/min**  
Carrier: **nitrogen, 25mL/min**  
Det.: **FID**  
Inj.: **1µL of 75ng/µL each component in carbon disulfide**




For more information, request Bulletins 775 and 865.

795-0184

### Figure 68. US EPA Method 607, Nitrosamines


Packing: **10% Carbowax 20M/2% KOH on 80/100 Chromosorb W AW**  
Cat. No.: **11739** (packing)  
Column: **2m x 4mm ID glass**  
Oven: **110°C and 220°C**  
Carrier: **helium, 40mL/min**  
Det.: **FID**  
Inj.: **2µL of 0.2mg/mL standard in methanol at concentrations indicated**



795-0185

### Figure 69. N-Nitrosamines

Packing: 10% Carbowax 20M/2% KOH on 80/100 Chromosorb W AW  
 Cat. No.: 11739 (packing)  
 Column: 2m x 4mm ID TightSpec glass  
 Oven: 150°C  
 Carrier: helium, 40mL/min  
 Det.: Hewlett-Packard NPD  
 Inj.: 1ng each component in 2µL methylene chloride




For more information, request Bulletin 775.

795-0186

### Figure 70. US EPA Method 608, Pesticides

Packing: 1.5% SP-2250/1.95% SP-2401 on 100/120 SUPELCOPORT  
 Cat. No.: 11947 (packing)  
 Column: 2m x 4mm ID glass (stock column available)  
 Oven: 200°C  
 Carrier: argon:methane (95:5), 6mL/min  
 Det.: EC Ni63  
 Inj.: 0.5µL of 0.02mg/mL of each pesticide




For more information, request Bulletin 775.

795-0187

## Figure 71. US EPA Method 609, Nitroaromatics and Isophorone

Packing: **1.5% SP-2250/1.95% SP-2401 on 100/120 SUPELCOPORT**  
 Cat. No.: **11947** (packing)  
 Column: 2m x 4mm ID glass  
 (stock column available)  
 Oven: 85°C and 145°C  
 Carrier: nitrogen, 44mL/min  
 Det.: FID for nitrobenzene &  
 isophorone  
 EC for dinitrotoluenes  
 Inj.: 2µL of 0.2mg/mL  
 standard on FID  
 0.5µL of 2ng/µL  
 standard on EC

1. Nitrobenzene
2. Isophorone
3. 2,6-Dinitrotoluene
4. 2,4-Dinitrotoluene




For more information, request Bulletin 775.

795-0188

## Figure 72. US EPA Method 610, Polynuclear Aromatic Hydrocarbons

Packing: **3% SP-2250 on 100/120 SUPELCOPORT**  
 Cat. No.: **11744** (packing)  
 Column: 2m x 2mm ID glass  
 Oven: 100°C (4 min) to 280°C at 8°C/min and hold  
 Carrier: nitrogen, 40mL/min  
 Det.: FID  
 Inj.: 0.2µg each component in 1µL methanol/methylene chloride

1. Naphthalene
2. Acenaphthylene
3. Acenaphthene
4. Fluorene
5. Phenanthrene & Anthracene
6. Fluoranthene
7. Pyrene
8. Benzo(a)anthracene
9. Chrysene
10. Benzo(b)fluoranthene & Benzo(k)fluoranthene
11. Benzo(a)pyrene
12. Dibenz(a,h)anthracene & Indeno(1,2,3-cd)pyrene
13. Benzo(ghi)perylene




For more information,  
request Bulletin 775.

795-0189

**Figure 73. US EPA Method 611, Haloethers**

Packing: **3% SP-1000 on 100/120 SUPELCOPORT**  
 Cat. No.: **11746** (packing)  
 Column: 2m x 2mm ID glass  
 Oven: 60°C (2 min) to 230°C at 8°C/min and hold  
 Carrier: helium, 40mL/min  
 Det.: FID  
 Inj.: 1µL of 0.2mg/mL of each haloether in methanol




For more information, request Bulletins 775 and 865.

795-0190

**Figure 74. US EPA Method 612, Chlorinated Hydrocarbons**

Packing: **1.5% OV-1/1.5% OV-225 on 80/100 SUPELCOPORT**  
 Cat. No.: **custom** (packing)  
 Column: 2m x 2mm ID glass  
 Oven: 75°C and 160°C  
 Carrier: argon:methane (95:5), 30mL/min  
 Det.: EC  
 Inj.: 2µL of synthetic mixture (0.1 to 10ng/µL) in hexane



For more information, request  
 Bulletins 775 and 865.

795-0191

## Figure 75. US EPA Method 625, Base/Neutrals

Packing: **3% SP-2250 on 100/120 SUPELCOPORT**

Cat. No.: **11756** (packing)

Column: **2m x 2mm ID glass**

Oven: **50°C (4 min) to 270°C at 8°C/min**


Carrier: **nitrogen, 30mL/min**

Det.: **FID, 300°C**

Inj.: **1µL methylene chloride, 0.125µg or 0.063µg (\*) each analyte**

1. 1,3-Dichlorobenzene
2. 1,4-Dichlorobenzene-d<sub>4</sub>
3. Bis(2-chloroethyl)ether
4. Hexachloroethane & 1,2-Dichlorobenzene
5. Bis(methyl-2-chloroethyl)ether
6. N-Nitroso-di-n-propylamine
7. Nitrobenzene & Hexachlorobutadiene
8. 1,2,4-Trichlorobenzene & Isophorone
9. Naphthalene & Bis(2-chloroethoxy)methane
10. Hexachlorocyclopentadiene
11. 2-Chloronaphthalene
12. Acenaphthylene
13. Acenaphthene
14. Dimethylphthalate
15. 2,6-Dinitrotoluene

16. Fluorene & 4-Chlorophenylphenyl ether
17. 2,4-Dinitrotoluene
18. Diethylphthalate & N-Nitrosodiphenylamine
19. 4-Bromophenylphenyl ether & Hexachlorobenzene
20. Phenanthrene & Anthracene
21. Di-n-butyl phthalate
22. Fluoranthene\*
23. Pyrene\*
24. Benzidine
25. Butyl benzyl phthalate
26. Bis(2-ethylhexyl)phthalate
27. Benzo(a)fluoranthene\*
28. Chrysene\*
29. 3,3'-Dichlorobenzidine
30. Di-n-octyl phthalate
31. Benzo(b)fluoranthene\* & Benzo(k)fluoranthene\*
32. Benzo(a)pyrene
33. Indeno(1,2,3-cd)pyrene\* & Dibenzo(a,h)anthracene\*
34. Benzo(ghi)perylene\*



For more information, request  
Bulletins 775 and 865.

713-1037

## Figure 76. US EPA Method 624, Volatiles

Packing: **60/80 Carboxack B/1% SP-1000**

Cat. No.: **11815** (packing, 15g/bottle)

Column: **8' x 1/8" stainless steel**

Cat. No.: **12543** (general configuration stock column; other stock columns available)

Oven: **45°C (3 min) to 220°C at 8°C/min (hold 15 min)**


Carrier: **helium, 40mL/min**

Det.: **FID, 250°C**

Inj.: **1µL of a synthetic mixture of volatile pollutants in dodecane, 200-500ng/µL each, 200°C**

1. Chloromethane
2. Methanol
3. Bromomethane
4. Vinyl chloride, Dichlorofluoromethane
5. Chloroethane
6. Methylene chloride
7. Trichlorofluoromethane
8. 1,1-Dichloroethylene
9. Bromochloromethane (IS)
10. 1,1-Dichloroethane
11. trans-1,2-Dichloroethylene
12. Chloroform
13. 1,2-Dichloroethane
14. 1,1,1-Trichloroethane
15. Carbon tetrachloride
16. Bromodichloromethane
17. 1,2-Dichloropropane
18. trans-1,3-Dichloropropene
19. Trichloroethylene
20. Benzene
21. cis-1,3-Dichloropropene, 1,1,2-Trichloroethane, Chlorodibromomethane
22. 1-Chloro-2-bromopropene (IS)
23. 2-Chloroethyl vinyl ether
24. Bromoform
25. 1,1,2,2-Tetrachloroethylene, 1,1,2,2-Tetrachloroethane
26. 1,4-Dichlorobutane (IS)
27. Toluene
28. Chlorobenzene
29. Ethylbenzene
30. 1,3-Dichlorobenzene
31. 1,2-Dichlorobenzene
32. 1,4-Dichlorobenzene


For more information, request  
Bulletins 775 and 865.



713-0704

### Figure 77. US EPA Method 625, Phenols

Packing: 1% SP-1240-DA on 100/120 SUPELCOPORT  
 Cat. No.: 11832 (packing)  
 Column: 2m x 2mm ID glass  
 Oven: 70°C (2 min) to 200°C at 8°C/min  
 Carrier: helium, 30mL/min  
 Det.: FID  
 Inj.: 0.5µL of a standard phenol mixture (Cat. No. 48859)




For more information, request Bulletins 775 and 865.

795-0194

### Figure 78. US EPA Method 625, Pesticides

Packing: 3% SP-2250 on 100/120 SUPELCOPORT  
 Cat. No.: 11878 (packing)  
 Column: 2m x 2mm ID glass  
 Oven: 50°C (4 min) to 270°C at 8°C/min  
 Carrier: nitrogen, 30mL/min  
 Det.: FID  
 Inj.: 0.1µL of a synthetic mixture of pesticides (1ng/µL each)




For more information, request Bulletins 775 and 865.

795-0193

## Figure 79. Ethylene Oxide and Ethylene Oxide Residues

Packing: **80/100 CarboPack C/0.8% THEED**  
Cat. No.: **11880-U** (packing)  
Column: 1m x 2mm ID glass  
(stock column available)  
Oven: 115°C  
Carrier: nitrogen, 20mL/min  
Det.: FID  
Inj.: 1µL water containing  
50ppm each component




For more information, request Application Note 109.

713-0732

## Figure 80. Trace Ethylene Oxide in Nitrogen

Packing: **HayeSep D**  
Cat. No.: **10292** (packing)  
Column: 10' x 1/8" stainless steel  
Oven: 130°C  
Carrier: helium, 30mL/min  
Det.: FID, 140°C  
Inj.: 250mL Valco valve, 100°C




For more information, request Application Note 109.

795-0195

**Figure 81. Ethylene Oxide by US Army Method**

Sample: HBr-coated Carboxen-564 spiked with 2-bromoethanol (1.2ng) and desorbed with toluene:acetonitrile (50:50)  
 Sampling Tube: **ORBO-78**  
 Cat. No.: **20355**  
 Packing: **10% DEGS on 80/100 SUPELCOPORT**  
 Cat. No.: **11999** (packing, 20g/bottle)  
 Column: 12' x 1/8" stainless steel  
 Oven: 155°C  
 Carrier: nitrogen, 30mL/min  
 Det.: ECD




For more information, request Application Note 109.

713-0696, 0697

**Figure 82. Fatty Acid Methyl Esters of Rapeseed Oil Mixture**

Packing: **5% DEGS-PS on 100/120 SUPELCOPORT**  
 Cat. No.: **11870-U** (packing, 20g/bottle)  
 Column: 6' x 1/8" stainless steel  
 Oven: 200°C  
 Carrier: nitrogen, 20mL/min  
 Det.: FID  
 Inj.: 0.5µL RM-3 Mix (Cat. No. O7256-1AMP), 10µg FAMES/µL in chloroform




For more information, request Bulletin 856.

713-0973



### Figure 83. Fatty Acid Methyl Esters of Rapeseed Oil Mixture

Packing: **GP 10% DEGS-PS on 80/100 SUPELCOPORT**  
Cat. No.: **11999** (packing, 20g/bottle)  
Column: 6' x 1/8" stainless steel  
Oven: 200°C  
Carrier: nitrogen, 20mL/min  
Det.: FID  
Inj.: 0.5µL RM-3 Mix (Cat. No. O7256-1AMP), 10µg FAMES/µL in chloroform




For more information, request Bulletin 856.

713-1072

### Figure 84. C18:3, C20:0, and C20:1 Methyl Esters

Packing: **GP 3% SP-2310/2% SP-2300 on 100/120 Chromosorb W AW**  
Cat. No.: **11833** (packing)  
Column: 6' x 1/8" stainless steel  
Oven: 190°C (2 min) to 220°C at 2°C/min  
Carrier: nitrogen, 20mL/min  
Det.: FID  
Inj.: 5µg mixed esters in 0.5µL methylene chloride




For more information, request Bulletin 856.

795-0260

### Figure 85. Methyl Elaidate and Oleate FAMES

Packing: **15% OV-275 on 100/120 Chromosorb P AW DMCS**  
 Cat. No.: **11844-U** (packing, 20g/bottle)  
 Column: 20' x 1/8" stainless steel  
 Oven: 220°C  
 Carrier: nitrogen, 10mL/min  
 Det.: FID  
 Inj.: 0.5µL isoctane containing 5µg each analyte




For more information, request Bulletin 856.

713-0974

### Figure 86. Whole Milk Free Acids

Packing: **10% SP-216-PS on 100/120 SUPELCOPORT**  
 Cat. No.: **11879** (packing)  
 Column: 3' x 2mm ID glass  
 Oven: 130°C to 200°C at 15°C/min  
 Carrier: nitrogen, 20mL/min




For more information,  
 request Bulletin 856.

795-0197

## Figure 87. C18:3, C20:0, and C20:1 Methyl Esters

Packing: **GP 3% SP-2310/2% SP-2300 on 100/120 Chromosorb W AW**  
Cat. No.: **11833** (packing)  
Column: **6' x 1/8" stainless steel**  
Oven: **190°C (2 min) to 220°C at 2°C/min**  
Carrier: **nitrogen, 20mL/min**  
Det.: **FID**  
Inj.: **5µg mixed esters in 0.5µL methylene chloride**




For more information, request Bulletin 856.

795-0198

## Figure 88. Polyunsaturated FAMES

Packing: **10% SP-2330 on 100/120 Chromosorb W AW**  
Cat. No.: **11851** (packing)  
Column: **6' x 1/8" stainless steel**  
Cat. No.: **13775** (general configuration stock column; other stock columns available)  
Oven: **210°C**  
Inj.: **15µg mixed esters in 1.5µL iso-octane**




For more information, request Bulletin 856.

795-0199

**Figure 89. Reference Standard: RM-3 Peanut Oil Type**

Packing: 10% SP-2340 on 100/120 Chromosorb W AW  
 Cat. No.: 11852 (packing)  
 Column: 2 meter x 2mm ID glass  
 Oven: 195°C  
 Carrier: nitrogen, 20mL/min  
 Inj.: 0.5µL, 10mg/cc isooctane




For more information, request Bulletin 856.

795-0200

**Figure 90. Grapefruit Oil**


Packing: 5% SP-2100/0.1% SP-401 on 100/120 SUPELCOPORT  
 Cat. No.: 11839 (packing)  
 Column: 6' x 2mm ID glass  
 Oven: 75°C (2 min) to 175°C at 4°C/min  
 Carrier: nitrogen, 20mL/min  
 Det.: FID  
 Inj.: 0.2µL



795-0201

## Figure 91. Lemon Oil

Packing: **5% SP-2100/0.1% SP-401 on 100/120 SUPELCOPORT**  
 Cat. No.: **11839** (packing)  
 Column: 6' x 2mm ID glass  
 Oven: 75°C (2 min) to 175°C at 4°C/min  
 Carrier: nitrogen, 20mL/min  
 Det.: FID  
 Inj.: 0.2µL




For more information, request Bulletin 856.

795-0203

## Figure 92. Orange Oil

Packing: **5% SP-2100/0.1% SP-401 on 100/120 SUPELCOPORT**  
 Cat. No.: **11839** (packing)  
 Column: 6' x 2mm ID glass  
 Oven: 75°C (2 min) to 175°C at 4°C/min  
 Carrier: nitrogen, 20mL/min  
 Det.: FID  
 Inj.: 0.2µL




795-0202

**Figure 93. Freon® and Chlorinated Hydrocarbons**

Packing: **80/100 Carboxpack B/1% SP-1000**  
 Cat. No.: **11815** (packing)  
 Column: 6' x 1/8" stainless steel  
 Cat. No.: **12485** (general configuration stock column;  
 other stock columns available)  
 Oven: 70°C  
 Carrier: nitrogen, 20mL/min  
 Det.: FID

1. Methyl chloride
2. Methyl bromide
3. Vinyl chloride & Freon 12
4. Chloroethane




For more information, request Bulletin 786.

795-0204

**Figure 94. Freon Gases**

Packing: **60/80 Carboxpack C/0.2% Carbowax 1500**  
 Cat. No.: **11826** (packing)  
 Column: 6' x 1/8" stainless steel  
 Cat. No.: **13860-U** (general configuration stock column;  
 other stock columns available)  
 Oven: 65°C to 150°C at 8°C/min  
 Carrier: nitrogen, 20mL/min  
 Det.: FID

1. Freon 23
2. Freon 13
3. Freon 22
4. Freon 115
5. Freon 12
6. Freon 21
7. Freon 114
8. Freon 11
9. Freon 113




For more information, request Bulletin 786.

795-0205

## Figure 95. Freon Gases

Packing: **60/80 Carboxpack B/5% Fluorcol®**  
Column: **10' x 1/8" SP alloy**  
Cat. No.: **12425**  
(general configuration stock column only)  
Oven: **50°C**  
Carrier: **nitrogen, 30mL/min**  
Det.: **FID**  
Inj.: **10µL synthetic mixture (appr. 20% each compound)**




For more information, request Bulletin 786.  
Figure provided by J.L. Glajch and W.G. Schindel, I.I. du Pont de Nemours and Co., Inc.,  
Experimental Station, Wilmington, DE, USA.

713-0705

## Figure 96. Formaldehyde

Packing: **100/120 HayeSep T**  
Cat. No.: **10311** (packing, 75cc/bottle)  
Column: **5' x 1/8" stainless steel**  
Oven: **132°C**  
Carrier: **helium, 30cc/min**  
Det.: **thermal conductivity (175mA)**  
Inj.: **0.2µL (165°C)**


1. Air
2. Formaldehyde
3. Water
4. Methanol



713-0690

### Figure 97. Formalin Solution

Packing: **45/60 Carboxen-1000**  
 Column: **2' x 1/8" stainless steel**  
 Cat. No.: **12370-U** (general configuration stock column; other stock columns available)  
 Oven: 150°C  
 Carrier: helium, 30mL/min  
 Det.: TCD  
 Inj.: 0.2µL formalin solution




795-0206

### Figure 98. Ethylene and Scott® Mix 216

Packing: **100/120 HayeSep DB**  
 Cat. No.: **custom**  
 Column: 30' x 1/8" stainless steel  
 Oven: 120°C  
 Carrier: helium, 30mL/min  
 Det.: P.E. 900 T.C., 300ma  
 Inj.: 250mL Valco valve, 120°C

1. Nitrogen
2. Carbon monoxide
3. Methane
4. Carbon dioxide
5. Acetylene (0.24ppm)
6. Ethylene
7. Water
8. Ethane




For more information, request Application Note 10.

795-0207



## Figure 99. Carbon Dioxide and C2 Hydrocarbons

Packing: **HayeSep T**  
Cat. No.: **10311** (packing)  
Column: 5' x 1/8" stainless steel  
Oven: 32°C  
Carrier: helium, 30mL/min  
Det.: 180°C  
Inj.: 50mL




For more information, request Application Note 10.

795-0208

## Figure 100. Mapp Gases, C1-C3 Unsaturated Hydrocarbons and CO<sub>2</sub>

Packing: **HayeSep R**  
Cat. No.: **10304** (packing)  
Column: 10' x 1/8" stainless steel  
Oven: 80°C  
Carrier: helium, 30mL/min  
Inj.: 15mL

1. Air
2. Methane
3. Carbon dioxide
4. Ethane
5. Propylene
6. Propane
7. Propadiene
8. Methyl acetylene




For more information, request Application Note 10.

795-0209

**Figure 101. Scott Mix 216, C1-C2, CO, and CO<sub>2</sub>**

Packing: **HayeSep D**  
 Cat. No.: **10293** (packing)  
 Column: 20' x 1/8" Ni  
 Oven: 40°C (2 min) to 110°C at 24°C/min  
 Carrier: helium, 30mL/min  
 Det.: P.E. 900 T.C., 225ma, 140°C  
 Inj.: 100mL, Valco valve, ambient

1. Nitrogen (balance)
2. Carbon monoxide, 1%
3. Methane, 1%
4. Carbon dioxide, 1%
5. Acetylene, 1%
6. Ethylene, 1%
7. Ethane, 1%




For more information, request Application Note 10.

795-0210

**Figure 102. C1-C2 Hydrocarbons**

Packing: **100/120 HayeSep D**  
 Cat. No.: **10293** (packing)  
 Column: 10' x 1/8" stainless steel  
 Oven: 80°C  
 Carrier: helium, 35mL/min  
 Det.: FID (3700 Varian)

1. Methane, 1%
2. Acetylene, 1%
3. Ethylene, 1%
4. Ethane, 1%




For more information, request Application Note 10.

795-0211

### Figure 103. Trace Acetylene in Ethylene

Packing: **45/60 Carboxen-1000**  
Column: **5' x 1/8" stainless steel**  
Cat. No.: **12380** (general configuration stock column; other stock columns available)  
Oven: 160°C  
Carrier: helium, 30mL/min  
Det.: FID  
Inj.: 0.5mL ethylene with trace impurities




For more information, request Application Note 10.

795-0212

### Figure 104. Trace Acetylene in Ethylene

Packing: **60/80 Carbosieve S-II custom**  
Cat. No.:  
Column: **5' x 1/8" stainless steel**  
Oven: 125°C  
Carrier: nitrogen, 30mL/min  
Det.: FID  
Inj.: 1mL, 100ppm acetylene in ethylene




For more information, request Application Note 10.

795-0213

**Figure 105. C1-C3 Hydrocarbons**

Packing: **45/60 Carboxen-1000**  
 Column: **2' x 1/8" stainless steel**  
 Cat. No.: **12370-U** (general configuration stock column; other stock columns available)  
 Oven: 165°C (1 min) to 220°C at 16°C/min  
 Carrier: helium, 60mL/min  
 Det.: FID  
 Inj.: 0.6mL Scott Gas Mix (Cat. No. 2-3470), 15ppm each component, appr. 10ng each component on column




For more information,  
 request Bulletin 743 and Application Note 10.

795-0214

**Figure 106. C1-C3 Hydrocarbons**

Packing: **60/80 Carbosieve G**  
 Cat. No.: **10198** (packing)  
 Column: **5' x 1/8" stainless steel**  
 Oven: 145°C to 195°C at 6°C/min  
 Carrier: nitrogen, 50mL/min  
 Det.: FID  
 Inj.: 1mL, appr. 15ppm each component in nitrogen




For more information,  
 request Bulletin 743 and Application Note 10.

795-0215

## Figure 107. Methanol in Propylene/Propane

Packing: **100/120 HaysSep T**  
Cat. No.: **10311** (packing)  
Column: **5' x 1/8" stainless steel**  
Oven: **120°C**  
Carrier: **helium, 30mL/min**  
Det.: **P.E. 900 T.C., 225ma, 150°C**  
Inj.: **132°C**

1. Air, Methane
2. Carbon dioxide
3. Ethane
4. Propylene, Propane
5. Butane
6. Butane
7. Methanol




For more information, request Bulletin 743  
and Application Note 10.

795-0216

## Figure 108. C3-C5 Saturated and Unsaturated Hydrocarbons

Packing: **80/100 Carboxack C/0.19% picric acid**  
Cat. No.: **11824** (packing)  
Column: **2m x 1/8" stainless steel**  
Cat. No.: **13867** (general configuration stock column only)  
Oven: **35°C to 70°C at 6°C/min**  
Carrier: **nitrogen, 30mL/min**  
Det.: **FID**

1. Propane
2. Propylene
3. Isobutane
4. 1-Butene
5. n-Butane
6. Isobutylene
7. cis-Butene-2
8. trans-Butene-2
9. 1,3-Butadiene
10. Isopentane
11. n-Pentane




For more information, request Bulletin 743  
and Application Note 142.

713-0757

### Figure 109. C3-C5 Saturated and Unsaturated Hydrocarbons

Packing: 10% EDO-1 on 100/120 Chromosorb P AW  
 Cat. No.: 11967 (packing)  
 Column: 20' x 1/8" stainless steel  
 Oven: 0°C  
 Carrier: nitrogen, 40mL/min  
 Det.: FID  
 Inj.: ASTM Hydrocarbon Mix  
 Section L, Blend 6




1. Propane
2. Propylene
3. Isobutane
4. n-Butane
5. 1-Butene
6. Isobutylene
7. trans-2-Butene
8. Isopentane
9. cis-2-Butene
10. n-Pentane
11. 1,3-Butadiene

For more information, request  
 Bulletin 743 and Application Note 142.

795-0217

### Figure 110. C1-C5 Saturated and Unsaturated Hydrocarbons

Packing: 15% EDO-1 on 60/80 Chromosorb P AW  
 Cat. No.: custom  
 Column: 20' x 3/16" stainless steel  
 Carrier: nitrogen, 40mL/min



1. Methane
2. Ethane
3. Ethylene
4. Propane
5. Propylene
6. Isobutane
7. n-Butane
8. 1-Butene
9. Isobutylene
10. trans-2-Butene
11. Isopentane
12. cis-2-Butene
13. n-Pentane
14. 3-Methyl-1-butene
15. 1,3-Butadiene
16. 1-Pentene
17. 2-Methyl-1-butene
18. trans-2-Pentene
19. cis-2-Pentene
20. 2-Methyl-2-butene


For more information, request Bulletin 743 and Application Note 10 and 142.

795-0218

### Figure 111. Trace C2-C6 Hydrocarbons

Packing: **100/120 HaysSep D**  
 Cat. No.: **10293** (packing)  
 Column: 10' x 1/8" stainless steel  
 Oven: 120°C to 200°C at 24°C/min  
 Carrier: helium, 35mL/min  
 Det.: FID (3700 Varian)

1.	Acetylene	16ppm
2.	Ethylene	15ppm
3.	Propylene	14.3ppm
4.	1-Butene	15ppm
5.	1-Pentene	14.75ppm
6.	1-Hexane	16ppm




For more information, request Bulletin 743 and Application Note 10 and 142.

795-0219

### Figure 112. C1-C5 Hydrocarbons

Packing: **100/120 HaysSep D**  
 Cat. No.: **10293** (packing)  
 Column: 10' x 1/8" stainless steel  
 Oven: 120°C to 200°C at 20°C/min  
 Carrier: helium, 35mL/min  
 Det.: FID (3700 Varian)

1.	Methane	0.1894%
2.	Ethane	0.0965%
3.	Propane	0.0989%
4.	Isobutane	0.1019%
5.	Butane	0.1019%
6.	n-Pentane	0.2002%




For more information, request Bulletin 743 and Application Note 10 and 142.

795-0220

**Figure 113. C1-C6 Saturated and Unsaturated Hydrocarbons**

Packing: **23% SP-1700 on 80/100 Chromosorb P AW**  
 Cat. No.: **11865** (packing)  
 Column: **30' x 1/8" stainless steel**  
 Cat. No.: **12809-U** (general configuration stock column only)  
 Oven: **70°C**  
 Carrier: **helium, 25mL/min**  
 Det.: **FID, 150°C**  
 Inj.: **0.6µL ASTM Section L, Blend 6,  
 plus C5 hydrocarbons, 100°C**




For more information, request Bulletin 743, and Application Note 142.

713-0758

**Figure 114. Propane and Propylene at High Concentrations**

Packing: **23% SP-1700 on 80/100 Chromosorb P AW**  
 Cat. No.: **11865** (packing)  
 Column: **30' x 1/8" stainless steel**  
 Cat. No.: **12809-U** (general configuration stock column only)  
 Oven: **70°C**  
 Carrier: **helium, 25mL/min**  
 Det.: **FID, 150°C**  
 Inj.: **0.5mL liquid propane gas, 100°C**

1. Pressure disturbance	
2. Ethane	0.42%
3. Propane	96.75%
4. Propylene	1.72%
5. Isobutane	0.93%
6. n-Butane	0.17%
7. 1-Butene, Isobutylene	trace
8. Isopentane	trace



For more information, request Bulletin 743 and Application Note 142.


795-0221



## Figure 115. Industrial Butane III

Packing: **23% SP-1700 on 80/100 Chromosorb P AW**  
Cat. No.: **11865** (packing)  
Column: **30' x 1/8" stainless steel**  
Cat. No.: **12809-U** (general configuration  
stock column only)  
Oven: **70°C**  
Carrier: **helium, 25mL/min**  
Det.: **FID, 110°C**  
Inj.: **0.75mL industrial butane III**

1. Pressure disturbance
2. Isobutane 2.02%
3. n-Butane 97.92%
4. Isopentane 0.05%
5. n-Pentane 0.01%




For more information, request Bulletin 743 and Application Note 142.

795-0222

## Figure 116. Permanent Gases in Air

Packing: **60/80 Molecular Sieve 5A**  
Cat. No.: **20302** (packing)  
Column: **3' x 1/8" stainless steel**  
Oven: **60°C**  
Carrier: **helium, 20mL/min**  
Det.: **thermal conductivity, 150ma**  
Inj.: **0.25cc, 1% concentrations of  
H<sub>2</sub>, O<sub>2</sub>, CH<sub>4</sub>, and CO in N<sub>2</sub>**




For more information, request Bulletin 743 and Application Note 142.

713-0739

### Figure 117. Permanent Gases in Air

Packing: **60/80 Chromosorb 102**  
 Cat. No.: **20200-U** (packing)  
 Column: 6' x 1/8" stainless steel  
 Oven: 60°C  
 Carrier: helium, 20mL/min  
 Det.: thermal conductivity, 150ma  
 Inj.: 0.2cc, 1% concentrations of CH<sub>4</sub>, and CO<sub>2</sub> from N<sub>2</sub>




For more information, request Application Note 10.

795-0223

### Figure 118. Scott Mix 237


Packing: **100/120 HayeSep D**  
 Cat. No.: **10293** (packing)  
 Column: 20' x 1/8" Ni  
 Oven: 25°C  
 Carrier: helium, 30cc/min  
 Det.: P.E. 900 T.C., 225ma, 140°C  
 Inj.: 50mL vapor, Valco valve, ambient



795-0224

### Figure 119. Nitrogen in Argon

Packing: **100/120 HayeSep D**  
Cat. No.: **10293** (packing)  
Column: 30' x 1/8" stainless steel  
Oven: 25°C  
Carrier: helium, 30cc/min  
Det.: P.E. 900 T.C, 300ma, 140°C  
Inj.: 250mL Valco valve, 25°C




795-0225

### Figure 120. Permanent Gases

Packing: **80/100 HayeSep A**  
Cat. No.: **10283** (packing)  
Column: 36' x 1/8"  
Oven: 25°C  
Carrier: helium, 23cc/min  
Det.: P.E. 900 T.C, 225ma  
Inj.: 25mL air plus hydrogen and carbon monoxide

1. Hydrogen 5%  
2. Nitrogen 48.5%  
3. Oxygen 13%  
4. Argon 0.5%  
5. Carbon monoxide 33%




795-0226

**Figure 121. Permanent Gases, Methane, and C2 Hydrocarbons**

Packing: **100/120 Carbosieve S-II**  
 Column: **10' x 1/8" stainless steel**  
 Cat. No.: **13821** (general configuration stock column;  
 other stock columns available)  
 Oven: 35°C (7 min) to 225°C  
 at 32°C/min  
 Carrier: helium, 30mL/min  
 Det.: thermal conductivity,  
 150ma  
 Inj.: 0.6mL

- |                    |       |
|--------------------|-------|
| 1. Hydrogen        | 2.5%  |
| 2. Oxygen          | 2.5%  |
| 3. Nitrogen        | 91.4% |
| 4. Carbon monoxide | 0.6%  |
| 5. Methane         | 0.6%  |
| 6. Carbon dioxide  | 0.6%  |
| 7. Acetylene       | 0.6%  |
| 8. Ethylene        | 0.6%  |
| 9. Ethane          | 0.6%  |



For more information, request Application Note 10.


795-0227

**Figure 122. Permanent Gases and C1-C2 Hydrocarbons**

Packing: **60/80 Carboxen-1000**  
 Column: **15' x 1/8" stainless steel**  
 Cat. No.: **12390-U** (general configuration stock column;  
 other stock columns available)  
 Oven: 35°C (5 min) to 225°C at 20°C/min  
 Carrier: helium, 30mL/min  
 Det.: TCD  
 Inj.: 0.6mL Scott Gas Mix (Cat. No. 23437) with oxygen added, 1% each  
 component in nitrogen, on column

- |                    |
|--------------------|
| 1. Oxygen          |
| 2. Nitrogen        |
| 3. Carbon monoxide |
| 4. Methane         |
| 5. Carbon dioxide  |
| 6. Acetylene       |
| 7. Ethylene        |
| 8. Ethane          |


For more information,  
 request Application Note 10.



795-0228

### Figure 123. Pesticides

Packing: **GP 1.5% SP-2250/1.95% SP-2401 on 100/120 SUPELCOPORT**  
 Cat. No.: **11947** (packing)  
 Column: **6' x 4mm ID glass**  
 (stock column available)  
 Oven: **200°C**  
 Carrier: **argon:methane (95:5), 70mL/min**  
 Det.: **EC NI63**  
 Inj.: **5µL CP Mix**




1. □-BHC
2. □BHC (Lindane)
3. □-BHC
4. Heptachlor
5. Aldrin
6. Heptachlor epoxide
7. 4,4'-DDE
8. Dieldrin
9. 2,4-DDD (TDE)
10. Endrin
11. 2,4'-DDT
12. 4,4'-DDD (TDE)
13. 4,4'-DDT

For more information, request Bulletins 758 and 775.

795-0229

### Figure 124. Pesticides

Packing: **GP 4% SE-30/6% SP-2401 on 100/120 SUPELCOPORT**  
 Cat. No.: **11948** (packing)  
 Column: **6' x 4mm ID glass**  
 Oven: **200°C**  
 Carrier: **argon:methane (95:5), 70mL/min**  
 Det.: **EC NI63**  
 Inj.: **5µL CP Mix**




1. □-BHC
2. □BHC (Lindane)
3. □-BHC
4. Heptachlor
5. Aldrin
6. Heptachlor epoxide
7. 4,4'-DDE
8. 2,4-DDD (TDE)
9. Dieldrin
10. 2,4'-DDT
11. Endrin
12. 4,4'-DDD (TDE)
13. 4,4'-DDT

For more information, request Bulletins 758 and 775.

795-0230

**Figure 125. Phenols**

Packing: **0.1% SP-1000 on 80/100 Carboxpack C**  
 Cat. No.: **11820** (packing)  
 Column: **6' x 2mm ID glass** (stock column available)  
 Oven: **225°C**  
 Carrier: **nitrogen, 20mL/min**  
 Det.: **FID**  
 Inj.: **0.5µL approx. 0.1% each in MEK**




For more information, request Bulletin 775.

795-0231

**Figure 126. Mixed Phenols**

Packing: **3% SP-1000 on 100/120 SUPELCOPORT**  
 Cat. No.: **11738** (packing)  
 Column: **6' x 2mm ID glass**  
 Oven: **130°C**  
 Carrier: **nitrogen, 20mL/min**  
 Det.: **FID, 175°C**  
 Inj.: **0.05µL**




For more information, request Bulletin 775.

795-0232

## Figure 127. Phenols

Packing: **60/80 Carboxpack F-TA**  
Column: 2m x 2mm ID glass  
Oven: 220°C  
Carrier: nitrogen, 15mL/min  
Det.: FID  
Inj.: 0.2µL approx. 0.1% each in MEK

- S Solvent  
1. Phenol  
2. o-Cresol  
3. m-Cresol  
4. p-Cresol  
5. o-Ethylphenol  
6. p-Ethylphenol  
7. p-Ethylphenol  
8. 2,6-Dimethylphenol  
9. 2,4-Dimethylphenol  
10. 2,5-Dimethylphenol  
11. 2,3-Dimethylphenol  
11. 3,4-Dimethylphenol




For more information,  
request Bulletin 775.

794-0139

## Figure 128. Butyl Phenols

Packing: **10% SP-1000 on 80/100 SUPELCOPORT**  
Cat. No.: **11872** (packing)  
Column: 6' x 1/8" stainless steel  
Oven: 180°C  
Carrier: nitrogen, 20mL/min  
Det.: FID  
Inj.: 0.5µL with 0.1% of each  
component in chloroform


1. 2,6-Di-tert-butylphenol  
2. 2,6-Di-tert-butyl-p-cresol



795-0233

**Figure 129. Phenols (Anisoles)**


Packing: **5% SP-1200/1.75% Bentone 34**  
 Cat. No.: **12134** (packing)  
 Column: **6' x 1/8" stainless steel**  
 Cat. No.: **12721** (general configuration stock column;  
 other stock columns available)  
 Oven: **150°C**  
 Carrier: **nitrogen, 20mL/min**



795-0234

**Figure 130. Pheromones**

Packing: **15% SP-2340/ on 100/120 Chromosorb P AW-DMCS**  
 Cat. No.: **custom**  
 Column: **20' x 1/8" stainless steel**  
 Oven: **225°C**  
 Carrier: **nitrogen, 10mL/min at 50psi**  
 Inj.: **E and Z-9-tetradecenyl acetate**




795-0238



### Figure 131. SIMDIS Calibration Blend (ASTM D3710)

Column: **Petrocol™ A, 20" x 1/8" stainless steel**  
 Cat. No.: **12445** (general configuration stock column only)  
 Oven: -20°C to 200°C at 20°C/min  
 Carrier: helium, 25mL/min  
 Det.: FID, 300°C  
 Inj.: 0.1µL of D3710 Qualitative Calibration Mix (Cat. No. 48884), 300°C




For more information, request Bulletin 864.

713-0772

### Figure 132. Gasoline by ASTM D3710

Column: **Petrocol A, 20" x 1/8" stainless steel**  
 Cat. No.: **12445** (general configuration stock column only)  
 Oven: -20°C to 250°C at 15°C/min  
 Carrier: helium, 25mL/min  
 Det.: FID, 300°C  
 Inj.: 0.1µL gasoline, 300°C




For more information, request Bulletin 864.

713-0771

**Figure 133. SIMDIS Calibration Blend (ASTM D2887)**

Column: **Petrocol B, 20" x 1/8" stainless steel**  
 Cat. No.: **12449** (general configuration stock column only)  
 Oven: **-25°C to 350°C at 15°C/min**  
 Carrier: **helium, 30mL/min**  
 Det.: **FID**  
 Inj.: **0.1µL Quantitative Calibration Mix (Cat. No. 48882)**




For more information, request Bulletin 864.

712-0427

**Figure 134. Reference Gas Oil Sample (ASTM D2887)**

Column: **Petrocol B, 20" x 1/8" stainless steel**  
 Cat. No.: **12449** (general configuration stock column only)  
 Oven: **-25°C to 350°C at 15°C/min**  
 Carrier: **helium, 30mL/min**  
 Det.: **FID**  
 Inj.: **0.1µL Reference Gas Oil Sample (Cat. No. 48873)**




For more information, request Bulletin 864.

713-0769

### Figure 135. PS-18-44D Wax Standard Mix (ASTM D2887)

Column: **Petrocol B, 20" x 1/8" stainless steel**  
Cat. No.: **12449** (general configuration stock column only)  
Oven: 80°C to 350°C at 10°C/min  
Carrier: nitrogen, 30mL/min  
Det.: FID, 380°C  
Inj.: 1µL PS-18-44D Mix (Cat. No. 48928), 350°C




For more information, request Bulletin 864.

794-0408

### Figure 136. PS-18-60D Wax Standard Mix (ASTM D2887)

Column: **Petrocol B, 20" x 1/8" stainless steel**  
Cat. No.: **12449** (general configuration stock column only)  
Oven: 50°C (1 min) to 350°C at 35°C/min  
Carrier: nitrogen, 60mL/min  
Det.: FID, 380°C  
Inj.: 0.5µL PS-12-60D Mix (Cat. No. 48932), 350°C




794-0165

**Figure 137. Crude Oil Qualitative Calibration Standard**

Column: **Petrocol C, 20" x 1/8" stainless steel**  
 Cat. No.: **12455** (general configuration stock column only)  
 Oven: -25°C to 350°C at 15°C/min  
 Carrier: helium, 30mL/min  
 Det.: FID, 360°C  
 Inj.: 1µL Crude Oil Qualitative Calibration Standard, 360°C

1. Propane (nC3), 1.0%
  2. n-Butane (nC4)
  3. n-Pentane (nC5), 3.4%
  4. n-Hexane (nC6), 5.2%
  5. n-Heptane (nC7)
  6. n-Octane (nC8), 6.2%
  7. n-Nonane (nC9), 4.7%
  8. n-Decane (nC10)
  9. n-Undecane (nC11)
  10. n-Dodecane (nC12)
  11. n-Tridecane (nC13)
  12. n-Tetradecane (nC14)
  13. n-Pentadecane (nC15)
  14. n-Hexadecane (nC16)
  15. n-Heptadecane (nC17)
  16. n-Octadecane (nC18)
  17. n-Eicosane (nC20)
  18. n-Tetracosane (nC24)
  19. n-Octacosane (nC28)
  20. n-Dotriacontane (nC32)
  21. n-Hexatriacontane (nC36)
  22. n-Tetracontane (nC40)
  23. n-Tetratetracontane (nC44)
- 4.4% each by weight, except as noted.




For more information, request Bulletin 864.

711-0072

**Figure 138. Crude Oil — North Sea, Oseberg**

Column: **Petrocol C, 20" x 1/8" stainless steel**  
 Cat. No.: **12445** (general configuration stock column only)  
 Oven: -25°C to 350°C at 15°C/min  
 Carrier: helium, 30mL/min  
 Det.: FID, 360°C  
 Inj.: 1µL oil:CS mix (1:1), 360°C

1. Tetradecane (nC14)
2. Pentadecane (nC15)
3. Hexadecane (nC16)
4. Heptadecane (nC17)




For more information, request Bulletin 864.

711-0076,0077

### Figure 139. SIMDIS, Regular Gasoline

Packing: **10% SP-2100 on 100/120 SUPELCOPORT**  
Cat. No.: **12140** (packing)  
Column: 10' x 1/8" stainless steel  
Cat. No. **13766** (general configuration stock column;  
other stock columns available)  
Oven: 75°C to 200°C at 8°C/min  
Carrier: 20mL/min  
Det.: FID  
Inj.: 0.1µL




For more information,  
request Bulletin 864.

711-0088

### Figure 140. Naphtha

Packing: **10% SP-2100 on 100/120 SUPELCOPORT**  
Cat. No.: **12140**  
Column: 10' x 1/8" stainless steel  
Cat. No.: **13766** (general configuration stock column;  
other stock columns available)  
Oven: 75°C to 200°C at 8°C/min  
Carrier: 20mL/min  
Det.: FID  
Inj.: 0.1µL




For more information, request Bulletin 864.

711-0087

### Figure 141. Kerosene (Jet A)

Packing: **10% SP-2100 on 100/120 SUPELCOPORT**  
 Cat. No.: **12140** (packing)  
 Cat. No.: **13766** (general configuration stock column;  
 other stock columns available)  
 Column: 10' x 1/8" stainless steel  
 (stock column available)  
 Oven: 75°C to 200°C at 8°C/min  
 Carrier: 20mL/min  
 Det.: FID  
 Inj.: 0.1µL

1. C6
2. C8
3. C10
4. C12
5. C14
6. C16




For more information, request Bulletin 864.

711-0089

### Figure 142. Microcrystalline Wax

Packing: **1% Dexsil® 300 on 100/120 SUPELCOPORT**  
 Cat. No.: **11972** (packing, 20g/bottle)  
 Column: 18' x 1/8" stainless steel  
 Oven: 175°C to 350°C at 8°C/min  
 Carrier: nitrogen, 20mL/min  
 Det.: 350°C  
 Inj.: 1µL chloroform containing 30µg wax (325°C)




For more information, request Bulletin 864.

713-0750

## Figure 143. Wax

Packing: **3% Dexsil® 300 on 100/120 SUPELCOPORT**  
Cat. No.: **11973 (packing)**  
Column: **6' x 1/8" stainless steel**  
Oven: **100°C to 360°C at 4°C/min**  
Carrier: **nitrogen, 20mL/min**

1. C10
2. C12
3. C14
4. C16
5. C18
6. C20
7. C22
8. C24
9. C26
10. C28
11. C30
12. C32
13. C34
14. C36
15. C38
16. C40
17. C44




For more information, request Bulletin 864.

795-0242

## Figure 144. Solvents — Alcohol Retarder Column

Packing: **15% THEED on 100/120 Chromosorb W AW**  
Cat. No.: **11823 (packing)**  
Column: **10' x 1/8" stainless steel**  
Oven: **80°C**  
Carrier: **nitro**  
Det.: **FID**  
Inj.: **0.5µLgen, 20mL/min**

1. Heptane
2. Octane
3. Ethyl acetate
4. n-Propyl acetate
5. Toluene
6. Ethyl benzene
7. Xylene
8. Isopropyl alcohol
9. Ethyl alcohol
10. n-Propyl alcohol




For more information, request Bulletins 816 and 824.

795-0243

## Figure 145. Solvents

Packing: **60/80 Carbowax C/0.2% Carbowax 1500**  
 Cat. No.: **11826** (packing)  
 Column: **6' x 1/8" stainless steel**  
 Cat. No.: **13860-U** (general configuration stock column,  
 other stock columns available)  
 Oven: **100°C**  
 Carrier: **nitrogen, 20mL/min**  
 Det.: **FID**  
 Inj.: **1.0µL**




1. Methanol
2. Ethanol
3. Acetone
4. Isopropyl alcohol
5. Ethyl ether
6. Ethyl acetate
7. Methyl ethyl ketone
8. Isobutyl alcohol
9. Benzene

For more information,  
 request Bulletins 816 and 824.

795-0244

## Figure 146. Solvents (Methanol through Pentanol)

Packing: **GP 80/100 Carbowax C/0.2% Carbowax 1500**  
 Cat. No.: **11826** (packing)  
 Column: **6' x 2mm stainless steel**  
 Cat. No.: **13860-U** (general configuration stock column;  
 other stock columns available)  
 Oven: **125°C**  
 Carrier: **nitrogen, 20mL/min, 31psi**  
 Det.: **FID**  
 Inj.: **0.02µL**



1. Methanol
2. Ethanol
3. 2-Propanol
4. 1-Propanol
5. 2-Methyl-2-propanol (1-butyl)
6. 2-Butanol (sec-butyl)
7. 2-Methyl-1-propanol (iso)
8. 1-Butanol
9. 2-Methyl-2-butanol (t-amyl)
10. 2,2-Dimethyl-1-propanol
11. 3-Methyl-2-butanol
12. 3-Pentanol
13. 2-Pentanol
14. 2-Methyl-1-butanol (active)
15. 3-Methyl-1-butanol
16. 1-Pentanol

For more information, request Bulletins 816 and 824.


795-0239



## Figure 147. Solvents (Ethyl Ether Through Styrene)

Packing: **10% SP-1000 on 80/100 SUPELCOPORT**  
 Cat. No.: **11872** (packing)  
 Column: 20' x 1/8" stainless steel  
 Cat. No.: **12719** (general configuration stock column;  
 other stock columns available)  
 Oven: 100°C (6 min) to 140°C at 20°C/min  
 Carrier: nitrogen, 30mL/min  
 Det.: FID  
 Inj.: 5µL CS<sub>2</sub> approx.  
 100ppm each solvent

1. Ethyl ether
2. Isooctane
3. Carbon disulfide
4. Acetone
5. Carbon tetrachloride,  
Methylchloroform
6. Methyl ethyl ketone
7. Methylene chloride
8. Benzene
9. Trichloroethylene
10. Chloroform
11. Tetrachloroethylene
12. Toluene
13. Ethylene dichloride
14. Dioxane
15. Ethylbenzene
16. p-Xylene
17. m-Xylene
18. o-Xylene
19. 1,1,2-Trichloroethane
20. Styrene




For more information, request Bulletins 816 and 824.

795-0240

## Figure 148. Solvents (Methanol Through Xylenes)

Packing: **20% SP-2100/0.1% Carbowax 1500 on 100/120 SUPELCOPORT**  
 Cat. No.: **11821** (packing)  
 Column: 10' x 1/8" stainless steel  
 Cat. No.: **12718** (general configuration stock column;  
 other stock columns available)  
 Oven: 100°C  
 Carrier: nitrogen, 20mL/min  
 Det.: FID  
 Inj.: 0.5µL, approx. equal proportions  
 of neat compounds

1. Methanol
2. Ethanol
3. Acetone
4. Methyl ethyl ketone
5. Ethyl acetate
6. Isobutyl alcohol
7. Methyl Cellosolve
8. Isopropyl acetate
9. 2-Nitropropane
10. MIBK
11. Isobutyl acetate
12. Toluene
13. n-Butyl acetate
14. Diacetone alcohol
15. Ethylbenzene
16. m- & p-Xylene
17. Cellosolve acetate
18. Butyl Cellosolve & o-Xylene




For more information, request Bulletins 816 and 824.

795-0241

### Figure 149. Solvents — Ketone-Retarding Column

Packing: **20% SP-2401/0.1% Carbowax 1500 on 100/120 SUPELCOPORT**  
 Cat. No.: **11822** (packing)  
 Column: 10' x 1/8" stainless steel  
 Oven: 100°C  
 Carrier: nitrogen, 20mL/min  
 Det.: FID  
 Inj.: 0.1µL, approx. equal proportions of neat components




For more information, request Bulletins 816 and 824.

795-0245

### Figure 150. Industrial Solvents

Packing: **80/120 Carbowax B/3% SP-1500**  
 Cat. No.: **11813-U** (packing)  
 Column: 10' x 1/8" stainless steel  
 Cat. No.: **12592** (general configuration stock column; other stock columns available)  
 Oven: 70°C to 225°C at 4°C/min




For more information, request Bulletin 824.

713-0787

## Figure 151. Industrial Solvents

Packing: **60/80 Carboxpack B/1% SP-1510**  
 Cat. No.: **11809** (packing)  
 Column: 10' x 1/8" stainless steel  
 Oven: 100°C to 225°C at 8°C/min  
 Carrier: nitrogen, 20mL/min  
 Det.: FID  
 Inj.: 0.3µL, approx. equal proportions of neat components




For more information, request Bulletin 824.

713-0788

## Figure 152. Industrial Solvents

Packing: **80/100 Carboxpack C/0.1% SP-1000**  
 Cat. No.: **11820** (packing)  
 Column: 10' x 1/8" stainless steel  
 Oven: 70°C to 225°C at 4°C/min  
 Carrier: nitrogen, 20mL/min  
 Det.: FID  
 Inj.: 0.3µL, 3% each component by volume (approx. 12µg each)




For more information, request Bulletin 824.

713-0789

**Figure 153. Industrial Solvents**

Packing: **60/80 Carbo-pack F-TA**  
 Column: 2m x 2mm ID glass  
 Oven: 75°C (2 min) to 225°C  
 at 4°C/min  
 Carrier: nitrogen, 20mL/min  
 Det.: FID  
 Inj.: 0.05µL solvents mix,  
 3% each solvent



1. Methanol
2. Methyl formate
3. Ethanol
4. Acetone
5. Methylene chloride
6. Isopropyl alcohol
7. Ethyl formate
8. Methyl acetate
9. n-Propyl alcohol
10. Tetrahydrofuran (THF)
11. Methyl ethyl ketone (MEK)
12. Methyl Cellosolve
13. sec-Butyl alcohol
14. Ethyl acetate
15. Isobutyl alcohol
16. 2-Nitropropane
17. n-Butyl alcohol
18. Ethyl Cellosolve
19. n-Propyl acetate
20. Cyclohexanone
21. sec-Butyl acetate
22. Isobutyl acetate
23. n-Butyl acetate
24. Toluene
25. Mesityl oxide
26. Cellosolve acetate
27. Butyl Cellosolve
28. Isoamyl acetate
29. Ethylbenzene
30. m-Xylene
31. p-Xylene
32. o-Xylene


For more information, request Bulletin 824.

795-0248

**Figure 154. Alcohols, C1-C5**

Packing: **80/100 Carbo-pack F-SL**  
 Column: 2m x 2mm ID TightSpec glass  
 Oven: 130°C  
 Carrier: nitrogen, 15mL/min

1. Methanol
2. Ethanol
3. i-Propanol
4. n-Propanol
5. t-Butanol
6. i-Butanol
7. 2-Methyl-1-propanol
8. n-Butanol
9. 2-Methyl-2-butanol
10. 2,2-Dimethyl-1-propanol
11. 3-Methyl-2-butanol
12. 3-Pentanol
13. 2-Pentanol
14. 2-Methyl-1-butanol
15. 3-Methyl-1-butanol
16. n-Pentanol




For more information, request Bulletin 824.

713-1069

### Figure 155. Sugars

Packing: **3% SP-2330 on 100/120 SUPELCOPORT**  
Cat. No.: **11802** (packing)  
Column: 6' x 2mm ID glass  
Oven: 225°C  
Carrier: nitrogen, 20mL/min  
Det.: FID  
Inj.: 5µL

1. Glyceraldehyde
2. Erythrose
3. Deoxyribose
4. Rhamnose
5. Fucose
6. Ribose
7. Arabinose
8. Xylose
9. Deoxyglucose
10. Mannose
11. Galactose
12. Glucose
13. Inositol




795-0249

### Figure 156. Sugars

Packing: **3% SP-2340 on 100/120 SUPELCOPORT**  
Cat. No.: **11863** (packing)  
Column: 6' x 2mm ID glass  
Oven: 225°C  
Carrier: nitrogen, 20mL/min  
Det.: FID  
Inj.: 4µL


1. Glyceraldehyde
2. Erythrose
3. Deoxyribose
4. Rhamnose
5. Fucose
6. Ribose
7. Arabinose
8. Xylose
9. Deoxyglucose
10. Mannose
11. Galactose
12. Glucose
13. Inositol



795-0247

**Figure 157. Sulfur Dioxide in Nitrogen**

Packing: **45/60 Carboxen-1000**  
 Column: **2' x 1/8" stainless steel**  
 Cat. No.: **12370-U** (general configuration stock column;  
 other stock columns available)  
 Oven: 190°C  
 Carrier: helium, 30mL/min  
 Det.: TCD  
 Inj.: 1µL of 1000ppm SO<sub>2</sub> in N<sub>2</sub>,  
 approx. 1µg SO<sub>2</sub> on column




For more information, request Bulletin 722.

795-0250

**Figure 158. Sulfur Gases and Hydrocarbons**

Packing: **HayeSep Q**  
 Cat. No.: **10301-U** (packing)  
 Column: 8' x 1/8"  
 Oven: 90°C  
 Carrier: helium, 30cc/min




For more information, request Bulletin 722.

795-0251

## Figure 159. Hydrogen Sulfide, Carbon Dioxide, and Water

Packing: **100/120 HaysSep D**  
Cat. No.: **10293** (packing)  
Column: 10' x 1/8" Ni  
Oven: 60°C  
Carrier: helium, 30cc/min  
Det.: P.E. 900 T.C., 225ma, 140°C  
Inj.: Valco valve, 50 microliters  
vapor (ambient), 100°C

1. Air (balance)
2. Carbon dioxide (approx. 0.1%)
3. Water (approx. 2.5%)
4. Hydrogen sulfide (approx. 2.5%)




For more information, request Bulletin 722.

795-0252

## Figure 160. Hydrogen Sulfide, Nitrogen, Carbon Dioxide, Nitrous Oxide and Water

Packing: **100/120 HaysSep D**  
Cat. No.: **10293** (packing)  
Column: 10' x 1/8" stainless steel  
Oven: 80°C  
Carrier: helium, 30cc/min  
Det.: P.E. 900 T.C., 225ma  
Inj.: Valco valve, 100 microliters, 140°C

1. Nitrogen (balance)
2. Carbon dioxide (2%)
3. Nitrous oxide (3%)
4. Water (0.5%)
5. Hydrogen sulfide (3%)



For more information, request Bulletin 722.

795-0253

## Figure 161. Light Sulfur Gases, Alkyl Sulfides, and Mercaptans

Packing: **40/60 Carbopack B HT 100**

Cat. No.: **20272** (packing)

Column: **1.4m x 1/8" OD Teflon® (FEP)**


Cat. No.: **11502-U** (general configuration stock column only)

Oven: **100°C**

Carrier: **nitrogen, 70mL/min**

Det.: **flame photometric**

Inj.: **0.5mL of synthetic mix**



For more information, request Bulletin 722.

795-0254

## Figure 162. Trace Sulfur Gases

Packing: **40/60 Carbopack B HT 100**

Cat. No.: **20272** (packing)

Column: **1.4m x 1/8" OD Teflon (FEP)**


Cat. No.: **11502-U** (general configuration stock column only)

Oven: **35°C**

Carrier: **nitrogen, 20mL/min**

Det.: **flame photometric**

Inj.: **0.3mL**




For more information, request Bulletin 722.

713-0806



### Figure 163. Kraft Pulp Mill Stack Gases

Packing: **Supelpak™-S**  
 Column: **30" (18" packed) x 1/8" OD  
 Teflon® (FEP)**  
 Cat. No.: **12255-U** (general configuration  
 stock column only)  
 Oven: 30°C (1 min) to 210°C at 40°C/min  
 Carrier: helium, 30mL/min  
 Det.: flame photometric  
 Inj.: 2mL




1. H<sub>2</sub>S
2. COS
3. SO<sub>2</sub>
4. CH<sub>3</sub>SH
5. (CH<sub>3</sub>)<sub>2</sub>S
6. (CH<sub>3</sub>)<sub>2</sub>S<sub>2</sub>

For more information, request Bulletin 722.

713-0805

### Figure 164. Trace Sulfur Gases with C<sub>2</sub> Isomeric Sulfur Compounds

Packing: **GP 60/80 Carboxack B/1.5% XE-60/1.0% H<sub>3</sub>PO<sub>4</sub>**  
 Cat. No.: **11828** (packing)  
 Column: 6' x 2mm ID glass  
 Oven: 50°C  
 Carrier: nitrogen, 35mL/min  
 Det.: flame photometric  
 Inj.: 0.25mL




1. H<sub>2</sub>S
2. SF<sub>6</sub>
3. COS
4. SO<sub>2</sub>
5. Methyl mercaptan
6. Ethyl mercaptan
7. Dimethyl sulfide
8. CS<sub>2</sub>

For more information, request Bulletin 722.

713-0807

**Figure 165. Alkyl Sulfides and Mercaptans**

Packing: **Chromosil 330**  
 Column: **8' (6' packed) x 1/8" OD Teflon (FEP)**  
 Cat. No.: **11496** (general configuration stock column only)  
 Oven: 65°C  
 Carrier: nitrogen, 20mL/min  
 Det.: flame photometric  
 Inj.: 0.5cc of synthetic mix




For more information, request Bulletin 722.

795-0255

**Figure 166. Trace Light Sulfur Gases and C1-C3 Mercaptans**

Packing: **Chromosil 330**  
 Column: **8' (6' packed) x 1/8" OD Teflon (FEP)**  
 Cat. No.: **11496** (general configuration stock column only)  
 Oven: 40°C  
 Carrier: nitrogen, 20mL/min  
 Det.: flame photometric  
 Inj.: 0.2mL, approx. 1ppm  
 each component  
 in nitrogen




For more information,  
 request Bulletin 722.

713-0804

### Figure 167. Light Sulfur Gases at Percent Concentrations

Packing: **Chromosil 310**  
Column: **6' x 4mm ID glass**  
Cat. No.: **11501** (general configuration stock column only)  
Oven: 40°C  
Carrier: helium, 50mL/min  
Det.: thermal conductivity  
Inj.: 0.3mL of synthetic mix




For more information, request Bulletin 722.

794-0844

### Figure 168. Trace Light Sulfur Gases

Packing: **Chromosil 310**  
Column: **8' (6' packed) x 1/8" OD Teflon (FEP)**  
Cat. No.: **11501** (general configuration stock column only)  
Oven: 50°C  
Carrier: nitrogen, 20mL/min  
Det.: flame photometric  
Inj.: 0.5mL approx. 1ppm each component in nitrogen




For more information, request Bulletin 722.

713-0803

**Figure 169. Mercaptans, Sulfides, and Disulfides**

Packing: 12% polyphenyl ether/0.5% H<sub>3</sub>PO<sub>4</sub> on 40/60 Chromosorb T  
 Column: 36' x 0.085" ID Teflon  
 Cat. No.: 11500 (general configuration stock column only)  
 Oven: 50°C  
 Carrier: nitrogen, 80mL/min  
 Det.: FPD  
 Inj.: 0.5cc synthetic mix




For more information, request Bulletin 722.

794-0850

**Figure 170. Sulfur Gases**

Packing: 12% polyphenyl ether/0.5% H<sub>3</sub>PO<sub>4</sub> on 40/60 Chromosorb T  
 Column: 36' x 0.085" ID Teflon  
 Cat. No.: 11500 (general configuration stock column only)  
 Oven: 50°C  
 Carrier: nitrogen, 100mL/min  
 Det.: flame photometric (sens. settings on figure)  
 Inj.: 0.2mL nitrogen approx. 1ppm each analyte




For more information, request Bulletin 722.

713-0808

### Figure 171. Cholesterol (Steroids)


Packing: 3% SP-2250 on 100/120 SUPELCOPORT  
Cat. No.: 11878 (packing)  
Column: 3' x 2mm ID glass  
(stock column available)  
Oven: 250°C  
Carrier: nitrogen, 20mL/min



795-0235

### Figure 172. Cholesterol (Steroids)


Packing: 3% SP-2100 on 80/100 SUPELCOPORT  
Cat. No.: 11987 (packing)  
Column: 3' x 2mm ID glass  
Oven: 250°C  
Carrier: nitrogen, 20mL/min



795-0236

**Figure 173. Cholesterol (Steroids)**


Packing: **3% SP-2401 on 100/120 SUPELCOPORT**  
Cat. No.: **11978** (packing)  
Column: **3' x 2mm ID glass**  
Oven: **250°C**  
Carrier: **nitrogen, 20mL/min**



795-0237

**Figure 174. Cholesterol (Steroids)**


Packing: **3% SP-2300 on 100/120 SUPELCOPORT**  
Cat. No.: **11888** (packing)  
Column: **3' x 2mm ID glass**  
(stock column available)  
Oven: **250°C**  
Carrier: **nitrogen, 20mL/min**



795-0283

## Figure 175. Plant Sterols


Packing: 3% SP-2250 on 100/120 SUPELCOPORT  
Cat. No.: 11878 (packing)  
Column: 3' x 2mm ID glass  
Oven: 260°C  
Carrier: nitrogen, 20mL/min



795-0257

## Figure 176. Ketosteroids


Packing: 3% OV-225 on 80/100 SUPELCOPORT  
Cat. No.: 11957-U (packing)  
Column: 6' x 4mm ID glass  
Oven: 250°C  
Carrier: nitrogen, 20mL/min



795-0258

**Figure 177. Triglycerides**


Packing: **1% Dexsil 300 on 100/120 SUPELCOPORT**  
Cat. No.: **11972** (packing, 20g/bottle)  
Column: **18" x 1/8" stainless steel**  
Oven: **275°C to 350°C at 8°C/min**  
Carrier: **nitrogen, 20mL/min**  
Det.: **FID, 350°C**  
Inj.: **1µL chloroform containing 1µg each triglyceride, 325°C**



713-0944

**Figure 178 Pentaerithritol Esters**

Packing: **1% Dexsil 300 on 100/120 SUPELCOPORT**  
Cat. No.: **11972** (packing)  
Column: **18" x 1/8" stainless steel**  
Oven: **125°C to 300°C at 8°C/min**  
Carrier: **nitrogen, 20mL/min**  
Det.: **350°C**




795-0259



## Figure 179. Cholesteryl Esters

Packing: **1% Dexsil 300 on 100/120 SUPELCOPORT**  
Cat. No.: **11972 (packing)**  
Column: **18" x 2mm ID glass**  
Oven: **300°C to 350°C at 6°C/min**  
Carrier: **nitrogen, 40mL/min**  
Det.: **350°C**



713-0943

For expert answers to your questions,  
contact our Technical Service Department  
Phone **800-359-3041** or **814-359-3041**  
FAX **800-359-3044** or **814-359-5468**

For more information, or current prices, contact your nearest Supelco subsidiary listed below. To obtain further contact information, visit our website ([www.sigma-aldrich.com](http://www.sigma-aldrich.com)), see the Supelco catalog, or contact Supelco, Bellefonte, PA 16823-0048 USA.

**ARGENTINA** · Sigma-Aldrich de Argentina, S.A. · Buenos Aires 1119

**AUSTRALIA** · Sigma-Aldrich Pty. Ltd. · Castle Hill NSW 2154

**AUSTRIA** · Sigma-Aldrich Handels GmbH · A-1110 Wien

**BELGIUM** · Sigma-Aldrich N.V./S.A. · B-2880 Bornem

**BRAZIL** · Sigma-Aldrich Quimica Brasil Ltda. · 01239-010 São Paulo, SP

**CANADA** · Sigma-Aldrich Canada, Ltd. · 2149 Winston Park Dr., Oakville, ON L6H 6J8

**CZECH REPUBLIC** · Sigma-Aldrich s.r.o. · 186 00 Praha 8

**DENMARK** · Sigma-Aldrich Denmark A/S · DK-2665 Vallensbaek Strand

**FINLAND** · Sigma-Aldrich Finland/YA-Kemia Oy · FIN-00700 Helsinki

**FRANCE** · Sigma-Aldrich Chimie · 38297 Saint-Quentin-Fallavier Cedex

**GERMANY** · Sigma-Aldrich Chemie GmbH · D-82041 Deisenhofen

**GREECE** · Sigma-Aldrich (o.m.) Ltd. · Ilioupoli 16346, Athens

**HUNGARY** · Sigma-Aldrich Kft. · H-1067 Budapest

**INDIA** · Sigma-Aldrich Co. · Bangalore 560 048

**IRELAND** · Sigma-Aldrich Ireland Ltd. · Dublin 24

**ISRAEL** · Sigma Israel Chemicals Ltd. · Rehovot 76100

**ITALY** · Sigma-Aldrich s.r.l. · 20151 Milano

**JAPAN** · Sigma-Aldrich Japan K.K. · Chuo-ku, Tokyo 103

**KOREA** · Sigma-Aldrich Korea · Seoul

**MALAYSIA** · Sigma-Aldrich (M) Sdn. Bhd. · 58200 Kuala Lumpur

**MEXICO** · Sigma-Aldrich Química S.A. de C.V. · 50200 Toluca

**NETHERLANDS** · Sigma-Aldrich Chemie BV · 3330 AA Zwijndrecht

**NORWAY** · Sigma-Aldrich Norway · Torshov · N-0401 Oslo

**POLAND** · Sigma-Aldrich Sp. z o.o. · 61-663 Poznań

**PORTUGAL** · Sigma-Aldrich Quimica, S.A. · Sintra 2710

**RUSSIA** · Sigma-Aldrich Russia · Moscow 103062

**SINGAPORE** · Sigma-Aldrich Pte. Ltd.

**SOUTH AFRICA** · Sigma-Aldrich (pty) Ltd. · Jet Park 1459

**SPAIN** · Sigma-Aldrich Quimica, S.A. · 28100 Alcobendas, Madrid

**SWEDEN** · Sigma-Aldrich Sweden AB · 135 70 Stockholm

**SWITZERLAND** · Supelco · CH-9471 Buchs

**UNITED KINGDOM** · Sigma-Aldrich Company Ltd. · Poole, Dorset BH12 4QH

**UNITED STATES** · Supelco · Supelco Park · Bellefonte, PA 16823-0048 · Phone 800-247-6628 or 814-359-3441 · Fax 800-447-3044 or 814-359-3044 · email: [supelco@sial.com](mailto:supelco@sial.com) **H**

Supelco is a member of the Sigma-Aldrich family. Supelco products are sold through Sigma-Aldrich, Inc. Sigma-Aldrich warrants that its products conform to the information contained in this and other Sigma-Aldrich publications. Purchaser must determine the suitability of the product for a particular use. Additional terms and conditions may apply. Please see the reverse side of the invoice or packing slip.

**To Order:**

**Phone 800-247-6628 or 814-359-3441**

**FAX 800-447-3044 or 814-359-3044**



**Supelco Park, Bellefonte, PA, 16823 USA**