



Ketones

Separation of several ketones on a wide-bore fused silica column

Application Note

Energy & Fuels

Authors

Agilent Technologies, Inc.

Introduction

Gas chromatography using an Agilent CP-Sil 5 CB column separates five ketones in seven minutes.




Agilent Technologies

Conditions

Technique : GC-capillary
Column : Agilent CP-Sil 5 CB, 0.53 mm x 10 m fused silica
WCOT CP-Sil 5 CB (5.0 μm) (Part no. CP7645)
Temperature : 50 °C \rightarrow 200 °C, 5 °C/min
Carrier Gas : N₂, 10 kPa (0.1 bar), 52 cm/s
Injector : direct
T = 250 °C
Detector : FID, 100 x 10⁻¹² Afs
T = 275 °C
Sample Size : 0.2 μL
Concentration Range : acetone

Peak identification

1. methanol
2. acetone
3. 2,3-butanedione (diacetyl)
4. acetylacetone
5. mesityloxyde (4-methyl-3-penten-2-one)
6. methylpyrrolidone



www.agilent.com/chem

This information is subject to change without notice.

© Agilent Technologies, Inc. 2011

Printed in the USA

31 October, 2011

First published prior to 11 May, 2010

A00015



Agilent Technologies