



Sulfur compounds

Trace analysis of sulfur compounds

Application Note

Energy & Fuels

Authors

Agilent Technologies, Inc.

Introduction

The inertness of the PoraBOND Q column allows the separation of volatile sulfur compounds at low levels with an excellent peak shape. Together with the high sensitivity of the Pulsed Discharge Detector, linear quantification down to the sub-ppm level is possible under these conditions. Also water can be detected.



Agilent Technologies


Conditions

Technique : GC-capillary
Column : Agilent PoraBOND Q, 0.32 mm x 15 m fused silica
PLOT (df = 5 μ m) (Part no. CP7351)
(as 25 m column)
Temperature : 35 °C (3 min) \rightarrow 250 °C, 20 °C/min
Carrier Gas : He, 50 kPa (0.5 bar, 7.2 psi)
Injector : Valve/Split
Detector : VICI Pulsed Discharge Detector,
T = 300 °C
Sample Size : 100 μ L
Concentration Range : 100 ppm

Courtesy : Jim Luong, Analytical Sciences,
Dow Western Canada Operations

Peak identification

1. water
2. hydrogen sulfide
3. carbonyl sulfide
4. methanethiol (methyl mercaptan)
5. ethanethiol (ethyl mercaptan)



www.agilent.com/chem

This information is subject to change without notice.

© Agilent Technologies, Inc. 2011

Printed in the USA

31 October, 2011

First published prior to 11 May, 2010

A01458



Agilent Technologies