

# New Application Specific Columns

Simon Jones  
GC columns Application Engineer

# What Columns are New?

DB-CLP1

DB-CLP2

DB-UI 8270D



Environmental

DB-Select 624 UI 467 — US Pharmacopeia

DB-Sulfur SCD

PLOT-PT

# What are the Benefits?

Address specific issues with certain methods

Response

Tailing factor

Resolution

Detector coking

Detector Spiking / Valve rotor damage

Named for easier identification

# DB-CLP1 & DB-CLP2

What are they?

Contract Laboratory Pesticides

Why 2 columns?

Dual column method ( $\mu$ ECD)

Primary and Confirmatory columns

Pesticides are active compounds!

Response, peak shape, tailing, breakdown

# DB-CLP1/CLP2


## CLP Pesticides

DB-CLP1 30m 0.32mm ID 0.25µm (cat.# 123-8232)  
DB-CLP2 30m 0.32mm ID 0.5µm (cat.# 123-8336)  
5m x 0.32mm ID deact. guard column  
Inert Tee CFT device (cat.# G3184-60065) 1:1 Split


Instrument: Agilent 7890 GC with dual µECD  
Sampler: Agilent 7693  
Sample: 50 ng/mL CLP Pesticides  
Inj. Vol.: 1 µL splitless  
Liner: Ultra Inert liner, single taper splitless  
(cat.#5190-2292)  
Inj. Temp.: 250°C  
Oven Temp: 150°C (hold 0.2 min),  
45°C/min to 250°C, 18°C/min to 300°C,  
30°C/min to 330°C, hold 2.5 min  
Carrier Gas: Helium, constant flow 3.5 mL/min  
Detector: µECD @ 340 °C

1. Tetrachloro-m-xylene*	9. γ-Chlordane	16. Endosulfan II
2. α-BHC	10. α-Chlordane	17. 4,4'-DDT
3. γ-BHC	11. Endosulfan I	18. Endrin aldehyde
4. β-BHC	12. 4,4'-DDE	19. Endosulfan sulfate
5. Heptachlor	13. Dieldrin	20. Methoxychlor
6. δ-BHC	14. Endrin	21. Endrin ketone
7. Aldrin	15. 4,4'-DDD	22. Decachlorobiphenyl*
8. Heptachlor epoxide		*surrogate standard

Agilent DB-CLP1



Agilent DB-CLP2



# DB-CLP1/DB-CLP2

## Additional Methods

**EPA 504.1** - 1,2-dibromoethane (EDB), 1,2-dibromo-3-chloropropane (DBCP), and 1,2,3-trichloropropane

**EPA 505** – organohalides

**EPA 508.1** – Chlorinated pesticides, herbicides and organohalides

**EPA 551** – chlorinated solvents, trihalomethanes (THMs), and disinfection byproducts (DBPs)

**EPA 552.3** – haloacetic acids and dalapon

**EPA 8081B** – (extended analyte list) – organochlorine pesticides

**EPA 8082A** – polychlorinated biphenyls (PCBs) and aroclors

**EPA 8151A** – chlorophenoxyacid herbicides

# EPA Method 504.1


## EDB, DBCP, and 123TCP

DB-CLP1 30m 0.32mm ID 0.25µm (cat.# 123-8232)  
DB-CLP2 30m 0.32mm ID 0.5µm (cat.# 123-8336)  
5m x 0.32mm ID deact. guard column  
Inert Tee CFT device (cat.# G3184-60065) 1:1 Split


Instrument: Agilent 7890 GC with dual µECD  
Sampler: Agilent 7693  
Sample: 100 ng/mL 504.1 analytes  
100 ng/mL Chlorinated Solvents + THMs  
Inj. Vol.: 2 µL splitless  
Liner: Ultra Inert liner, single taper splitless  
(cat.#5190-2292)  
Inj. Temp.: 200°C  
Oven Temp: 50°C (hold 1.5 min) to 95°C at 20°C/min,  
40°C/min to 175°C, hold 1.25 min  
Carrier Gas: Helium, constant flow 3.75 mL/min  
Detector: µECD @ 300 °C

- |                          |  |
|--------------------------|--|
| 1. Chloroform            | 7. 1,1,2-Trichloroethane               |
| 2. 1,1,1-Trichloroethane | 8. Dibromochloromethane                |
| 3. Carbon tetrachloride  | 9. 1,2-Dibromoethane (EDB)             |
| 4. Trichloroethene       | 10. Bromoform                          |
| 5. Bromodichloromethane  | 11. 1,2,3-Trichloropropane (123TCP)    |
| 6. Tetrachloroethene     | 12. 1,2-Dibromo-3-chloropropane (DBCP) |

Agilent DB-CLP1



Agilent DB-CLP2



# EPA Method 505


## Organohalide Pesticides

DB-CLP1 30m 0.32mm ID 0.25µm (cat.# 123-8232)  
DB-CLP2 30m 0.32mm ID 0.5µm (cat.# 123-8336)  
5m x 0.32mm ID deact. guard column  
Inert Tee CFT device (cat.# G3184-60065) 1:1 Split


Instrument: Agilent 7890 GC with dual µECD  
Sampler: Agilent 7873B  
Sample: 100 ng/mL 505 analytes  
Inj. Vol.: 2 µL splitless  
Liner: Ultra Inert Liner double taper, splitless  
(cat.# 5190-3983)  
Inj. Temp.: 250°C  
Oven Temp: 90°C (hold 0.5 min), 35°C/min to 175°C,  
12°C/min to 300°C, hold 1.75 min  
Carrier Gas: Helium, constant flow 2.5 mL/min  
Detector: µECD @ 325 °C

- |                              |                              |
|------------------------------|------------------------------|
| 1. Hexachlorocyclopentadiene | 9. Heptachlor epoxide        |
| 2. Hexachlorobenzene         | 10. γ-Chlordane              |
| 3. Atrazine (25 µg/mL)       | 11. trans-Nonachlor          |
| 4. Simazine (25 µg/mL)       | 12. α-Chlordane              |
| 5. γ-BHC                     | 13. Dieldrin                 |
| 6. Heptachlor                | 14. Endrin                   |
| 7. Alachlor (1 µg/mL)        | 15. cis-Nonachlor            |
| 8. Aldrin                    | 16. Methoxychlor (0.5 µg/mL) |

Agilent DB-CLP1



Agilent DB-CLP2





# EPA Method 508.1


## Chlorinated Pesticides and Herbicides

DB-CLP1 30m 0.32mm ID 0.25µm (cat.# 123-8232)  
 DB-CLP2 30m 0.32mm ID 0.5µm (cat.# 123-8336)  
 5m x 0.32mm ID deact. guard column  
 Inert Tee CFT device (cat.# G3184-60065) 1:1 Split


Instrument: Agilent 7890 GC with dual µECD  
 Sampler: Agilent 7873B  
 Sample: 100 ng/mL 508.1 analytes  
 100 ng/mL Pesticide Surrogate Mix  
 Inj. Vol.: 2 µL splitless  
 Liner: Ultra Inert liner, single taper splitless  
 (cat.#5190-2292)  
 Inj. Temp.: 250°C  
 Oven Temp: 80°C (hold 0.5 min) to 175°C at 26°C/min,  
 6.5°C/min to 235°C, 15°C/min to 300°C,  
 hold 6 min  
 Carrier Gas: Helium at 35 cm/sec (constant flow)  
 Detector: µECD @ 340 °C

1. Hexachlorocyclopentadiene	11. γ-BHC	21. Heptachlor epoxide	31. 4,4'-DDD
2. Etriazole	12. β-BHC	22. Cyanazine	32. Endosulfan II
3. Chloroneb	13. Heptachlor	23. Butachlor	33. 4,4'-DDT
4. Trifluralin	14. Alachlor	24. γ-Chlordane	34. Endrin aldehyde
5. TCMX*	15. δ-BHC	25. α-Chlordane	35. Endosulfan sulfate
6. Propachlor	16. Chlorothalonil	26. Endosulfan I	36. Methoxychlor
7. Hexachlorobenzene	17. Aldrin	27. 4,4'-DDE	37. cis-Permethrin
8. α-BHC	18. Metribuzin	28. Dieldrin	38. trans-Permethrin
9. Atrazine	19. Metolachlor	29. Chlorobenzilate	39. Decachlorobiphenyl*
10. Simazine	20. DCPA	30. Endrin	<i>*surrogate standard</i>

Agilent DB-CLP1



Agilent DB-CLP2




# EPA Method 551

## Chlorinated Solvents, THMs, and DBPs


DB-CLP1 30m 0.32mm ID 0.25µm (cat.# 123-8232)  
DB-CLP2 30m 0.32mm ID 0.5µm (cat.# 123-8336)  
5m x 0.32mm ID deact. guard column  
Inert Tee CFT device (cat.# G3184-60065) 1:1 Split

Instrument: Agilent 7890 GC with dual µECD  
Sampler: Agilent 7693  
Sample: 100 ng/mL Chlorinated Solvents, THMs, and DBPs  
Inj. Vol.: 2 µL splitless  
Liner: Ultra Inert liner, single taper splitless  
(cat.#5190-2292)  
Inj. Temp.: 200°C  
Oven Temp: 35°C (hold 5.75 min) to 95°C at 20°C/min,  
40°C/min to 200°C, hold 1.25 min  
Carrier Gas: Helium at 45 cm/sec (constant flow)  
Detector: µECD @ 300 °C

- |                             |                                 |
|-----------------------------|---------------------------------|
| 1. Chloroform               | 11. Tetrachloroethene           |
| 2. 1,1,1-Trichloroethane    | 12. 1,1,2-Trichloroethane       |
| 3. Carbon tetrachloride     | 13. Dibromochloromethane        |
| 4. Trichloroacetonitrile    | 14. 1,2-Dibromoethane           |
| 5. Trichloroethene          | 15. 1,1,1-Trichloro-2-propanone |
| 6. Chloral hydrate          | 16. Bromochloroacetonitrile     |
| 7. Bromodichloromethane     | 17. Bromoform                   |
| 8. 1,1-Dichloro-2-propanone | 18. 1,2,3-Trichloropropane      |
| 9. Dichloroacetonitrile     | 19. Dibromoacetonitrile         |
| 10. Chloropicrin            | 20. 1,2-Dibromo-3-chloropropane |



Agilent DB-CLP1



Agilent DB-CLP2

# EPA Method 552.3


## Haloacetic Acids and Dalapon

DB-CLP1 30m 0.32mm ID 0.25µm (cat.# 123-8232)  
 DB-CLP2 30m 0.32mm ID 0.5µm (cat.# 123-8336)  
 5m x 0.32mm ID deact. guard column  
 Inert Tee CFT device (cat.# G3184-60065) 1:1 Split


Instrument: Agilent 7890 GC with dual µECD  
 Sampler: Agilent 7693  
 Sample: 10-100 ng/mL Haloacetic acids and Dalapon (methyl esters)  
 Inj. Vol.: 1 µL splitless  
 Liner: Ultra Inert liner, single taper splitless (cat.#5190-2292)  
 Inj. Temp.: 180°C  
 Oven Temp: 40°C (hold 0.5 min) to 95°C at 10°C/min, 30°C/min to 200°C, hold 1 min  
 Carrier Gas: Helium at 54.79 cm/sec (constant flow)  
 Detector: µECD @ 340 °C

- |  |   |
|--|---|
| 1. Methyl chloroacetate (30ng/mL)        | 7. Methyl bromochloroacetate (20ng/mL)    |
| 2. Methyl bromoacetate (20ng/mL)         | 8. Methyl bromodichloroacetate (20ng/mL)  |
| 3. Methyl dichloroacetate (30ng/mL)      | 9. Methyl dibromoacetate (10ng/mL)        |
| 4. Dalapon methyl ester (20ng/mL)        | 10. Methyl dibromochloroacetate (50ng/mL) |
| 5. Methyl trichloroacetate (10ng/mL)     | 11. Methyl tribromoacetate (100ng/mL)     |
| 6. 1,2,3-Trichloropropane (IS) (50ng/mL) |   |

Agilent DB-CLP1



Agilent DB-CLP2



# EPA Method 8081B (extended)


## Organochlorine Pesticides

DB-CLP1 30m 0.32mm ID 0.25µm (cat.# 123-8232)  
DB-CLP2 30m 0.32mm ID 0.5µm (cat.# 123-8336)  
5m x 0.32mm ID deact. guard column  
Inert Tee CFT device (cat.# G3184-60065) 1:1 Split


Instrument: Agilent 7890 GC with dual µECD  
Sampler: Agilent 7693  
Sample: 50 ng/mL 8081B analytes  
Inj. Vol.: 2 µL splitless  
Liner: Ultra Inert liner, single taper splitless (cat.#5190-2292)  
Inj. Temp.: 250°C  
Oven Temp: 80°C (hold 0.5 min) to 150°C at 20°C/min, 5°C/min to 235°C, 15°C/min to 300°C, hold 5 min  
Carrier Gas: Helium at 43.5 cm/sec (constant flow)  
Detector: µECD @ 325 °C

1. 1,2-Dibromo-3-chloropropane	11. α-BHC	21. DCPA	31. Chlorobenzilate (250ng/mL)	41. Captafol
2. Hexachlorocyclopentadiene	12. Pentachloronitrobenzene	22. Isodrin	32. Perthane (250ng/mL)	42. Methoxychlor
3. 1-Bromo-2-nitrobenzene	13. γ-BHC	23. Kelthane	33. Chloropropylate (250ng/mL)	43. Endrin ketone
4. Etriazazole	14. β-BHC	24. Heptachlor epoxide	34. Endrin	44. Mirex
5. Chloroneb	15. Heptachlor	25. γ-Chlordane	35. Nitrofen	45. cis-Permethrin
6. Trifluralin	16. Diclhone	26. trans-Nonachlor	36. 4,4'-DDD	46. trans-Permethrin
7. TCMX*	17. Alachlor	27. α-Chlordane	37. Endosulfan II	47. Decachlorobiphenyl*
8. Propachlor	18. δ-BHC	28. Endosulfan I	38. 4,4'-DDT	*surrogate standard
9. Diallate isomers (250ng/mL)	19. Chlorothalonil	29. 4,4'-DDE	39. Endrin aldehyde	
10. Hexachlorobenzene	20. Aldrin	30. Dieldrin	40. Endosulfan sulfate	

Agilent DB-CLP1



Agilent DB-CLP2



# EPA Method 8082A


## PCB Congeners

DB-CLP1 30m 0.32mm ID 0.25µm (cat.# 123-8232)  
DB-CLP2 30m 0.32mm ID 0.5µm (cat.# 123-8336)  
5m x 0.32mm ID deact. guard column  
Inert Tee CFT device (cat.# G3184-60065) 1:1 Split


Instrument: Agilent 7890 GC with  $\mu$ ECD  
Sampler: Agilent 7693  
Sample: 100 ng/mL PCB congeners  
100 ng/mL Pesticide Surrogate Mix  
Inj. Vol.: 2 µL splitless  
Liner: Ultra Inert liner, single taper splitless  
(cat.#5190-2292)  
Inj. Temp.: 250°C  
Oven Temp: 125°C (hold 0.25 min),  
20°C/min to 210°C, hold 0.5 min,  
7°C/min to 235°C, hold 0.75 min,  
25°C/min to 325°C, hold 2 min,  
Carrier Gas: Helium, constant flow 3 mL/min  
Detector:  $\mu$ ECD @ 340 °C

1. BZ #1	8. BZ #66	15. BZ #138
2. TCMX*	9. BZ #101	16. BZ #187
3. BZ #5	10. BZ #87	17. BZ #183
4. BZ #18	11. BZ #110	18. BZ #180
5. BZ #31	12. BZ #151	19. BZ #170
6. BZ #52	13. BZ #153	20. BZ #206
7. BZ #44	14. BZ #141	21. Decachlorobiphenyl* *surrogate standard

Agilent DB-CLP1



Agilent DB-CLP2





# EPA Method 8082A

## Aroclors

DB-CLP1 30m 0.32mm ID 0.25 $\mu$ m (cat.# 123-8232)  
DB-CLP2 30m 0.32mm ID 0.5 $\mu$ m (cat.# 123-8336)  
5m x 0.32mm ID deact. guard column  
Inert Tee CFT device (cat.# G3184-60065) 1:1 Split

Instrument: Agilent 7890 GC with dual  $\mu$ ECD  
Sampler: Agilent 7693  
Sample: 100 ng/mL Individual Aroclor  
100 ng/mL Decachlorobiphenyl  
Inj. Vol.: 2  $\mu$ L splitless  
Liner: Ultra Inert liner, single taper splitless  
(cat.#5190-2292)  
Inj. Temp.: 250°C  
Oven Temp: 160°C (hold 0.25 min),  
35°C/min to 330°C, hold 3 min  
Carrier Gas: Helium, constant flow 3.75 mL/min  
Detector:  $\mu$ ECD @ 340 °C



# DB-8270D UI

EPA Semivolatiles method


Dirty Samples

Issues with Endrin and DDT breakdown \*

Response of 2,4-Dinitrophenol

Resolution

# 78 Semi-Volatile Components on a DB-UI 8270D 30m x 0.25 x 0.25




## 25 minute Semi-Volatile Analysis

Application note 5991-0250EN




# Resolution of benzo-b & k fluoranthene isomers



## Positional isomers

Column: Agilent DB-UI 8270D, 30 m x 0.25 mm, 0.25  $\mu$ m (p/n 122-9732)  
Liner: Dual taper direct connect liner (p/n G1544-80700)  
Inlet: MMI in non-pulsed splitless mode 1  $\mu$ L at 275  $^{\circ}$ C  
Carrier: He, 1.2 mL/min, constant flow  
Septum purge flow: 3 mL/min, purge time on 0.7 min 50 mL/min  
Oven program: 30  $^{\circ}$ C (1.0 min), 15  $^{\circ}$ C/min to 100  $^{\circ}$ C, 20  $^{\circ}$ C/min to 240  $^{\circ}$ C (0.5 min),  
15  $^{\circ}$ C/min to 325  $^{\circ}$ C (6.7 min) Gas saver Off  
GC/MSD: Agilent 7890/5975C, 325  $^{\circ}$ C transfer line, 280  $^{\circ}$ C source,  
150  $^{\circ}$ C quad, 35-500 AMU range  
Sampler: Agilent 7693, 10.0  $\mu$ L syringe (p/n G4513-80216)

# 10 ng/ul Semivolatile Checkout Standard on 20m x 0.18mm x 0.36um DB-8270D Capillary GC Column using an Ultra Inert Liner with Wool




1	N-Nitrosodimethylamine	11	4-nitrophenol	21	Phenanthrene-d10
2	Aniline	12	2,4-dinitrotoluene	22	Aldrin
3	1,4-Dichlorobenzene-d4	13	Flourene	23	Heptachlor epoxide
4	Isophorone	14	4,6-dinitro-2-methyl phenol	24	Endrin
5	1,3-dimethyl-2-nitrobenzene	15	Trifluralin	25	4,4'-DDT
6	Naphthalene	16	Simazine	26	3,3'-dichlorobenzidine
7	hexachlorocyclopentadiene	17	Atrazine	27	Chrysene d-12
8	Mevinphos	18	pentachlorophenol	28	benzo[b]flouranthene
9	Acenaphthene-d10	19	Terbufos	29	Perylene-d12
10	2,4-dinitrophenol	20	Chlorothanlonil		

**16 minute Semi-Volatile Analysis**

Checkout Standard mix Agilent part # 5190-0473

# Excellent Peak Shape and Response for 2,4 Dinitrophenol



Peak #	Compound
9	Acenaphthalene D10
10	2,4 Dinitrophenol
11	4 Nitro phenol
12	2,4 Dinitro toluene

Column 1: Agilent J&W 8270-D 20 m x 0.18mm x 0.36  $\mu$ m

Column 2: 1 m 0.15 mm ID deactivated 2 psi with 5ml/min bleeder

Carrier: 1 ml/min constant flow at 40°C, AUX EPC 2 PSI constant pressure

Backflush: postrun 75 psi AUX EPC 2 PSI inlet

Oven: 40°C (2.5min) to 320°C (25°C/min) 4.8 min hold

Injection: 280°C Pulsed splitless 44 psi until 0.3 min, purge flow 60ml/min on at 0.35 min

Liner: UI single taper

MSD: Transfer Line 290°C, Source 300°C, Quad 150°C 50 -550 AMU

# Additional Enhancements

## Ultra Inert Liners



### Certificate of Performance

5190 -2293 Ultra Inert Liner

Splitless , Sngl taper, Glass Wool

Liner Body Lot: 0023A

Deactivation Lot: B11002

Tested for: 2ng 4-Aminopyridine  
2ng 2,4-Dinitrophenol

## Ultimetal Plus



Inlet Weldment

## Ultra Inert Gold Seal



Flexible Metal Ferrules



# DB-Select 624 UI 467

Residual Solvents Method USP<467>

Procedure A uses “G43” stationary phase

6% cyanopropyl phenyl -94% dimethylpolysiloxane

30 m X 0.53 mm X 3.0  $\mu\text{m}$

30 m X 0.32 mm X 1.8  $\mu\text{m}$

# USP 467 Application Development

## ➤ **Column Performance Evaluation**

- USP Method 467 Residual Solvents

- Standards Criteria

- Class 1

- Carbon Tetrachloride Signal to Noise  $\geq 3$

- Class 2A

- Resolution of Acetonitrile:Dichloromethane

- Resolution of Methylcyclohexane:1,4-Dioxane


- Tailing of Methylcyclohexane

- Class 2B

- Tailing for Pyridine


# USP 467 Class 1

DB-Select 467UI Serial #USC9275175



# USP 467 Class 1


DB-Select 467UI Serial #USC284013H






# USP 467 Class 2A

DB-624 UI 467 Select Serial #USC9261831




# USP 467 Class 2B

DB-624 UI 467 Select Serial #USC9261831



# Impact of Data Acquisition Rate on S/N Determination




# DB-Sulfur SCD

Thick film PDMS phase generally coat ceramics in SCD

Must have good inertness for H<sub>2</sub>S


# Traditional PDMS column- Coating of Reactor Tubes

Overlay of before (green) and after 2 x 2 uL neat toluene Injection (red)



Data courtesy of Jim Luong, Ronda Gras, Myron Hawryluk of Dow Chemical Canada


# Configuration to test SCD Quenching Issue



*Data courtesy of Jim Luong, Ronda Gras, Myron Hawryluk of Dow Chemical Canada*

# New DB-Sulfur SCD column

Last Three Runs of the day (n=20, 100 ppmv std)



Data courtesy of Jim Luong, Ronda Gras, Myron Hawryluk of Dow Chemical Canada

# PLOT-PT columns

Solid stationary phase 'Particles' can cause:


Detector spiking

Valve scoring


Clogged Unions



# No Detector Spikes Observed on PT Columns with Repeated Temperature and Pressure Cycling




- Temperature: 150°C + 20°C/min → 250°C; 15 times
- Pressure 3x optimum
- Each run switch off/on carrier gas 10 times



*The unusual “chromatogram” shows the detector signal profile of the temp and pressure cycling*

# Integrated Particle Trap PLOT Column


2.5 meter integrated particle traps on both ends




Compatible with GC/MS and valve switching systems,  
and Capillary Flow Technology

Very similar selectivity, plates and peak shape to  
existing Agilent PLOT columns

# The PLOT Column and Integrated Particle Trap



# No spikes at Fixed Gases Analysis on CP-Molsieve 5Å PLOT PT column



# Halocarbons by GCMS

**Column: PoraPLOT Q PT, 25m×0.32mm,10um (P/N CP7551PT)  
(30m total length)**

Carrier : Helium, 42cm/s @55 °C

Oven: : 55 °C for 5min

55 °C - 200 °C at 12°C/min


200 °C for 10min

Injection: 250 °C, splitless, 0.2min purge activation time

Detector: MSD, 280 °C Transfer line, full Scan at m/z 45-185

Sample: 1uL

1. Fluoroform (Freon-23)
2. 1,1,1-trifluoroethane (Freon-143a)
3. Pentafluoroethane (Freon-125)
4. Bromotrifluoromethane (Freon-13b1)
5. 1,1,1,2-Tetrafluoroethane (Freon-134a)
6. 1,1-difluoroethane (Freon-152a)
7. Difluorochloromethane (Freon-22)
8. 1,1,2,2-tetrafluoroethane (Freon-134)
9. 1-chloro-1,1-difluoroethane (Freon-142)
10. Bromochlorodifluoromethane (Freon-12b1)
11. Ethyl Chloride (Freon-160)
12. Fluorodichloromethane (Freon-21)
13. Trichloromonofluoromethane (Freon-11)
14. 1,1-Dichloro-1-fluoroethane (Freon-141)
15. 2,2-dichloro-1,1,1-trifluoroethane (Freon-123)
16. 1,1,2-trichloro-1,2,2-trifluoroethane (Freon-113)
17. 1,2-dibromo-1,1,2,2-tetrafluoroethane (Freon-114b2)
18. Trichloromethane (Freon-20)
19. 1,2-dichloroethane
20. 1,1,1-trichloro-ethane
21. Trichloroethylene
22. 1,1,2-trichloroethane



# Agilent J&W

## PLOT PT Columns

### Available PLOT PT columns:

- Porous polymers:
  - PoraPLOT Q
  - PoraBOND Q
  - PoraBOND Q HT
  - HP-PLOT Q
  - GS-Q
  - PoraPLOT U
  - HP-PLOT U
- **NEW!** Aluminum oxide
  - HP-PLOT Al<sub>2</sub>O<sub>3</sub> S
  - HP-PLOT Al<sub>2</sub>O<sub>3</sub> M
  - HP-PLOT Al<sub>2</sub>O<sub>3</sub> KCl
  - GS-Alumina
  - GS-Alumina/KCl
  - CP-Al<sub>2</sub>O<sub>3</sub> KCl
  - CP-Al<sub>2</sub>O<sub>3</sub> Na<sub>2</sub>SO<sub>4</sub>
- **NEW!** Molsieve
  - CP-Molsieve 5A

Custom PLOT PT columns are available for these phases

Phase type	Part number	Description	Dimensions
PLOT Q	CP7348PT	PoraBOND Q PT	25m x 0.25mm x 3µm
	CP7351PT	PoraBOND Q PT	25m x 0.32mm x 5µm
	CP7352PT	PoraBOND Q PT	50m x 0.32mm x 5µm
	CP7353PT	PoraBOND Q PT	10m x 0.53mm x 10µm
	CP7354PT	PoraBOND Q PT	25m x 0.53mm x 10µm
	CP7550PT	PoraPLOT Q PT	10m x 0.32mm x 10µm
	CP7551PT	PoraPLOT Q PT	25m x 0.32mm x 10µm
	CP7554PT	PoraPLOT Q PT	25m x 0.53mm x 20µm
	CP7557PT	PoraPLOT Q-HT PT	25m x 0.32mm x 10µm
	115-3432PT	GS-Q PT	30m x 0.53mm
	19091P-QO3PT	HP-PLOT Q PT	15m x 0.32mm x 20µm
	19091P-QO4PT	HP-PLOT Q PT	30m x 0.32mm x 20µm
	19095P-QO3PT	HP-PLOT Q PT	15m x 0.53mm x 40µm
19095P-QO4PT	HP-PLOT Q PT	30m x 0.53mm x 40µm	
PLOT U	CP7584PT	PoraPLOT U PT	25m x 0.53mm x 20µm
	19095P-UO4PT	HP-PLOT U PT	30m x 0.53mm x 20µm
Al <sub>2</sub> O <sub>3</sub> KCl deactivated	CP7515PT	CP-Al <sub>2</sub> O <sub>3</sub> /KCl PT	50m x 0.32mm x 5µm
	CP7517PT	CP-Al <sub>2</sub> O <sub>3</sub> /KCl PT	25m x 0.53mm x 10µm
	CP7518PT	CP-Al <sub>2</sub> O <sub>3</sub> /KCl PT	50m x 0.53mm x 10µm
	19091P-K15PT	HP-PLOT Al <sub>2</sub> O <sub>3</sub> KCl PT	50m x 0.32mm x 8µm
	19095P-K23PT	HP-PLOT Al <sub>2</sub> O <sub>3</sub> KCl PT	30m x 0.53mm x 15µm
	19095P-K25PT	HP-PLOT Al <sub>2</sub> O <sub>3</sub> KCl PT	50m x 0.53mm x 15µm
	115-3352PT	GS-Alumina/KCl PT	50m x 0.53mm
Al <sub>2</sub> O <sub>3</sub> Na <sub>2</sub> SO <sub>4</sub> deactivated	CP7565PT	CP-Al <sub>2</sub> O <sub>3</sub> /Na <sub>2</sub> SO <sub>4</sub> PT	50m x 0.32mm x 5µm
	CP7568PT	CP-Al <sub>2</sub> O <sub>3</sub> /Na <sub>2</sub> SO <sub>4</sub> PT	50m x 0.53mm x 10µm
	19091P-S12PT	HP-PLOT Al <sub>2</sub> O <sub>3</sub> S PT	25m x 0.32mm x 8µm
	19091P-S15PT	HP-PLOT Al <sub>2</sub> O <sub>3</sub> S PT	50m x 0.32mm x 8µm
	19095P-S23PT	HP-PLOT Al <sub>2</sub> O <sub>3</sub> S PT	30m x 0.53mm x 15µm
	19095P-S25PT	HP-PLOT Al <sub>2</sub> O <sub>3</sub> S PT	50m x 0.53mm x 15µm
Al <sub>2</sub> O <sub>3</sub> with proprietary deactivation	115-3532PT	GS-Alumina PT	30m x 0.53mm
	115-3552PT	GS-Alumina PT	50m x 0.53mm
	19095P-M25PT	HP-PLOT Al <sub>2</sub> O <sub>3</sub> M PT	50m x 0.53mm x 15µm
Molsieve	CP7534PT	CP-Molsieve 5A PT	30m x 0.32mm x 10µm
	CP7536PT	CP-Molsieve 5A PT	25m x 0.32mm x 30µm
	CP7538PT	CP-Molsieve 5A PT	25m x 0.53mm x 50µm
	CP7539PT	CP-Molsieve 5A PT	50m x 0.53mm x 50µm

# Conclusions

The Application Specific columns are designed to help address some of the known issues with certain methods.

Can be used for other methods as well

**\*\*Samples still kill the columns\*\***

Activity Issues don't reside with the column only!

Depending upon the method, UI liners and gold seals can help reduce or eliminate activity issues

# Column Part Numbers

DB-CLP1 123-8232 30 m X 0.32 mm X 0.25  $\mu$ m

DB-CLP2 123-8336 30 m X 0.32 mm X 0.50  $\mu$ m

## DB-Select 624 UI for <467>

122-0334UI 30 m X 0.25 mm X 1.4  $\mu$ m

122-0364UI 60 m X 0.25 mm X 1.4  $\mu$ m

123-0334UI 30 m X 0.32 mm X 1.8  $\mu$ m

123-0364UI 60 m X 0.32 mm X 1.8  $\mu$ m

125-0334UI 30 m X 0.53 mm X 3.0  $\mu$ m

DB-UI 8270D 121-9723 20 m X 0.18 mm X 0.36  $\mu$ m

621-9723 Pack of 6

122-9732 30 m X 0.25 mm X 0.25  $\mu$ m

622-9732 Pack of 6

122-9736 30 m X 0.25 mm X 0.50  $\mu$ m



# Column Part Numbers

## DB-Sulfur SCD

G3903-63001 60 m X 0.32 mm X 4.2  $\mu$ m

G3903-63002 40 m X 0.32 mm X 0.75  $\mu$ m

G3903-63003 70 m X 0.53 mm X 4.3  $\mu$ m

G3903-63004 40 m X 0.32 mm X 3  $\mu$ m

# Application Notes

Residual Solvent Analysis with a Specifically Designed and Tested Agilent J&W DB-Select 624UI for USP <467> Column

**5991-0616EN**

Evaluating CLP and EPA Methods for Pesticides in Water Using Agilent J&W DB-CLP1/DB-CLP2 GC Columns

**5991-0615EN**

Semivolatile Analysis with Specially Designed Agilent J&W DB-UI 8270D Columns

**5991-0250EN**

Analysis of Sulfur Compounds in a Petroleum Fraction using Agilent J&W DB-Sulfur SCD GC Column and Inert Flow Path

**5991-3108EN**

# Agilent J&W Scientific Technical Support

**800-227-9770 (phone: US & Canada)\***

*\* Select option 3, then 3, then 1.*

**866-422-5571 (fax)**

**GC-Column-support@agilent.com**



**www.chem.agilent.com**

