

# Separation Superhero Seminar Tour

Battle the villains of bad chromatography



Solutions for light hydrocarbons and gasses:  
PLOT columns

# Solutions for Light Hydrocarbons and Gasses: PLOT Columns



Allen K. Vickers

Agilent Technologies


Folsom Technology Center

California, USA


April 18, 2019

# Table of Boiling Point Fractions

Carbon No.	C <sub>1</sub>	C <sub>2</sub>	C <sub>3</sub>	C <sub>4</sub>	C <sub>5</sub>	C <sub>6</sub>	C <sub>7</sub>	C <sub>8</sub>	C <sub>9</sub>	C <sub>10</sub>	C <sub>11</sub>	C <sub>12</sub>	C <sub>13</sub>	C <sub>14</sub>	C <sub>15</sub>	C <sub>16</sub>
Bpt of n-Paraffin at 760 mm Hg																
Centigrade	-161	-89	-42	-0.5	+36	69	98	126	151	174	196	216	235	253	270	287
Fahrenheit	-259	-127	-44	+31	97	156	209	258	303	345	384	421	421	488	519	548




# Wall Coated Open Tubular (WCOT) Columns



$K_c \Rightarrow \text{Large}$

$K_c \Rightarrow \text{Small}$

# WCOT Ethylene Analysis




# PLOT Ethylene Analysis

## Ethylene

**Column:** GS-Alumina  
50 m x 0.53 mm I.D.  
**J&W P/N:** 115-3552  
**Carrier:** Helium at 11 mL/min, measured at 35°C  
**Oven:** 35°C for 2 min  
35-190°C at 6°/min  
190°C for 3 min  
**Injector:** Split 1:30, 200°C  
0.2 mL of trace hydrocarbons in ethylene  
**Detector:** FID, 200°C  
Nitrogen makeup gas at 20 mL/min



1. Methane
2. Ethane
3. Ethylene
4. Propane
5. Cyclopropane
6. Propylene
7. Isobutane
8. Acetylene
9. *n*-Butane
10. Propadiene
11. *trans*-2-Butene
12. 1-Butene
13. Isobutylene
14. *cis*-2-Butene
15. Isopentane
16. *n*-Pentane
17. Propyne
18. 1,3-Butadiene
19. 1-Pentene
20. *n*-Hexane



C585

# What Is a PLOT Column?

- Porous layer (surface) on the inner wall of the capillary
- Porosity achieved by the deposition of porous particles on the wall from a suspension
- Porous layer serves as stationary phase
- Separation mechanism (gas-solid chromatography) differences in analyte distribution between carrier gas and the surface of the adsorbent (stationary phase)
- Separation mechanism for gas-liquid chromatography differences in analyte solubility in liquid phase (stationary phase)
- Benefit is higher selectivity ( $\alpha$ ) and retention ( $k$ ) of highly volatile solutes



# Selectivity Interactions in PLOT Phases





# Porous Layer Open Tubular (PLOT) Columns

## Advantages

- Retention for high vapor pressure solutes
  - No cryogenics needed
  - Useful for multicolumn and valve switching, heart cutting techniques.
- Variety adsorbents
  - Porous polymers, alumina, molecular sieve, and so on
  - Selectivity for isomeric compounds
  - Gases and solvents separations
- Efficiency capillary PLOT
  - Preferred over packed
  - Packed benefits in specific applications

# Agilent PLOT Columns Portfolio

## Porous Polymers

HP-PLOT Q, U  
PoraBOND Q, U  
PoraPLOT Q, U, S  
GS-Q  
PoraPLOT Q-HT

## Zeolites

HP-Molesieve  
CP- Molsieve 5A

## Oxygenates

Lowox  
GS-OxyPLOT

## Alumina

Al<sub>2</sub>O<sub>3</sub>/KCl, HP-PLOT Al<sub>2</sub>O<sub>3</sub>/KCl  
Al<sub>2</sub>O<sub>3</sub>/Na<sub>2</sub>SO<sub>4</sub>, HP-PLOT Al<sub>2</sub>O<sub>3</sub> S  
GS-Alumina

## Select MAPD

## Silica Porous

GS-GasPro  
SilicaPLOT

## Graphitized Carbon

CarboPLOT P7  
CarboBOND  
GS-CarbonPLOT

# Compounds and Columns

- Permanent gases
  - Noble gases, O<sub>2</sub>, N<sub>2</sub>, H<sub>2</sub>, CO, CO<sub>2</sub>, CH<sub>4</sub>
- Light hydrocarbons
  - C1 – C8, saturated, unsaturated
  - Aromatics, natural gas, C2, C3, C4 streams
- Volatile sulfur compounds
  - H<sub>2</sub>S, COS, mercaptans
  - Sulfides, disulfides
- Oxygenates
  - Alcohols, ketones, ethers
- Solvents
  - Oxygenates, aromatics, alkanes
  - Chlorinated hydrocarbons
- Chlorinated and fluorinated hydrocarbons

Molsieve

Select Perm. Gases

Alumina

Silica

Carbon

Porous polymers

Select Low Sulfur

Silica

Porous polymers

Lowox

OxyPLOT

Porous polymers

Silica

Porous polymers

# Porous Polymer PLOT Columns


- Analyses of polar and nonpolar volatile compounds
  - Oxygenates, gases, halogenated compounds, hydrocarbons C1 – C6, ketones, solvents
- Most versatile adsorbent materials
- Suitable for aqueous injections
- Elution of water as a sharp peak + quantification of water
- Recommended for column switching systems
- Divinyl benzene copolymer = nonpolar Q type
- Styrene-glycol methacrylate copolymer = polar U type

# Porous Polymer PLOT Columns

- **HP-PLOT Q, GS-Q and PoraPLOT Q**
  - Separation ethane/ethylene, propane/propylene
  - % levels, not ppm impurities in C2, C3 matrix
- **H<sub>2</sub>S and COS**
  - Refinery gas analysis
  - 100 ppm – %, not low ppm
- **Precolumn with CP-Molsieve columns for O<sub>2</sub>, N<sub>2</sub>, CH<sub>4</sub>, CO**
- **Poor C4 isomer separation**
  - Alumina preferred
- **Not sensitive to moisture**
- **Good choice for solvents and CFCs**


# Refinery Fuel Gas and Porous Polymer Separation

## HP- PLOT Q




# Porous Polymer PLOT, Aqueous Injections

Before




After

5 x 5  $\mu$ L water, **splitless**, at 80 °C



Retention times are the same for all compounds




Samples containing water can be analyzed

Isothermal and short cycle time


# Agilent "BOND" Technology

## Particle PLOT Q




**Particle size: 0.1 - 2  $\mu\text{m}$**   
**Number:  $>10^{12}$**

## PoraBOND Q



**Bonded polymer layer**  
**No particle traps**






# Solvent Analysis on PoraBOND Q

Column : CP-PoraBOND Q fused silica PLOT  
25 m x 0.53 mm, df = 10 µm, Cat. no. 7354  
Temperature : 100 °C (2 min) → 300 °C, 5 °C/min  
Carrier gas : He, 25 kPa (0.25 bar, 3.5 psi)  
Injector : Split, T = 250 °C  
Detector : FID, T = 250 °C  
Sample size : 0.5 µl  
Concentration range : 0.1% per compound  
Solvent sample : DMSO

## Peak identification:

- |                             |                           |
|-----------------------------|---------------------------|
| 1. methane                  | 8. dimethyl sulfide       |
| 2. methanol                 | 9. diethyl ether          |
| 3. ethanol                  | 10. 1-propanol            |
| 4. acetonitrile             | 11. pentane               |
| 5. acetone                  | 12. 2-butanone            |
| 6. dichloromethane          | 13. trichloromethane      |
| 7. 2-propanol (isopropanol) | 14. tetrahydrofuran       |
|                             | 15. ethyl acetate         |
|                             | 16. 2-methoxyethanol      |
|                             | 17. isobutanol            |
|                             | 18. butanol               |
|                             | 19. hexane                |
|                             | 20. benzene               |
|                             | 21. trichloroethylene     |
|                             | 22. cyclohexane           |
|                             | 23. 1,4-dioxane           |
|                             | 24. pyridine              |
|                             | 25. N,N-dimethylformamide |
|                             | 26. dimethyl sulfoxide    |



# Alumina Adsorbent and PLOT Columns

## **Best selectivity for hydrocarbon separations**

- **General C1 – C6 (C9) hydrocarbons**
- **Natural gas**
- **Ethylene streams, impurities**
- **Impurities in propylene**
- **Butylene streams, impurities and complex C4 composition**
- **Environmental hydrocarbons distributions**


# Alumina Adsorbent and Columns

Separation depends on:


- **Degree of hydrocarbon saturation**
  - Elution order: alkane, alkene, alkyne, (dialkenes)
- **Types of deactivation**
  - KCl, Na<sub>2</sub>SO<sub>4</sub> and proprietary
- **Operating conditions**
  - Column flow and oven temperature
- **Presence of water**
  - Al<sub>2</sub>O<sub>3</sub> columns are sensitive to moisture in carrier gas: Gas Clean filters

# Selectivity Difference KCl and Na<sub>2</sub>SO<sub>4</sub>

KCl



Na<sub>2</sub>SO<sub>4</sub>




1. Methane
2. Ethane
3. Ethene
4. Propane
5. Propene
6. Cyclopropane
7. Ethyne
8. Iso-butane
9. Propadiene
10. n-Butane
11. t-2-Butene
12. 1-Butene
13. Iso-butene
14. c-2-Butene
15. 1,3-Butadiene
16. Propyne


# Impurities in 1,3 Butadiene

50 m x 0.32 mm Al<sub>2</sub>O<sub>3</sub>/KCl, 5 μm


100 °C – 200 °C, 6 °C/min




# Select Alumina Column



# Response of Hydrocarbons



# Improved Response on Select for MAPD



Courtesy: J. Luong, Dow Chemical Canada



# Silica

- **Light hydrocarbon separation, C1 – C4**
- **Extended hydrocarbon range compared to other PLOT substrates (+C10)**
- **Inert enough for light sulfurs, H<sub>2</sub>S, COS, mercaptans**
- **CFCs**
- **GS-GasPro**
- **CP-SilicaPLOT**
- **Not sensitive to moisture in carrier gas**
- **MSD compatible**


# GS-GasPro: Inorganic Gases

## Inorganic Gases

**Column:** GS-GasPro  
30 m x 0.32 mm I.D.  
**J&W P/N:** 113-4332

**Carrier:** Helium at 53 cm/sec  
**Oven:** 25°C for 3 min  
25-200°C at 10°/min  
200°C Hold

**Injector:** Split 1:50, 200°C, 50 µL  
**Detector:** TCD, 250°C




# GS-GasPro: Halocarbons

## Halocarbons

**Column:** GS-GasPro  
30 m x 0.32 mm I.D.  
**J&W P/N:** 113-4332  
**Carrier:** Helium at 30 cm/sec  
**Oven:** 130°C for 4 min  
130-225°C at 10°/min  
225°C Hold  
**Injector:** Split 1:67, 1 µL, 250°C  
**Detector:** FID, 250°C

1. CH<sub>4</sub>
2. CHClF<sub>2</sub> (Freon 22)
3. CCl<sub>2</sub>F<sub>2</sub> (Freon 12)
4. ClCF<sub>2</sub>CF<sub>2</sub>Cl (Freon 114)
5. CHCl<sub>2</sub>F (Freon 21)
6. CCl<sub>3</sub>F (Freon 11)
7. CF<sub>2</sub>Br<sub>2</sub> (Freon 12B2)
8. CH<sub>3</sub>I
9. CH<sub>2</sub>Cl<sub>2</sub>
10. *trans*-ClCH=CHCl
11. CF<sub>3</sub>CCl<sub>3</sub> (Freon 113)
12. *cis*-ClCH=CHCl
13. CHCl<sub>3</sub>
14. ? from CCl<sub>4</sub>
15. CCl<sub>4</sub>
16. ? from CCl<sub>4</sub>
17. CH<sub>3</sub>CH<sub>2</sub>I
18. CH<sub>2</sub>Br<sub>2</sub>
19. CHCl<sub>2</sub>Br
20. C<sub>4</sub>F<sub>9</sub>I
21. CHClBr<sub>2</sub>
22. CH<sub>3</sub>CH<sub>2</sub>CH<sub>2</sub>I



# Molecular Sieves

- **Noble gases, Kr, He, Ar, Xe**
- **N<sub>2</sub>, O<sub>2</sub>, H<sub>2</sub>, CH<sub>4</sub>, CO (not CO<sub>2</sub>)**
  
- **CP-Molsieve 5 A**
- **HP-Molesieve**
  
- **Very sensitive to moisture/CO<sub>2</sub> in carrier gas/sample**
  - **Conditioning needed**
  - **Gas Clean moisture filters in carrier gas essential**

# PLOT Columns for Oxygenate Analysis


- Light oxygenates in hydrocarbon streams (ppm)
- C1 – C5 alcohols, ethers, ketones
- One column solution for ethylene, propylene, butylene feeds
  
- CP-Lowox
- GS-OxyPLOT
  
- Highly sensitive to moisture in carrier gas
  - Conditioning needed
  - Gas Clean moisture filters essential

# Oxygenates in Hydrocarbons

- Oxygenates blended in gasoline
  - Additives to boost octane content, prevent engine “knocking”
  - MTBE, ETBE, ethanol
  - % level oxygenates
  - GC analyses on WAX or TCEP polar liquid phase columns
- Oxygenate in intermediates (monomers, naphthas)
- Lower catalyst effectiveness, lower yield
  - Higher catalyst costs
  - More refinery downtime
  - ppm level oxygenates
  - GC analyses on Lowox, OxyPLOT columns, Agilent exclusives


# Analyzing Oxygenates in Hydrocarbon Matrix

- Low ppm concentration level oxygenates
- FID detection (MS rarely)
- High selectivity columns hydrocarbons/oxygenates




# Rotary Valve Configuration for ASTM Trace Oxygenate Analysis Methods

**Transfer of Oxygenates  
Valve Off**




**Venting Hydrocarbons  
Valve On**





# Hydrocarbons and Oxygenates Separation using DB-1 Stripper Column and GS-OxyPLOT Separation Column



# GS-CarbonPLOT

Monolithic carbon molecular sieve

Phase formed *in situ*

Extended temperature limit of 360 °C

Unique selectivity

# C1 to C3 Hydrocarbon Split Injector

**GS-CarbonPLOT**

**30 m x 0.32 mm id, 1.5  $\mu$ m**

**Carrier: He at 30 cm/s**


**Oven: 125 °C for 2 min**

**125 °C to 225 °C at 20 °C/min**

**Injector: Split 1:20, 250 °C**

**Detector: MSD, SIM**

- 1. Methane**
- 2. Acetylene**
- 3. Ethylene**
- 4. Ethane**
- 5. Propyne**
- 6. Propylene**
- 7. Propane**





# PLOT Column with Integrated Particle Trap

- **Zero particle shedding when using steep temperature gradients and pressure ranges**
- **Integrated to column – no unions and fittings**
- **Compatible with multicolumn valve switching systems and systems with CFT technology**
  - Particle traps integrated on both ends – supports backflush apps
- **Similar selectivity, plates, and peak shape performance to existing Agilent porous polymer PLOT columns**
  - Minimum method redevelopment required
  - Available in a wide variety of our most popular porous polymer PLOT columns configurations

# Proof of Agilent PLOT PT Column Performance


- 150 °C to 250 °C at 20 °C/min
- Pressure 3x higher than optimal pressure
- 15 cycles
- Carrier gas off and on 10 times



# PLOT Column with Integrated Particle Trap

**Integrated particle trap at front and back of the column**

**PLOT column part**



# Example of Chromatographic performance

Column: PoraBond Q PT, 30 m × 0.25 mm, 3 μm

Carrier: Helium, 43 cm/s at 90 °C

Oven: 90 °C – 140 °C at 10 °C/min

140 °C for 6 min

140 °C – 200 °C at 5 °C/min


200 °C for 10 min

Injection: Split, 250 °C, split ratio 1:160

Detector: MSD, 280 °C

Transfer line, full scan at m/z 30-350

- |                       |                           |                               |
|-----------------------|---------------------------|-------------------------------|
| 1. Methyl Alcohol     | 14. Trichloromethane      | 27. Pyridine                  |
| 2. Acetaldehyde       | 15. 2-Butanone (MEK)      | 28. Dimethyl Formamide (DMF)  |
| 3. Ethanol            | 16. Ethyl Acetate         | 29. Isoamyl Alcohol           |
| 4. Acetonitrile       | 17. sec-Butyl alcohol     | 30. Dimethyl Sulfoxide (DMSO) |
| 5. Acetone            | 18. MTBE                  | 31. Toluene                   |
| 6. Methylene Chloride | 19. 2-chlorobutane        | 32. Heptane                   |
| 7. Isopropyl Alcohol  | 20. 1-Butanol             | 33. Paraldehyde               |
| 8. 2-Propanamine      | 21. Benzene               | 34. Chlorobenzene             |
| 9. Ethyl Formate      | 22. 1,1,1-Trichloroethane | 35. Ethylbenzene              |
| 10. 1-Propanol        | 23. 1-chlorobutane        | 36. m-Xylene                  |
| 11. Ethyl ether       | 24. Carbon Tetrachloride  | 37. p-Xylene                  |
| 12. t-Butyl alcohol   | 25. Hexane                | 38. o-xylene                  |
| 13. 1,2-Ethanediol    | 26. 1,4-Dioxane           |                               |



# Considerations for PLOT Column Analysis

- Inlet issues
  - Split versus direct injection
  - Gas sampling valves
  - Low dead volume
  - Column id and flow rate
- Detector issues
  - Particle generation or “spiking”; particle traps
  - Column id and flow rate



# Considerations for PLOT Column Analysis

- Column issues
  - Selectivity
  - Capacity; overloaded peaks
  - Inertness
  - Temperature limits
- Column contamination
  - Efficiency loss; “ghost peaks”; increase in bleed
  - Water, CO<sub>2</sub>, high molecular weight hydrocarbons?
  - Carrier gas purifiers

# Gas Clean Offer More Reliable Data

## Bottom Line:

**By removing impurities from the carrier gas stream, Gas Clean filters improve separation and repeatability while using PLOT columns.**

## We suggest:

**For a 7890 GC:**

**CP17988 Gas Clean carrier gas kit (1/8 in fittings)**

**For a 6890 GC:**

**CP17974 Gas Clean carrier gas kit (1/8 in fittings)**



**Agilent 7890 Gas  
Clean carrier gas  
kit installed on a  
7890B GC**

# Summary

- Agilent supplies the largest selection of PLOT columns in the market for all gases and volatiles applications. Agilent has dedicated columns for challenging analyses in the petrochemical industry.
- Fully QC tested to assure column-to-column reproducibility with excellent peak shape performance and separation for the best data accuracy.
- Agilent PLOT columns come with the lowest levels of particle shedding for better baseline stability and trouble free analyses.
- The PLOT-PT columns with integrated particle traps to assure “spike free” detection, mass spectrometer compatibility, and improved system performance with complex valve applications.

# Contact Agilent Chemistries and Supplies Technical Support



**1-800-227-9770 option 3, option 3:**

**Option 1 for GC and GC/MS columns and supplies**

**Option 2 for LC and LC/MS columns and supplies**

**Option 3 for sample preparation, filtration and QuEChERS**

**Option 4 for spectroscopy supplies**

**Option 5 for chemical standards**

**Available in the USA and Canada 8-5 all time zones**

**[gc-column-support@agilent.com](mailto:gc-column-support@agilent.com)**

**[lc-column-support@agilent.com](mailto:lc-column-support@agilent.com)**

**[spp-support@agilent.com](mailto:spp-support@agilent.com)**

**[spectro-supplies-support@agilent.com](mailto:spectro-supplies-support@agilent.com)**

**[chem-standards-support@agilent.com](mailto:chem-standards-support@agilent.com)**