



Agilent MassHunter Workstation Software – 7250 Accurate-Mass Quadrupole Time of Flight GC/MS

Familiarization Guide

Before you begin	3
Prepare your system	3
Prepare the samples required for data acquisition	4
Exercise – Develop an acquisition method for the 7250	5
Task 1. Set the inlet and injection parameters	5
Task 2. Check the GC Configuration	7
Task 3. Optimize Base Ion Abundance and Perform a Mass Calibration	9
Task 4. Enter GC acquisition parameters	13
Task 5. Create a Qual acquisition method for scanning ions	17
Task 6. Acquire MS scan data (Optional)	19
Task 7. Using a Sequence to Schedule Mass Calibrations	21
Exercise – Analyze data	23
Task 1. Open a data file in the Qualitative Navigator program	23
Task 2. Configure the qualitative analysis user interface	25
Task 3. Identify Peaks in Qualitative Navigator	28
Task 4. Identify compounds using Qualitative Analysis Workflows	30
Task 5. Configure Method Automation Reports	32
Task 6. Generate the Method Automation workflow	34
Task 7. Review the results	36
Task 8. Identify compounds using Unknowns Analysis	40



Reference material

This guide shows how to use the Agilent 7250 Q-TOF GC/MS System to acquire and analyze sample data. If you want to skip the data acquisition steps in this guide, use the demo data files shipped with MassHunter (See the Reference Material section below).

In this guide, you learn how to determine the best acquisition settings for analyzing your compounds of interest. These instructions help you understand not only how to set up a method to optimize instrument parameters for best sensitivity in acquisition, but also how to use the Qualitative Analysis program to identify parameter values producing optimum signal response.

Reference material

Available on the 7250 QTOF installation of the MassHunter Software User Information. This documentation application is delivered to you as part of your MassHunter Software.

- *MassHunter Qualitative Analysis Familiarization Guide for GC/MS* for an introduction to many Qualitative Analysis program features than are not covered in this guide including the Qualitative Analysis Navigator program.
- *Qualitative Analysis Training Videos* for those seeking visual and audio lessons presenting a comprehensive use of the MassHunter Qualitative Analysis Navigator and MassHunter Qualitative Workflows.
- *On-line Help* for detailed information on how Qualitative Analysis works.
- *Demo data files and accurate mass library* that allow you to perform all the analysis steps demonstrated here using you own installation of Qualitative Analysis without acquiring compound data or owning a library license.
- *Quick Start Guide* explains what documentation is included in the applications and what information is on each document.

Available on the 7250 QTOF installation of the User Manuals and Tools DVD. This documentation application is delivered to you with your 7250 QTOF instrument.

- *Concepts Guide* to learn more about how the 7250 Q-TOF GC/MS System works.
- *Quick Start Guide* explains what documentation is included in the applications and what information is on each document.
- *Hardware manuals* to learn how to operate and perform maintenance on the 7250 Q-TOF.

Before you begin

In this manual, each task is presented in a table with three columns:

- Steps – Use these general instructions to proceed on your own to explore the program.
- Detailed instructions – Use these if you need help or prefer to use a step-by-step learning process.
- Comments – Read these to learn tips and additional information about each step in the exercise.

Before you begin, you need to check that your system is ready. If you plan to acquire data, you also need to set up the instrument.

Prepare your system

- 1 Check that:
 - MassHunter Acquisition, MassHunter Qualitative Analysis, and MassHunter Quantitative Analysis are installed.
 - Your system uses an Agilent 7890 GC with split/splitless or MultiMode (MMI) inlet and automatic liquid sampler.
 - The acquisition uses a 10 uL ALS syringe tapered, fixed, with 23-26s needle. A suitable syringe may be substituted.
 - The 7250 Q-TOF GC/MS System is configured and has a valid tune.
 - The performance is verified.
 - The system is turned on.
 - A suitable column is installed. The J&W model 122-3832 DB-35MS: 30 m x 250 μ m, 0.25 μ m column is used for the examples in this guide.
- 2 Configure the GC for the installed column.
- 3 If needed, copy the *Demo data files and accurate mass library* noted in “Reference material” on page 2 to any location on your hard disk. This data file and accurate mass library file are needed for this exercise if you are not acquiring data and do not have an accurate mass library of the compounds shown in “Sample Compound list” on page 4.

Before you begin

Prepare the samples required for data acquisition

Prepare the samples required for data acquisition

If you do not intend to acquire data but want to learn how to use the Qualitative Analysis program you can skip the sample preparation and actual acquisition and use the data file shipped with this guide. It is recommended that you read the “Exercise – Develop an acquisition method for the 7250” on page 5 to understand settings unique to the Agilent instrument.

Materials required for sample preparation:

- Sample (p/n 05970-60045 or p/n 5074-3025 Japan only)
- Isooctane for sample dilution
- Sample vials

The sample compounds are in an isooctane solvent contained in 1 mL ampules of 10 ng/μL, 100 ng/μL, and 100 pg/μL concentrations and are shown in Table 1.

Table 1 Sample Compound list

Compound	<i>m/z</i>	Formula
Dodecane	170.2029	C ₁₂ H ₂₆
Biphenyl	154.0777	C ₁₂ H ₁₀
4-Chlorobiphenyl (p/n 05970-60045 only)	188.0387	C ₁₂ H ₉ Cl
Methyl palmitate	270.2553	C ₁₇ H ₃₄ O ₂

Prepare the Qualitative Analysis sample by emptying the contents of the 10 ng/μL ampoule into an ALS sample vial and cap the vial.

Fill an ALS wash vial with isooctane.

Exercise – Develop an acquisition method for the 7250

Task 1. Set the inlet and injection parameters

Steps	Detailed instructions	Comments
1 Set up the inlet, injection source, and enable the 7250.	<ol style="list-style-type: none"> Double-click the Data Acquisition icon on the windows desktop. Click the Inlet and Injection Parameters icon. Select GC for the sample inlet and the installed ALS for the injection source. Select the Use MS check box. 	<ul style="list-style-type: none"> The Data Acquisition window shown in Figure 1 is displayed. Hover over an icon to display a tag identifying the icon. The Inlet and Injection Parameters dialog box shown in Figure 2 on page 6 is displayed.




Figure 1 Agilent MassHunter Workstation Software – Data Acquisition window

Exercise – Develop an acquisition method for the 7250

Task 1. Set the inlet and injection parameters

Inlet and Injection Parameters

Sample Inlet: GC

Injection Source: GC ALS

Use MS

Inlet Location: Front Rear Dual

MS Connected to: Front Rear

OK Cancel Help

Figure 2 Inlet and Injection Parameters

Task 2. Check the GC Configuration

In this exercise, you review the GC hardware setup for the analysis.

Steps	Detailed instructions	Comments
1 Check that the GC hardware configuration is suitable for the analysis.	<p>a Click the GC Edit Parameters icon.</p> <p>b In the navigation menu select Configuration > Miscellaneous.</p> <p>c Set the Pressure Units to psi.</p> <p>d In the Oven area the Slow Fan mode is unchecked.</p> <p>e In the navigation menu select Configuration > Columns and set Column 1 to a J&W 122-3832 column or one that is similar. Set the Inlet to Front (or Rear) Inlet and the Outlet to MSD. Heated By is set to Oven.</p> <p>f In the navigation menu select Configuration > Modules and set the SS inlet gas to He and the Collision Cell EPC gas to N2.</p> <p>g In the navigation menu select Configuration > ALS and set the Syringe Size to 10 uL and the Solvent Wash Mode to A, B.</p> <p>h Select the OK button.</p>	<ul style="list-style-type: none"> • See Figure 1. The GC edit parameters window shown in Figure 3 is displayed. • If using a different column you must adjust your GC parameter settings accordingly for acceptable chromatography. • 10 uL ALS syringe tapered, fixed, with 23-26s needle. A suitable syringe may be substituted. • The GC parameters are downloaded to the GC and the window closes.

Exercise – Develop an acquisition method for the 7250

Task 2. Check the GC Configuration

The screenshot shows the 'GC Edit Parameters' dialog box. On the left is a tree view with 'Miscellaneous' selected. The main area is divided into 'Pressure Units' and 'Valve Configuration'.

Pressure Units: The unit is set to 'psi'. There is an 'Oven' section with a 'Slow Fan' checkbox. Below is a table for 'Thermal Aux Type':

	Thermal Aux Type
1	Not Installed
2	MSD Transfer Line
3	Not Installed

Valve Configuration: A table with columns for Valve #, Valve Type, Name, and Parameters.

	Valve Type	Name	Parameters
1	Not Installed	(Valve #1)	
2	Not Installed	(Valve #2)	
3	Not Installed	(Valve #3)	
4	Not Installed	(Valve #4)	
5	Not Installed	(Valve #5)	
6	Not Installed	(Valve #6)	
7	Not Installed	(Valve #7)	
8	Not Installed	(Valve #8)	

Figure 3 The Configuration Settings

Task 3. Optimize Base Ion Abundance and Perform a Mass Calibration

In this exercise you optimize the abundance of the base ion and perform a mass calibration. A mass calibration is completed in less than two minutes and it is good practice to calibrate the instrument daily or even every couple of hours. A sequence table keyword allows automatic mass calibration between samples in a sequence. See the on-line help for more information.

Steps	Detailed instructions	Comments
1	<p>Set the m/z range of acquired data and the range of that data to store for analysis.</p> <p>a In MassHunter Instrument Control view click the MS tune icon.</p> <p>b Click the Manual Tune tab, then click the Acquisition tab.</p> <p>c Select 1 Hz from the Acquisition Rate dropdown.</p> <p>d Select Low from the Maximum Mass Range dropdown.</p> <p>e Enter 25 for the Low end of the range and 650 for the High.</p>	<ul style="list-style-type: none"> The GC/Q-TOF Tune window is displayed. See Figure 4 below. This is the rate used during calibration. Data will be scanned from 20 to 650 m/z. There are 2 other mass ranges available for scanning data. The Standard range for 20 to 1200 m/z and the Extended range for 20 to 3000 m/z. Here we select the Low range to get the highest sensitivity for our data. Acquired data between 25 and 650 m/z is displayed in the tune spectrum window.




Figure 4 Selecting the m/z range of acquired data

Exercise – Develop an acquisition method for the 7250

Task 3. Optimize Base Ion Abundance and Perform a Mass Calibration

Steps	Detailed instructions	Comments
2 Optimize the base ion abundance. This step is usually done before calibration is performed.	<p>a Click the Ion Source tab and in the Tune Masses area select Enabled.</p> <p>b Enable the Emission and EI Cal Valve.</p> <p>c Adjust the Emission current so that the abundance of the ion of interest is ideally between 0.8×10^6 and 1.2×10^6 counts.</p>	<ul style="list-style-type: none"> • See Figure 5 on page 10. • To enable calibrant flow ionization. • Higher values will saturate the signal and lower values will not provide sufficient ion statistics for optimal mass accuracy.




Figure 5 Optimizing base ion abundance

Exercise – Develop an acquisition method for the 7250
Task 3. Optimize Base Ion Abundance and Perform a Mass Calibration

Steps	Detailed instructions	Comments
3 Perform a mass calibration.	<p>a Select the TOF Mass Calibration tab from the GC/Q-TOF Tune window.</p> <p>b Click the Run Calibration button.</p> <p>c Select the Close button.</p> <p>d Click the File and Reports tab and save the tune file.</p> <p>e Select the Close button.</p>	<ul style="list-style-type: none"> • See Figure 6 on page 11. Also see “Task 6. Acquire MS scan data (Optional)” on page 19. • When the calibration completes the TOF Mass Calibration Results window displays. Mass Accuracy (PPM) should typically be below 1PPM for all ions used in calibration. See Figure 7 on page 12. • The TOF Mass Calibration Results closes. • Save the tune file as atune-lowmz_ <i>date</i>.ei.tune.xml. Where <i>date</i> is today’s date. • The GC/Q-TOF Tune window closes and you are returned to Instrument Control view.





Figure 6 TOF Mass Calibration tab

Exercise – Develop an acquisition method for the 7250

Task 3. Optimize Base Ion Abundance and Perform a Mass Calibration

TOF Mass Calibration Results



A = 0.000345801089478936, T0 = 1240.52592904741

✓	Target Mass	Actual Mass	Accuracy (PPM)	Previous Mass	Previous
✓	68.9947	68.9947	0.01	68.9949	3.27
✓	130.9915	130.9915	-0.03	130.9918	2.29
✓	218.9851	218.9851	0.09	218.9856	2.25
✓	263.9866	263.9866	-0.07	263.9871	2.03
✓	413.9770	413.9769	-0.12	413.9776	1.62
✓	463.9738	463.9739	0.27	463.9747	2.06
✓	501.9705	501.9705	-0.17	501.9715	1.82
✓	613.9642	613.9642	0.01	613.9666	3.99

For enabled m/z values over 50, average PPM error 0.00, maximum PPM error 0.3
limits for average PPM error 3.0, maximum PPM error 8.0

Show Detailed Chart

Figure 7 TOF Mass Calibration Results

Task 4. Enter GC acquisition parameters

In this exercise, you enter the GC conditions for the analysis.

Steps	Detailed instructions	Comments
1	<p>Enter GC parameters appropriate for the sample. See Table 2.</p> <p>a Click the GC Edit Parameters icon (Figure 1).</p> <p>b In the navigation menu select Columns then select column 1 in the Selection column.</p> <p>c Select control mode On and then select Constant Flow mode. Enter 1.1 mL/min for the initial Flow.</p> <p>d Select the Collision Cell EPC in the Selection column and then in the Collision Cell EPC area, set the N2 Collision Gas on at 1.5 mL/min.</p> <p>e In the Collision Cell EPC area, uncheck the He Quench Gas.</p> <p>f In the navigation menu select Inlets > SSL and enter the inlet parameters listed in Table 2.</p> <p>g Select the Oven icon and enter the oven parameters listed in Table 2.</p> <p>h In the navigation menu select ALS > Front Injector and enter the injector parameters listed in Table 2.</p> <p>i In the navigation menu select Aux Heaters, enable, and set the temperature to 280 °C.</p> <p>j Select the OK button.</p>	<ul style="list-style-type: none"> • The GC edit parameters window shown in Figure 8 on page 14 is displayed. • With the window selected, mouse over the icons to identify the icon from the tool tip. • If the current flow value of the collision cell N2 gas is not 1.5 mL/min and you change it to this value, an autotune will be required. • If your ALS is attached to the Back Inlet select the Back Injector tab. • This is the MSD transfer line heater. • The GC parameters are downloaded to the GC and the window closes.

Exercise – Develop an acquisition method for the 7250

Task 4. Enter GC acquisition parameters



Figure 8 GC Edit Parameters window with **Oven** icon selected

Table 2 GC parameters for data acquisition method

Parameter	Value
Oven	
Equilibration Time	0.1 min
Oven Program	80 °C for 3 min, 25 °C/min to 250 °C, hold for 2.2 min
Run Time	12 min
Front SS Inlet	
Mode	Split
Heater	On 250 °C
Pressure	On Value automatically set with column flow
Septum Purge Flow	On 3 mL/min
Gas Saver	On 20 mL/min after 3 min
Split Flow	220 mL/min
Split Ratio	200:1
Thermal Aux 2 {MSD Transfer Line}	
Heater	On
Temperature	280 °C
Column # 1	
In	Front SS Inlet He
Out	Vacuum
(Initial)	80 °C
Flow	1.1 mL/min
Flow Program	Off
Front Injector	
Syringe Size	10 µL
Injection Volume	1 µL
Solvent A Washes (Prelinj)	2


Exercise – Develop an acquisition method for the 7250

Task 4. Enter GC acquisition parameters

Parameter	Value
Solvent A Washes (PostInj)	2
Solvent A Volume	8 μ L
Solvent B Washes (PreInj)	2
Solvent B Washes (PostInj)	2
Solvent B Volume	8
Sample Washes	0
Sample Wash Volume	8 μ L
Sample Pumps	4
Dwell Time (PreInj)	0 min
Dwell Time (PostInj)	0 min
Solvent Wash Draw Speed	300 μ L/min
Solvent Wash Dispense Speed	6000 μ L/min
Sample Wash Draw Speed	300 μ L/min
Sample Wash Dispense Speed	6000 μ L/min
Injection Dispense Speed	6000 μ L/min
Viscosity Delay	0 sec
Sample Depth	Disabled
Collision cell EPC Module	
Nitrogen	On 1.5 mL/min
Helium	Off

Task 5. Create a Qual acquisition method for scanning ions

This exercise starts with the GC parameters entered in the method from Task 4. In this task you will enter the 7250 parameters for ion scanning and save to the method.

Steps	Detailed instructions	Comments
1	<p>Enter MS parameters appropriate for the sample and save the method as OFN EI 70eV.M.</p> <p>a Click the QTOF Method Editor icon (Figure 1).</p> <p>b In the Tune file area, click the  icon.</p> <p>c In the Ion Source area, set the Source temperature to 200 °C, set the Emission to Fixed with a value of 5.0 entered, and set the Electron energy to Tune Setting.</p> <p>d Set the Solvent delay to 5 minutes.</p> <p>e In the Time Filtering area select Peak width and set it to 0.7 seconds.</p> <p>f In the Data Threshold area enter 10 for counts.</p> <p>g Select Apply to Profile data.</p> <p>h In the Time segment area, select a Scan Type of MS from the Acq mode drop-down list.</p> <p>i Select Both for Data stored.</p> <p>j In the MS mode section, for the Mass range enter 40 for the start mass, 600 for the end mass, and 5.00 spectra/s for the Acq rate.</p> <p>k Click OK to close the window.</p> <p>l From the main window select Method > Save Method As and save the method as OFN EI 70eV.M.</p>	<ul style="list-style-type: none"> The QTOF Method Editor window shown in Figure 9 on page 18 opens. Select the tune file created in Select the tune file created at the end of Task 3. The 7250 starts collecting data at 5 minutes due to the Solvent delay setting. Filter out unwanted peaks to reduce data storage. Filter out unwanted noise to reduce data storage. To apply the Data Threshold filter to profile data to reduce storage. If we were doing an MS/MS acquisition we would enter the counts here to reduce data storage. Selecting Both stores both a peak's profile data and centroid data for data analysis. All data up to 650 m/z is always acquired but only the data selected here (40 to 600) is saved to disk.

Exercise – Develop an acquisition method for the 7250
Task 5. Create a Qual acquisition method for scanning ions

QTOF Method Editor

Ion source

Ion source: EI

Source temp.

Tune Setting 200 °C

Fixed

EI Mode

Standard Low Energy

Emission

Tune Setting 1.3 0.5 μA

Fixed 5.0 μA

By time segment

Electron energy

Tune Setting 70.0 15.0 eV

Fixed eV

By time segment eV

Tune file

atunes.eihs.tune.xml

Run time 5 min

Solvent delay 6 min

Time filtering

Peak width 0.700 sec

Data Threshold

MS 10 counts 0 %

MS/MS 0 counts 0 %

Apply to Profile data

Time segments

Time	Acq mode	EI Mode	Emission	Electron energy	Data storage	Data stored
1	0.00 MS	Standard	5	70	Both	<input checked="" type="checkbox"/>

Timed events

Time	Type	Address	Value

Display Timed Events

Acquisition Reference Mass Instrument Chromatogram

MS mode

Quad TTI cutoff mass

Default 50 amu

Override amu

Mass range 40 to 600 amu

Acq rate 5.00 spectra/s

Acq time 200 ms/spectrum

Transients/spectrum 1941

OK Apply

Figure 9 QTOF Method Editor

Task 6. Acquire MS scan data (Optional)

In this task, you acquire the scan data using the method developed in the previous tasks. This task is optional because you can perform the next task with an example data file that is provided with MassHunter in the location shown in “Reference material” on page 2. However, if you prefer, you can acquire your own data file as described in this task.

Steps	Detailed instructions	Comments
1 Acquire data (optional). <ul style="list-style-type: none"> Name the data file Sol_A.D. Designate a directory path to hold your data files and method. 	<p>a Click the Start Run (green arrow) icon.</p> <p>b In the Data Path enter the directory to save the data file that is acquired by this run.</p> <p>c In the Inlet section you are using, enter Sol_A.D for the Data File Name.</p> <p>d Enter the Vial Number location in the auto sampler tray.</p> <p>e Select Current Method for the Injection Volume.</p> <p>f In the Method Sections to Run section, select Data Acquisition.</p> <p>g Click the OK and Run Method button.</p>	<ul style="list-style-type: none"> The Start Run dialog box shown in Figure 10 on page 20 is displayed. The Injection volume that you entered in Task 4 is used. The method is sent to the GC and the Q-TOF. When the instruments are ready the sample is injected and the data is collected and sent to the data directory specified.

Exercise – Develop an acquisition method for the 7250
Task 6. Acquire MS scan data (Optional)

Start Run

Basic Advanced

Current Method Injection Style: GC ALS

Inlet Location: Front Rear Dual MS Connected to: Front Inlet Rear Inlet

Operator Name:

Data Path: D:\MassHunter\GCMS\2\DATA

Front Inlet

Data File Name: SoL_A.D

Sample Name:

Misc. Info:

Expected Barcode:

Sample Amount: 0

Multiplier: 1

Vial Number: 1

Tray Name: Agilent ALS

Injection Volume: Current Method 1 µL Override using 0 µL

Rear Inlet

Data File Name: EVALDEMO.D

Sample Name:

Misc. Info:

Expected Barcode:

Sample Amount: 0

Multiplier: 1

Vial Number:

Tray Name: Agilent ALS

Injection Volume: Current Method 0 µL Override using µL

9,854,193,664 bytes free on drive D:

Method Sections to Run:

Data Acquisition

Data Analysis (MassHunter DA)

Figure 10 Start Run dialog box

Task 7. Using a Sequence to Schedule Mass Calibrations

This automated procedure is used to schedule mass calibrations at the start of a sequence of sample runs and at timed intervals during those runs. It is recommended to do a mass calibration about every 2 hours when continuously running samples. This mass calibration only takes a couple of minutes but allows you to maintain higher mass accuracy and immunity to drift.

Steps	Detailed instructions	Comments
1 Add a mass calibration at the start of sequence.	<ol style="list-style-type: none"> a Insert an entry for the running of the mass calibration. b Add the MASSCAL keyword to this entry. c Enter the Method used for processing samples that follow and select CAL from the Type dropdown. d Save the sequence. 	<ul style="list-style-type: none"> • This entry might follow the running of a sample blank. • See Figure 11 below. • This entry will also be used after running samples every 2 hours in our sequence.

	Name	Vial	Type	Keyword	Method Path	Method File	Data Path
1	Hexane	1	DoubleBlank		D:\MassHunter\GCMS\1\methods	OFN EI 70eV.m	D:\MassHunter\GCMS\1\data
2			Cal	MassCal	D:\MassHunter\GCMS\1\methods	OFN EI 70eV.m	D:\MassHunter\GCMS\1\data
3	OFN 100fg-1	5	Sample		D:\MassHunter\GCMS\1\methods	OFN EI 70eV.m	D:\MassHunter\GCMS\1\data
4	OFN 100fg-2	6	Sample			OFN EI 70eV.m	D:\MassHunter\GCMS\1\data
5	OFN 100fg-3	7	Sample			OFN EI 70eV.m	D:\MassHunter\GCMS\1\data
6	OFN 100fg-4	8	Sample			OFN EI 70eV.m	D:\MassHunter\GCMS\1\data

Figure 11 Sequence table setup for automated mass calibration

Exercise – Develop an acquisition method for the 7250
Task 7. Using a Sequence to Schedule Mass Calibrations

Steps	Detailed instructions	Comments
<p>2 You can skip this step and use the run time procedure shown in the step below instead.</p> <p>Use this step to add a mass calibration entry at 2 hour intervals during sample runs by creating mass calibration entries in the sequence table.</p>	<p>a Calculate the number of samples to run before performing the mass calibration.</p> <p>b Copy the mass calibration entry added in the above step.</p> <p>c Select the sample entry location for the mass calibration and select Insert sample from the context menu.</p> <p>d Click paste with this entry selected.</p> <p>e Repeat step 2 as required.</p> <p>f Save the sequence.</p>	<ul style="list-style-type: none"> • (120 min / <i>sample run time</i>) For a 10 minute sample run this means 12 samples can be processed before performing a mass calibration. • $2 + 12 + 1 =$ entry 15 for our example. An empty entry is created. • The copied calibration is entered.
<p>3 The same result for the mass calibration interval created in the above step can be obtained at run time as follows.</p>	<p>a Select Sequence > Run Sequence.</p> <p>b Fill out the required information in the Basic tab.</p> <p>c Click the Recurring Keyword tab.</p> <p>d Select Hours as the Recurring Type.</p> <p>e Enter 2 hours for the Recurring Interval.</p> <p>f Select MassCal as the Recurring Keyword string.</p> <p>g Click Run Sequence.</p>	<ul style="list-style-type: none"> • The Start Sequence dialog opens. • See online help for additional information. • You could also select Method change or Hours if you want to recalibrate when the method changes or at a timed interval. • The mass calibration runs at the specified interval. • The loaded sequence runs.





Figure 12 Using the Recurring Keyword feature

Exercise – Analyze data

In this exercise, you analyze data acquired from the previous exercises in this manual. If you did not acquire this data you can use the example data file **Sol_A.D** provided in the location shown in “Reference material” on page 2.

For additional details on using this program, see the *MassHunter Workstation Qualitative Analysis Familiarization Guide for GC/MS* and the *Qualitative Analysis Training Videos* that are provided with MassHunter in the location shown in “Reference material” on page 2.

Task 1. Open a data file in the Qualitative Navigator program

Steps	Detailed instructions	Comments
1 Start the Qualitative Analysis program.	<p>a Double-click Qualitative Navigator shortcut on your desktop.</p> 	<ul style="list-style-type: none"> The system displays the Open Data File dialog. See Figure 13 on page 24. You can get help by: <ul style="list-style-type: none"> Pressing the F1 key when a window is active Selecting Help > Contents in the main menu Selecting the Help icon in the active window
2 Load a Data file.	<p>b Navigate to the location where your data file is located and then select your acquired data file or the demo file provided for this exercise Sol_A.D.</p> <p>c Under Options, select Use current method and clear Load result data.</p> <p>d Click Open.</p>	<ul style="list-style-type: none"> See Figure 13 on page 24. The data file is loaded and a TIC of the data is displayed. See Figure 14 on page 24.

Exercise – Analyze data

Task 1. Open a data file in the Qualitative Navigator program




Figure 13 Opening the data file




Figure 14 Sample Chromatogram for Sol_A.d loaded

Task 2. Configure the qualitative analysis user interface

Steps	Detailed instructions	Comments
1 Create a Qualitative Analysis method.	<p>a From the main menu, select Method > Open.</p> <p>b Select default-GCMS.m and click Open.</p> <p>c From the main menu, select Configuration > Show Advanced Settings.</p> <p>d From the main menu, select Configuration > User Interface Configuration.</p> <p>e Click OK.</p> <p>f From the main menu, select Method > Save As.</p> <p>g In the File Name enter QTOF-GCMS.m and click OK.</p>	<ul style="list-style-type: none"> The Open Method dialog is displayed. The default method for a GC QTOF is loaded. Advanced menu items are now available. The User Configuration dialog box opens. See Figure 16 on page 26. The settings selected here are based on the data file loaded and the default method previously selected. You may edit these if required. To close the dialog. The Save Method dialog is displayed. Rename the method since you can not overwrite the default method.
2 Configure the MS and MS/MS Spectra Display options.	<p>a From the main menu, select Configuration > MS and MS/MS Spectra Display Options.</p> <p>b Click the Spectrum Peak Label Options tab and set the values for the peak labels.</p> <p>c Save the method.</p>	<ul style="list-style-type: none"> The MS and MS/MS Spectra Display Options dialog is displayed. Use the values displayed in Figure 17 on page 27. These values result in the horizontal display of m/z and Formula & Ion Species value above identified spectra.
3 Assign a library for identifying spectra.	<p>a Click the Method Editor icon to open the method editor and click on Identify Spectra > Identification Workflow.</p> <p>b Set the values for the library used to identify spectra. The demo.L library is included with MassHunter and includes the 4 compounds used in our sample</p>	<ul style="list-style-type: none"> For easier viewing, float the editor pane outside the Qualitative Analysis Navigator window. Use the values displayed in Figure 18 on page 27.

Exercise – Analyze data

Task 2. Configure the qualitative analysis user interface




Figure 15 Default layout for a GC QTOF




Figure 16 Configuring the User Interface

Task 2. Configure the qualitative analysis user interface




Figure 17 Configuring spectrum peak labels







Figure 18 Assign a library for identifying spectra

Task 3. Identify Peaks in Qualitative Navigator

Steps	Detailed instructions	Comments
1	<p>Identify the first major peak in the Chromatogram.</p> <p>a In the Chromatogram Results pane, right-click and drag around the peak at 4.676 RT.</p> <p>b With Range Select  selected, click and drag the mouse to select background from RT 4.7 to 4.71.</p> <p>c Right-click inside the shaded area and select Extract MS Spectrum to Background from the context menu.</p> <p>d Click and drag the mouse to select an area containing the peak apex from RT 4.675 to 4.68.</p> <p>e Right-click inside the shaded area and select Extract MS Spectrum from the context menu.</p> <p>f In the MS Spectrum Results tab, with Autoscale Y-axis during Zoom  and Show Predicted Isotope Distribution  selected, from the main menu select Identify > Search Library/DB for Spectra.</p> <p>g Zoom the <i>m/z</i> axis to display the range of values between 168 through 179.</p> <p>h If desired, repeat this procedure to identify other peaks in this sample.</p> <p>i In the Method menu select Save.</p>	<ul style="list-style-type: none"> The single zoomed peak is displayed. This area selected includes background spectra. The extracted spectrum displays in the MS Spectrum Results pane, and is also selected under Background Spectra in the Data Navigator. It will be automatically subtracted from any extracted MS spectrum that follows. This spectra is extracted and displayed in the MS Spectrum Results top pane and is also selected under Spectra in the Data Navigator. The background spectra was subtracted from it and noted by the Subtract label in the MS Spectrum Results top pane. The Spectrum Identification Results tab for this scan displays Dodecane as the compound found by the Library. The predicted isotopes for Dodecane are displayed surrounded by a red outline indicating an isotope of Dodecane. See Figure 19 on page 29.

Exercise – Analyze data
Task 3. Identify Peaks in Qualitative Navigator






Figure 19 MS Spectrum Results showing predicted isotopes

Task 4. Identify compounds using Qualitative Analysis Workflows

Steps	Detailed instructions	Comments
1	<p>Edit the Qualitative Analysis method settings in the Compound Identification workflow.</p> <p>a Double-click Qualitative Workflows shortcut on your desktop.</p>  <p>b Select the SolA_70eV_2.D data file.</p> <p>c In the Method Editor, click Method Automation/Workflow.</p> <p>d In the Options tab select Target/Suspect Screening for the Workflow.</p> <p>e In the Method Editor, click Target/Suspect Screening > Find by Fragments and click the Target Source tab.</p> <p>f Click <input type="text" value="..."/> and select SolA.cdb as the target source.</p>	<ul style="list-style-type: none"> The Open Data File dialog displays. If you are in the Navigator you can also select Launch > Qualitative Workflows from the main menu. This is the same data file loaded in the Qualitative Navigator in the previous task. The Options tab displays. Find by Fragments is automatically selected for Compound mining and is the only type available here. You can find this file in the location shown in “Reference material” on page 2.
2	<p>Configure the Find by Fragments algorithm.</p> <p>a In the Method Editor click Fragment Options and set its values.</p> <p>b Click the Find by Fragments icon in the Method Editor toolbar.</p>	<ul style="list-style-type: none"> Use the values shown in Figure 20 on page 31. It is important to have these values set correctly. The default values in the other tabs are OK for processing this sample. The results are listed in the Sample Table, Compound List, Sample Chromatogram Results, Compound Chromatogram Results, and Compound MS Spectrum Results panes.
3	<p>Save the method.</p> <p>a In the Method Editor toolbar click Save Method .</p>	

Task 4. Identify compounds using Qualitative Analysis Workflows

The screenshot shows the 'Method Editor: Find by Fragments' dialog box. The 'Find by Fragments' tab is selected in the left-hand navigation pane. The main area is divided into several sections:

- Fragment Annotation**: Includes 'Fragment Peak Filter' and 'Target Source'.
- EIC Integration**: Includes 'Calculate Signal-to-Noise' and 'Fragment Options'.
- EIC Peak Filters**: Includes 'Results' and 'Match Tolerance'.


The 'Fragment Options' section is expanded and contains the following settings:

- Fragment ion source**:
 - Use spectral library only
 - Use average fragment spectrum if spectral library not available
 - Number of most specific ions from spectral library:
 - Number of most specific ions from average fragment spectrum:
- Fragment ion EIC qualification settings**:
 - RT difference +/-: min. of expected RT
 - S/N ratio >=
 - Coelution score >=
- Fragment ion confirmation criteria**:
 - Minimum number of qualified fragments:
 - Minimum percent of qualified fragments:

Figure 20 Setting values in the Fragments Options tab

Task 5. Configure Method Automation Reports

You can print an analysis report interactively or generate it as part of a **Method Automation** workflow as we are doing here. An analysis report can contain the results from extracting and integrating chromatograms, extracting spectra, finding compounds, searching the database for peak spectra or generating formulas from peak spectra.

Steps	Detailed instructions	Comments
1 Edit the Qualitative Analysis method to include the Reports workflow.	a In the Method Editor, click Method Automation and select Reports .	<ul style="list-style-type: none"> • See Figure 21 on page 33.
	b Set the entries for the Destination tab to the values shown in Figure 21 on page 33.	<ul style="list-style-type: none"> • Select Save Report and enter a folder location for the PDF-based report.
	c Set the entries for the Templates tab to the values shown in Figure 21 on page 33.	<ul style="list-style-type: none"> • For this exercise, based on the c workflow selection, we are using the Target screening report template.
	d Set the entries for the Layout tab to the values shown in Figure 21 on page 33.	<ul style="list-style-type: none"> • Select the compound tables, chromatograms, spectrum types, peak tables and library search results spectrum types to include.
	e Set the entries for the Contents tab to the values shown in Figure 21 on page 33.	
2 Save the method.	f In the Method Editor toolbar click Save Method  .	
	g	

Exercise – Analyze data
Task 5. Configure Method Automation Reports

The Destination tab contains the following settings:

- Print report:** Print report. Printer name: <Default>
- Save report:** Save report.
 - Inside data file's reports subdirectory
 - At specified directory: D:\MassHunter\reports
- If report file already exists:**
 - Overwrite existing report
 - Auto-generate new report file name

Destination tab

The Templates tab contains the following settings:

- Use Microsoft Excel® for reporting
- Use PDF-based reporting
- Report template folder:** D:\MassHunter\Report Templates\Qual\B.08.00\en-U
- Report templates:**
 - Target screening report template: TargetCompoundScreeningReport.xslt
 - Compound Discovery report template: CompoundReportWithIdentificationHits.xslt
 - Sample purity report template: FormulaConfirmationReport.xslt
 - Compound report template: CompoundReport.xslt

Templates tab

The Layout tab contains the following settings:

- General:**
 - Page orientation: Portrait, Landscape
 - Page size: A4, Letter
 - Show sample information
- Peak table limits (include all peaks unless limits are specified):**
 - Chromatogram peaks: 10 largest peaks per table
 - Mass spectrum peaks: 10 largest peaks per table
- Custom plot limits (autoscaled unless limits are specified):**
 - Chromatograms: 2,000-10,000 min
 - Compound chromatograms: min
 - MS spectra: 100,0000-1200,0000 m/z
 - Deconvoluted spectra: 15000,00-35000,00 Da
 - UV spectra: 240-330 nm

Layout tab

The Contents tab contains the following settings:


- Compounds:**
 - Show compound table. Sort by: Retention time. Sort order: Increasing.
 - Exclude details for unidentified compounds
- Chromatograms:**
 - Show user chromatogram(s)
 - Show compound chromatogram(s)
 - Overlay compound chromatogram(s)
- Compound spectrum (MS):**
 - Show MS spectrum. Show MS peak table
 - Show predicted isotope match table
 - Show MS spectrum (zoomed in on special peaks). Zoom padding: - 30.0 + 30.0 m/z
 - Overlay predicted isotope distribution
- Library search results:**
 - Show library spectrum. Show difference spectrum

Contents tab

Figure 21 Reports tab settings

Task 6. Generate the Method Automation workflow

This task automates the generation of the Find by Fragments and Reports workflows in a single saved method. Here we will perform the analysis based on the method created in the previous tasks. You can also use a saved data analysis method containing several workflows to generate the analysis automatically at the end of your sample data acquisition run.

Steps	Detailed instructions	Comments
1 Review the Method Automation Workflow for finding compounds.	a In the Method Editor select Method Automation > Workflow .	<ul style="list-style-type: none"> See Figure 22 on page 35. This workflow was previously configured. See “Task 4. Identify compounds using Qualitative Analysis Workflows” on page 30.
2 Review the Method Automation Reports for generating a report of the results.	a In the Method Editor select Method Automation > Reports .	<ul style="list-style-type: none"> This workflow was previously configured. See “Task 5. Configure Method Automation Reports” on page 32.
3 Run the Automation Workflow.	a From the dropdown menu in the toolbar select Run Method Automation 	<ul style="list-style-type: none"> The compounds are found, identified, and a report is generated and saved as a pdf in the location specified. See Figure 29 on page 39.




Figure 22 Setting values in Method Automation Workflow







Figure 23 Setting values in Method Automation Reports

Task 7. Review the results

This task takes a brief look at the results shown in various windows of the Quantitative Analysis Workflow program. The first thing you will notice is that all windows are now populated with various results data.

Steps	Detailed instructions	Comments
1 Determine the number of compounds found and how many of these were identified.	<p>a In the Sample Table section, the Result Summary column shows there were 4 qualified targets found.</p>	<ul style="list-style-type: none"> • (See Figure 24 on page 36). • Scroll to the right to review the data in the other columns.
2 View the properties of the identified compounds in the Compound List .	<p>a In the Compound List toolbar click Hide any current empty columns .</p> <p>b In the Compound List toolbar click Auto Size All Columns .</p> <p>c Scroll to the right and review the results for the four identified compounds.</p>	<ul style="list-style-type: none"> • (See Figure 25 on page 37) Saves time reviewing the parameters • Eliminates adjusting columns.
3 View the Compound Identification Results .	<p>a Click on the Compound Identification Results tab at the bottom of the window.</p> <p>b In the Compound List scroll right to view the Compound Identification section and select Dodecane.</p> <p>c Scroll to the right and review the results for Dodecane.</p> <p>d Select the other three compounds in this list and see the view their plots in the results panes.</p>	<ul style="list-style-type: none"> • The Compound Identification Results pane is displayed. • Dodecane is the compound now displayed in the Compound Identification Results window. (See Figure 26 on page 37).



Result Summary	Saved Results Method	Last Run Method	Workflow	Target Source	Sample Name	FileName	Sample Position	Acquisition
▶ 4 qualified (4 targets)	QTOF-GCMS.m	QTOF-GCMS.m	Target/Suspect Screening	D:\MassHunter\Library\SolA.cdb	Sol_A.D	1	DFN EI 70eV.f	

Figure 24 Sample Table results

Exercise – Analyze data Task 7. Review the results

Compound List: 4 found, 4 shown, filtered on Formula

General													Compound Identification			
Formula	m/z	Mass	RT	Width	Height	Area	Score	Base Peak	Ions	Mining Algorithm	ID Source	File	Cpd	Name	Hits	CAS
C12 H10	154.0773	154.0779	5.445	0.014	6014724	5358926	98.26	154.0773	3	Find By Fragment	FBF-FragConfim	SolA_70eV_2.D	1	Biphenyl	1	92-52-4
C12 H9 Cl	188.0382	188.0388	6.111	0.013	7524320	6279273	95.29	188.0382	3	Find By Fragment	FBF-FragConfim	SolA_70eV_2.D	2	4-Chlorobiphenyl	1	2051-62-
C17 H34 O2	270.2544	270.255	6.96	0.01	355478	250343	96.08	74.0361	4	Find By Fragment	FBF-FragConfim	SolA_70eV_2.D	3	Methyl palmitate	1	112-39-0
C12 H26	170.203	170.2035	4.676	0.01	34553	21368	99.99	57.0698	3	Find By Fragment	FBF-FragConfim	SolA_70eV_2.D	4	Dodecane	1	112-40-3

Figure 25 Compound List results


Compound Identification Results: Cpd 4: Dodecane; C12 H26; 4.676

ID Techniques Applied

- FBF-FragConfim

Coelution Score	FragMassDiff(ppm)	Flags(Fis)	Height	Abundance(Lib)	mz(Lib)	m/z	ObsPkHeight(MS)	Compound Name	RT	RT Diff
99.9	1.2	Qualified	356014.1	19.3	55.0542	55.0543	80604.8	Dodecane	4.676	0
99.69	4	Qualified	116183.1	6.3	98.109	98.1094	23855.2	Dodecane	4.676	0
99.85	0	Qualified	281363.3	15.2	70.0777	70.0777	60018.7	Dodecane	4.676	0
99.94	0.3	Qualified	645282.6	35	85.1012	85.1012	135863.5	Dodecane	4.676	0
99.95	1.1	Qualified	1844515	100	57.0699	57.0698	412740.6	Dodecane	4.676	0
100	0.5	Reference ion	1348610.6	73.1	71.0855	71.0855	292281.5	Dodecane	4.676	0

Figure 26 Compound Identification Results for Dodecane

Steps	Detailed instructions	Comments
4 Review the Sample Chromatogram Results window	a Click Compound Overlaid mode .	<ul style="list-style-type: none"> The Dodecane compound is prominently displayed in the TIC at a RT of 4.676 min. (See Figure 27 on page 38).
5 Review the Compound Chromatogram Results window and display its Overlaid mode.	a On the Compound Chromatogram Results window click on  to change from List mode to Overlaid mode.	<ul style="list-style-type: none"> Overlaid mode shows the EIC as a peak outline and the ECC as a filled peak. In our example they are aligned. (See Figure 28 on page 38).

Exercise – Analyze data

Task 7. Review the results




Figure 27 Sample Chromatogram Results for Dodecane in compound overlaid mode




Figure 28 Compound Chromatogram and MS Spectrum Results for Dodecane

Steps	Detailed instructions	Comments
6 Review the compound identification report.	a In file explorer navigate to the folder containing the compound identification report pdf file and Open it.	• For our demo we specified to store it in D:\MassHunter\reports. (See Figure 29 on page 39)

Target Compound Screening Report

Data File	SciA_70eV_2.D	Sample Name	4stds_70eV_10ng.sp50
Sample Type		Position	1
Instrument Name	7250A Marketing	User Name	
Acq Method	4_std5-70eV_SSI_033017.M	Acquired Time	3/30/2017 1:41:10 PM (UTC-07:00)
IRM Calibration Status	Success	DA Method	QTOF-GCMS.m
Comment		Sample Amount	
Expected Barcode		TuneName	PP4-03-14-2017_ehls.tune.xml
Dual Inj Vol	1	TuneDateStamp	2017-03-30T19:40:19-07:00
TimePath	D:\MassHunter\GCMS\1\7250\	OperatorName	
MSFirmwareVersion	G.7250.02.01E	Acquisition Time (Local)	3/30/2017 4:41:10 PM (UTC-04:00)
RunCompletedFlag	True	QuadrupoleTimeOff	
Acquisition SW Version	MassHunter GC/MS Acquisition 8.07.06.2628 28-Mar-2017 Copyright © 1989-2016 Agilent Technologies, Inc.	Light Driver Version	MSQTOFDriver 7.6.0.0
QuadrupoleTimeOff	G.7250.02.01E		
Light Firmware Version			

Compound Table

Label	Tgt Name	Tgt Score	RT Diff	Mass Error (ppm)	Tgt Formula	Tgt RT	Obs. RT	Ref. Mass	Obs. Mass
Cpd 4: Dodecane; C12 H26; 4.676	Dodecane	99.99	0	0.29	C12 H26	4.676	4.676	170.2035	170.2035
Cpd 1: Biphenyl; C12 H10; 5.445	Biphenyl	98.26	-0.003	-2.03	C12 H10	5.449	5.445	154.0783	154.0779
Cpd 2: 4-Chlorobiphenyl; C12 H9 Cl; 6.111	4-Chlorobiphenyl	95.29	-0.008	-2.46	C12 H9 Cl	6.119	6.111	188.0393	188.0388
Cpd 3: Methyl palmitate; C17 H34 O2; 6.960	Methyl palmitate	96.08	-0.002	-3.22	C17 H34 O2	6.962	6.96	270.2559	270.255

Name	Obs. m/z	Obs. RT	Obs. Mass	Tgt RT	Tgt Formula	Tgt Mass	Tgt Mass Error (ppm)	RT Diff.	Find Cpds Algorithm
Dodecane	170.203	4.676	170.2035	4.676	C12 H26	170.2035	0.29	0	Find By Fragment

Compound Chromatograms





Figure 29 Target Compound Screening Report

Exercise – Analyze data

Task 8. Identify compounds using Unknowns Analysis

Task 8. Identify compounds using Unknowns Analysis

This task uses the Unknown Analysis application supplied with MassHunter Quantitative Analysis to find unknown compounds using SureMass.

Steps	Detailed instructions	Comments
1 Start the Unknown Analysis application.	<p>a Double click the Unknown Analysis application shortcut on your desktop.</p>  <p>b Select File > New Analysis from the main menu.</p> <p>c Enter UnknownsDemo.uaf for the File name.</p>	<ul style="list-style-type: none">• The application opens.• The New Analysis dialog opens.• The application title bar displays this name for the analysis.
2 Add samples to the new analysis.	<p>a Select File > Add samples from the main menu.</p> <p>b Select the Sol_A.D file from the Quant batch folder and click OK.</p>	<ul style="list-style-type: none">• The Add Samples dialog opens.• If the data file does not exist in the Quant batch folder, click Browse to Copy Samples and select the sample from another location. The sample is copied to the Quant batch folder.
3 Convert the sample data to SureMass format.	<p>a Select Tools > Convert Accurate Mass Samples from the main menu.</p> <p>b Browse to the Batch folder that contains the sample loaded above.</p> <p>c Select (highlight) the Sol_A.D file from the batch folder and click OK.</p> <p>d In the Convert section of the dialog select Convert to SureMass format and click Convert.</p>	<ul style="list-style-type: none">• The Convert Samples dialog opens.• See step 2 above.• The data file is converted to the SureMass format. Click Close when it finishes.

Steps	Detailed instructions	Comments
4 Edit the analysis method.	<p>a Select Method > Edit from the main menu.</p> <p>b On the Peak Detection tab select SureMass under Peak Detection.</p> <p>c Click on the Deconvolution tab and under Resolution change the RT windows size factor to 100.</p> <p>d Under Component shape deselect Use base peak shape.</p> <p>e Click on the Library Search tab and click on Change Library.</p> <p>f Select the SolA.cdb file.</p> <p>g Under Search criteria select None for the Pre-search type and select Remove Duplicate Hits.</p> <p>h Under Match Factor select Use RT Match.</p> <p>i Under RT penalty function select Trapezoidal, set the RT range to 20 sec and the Penalty-free RT range to 15 sec.</p> <p>j Click on the Compound Identification tab and set the Min MZ to 60.</p> <p>k Click Apply to All Samples.</p> <p>l Click Advanced and then click on the Library Search tab.</p> <p>m Set the Accurate Mass Tolerance to 50 and click OK.</p>	<ul style="list-style-type: none"> • The Method editor dialog opens. • See Figure 30 on page 42. • See Figure 31 on page 42. • See Figure 32 on page 43. • The Open dialog opens. • Browse to the location where the SolA.cdb file is located. • See Figure 33 on page 44. • Scroll to the last column in the table. The Method Editor closes. See Figure 34 on page 44.

Exercise – Analyze data

Task 8. Identify compounds using Unknowns Analysis

The screenshot shows the 'Method' dialog box with the 'Peak Detector' tab selected. The 'Peak detection' section has a dropdown menu set to 'SureMass' and a button labeled 'Specify scan/signal for TIC Analysis...'. The 'Peak filter' section includes an 'Excluded m/z' field with the value '28' and an example 'Example: 28,91,149'. The 'SNR threshold' field is set to '0'.

Figure 30 Method Editor Peak Detector tab parameter changes

The screenshot shows the 'Method' dialog box with the 'Deconvolution' tab selected. The 'Resolution' section has an 'RT window size factor' field set to '100'. The 'Extraction window' section includes 'Left m/z delta' (0.3), 'Right m/z delta' (0.7), and 'm/z delta units' (AMU). There is an unchecked checkbox for 'Use integer m/z values'. The 'Component shape' section has an unchecked checkbox for 'Use base peak shape' and a 'Sharpness threshold' field set to '25 %'.

Figure 31 Method Editor Deconvolution tab parameter changes




Figure 32 Method Editor Library Search tab parameter changes

Exercise – Analyze data

Task 8. Identify compounds using Unknowns Analysis




Figure 33 Method Editor Compound Identification tab parameter changes




Figure 34 Method Editor Advanced Library Search parameter changes

Steps	Detailed instructions	Comments
5 Run the analysis.	<p>a In the Unknowns Analysis toolbar click the Analyze All icon.</p> <p>b Click on Non-Target in the toolbar to display the identified compounds.</p> <p>c Right-click inside the Spectrum window and select Header to Tail if not already selected.</p>	<ul style="list-style-type: none">• The sample file is analyzed according to the parameters set in the method.• Results are displayed in Figure 35 on page 45.• This shows the Library spectra compared against the sample data.

Exercise – Analyze data
Task 8. Identify compounds using Unknowns Analysis




Figure 35 Unknown Analysis results

www.agilent.com

© Agilent Technologies, Inc. 2017

First Edition, July 2017



G7250-90005



Agilent Technologies