



Automatic Determination of Greenhouse Gases by GC

Massimo Santoro, *Thermo Fisher Scientific*

Cristiane de Oliveira Silva, Henrique Franciscato Melo

Danilo Vinicius Pierone, *NovaAnalitica, Brazil*

Introduction

- Greenhouse gases (GHG) are defined as trace components of the atmosphere that absorb infrared radiation emitted by the Earth's surface.
- Increases in GHG have the effect of reducing atmospheric heat loss into space and keeping the Earth warmed.^{1,2}

Introduction

- In 1997, the Kyoto Protocol established commitments for reducing or limiting GHG emissions to be met by industrialized countries between 2008 and 2012. In 2011 at Durban, South Africa, the deadline for the implementation of those commitments was extended for another five to eight years.³
- At COP 194, held in Warsaw (PL) in November 2013, it was decided that a new global agreement to reduce emissions will have to be approved by the first quarter of 2015.⁴
- Brazil is included in the group of developing countries and so has no quantified targets to meet, however, it established the National Policy on Climate Change which sets a national voluntary commitment to reduce its GHG emissions.⁵

Introduction


- The Greenhouse Gas compounds usually considered in the emission estimations are:
 - Carbon Dioxide (CO₂)
 - Nitrous Oxide (N₂O)
 - Methane (CH₄)
 - Hydrofluorocarbons (HFCs)
 - Perfluorocarbons (PFCs)
 - Sulfur Hexafluoride (SF₆)
- Changes in the atmospheric GHG concentration are usually determined by gas chromatography and used for calculating the rates of emission or absorption.

Greenhouse Gas Analysis Instrumentation


- The analysis is performed using the Thermo Scientific™ TriPlus™ RSH autosampler with the Thermo Scientific™ TRACE™ 1310 Gas Chromatograph controlled by the Thermo Scientific™ Chromeleon™ 7.2 Chromatography Data System software.



Two-Detector Configuration, Load Position



Two-Detector Configuration, Inject Position




Chromatographic Parameters of the GHG Analysis


- Loop purge: 20 mL/min of He for 12 s
- GC oven temperature: 50 °C
- Auxiliary oven temperature : 50 °C isothermal
- Carrier gas: He, flow rate: 18 mL/min

• Detector	Temperature (°C)	Flow rate (mL/min)
• FID	250	350 (synthetic air) 25 (H ₂) 12 (N ₂ make-up)
• μECD	350	30 (N ₂ make-up)*
• TCD	110	1.0 (N ₂ make-up)


Three-Detector Configuration, for Oxygen or High CO2 Concentrations



TRACE 1300 Series GC: Tailor Instrument Configuration



- Proprietary, patent-pending Thermo Scientific Instant Connect modules
- Modules are user-installable in less than two minutes
 - Just remove three screws and put the new module in place
 - No special training, dedicated tools or on-site service engineers required
- Every injector and detector is self-sufficient
 - Contains the Integrated Electronic Control (IEC) modules
 - Storing module calibration





Standards and Sample

- Four standard mixtures of CO₂, CH₄, N₂O and SF₆ in different concentrations in helium were analyzed to plot the calibration curves.
- A sample of rumen gases from cattle was analyzed by external standard calibration.

GHG Concentrations in the Standard Mixtures

Mixture	1	2	3	4
CO ₂ (ppm)	252.7	502.2	1027	1998
CH ₄ (ppm)	0.514	1.027	3.150	5.117
N ₂ O (ppb)	253.7	506.2	1000	2096
SF ₆ (ppt)	34	100	1009	n.c.

Typical Chromatograms Obtained for One Mixture



Repeatability

- The repeatability was evaluated through the relative standard deviation (RSD) of the peak area average. Mixtures of GHG were analyzed in triplicate and the peak areas were used for the calculation of the RSD. The low RSD values obtained indicate excellent repeatability.

GHG Concentrations in the Standard Mixtures

GHG	Concentration	Area average (n=3)	Standard deviation	RSD %
CO ₂	252.7 ppm	22.6962	0.1738	0.8
CH ₄	1.027 ppm	0.092	0.0008	0.8
N ₂ O	506.2 ppb	0.019	0.00016	0.8
SF ₆	34 ppt	0.014	0.0002	1.4


Linearity

- Linearity was evaluated by the correlation coefficients (r) of the calibration curves. These calibration curves were obtained using the normal method of quadratic least squares fit. The (r) values greater than 0.99 indicate a good linear correlation achieved between the peak areas and the GHG concentrations determined with FID and ECD detectors.

CH₄ $r = 0.9990$

Conc (ppm) Area average (n=3)

0.514	0.070
1.027	0.123
3.150	0.332
5.117	0.525



Linearity

CO₂ r = 0.9996

Conc (ppm)

Area average (n=3)

252.7

22.696

502.2


54.453

1027

103.725

1998

197.388



N₂O r = 0.9979

Conc (ppb)

Area average (n=3)

253.7

0.004

506.2


0.007

1000.


0.013

2096

0.025



Linearity



SF₆ r = 0.9997

Conc (ppt)

Area average (n=3)

34

0.014

100

0.022

1009

0.206



Limits of Detection and Quantification

- The limits of detection (LOD) of the developed analytical method were determined in successive chromatograms of GHG mixtures with decreasing concentrations. The lowest concentrations that generate analytical signals were considered as LOD. The limits of quantification (LOQ) were calculated based on the 10:1 ratio, i.e. 10LOQ: 1LOD.

Limits of Detection and Quantification



GHG	LOD	LOQ
CO ₂ (FID)	7.432 ppm	74.32 ppm
CH ₄	0.056 ppm	0.56 ppm
N ₂ O	32.76 ppb	327.6 ppb
SF ₆	4.35 ppt	43.5 ppt

Low Concentration Chromatograms



32.76 ppb of N₂O detected by ECD; (b) 0.056 ppm of CH₄ and 7.432 ppm of CO₂ detected by FID/Methanizer

Low Concentration Chromatograms



Low concentrations of SF₆: (a) 100 ppt (b) 4.35 ppt

Rumen Sample from Cattle

- Two chromatograms of a ruminal gas sample and the concentrations of N_2O , CH_4 and CO_2 quantified by the external standard method



ECD Concentration

N_2O 337.4 ppb

FID

CH_4 2.82 ppm

CO_2 981.1 ppm



Conclusions

- The results obtained for area repeatability, linearity, efficiency of the analyte separation, limits of detection, and quantification show that the TRACE 1300 Series GC, in the configuration presented in this study, is a system perfectly suited for the analysis of greenhouse gases.
- The approach is very simple and easily automated and applicable to low- and high-level calibrations.
- Samples are completed in less than 10 minutes, giving high productivity.

References

- 1. *Climate Change 2007: Working Group I: The Physical Science Basis*; IPCC, Intergovernmental Panel on Climate Change, Fourth Assessment Report 2007, http://www.ipcc.ch/publications_and_data/ar4/wg1/en/faq-1-3.html (Accessed March 10, 2014).
- 2. *Climate Change 2013: Working Group I: The Physical Science Basis*; Summary for Policymakers; IPCC, Intergovernmental Panel on Climate Change, Summary for Policymakers Report 2013, http://www.ipcc.ch/pdf/assessment-report/ar5/wg1/WG1AR5_SPM_FINAL.pdf (Accessed March 20, 2014).
- 3. The United Nations. <http://www.un.org/wcm/content/site/climatechange/pages/gateway/the-negotiations/durban> (Accessed March 10, 2014).
- 4. United Nations Framework Convention on Climate Change. Warsaw Climate Change Conference, November 2013. https://unfccc.int/meetings/warsaw_nov_2013/meeting/7649.php (Accessed March 10, 2014).
- 5. *Annual emission estimates of greenhouse gases in Brazil*. Ministry of Science, Technology and Innovation. General Coordination on Global Climate Change. <http://www.mct.gov.br/> (Accessed March 10, 2014).
- 6. Interdisciplinary Research Group on Global Climate Change. <http://www.cpa.unicamp.br/alcscons/abc.php?idAbc=18> (Accessed December 13, 2013).
- 7. Nicoloso, R. S.; Bayer, C.; Denega, G. L.; de Oliveira, P. A. V.; Higarashi, M. M.; Correa, J. C.; Lopes, L. S. Gas chromatography and photoacoustic spectroscopy for the assessment of soil greenhouse gases emissions. *Ciência Rural*, **2013**, 43 (2), 262-269.

Acknowledgements

- The authors thank Prof. Dr. Paulo Henrique Mazza Rodrigues from the Department of Nutrition and Animal Production, University of São Paulo, Pirassununga, SP, BR, for contributing to the improvement of the chromatographic system configuration, the Embrapa São Carlos, SP, BR, by giving the mixtures of GHG standards and Silvana Odete Pisani, Ph.D., from Nova Analítica, São Paulo, SP, BR, by the revision of this writing.

HayeSep is a trademark of Hayes Separations. All other trademarks are the property of Thermo Fisher Scientific and its subsidiaries.

Thank You for Your Attention!



Questions?

Stay connected with us



Twitter

[@ChromSolutions](https://twitter.com/ChromSolutions)



Chromatography Solutions Blog

<http://chromblog.thermoscientific.com/blog>



YouTube

<http://www.youtube.com/ChromSolutions>



Facebook

<http://www.facebook.com/ChromatographySolutions>



Pinterest

<http://pinterest.com/chromsolutions/>