

# A Holiday Gift Buying Guide for your GC: Latest Supply Innovations

Alexander Ucci  
Online Application Engineer  
November 12, 2019



# Exclusive Offer for Attending Agilent Chemistries Webinars

Receive a 25% discount off any Agilent J&W GC columns, Agilent LC columns, sample preparation products, chemical standards, and your most often used GC and LC supplies.\*


- Discount applies on future purchase of up to \$3,000 (USD)
- Offer valid for up to 30 days after the Agilent webinar
- Use promotion code 9969

*\*Some restrictions apply on supplies*




# An Overview of the Latest Supply Innovations

Gas traps  
Crimpers  
Sample cleanup  
ADM Flowmeter  
Standards



# Let's Talk About Gas Quality and Filters

- Oxygen in carrier gas is detrimental to GC/MS
  - Reduced response
  - Elevated background
  - Irreversible column damage
  - Impaired electron multiplier function
  - Premature filament, liner lifetime
- Use UHP carrier gases
  - 99.9995% or greater
- Use Gas Clean carrier gas filters




GC/MS filter  
Agilent p/n  
CP17973



# Let's Talk About Gas Quality and Filters

If you used lower quality gases, how much O<sub>2</sub> could the filter clean up?




Installing (and properly purging) the Gas Clean carrier gas filter lowered the O<sub>2</sub> signal by a **factor of 5**.



# Let's Talk About Gas Quality and Filters

If you used lower quality gases, how much H<sub>2</sub>O could the filter clean up?




Gas Clean filter lowered the H<sub>2</sub>O signal by a **factor >10**

You can further increase Gas Clean and column lifetime by installing the universal filter before the Gas Clean filter.



# Let's Talk About Gas Quality and Filters

If lower quality gases were used, how much background could the filter clean up?



- Install the Universal trap vertically – use the mounting bracket(s)
- Extend the lifetime of your Gas Clean (indicating) filter **and** your column

# New Agilent A-Line Electron Crimpers and Decappers




- More vials crimped per battery charge
- Increased crimping speed when compared to previous models
- Lighter weight means less hand strain and effort
- OLED screen for easier viewing interface
- More efficient charging
- New motor life is significantly longer leading to extended productivity

<https://www.agilent.com/en/products/lab-supplies/chromatography-spectroscopy/vials-closures/crimpers-decappers-accessories>



# Vials

- Choose high-quality vials and caps
  - Poorly constructed vial septa → siloxanes → bleed peaks
  - Low-quality vial → leach contaminants into sample
  - Choose the right cap/septa for your solvent



|                             | High performance septa                  | Thin PTFE   | PTFE/Silicone*  | PTFE/Silicone/PTFE*  | PTFE/Red rubber   | Flouoroelastomer                          | Butyl  |
|-----------------------------|---|---|---|--|---|---|--|
| Temperature range           | 40 °C to 300 °C**                       | Up to 260 °C  | -40 °C to 200 °C  | -40 °C to 200 °C   | -40 °C to 90 °C   | -40 °C to 260 °C                          | -50 °C to 150 °C                                     |
| Use for multiple injections | No                                      | No  | Yes   | Yes  | No  | No  | No   |
| Price                       | More expensive                          | Very economical   | Economical  | Most expensive   | Very economical   | Economical                                | Economical   |
| Resistance to coring        | Excellent                               | None  | Excellent   | Excellent  | None  | None                                      | None   |
| Recommended for storage     | No                                      | No  | Yes   | Yes  | No  | No  | No   |
| Best for                    | High temperature headspace applications | Superior chemical inertness, short cycle times, and single injections | Most common HPLC and GC analyses, not as resistant to coring as P/S/P | Superior performance for ultra trace analysis, repeat injections, and internal standards | Chlorosilanes, more economical option for single injections | Chlorinated solvents, higher temperatures | Organic solvents, acetic acids, impermeable to gases |

\* Agilent silicone is platinum cured (versus peroxide cured), making it more inert and less likely to interact with samples.

\*\* For up to 1 hour.

# ADM Flow Meter: Easy No-Return Recalibration





- Replaceable calibration cartridge
- Automatic notification of cartridge replacement
- Ergonomic and robust design
- Universal 3AA or USB power
- USB connects to Web interface for added functionality
- Easy to view OLED Screen
- Kickstand

Technical overview 5991-7685EN

# The Matrix

If your target ions are buried beneath matrix peaks, it might be time to trim the column or do sample cleanup



(or improve your sample cleanup)




# How Did it Become Contaminated in the First Place?





# Offline Options for Sample Matrix Removal



QuEChERS



Captiva EMR-Lipid filtration cartridges and plates

Bond Elut Solid Phase Extraction cartridges and plates



Captiva syringe filters



Filter vials




Synthetic Chem Elut S

# SPME Fiber and Arrow Offering from Agilent

## Solid Phase Microextraction (SPME)

- Environmental analyses of water samples
  - Odor analyses (ppt)
  - Flavor analyses of food products
  - Forensic analyses of arson/explosives samples
  - Toxicology analyses: blood alcohol or drugs in urine/serum
  - Surfactants, other industrial applications
- Trace analysis in food
  - Drugs and pharmaceuticals
  - Herbicides/pesticides
  - Medical diagnostics
  - Trace impurities in polymers and solid samples
  - Solvent residues in raw materials
  - Explosives



SPME fibers




SPME arrows


# Examination of Lower Molecular Weight PAHs in Drinking Water Using Agilent PDMS SPME Fibers

Polycyclic aromatic hydrocarbons (PAHs) are a large class of organic compounds containing two or more fused aromatic rings. PAHs are considered compounds of concern by every environmental organization; their concentration in water is strictly regulated.

5994-1301EN




SIM chromatogram of naphthalenes with PDMS fibers (black trace = 100  $\mu\text{m}$ ; green trace = 30  $\mu\text{m}$ ; blue trace = 7  $\mu\text{m}$ )



SIM chromatogram of fluoranthene and pyrene with PDMS fibers (black trace = 100  $\mu\text{m}$ ; green trace = 30  $\mu\text{m}$ ; blue trace = 7  $\mu\text{m}$ )

# Liquid/Liquid Extraction (LLE)



- LLE has been successfully used as a method of sample preparation for many years.
- It separates the more organic solvent soluble compounds from the more water-soluble compounds using water immiscible organic solvents.
- It can remove many interfering substances like salts
- Modulating pH can selectively extract or eliminate specific compound types.






# Drawbacks of Liquid/Liquid Extraction

- LLE does have drawbacks
  - Inconsistent results from one analyst to another
    - Shaking time
    - Shaking motion
    - Determination of where to cut between layers
  - Emulsions
  - Labor intensive
  - Quite tedious with small sample sizes (<5 mL)
  - Challenging with large numbers of samples
  - Difficult to automate for large numbers of samples



How many of these problems can be fixed with Solid Supported Liquid Extraction?


# How Does SLE Work?



1. Add aqueous sample to SLE sorbent
2. Wait 5 minutes
3. Elute with water immiscible solvent

# What is SLE Sorbent?


- There are two types of SLE media
  - Diatomaceous earth (DE) based products like our Chem Elut brand of SLE products
    - A mined fossil diatom material, which is heterogeneous and inconsistent from one mine to the next



## Diatomaceous Earth in Chem Elut

- ✗ Naturally occurring; mined
- ✗ Broad particle size distribution
- ✗ Supplier reliability issues
- ✗ Poor lot-to-lot consistency

- Synthetic media we use in Chem Elut S
  - Controlled synthesis to be consistent batch after batch



## Synthetic SLE sorbent


- ✓ Large scale synthesis
- ✓ Narrow particle size distribution
- ✓ Reliable supplier
- ✓ Controlled manufacturing

# Chem Elut S – 15 Minute Hold Time


## Large scale format comparison with aromatic amines

Recovery (30 µg/mL; n = 3), 87-119%

Chem Elut S 5 mL Chem Elut S 10 mL Chem Elut S 20 mL



Reproducibility (30 µg/mL; n = 3), RSD 0.4-8.2%







5994-0951EN



# Enhanced Matrix Removal: Agilent Captiva EMR-Lipid

EMR-Lipid sorbent technology effectively traps lipids through two mechanisms:


- Size exclusion – unbranched hydrocarbon chains (lipids) enter the sorbent; bulky analytes do not
- Sorbent chemistry – lipid chains that enter the sorbent are trapped by hydrophobic interactions



# Determination of 19 Polycyclic Aromatic Hydrocarbon Compounds in Salmon and Beef

Using Captiva EMR-Lipid cleanup by GC/MS/MS


Salmon dried residue



No cleanup

Captiva EMR-Lipid  
cleanup

Beef dried residue




No cleanup

Captiva EMR-Lipid  
cleanup

5994-0553EN

# Captiva EMR-Lipid Cleanup Improves Analytes S/N Ratio and Integration Accuracy on GC/MS(/MS)



5994-0405EN

# ULTRA Scientific is Now Part of Agilent Technologies

Agilent ULTRA Chemical Standards have:


- Best in class online search, compare, and ordering capabilities
- Rapid shipping: 99.9% of orders dispatched within 24 to 48 hours
- Custom solutions
- Sample preparation materials, columns, supplies, instrumentation, and reference materials from a single source
- Rigorously tested and manufactured under ISO 9001, ISO 17025, and Guide 34 certifications



You can calibrate with confidence and maximize accuracy

# An Overview of the Latest Supply Innovations


Gas traps  
Crimpers  
Sample cleanup  
ADM Flowmeter  
Standards





# Inlet

- Injection efficiency:
  - Main function of the inlet is to produce a narrow sample band at the head of the column
  - One of the most important aspects to any high-resolution GC method
- Must be reproducible
- The liner volume must be large enough to accommodate the solvent's phase transformation into a vapor (backflash)
- Most chromatography problems are “front-end” related
- Many consumables to replace: septa, liner, gold seal
- Inlet body must be cleaned/solvent rinsed periodically (**no steel brushes, please**)




# 30+ Years of R&D Focused on Surface Deactivation

Two unique chemistries were developed to treat the surfaces in a GC flow path.

## UltiMetal Plus

- Inorganic vapor deposition in Middelburg
- CFT device
- Ultimate union
- Flexible metal ferrule
- Inlet welding
- Detector parts
- GC tubing



## Ultra Inert

- Organic vapor deposition
  - Ultra Inert inlet liners
  - Ultra Inert gold seals
- Unique proprietary deactivation process
  - UI columns




# Septa

- Typical cost of one premium septum (list), \$1.25
- Typical cost of one GC column, 30 m x 0.25 mm id, \$450
- Proactively change inlet septa
- Agilent packing eliminates contamination of septa
- “CenterGuide septa” puts less strain on syringe compared to solid septa
- Do not overtighten septum nut; septum can begin to “bulge” out
- Should tighten nut until c-clamp on top stops turning, then  $\frac{1}{2}$  to  $\frac{3}{4}$  turn more




Septum nut




# Ultra Inert Inlet Liners

1. Ultra Inert deactivated inlet liners provide higher response for sensitive compounds
2. Ultra Inert **glass wool liners** deliver benefits of glass wool without loss of active compounds
3. QC tested and certified for consistent performance



## Productivity

Touchless packaging with pre-installed O-ring: quick and easy hassle free installation



# Agilent Assurance

UI liner lot QC with demanding test probes (dinitrophenol and others)

Assurance of consistent performance

Label for p/n, batch, and lot testing

Traceability

Deactivation lot number is on certificate

- Liner lot number (and part number) etched on glass

## *Certificate of Performance*

**5190-2293 Ultra Inert Liner**

Splitless, Sngl taper, Glass Wool

Liner Body Lot: 0023A

Deactivation Lot: B11002

**Tested for: 2ug 4-Aminopyridine  
2ug 2,4-Dinitrophenol**



# Packaging, Ease, and Productivity



Individual liner



Convenient 5/pk in  
Touchless packaging




25/pk in Touchless  
dispenser  
High throughput



100/pk, bulk, blister  
(no O-ring)  
High throughput

# What's New? Glass Wool Alternative Liners

## Ultra Inert liners with sintered glass frits





# Modern Day Solution: The Glass-Frit



Substituting glass wool for a sintered glass frit provides equivalent vaporization properties while enhancing performance consistency and quality


# The Benefits of the Glass Frit

Dislodged glass wool fibers expose active sites that interact with sensitive analytes



# 8270 Chromatogram Overlay – Phenols/Low Frit Liner

Low frit liner gives excellent response and peak shape for tricky analytes like phenols




Blue = Initial injection  
Red = After 100 matrix injections  
Green = After liner replacement

5 ppm, 1  $\mu$ L splitless injection  
FID detector

# Pesticides by FPD



Results for fritted liner indicate highly inert, consistent response






# Agilent UI Gold Seal: Deactivated Gold Surface

- Soft gold plating is essential for proper sealing
- Ultra Inert chemistry blocks active sites (gold is **not** inert)
- Smooth surface doesn't leak (injection molded)
- Part numbers 5190-6144 ea. 5190-6145, 10/pk.  
5190-6149, 50/pk.





Agilent MIM seal



Competitor's  
machined seal

Reliable ppb and ppt  
measurements require  
attention to the little things

# Column Installation: Self-Tightening Column Nut



- Spring driven piston continuously presses against ferrule
- Automatically retightens when ferrule shrinks
- No leaks, no downtime, no frustration
- Wing design for finger tightening
- No tools needed
- No polymer materials for durability
- Compatible with **only** short graphite
- Vespel ferrules

# Increasing Ease of Use Through Continued Innovation: Self-Tightening Nuts



For GC inlet or detector

Innovating for the next generation



For mass spectrometry transfer line

Innovating for the next generation



- Easier and faster to install
- Collar holds column in place
- Single-hand installation into inlet
- No tools needed




# Self-Tightening Nuts: No leaks, No Downtime, No Frustration

- Spring-driven piston continuously presses against ferrule
- Automatically retightens when ferrule shrinks
- Wing design for finger tightening
- No tools needed


| Part Number | Description                                       |
|-------------|---|
| G3440-81013 | Column Nut, Collared Self-Tightening MSD          |
| G3440-81011 | Column nut, Collared Self Tightening Inlet/Detect |
| G3440-81012 | Collar for Self Tightening Nut                    |




Standard column nuts new fitting



Standard column nuts after 25 injections




Agilent Self Tightening Column Nuts after 400 injections



400 injections

# An Overview of the Latest Supply Innovations

Gas traps  
Crimpers  
Sample cleanup  
ADM Flowmeter  
Standards





# Agilent J&W Column Portfolio- DB, HP, CP, VF

| Low Polarity      |                    |                  | Mid Polarity     |                  |           | High Polarity |                   |              |
|-------------------|--------------------|------------------|------------------|------------------|-----------|---------------|-------------------|--------------|
| CP-Sil 2          | DB & HP-1MS UI     | DB & HP-5MS UI   | DB-XLB           | DB-225MS         | DB-ALC1   | HP-88         | DB-WAX            | DB-WAX UI    |
| DB-MTBE           | DB & HP1-MS        | DB & HP5-MS      | VF-XMS           | DB-225           | DB-Dioxin | CP-Sil 88     | DB-WAX ETR        | DB-HeavyWAX  |
| CP-Select CB MTBE | VF-1MS             | VF-5MS           | DB-35MS UI       | CP-Sil 43 CB     | DB-200    | DB-23         | HP-INNOWax        | DB-FATWAX UI |
|                   | DB & HP-1          | DB & HP-5        | DB & VF-35MS     | VF-1701 MS       | VF-200MS  | VF-23 MS      | VF-WAX MS         |              |
|                   | CP-Sil 5 CB        | CP-Sil 8 CB      | DB & HP-35       | DB-1701          | DB-210    |               | CP-WAX 57 CB      |              |
|                   | Ultra 1            | Ultra 2          | DB & VF-17MS     | CP-Sil 19 CB     | DX-4      |               | DB & HP-FFAP      |              |
|                   | DB-1HT             | VF-DA            | DB-17            | HP-Blood Alcohol |           |               | DB-WAX FF         |              |
|                   | DB-2887            | DB-5.625         | HP-50+           | DB-ALC2          |           |               | CP-FFAP CB        |              |
|                   | DB-Petro/PONA      | DB & VF-5HT      | DB-17HT          | DX-1             |           |               | CP-WAX 58 FFAP CB |              |
|                   | CP-Sil PONA CB     | CP-Sil PAH CB    | DB-608           |                  |           |               | CP-WAX 52 CB      |              |
|                   | DB-HT SimDis       | Select Biodiesel | DB-TPH           |                  |           |               | CP-WAX 51         |              |
|                   | CP-SimDis          | SE-54            | DB-502.2         |                  |           |               | CP-Carbowax 400   |              |
|                   | CP-Volamine        |                  | HP-VOC           |                  |           |               | Carbowax 20M      |              |
|                   | Select Mineral Oil |                  | DB-VRX           |                  |           |               | HP-20M            |              |
|                   | HP-101             |                  | DB-624           |                  |           |               | CAM               |              |
|                   | SE-30              |                  | VF-624MS         |                  |           |               | CP-TCEP           |              |
|                   |                    |                  | CP-Select 624 CB |                  |           |               |                   |              |
|                   |                    |                  | DB-1301          |                  |           |               |                   |              |
|                   |                    |                  | VF-1301MS        |                  |           |               |                   |              |
|                   |                    |                  | CP-Sil 13 CB     |                  |           |               |                   |              |

**Agilent J&W columns have over 50 different stationary phase offerings**

# Introducing DB-HeavyWAX

- WAX column with increased MAOT compared to existing columns on the market
  - 280 °C isothermal and 290 °C programmed
- Provides increased thermal stability
- Has a low bleed level

## Advantages

### General GC

- Shorter runtimes when late eluters are present
- Better S/N ratio, improved detection
- Better thermal stability
- Faster column bake-out


### GC/MS

- Desire for “zero” bleed
- Avoid MS contamination by column bleed for longer system uptime and column lifetime
- Improve detection limit


### GC×GC

- Extended scope of compounds

# Bleed Summary at 280 °C Over 100 Hours




# Thermal Stability and Retention Time Shifting for a Traditional WAX Phase



# DB-HeavyWAX

Application note 5991-9035EN


| Peak | Compound     |
|------|--------------|
| 1    | Methanol     |
| 2    | Benzene      |
| 3    | Toluene      |
| 4    | Ethylbenzene |
| 5    | p-Xylene     |
| 6    | m-Xylene     |
| 7    | o-Xylene     |







# Higher Final Temperatures Mean Earlier Elution of Heavy Compounds, Sharper Late Eluting Peaks, and Reduced Carryover

Pink grapefruit essential oil (cold pressed)




Application note 5991-9078EN

# Better Peak Shape and Decreased Runtime



# Bleed Reductions for FAMEs in Biodiesel (IP-585) Using a DB-HeavyWAX, 60 m x 0.25 mm x 0.50 μm




# DB-FATWAX UI Provides the Desired Inertness and Thermal Stability to Separate Challenging Organic Acids and Fatty Acids

- Great need in the food, forensic, and cosmetic industries to monitor the content of free and natural occurring fatty acids.
- Analysis of underivatized organic acids and free acids is desirable to eliminate the problems associated with derivatization, extraction, and cleanup procedures
- Volatile organic acids and fatty acids are difficult to quantify accurately by standard WAX columns. These acids often elute as tailing or poorly resolved peaks. For some acids, adsorption can become irreversible.
- Normally:
  - Fatty acids are derivatized to the methyl ester (FAME)
  - Free fatty acids are analyzed using acid-deactivated WAX columns. The acid modifier, nitroterephthalic acid, however, reduces thermal stability, operating temperature and reacts with humidity, reducing column life.

# Comparison Between Acid-Modified WAX and DB-FATWAX UI After 50 h at 250 °C

DB-FATWAX Ultra Inert shows superior inertness and thermal stability compared to acid-modified WAX




# DB-FastFAME


## 20 m x 0.18 mm x 0.20 $\mu\text{m}$

|                  |   |
|------------------|---|
| <b>Column</b>    | Agilent J&W DB-FastFAME, 20 m x 0.18 mm, 0.20 $\mu\text{m}$                             |
| <b>Gas</b>       | <b>Hydrogen</b> , 28 psi, constant pressure mode  |
| <b>Inlet</b>     | Split/splitless, 250 °C, split ratio 50:1   |
| <b>Oven</b>      | 80 °C (0.5 min), 65 °C/min to 175 °C, 10 °C/min to 185 °C (0.5 min), 7 °C/min to 230 °C |
| <b>FID</b>       | 280 °C, hydrogen: 40 mL/min; air: 400 mL/min; make-up gas: 25 mL/min.                   |
| <b>Injection</b> | 1 $\mu\text{L}$   |

Strong interaction between cis isomers and the dipoles of the cyano propyl ligands. That allows the trans to elute after the cis isomers.



$R_s \geq 1.95$  for cis/trans isomers



7 min analysis time



# Introducing the Agilent 8890 GC System

Flexible and expandable to meet your needs today and tomorrow



## Future-proof: Ready for anything

- Powerful next generation electronic architecture
- Expanded smart-connected functionality
- Full suite of inlets, detectors, and accessories, CFT, Deans switch, backflush, GC x GC, dual simultaneous injection
- Six valves, eight heated zones, plus LVO
- Generation 6 precision EPC
- Smart keys
- 7-inch color touch display



# Agilent 8890 GC System

## Smart-connected GC

### Modern intuitive interface

#### 7-inch color touch screen

- Configuration
- Status
- Methods
- Sequence info
- Troubleshooting, diagnostics, and help

#### Real-time chromatographic evaluation

- Blank evaluation
- Detector evaluation

# GC Columns with Smart Key (for the Agilent 8890 GC only)

For immediate identification and usage monitoring of your GC column

- Available with the Agilent 8890 GC only
- Can track use of a GC column
- Smart key contains GC column information, including:
  - Part and serial numbers
  - Number of injections/runs
  - Time at/above temperature limits
  - Date installed
  - Temperature limits – GC columns
    - If more than one column is installed, the temperature is determined by lowest column smart key installed (DB-WAX vs DB-5)
  - Column length/trimming done in “column maintenance mode” in software and rewritten to the smart key
  - S/N of last instrument installed if it was in an 8890 GC




# GC Columns with Smart Key

## Feature, advantages, and benefits

| Feature   | Advantage  | Benefit  | Economic Benefit   |
|---|--|--|--|
| <p>Smart key included with the GC column p/n must include "KEY"</p> <p>Designation example *122-5532UI-KEY*</p> | <ol style="list-style-type: none"> <li>Reduce possible errors of manual input of method parameters</li> <li>Optimized maintenance schedules with usage tracking</li> </ol> | <ol style="list-style-type: none"> <li>Better data quality as the system is aware of the configuration<br/>Ease-of use</li> <li>Ability to better plan preventive maintenance before issues occur</li> </ol> | <ol style="list-style-type: none"> <li>Ease-of-use, no risk of faulty information in GC, fewer reruns</li> <li>Limited downtime as preventive maintenance is easier to plan</li> </ol> |

|   | Literature   | Web page  |
|---|--|---|
| 1 | Agilent 8890 GC brochure                             | <a href="https://www.agilent.com/cs/library/brochures/brochure-gc-8890-5994-0476en-agilent.pdf">https://www.agilent.com/cs/library/brochures/brochure-gc-8890-5994-0476en-agilent.pdf</a>   |
| 2 | Smart key product page (not for ordering smart keys) | <a href="http://www.Agilent.com/chem/smartkey8890">www.Agilent.com/chem/smartkey8890</a>  |
| 3 | Instruction sheet                                    | <a href="https://www.agilent.com/cs/library/instructionsheet/public/insert-smart%20key-8890-5994-0700en-agilent.pdf">https://www.agilent.com/cs/library/instructionsheet/public/insert-smart%20key-8890-5994-0700en-agilent.pdf</a> |






# GC Columns with Smart Key

## Smart key clarifications

- It is not a USB flash drive. It comes already programmed from the factory with specific information about the GC column that it belongs to
- It is not tethered to the GC column, neither is connected through it with WIFI or Bluetooth.



# Expanding the Intuvo Column Family

## Key features

- VF-, CP-, and Select- applications specific: 24 new phases and 39 new part numbers
- New Intuvo column dimension: 0.15 mm id to add to our 0.18 mm id for a greater range of high-efficiency columns
- We now have PLOT columns
- Custom column options



| PDMS        | 5% Phenyl                 | WAX                    | Specialty             | Siloxane/<br>Cyanopropyl | PLOT               |
|-------------|---------------------------|------------------------|-----------------------|--------------------------|--------------------|
| VF-1ms      | VF-5ms                    | CP-WAX 52<br>CB        | CP-Volamine           | VF-17ms                  | HP-PLOT Q          |
| CP-Sil 5 CB | CP-Sil 8 CB               | VF-WAXms               | Select Mineral<br>Oil | VF-1701ms                | HP-PLOT<br>Al2O3 S |
|             | CP-Sil 8 CB for<br>Amines | CP-WAX 57<br>CB        | Select PAH            | VF-Xms                   |                    |
|             |                           | CP-WAX 58<br>FFAP      | Select Low<br>Sulfur  | VF-624ms                 |                    |
|             |                           | CP-Wax for<br>Volamine | Select-FAME           | CP-Select 624<br>CB      |                    |
|             |                           |                        | CP-Sil 88 for<br>FAME | CP-Sil 88                |                    |
|             |                           |                        | CP-TAP                | VF-200ms                 |                    |






# An Overview of the Latest Supply Innovations

Gas traps  
Crimpers  
Sample cleanup  
ADM Flowmeter  
Standards

Septa  
Liners  
Gold seals  
Column  
nuts

A-line vials



# New Agilent Universal Fit GC Detector Jets

- Easier column installation and jet replacement reducing the risk of column damage
- Lubricant free threads reducing the risk of contamination
- Made of strong material reducing the risk of deforming
- Universally fits in both capillary column and packed column (adaptable) FID detectors



| Previous Jets   |                          |                      |                              | New Universal Fit Jets        |                          |                      |                              |
|-----------------|--------------------------|----------------------|------------------------------|-------------------------------|--------------------------|----------------------|------------------------------|
| Previous Jet PN | Jet Orifice ID (inch/mm) | Jet Length (inch/mm) | Fit of Detector Fitting Type | New Jet PN (use for re-order) | Jet Orifice ID (inch/mm) | Jet Length (inch/mm) | Fit of Detector Fitting Type |
| 19244-80560     | 0.011 / 0.29             | 2.4 / 62             | FID, Adaptable               | 5200-0176                     | 0.011 / 0.29             | 1.2 / 31             | FID, Capillary & Adaptable   |
| G1531-80560     | 0.011 / 0.29             | 1.7 / 43             | FID, Capillary               |                               |                          |                      |                              |
| 18710-20119     | 0.018 / 0.47             | 2.5 / 64             | FID, Adaptable               | 5200-0177                     | 0.018 / 0.47             | 1.2 / 31             | FID, Capillary & Adaptable   |
| 19244-80620     | 0.018 / 0.47             | 2.4 / 62             | FID, Adaptable               |                               |                          |                      |                              |
| G1531-80620     | 0.018 / 0.47             | 1.7 / 43             | FID, Capillary               |                               |                          |                      |                              |
| 18789-80070     | 0.030 / 0.76             | 2.5 / 64             | FID, Adaptable               | 5200-0178                     | 0.030 / 0.76             | 1.2 / 31             | FID, Capillary & Adaptable   |
| G1534-80580     | 0.011 / 0.29             | 2.0 / 52             | NPD, Capillary               | 5200-0179                     | 0.011 / 0.29             | 1.6 / 40             | NPD, Capillary & Adaptable   |
| G1534-80590     | 0.011 / 0.29             | 2.8 / 71             | NPD, Adaptable               |                               |                          |                      |                              |

# Upcoming GC Seminars

|  |                   |                      |
|--|-------------------|----------------------|
| <a href="#">How to combat the Helium Shortage: Making the Switch from Helium to Hydrogen or Nitrogen Carrier Gas</a> | December 5, 2019  | 11:00am ET/8:00am PT |
| <a href="#">Optimizing Conditions for GC/MS Analyses</a>   | January 23, 2020  | 11:00am ET/8:00am PT |
| <a href="#">Plunge Deep into the World of Syringes: Simplifying GC Syringes and Sample Introduction</a>              | February 27, 2020 | 11:00am ET/8:00am PT |
| <a href="#">An Agilent Guide to SPME Analysis</a>  | March 26, 2020    | 11:00am ET/8:00am PT |

Have an idea for a topic? Let us know.



# Contact Agilent Chemistries and Supplies Technical Support



1-800-227-9770 option 3, option 3:

**Option 1 for GC and GC/MS columns and supplies**

Option 2 for LC and LC/MS columns and supplies

Option 3 for sample preparation, filtration, and QuEChERS

Option 4 for spectroscopy supplies

Option 5 for chemical standards

**Available in the USA and Canada 8–5, all time zones**



[gc-column-support@agilent.com](mailto:gc-column-support@agilent.com)

[lc-column-support@agilent.com](mailto:lc-column-support@agilent.com)

[spp-support@agilent.com](mailto:spp-support@agilent.com)

[spectro-supplies-support@agilent.com](mailto:spectro-supplies-support@agilent.com)

[chem-standards-support@agilent.com](mailto:chem-standards-support@agilent.com)