

APPLICATIONS

ASTM D6730 Detailed Hydrocarbon Analysis of Spark Ignition Fuel using Zebtron™ ZB-DHA-PONA Gas Chromatography Column

Dr. Ramkumar Dhandapani, Zeshan Aqeel, and Tim Nelson
Phenomenex, Inc., 411 Madrid Ave., Torrance, CA 90501 USA



Tim Nelson
Global Marketing Manager
Fuels and Chemicals

Tim is an avid outdoorsman who loves to hike and ski. His most recent exploration is tall ship sailing in our local Pacific Ocean. Tim loves history and everything about the stars and space.

Introduction

The motor fuel industry continues to make improvements in their formulas to lower emission and to provide higher operational efficiency. In addition, for those same reasons there are new combustion engines that are being designed around alternative formulated fuels. Fuels are composed of very complex matrices that have several hundred different components and there are many different refinery process stream feeds that can be used to blend with the straight run material. The hundreds of different components in the straight run spark ignition fuel are then further complicated exponentially with the numerous different sources of refinery process feed and their hundreds of components, and so it is mandated that the fuel industry use regulated test methods to measure the group type chemical components to make sure that formula guidelines are met. Not only is the component composition needed for the regulatory specifications, but also for the refinery process quality control.


The petroleum fractions in the spark ignition fuel are in the carbon number range of C1 to C12, and the refinery process feedstock are virgin straight run naphtha, along with other blending streams from the FCC, coker, steam cracking, alkylation, and synthetic sources, along with the chemical additives such as oxygenates. Detailed Hydrocarbon Analysis (DHA) are test methods developed specifically to measure the primary chemical group types, such as paraffin, isoparaffin, olefin, naphthenes, and aromatics.

ASTM-D6730-01 is a specific DHA test method for the analysis of hydrocarbon components along with the oxygenated additives that are in spark-ignition engine fuels. This test method covers the determination of individual hydrocarbon components of spark-ignition engine fuels with boiling ranges up to 225 °C. The accuracy for total PONA analysis decreases with higher olefin content and carbon numbers, though the aromatics have reasonable accuracy. A precolumn is used to improve resolution of some key components such as benzene and toluene from paraffins and cycloparaffins.

Table 1.
Zebtron ZB-DHA-PONA Column Dimensional Options for Detailed Hydrocarbon Analysis (ASTM and alternative methods)

Part No.	Dimensions	Stationary Phase	Benefits
7JE-G042-17	50 m x 0.20 mm x 0.50 μm	100 % Dimethylpolysiloxane	50 m column provides shorter run time while the 0.20 mm tighter ID provides higher efficiency
7MG-G042-17	100 m x 0.25 mm x 0.50 μm	100 % Dimethylpolysiloxane	100 m length provides high plate count / efficiency
7QG-G042-22	150 m x 0.25 mm x 1.0 μm	100 % Dimethylpolysiloxane	150 m with a 1.0 μm thicker film provides better separation of lower boiling fractions and maintains high efficiency from the column length
7AG-G042-22	5 m x 0.25 mm x 1.0 μm	5 % Phenyl 95 % Dimethylpolysiloxane	Optional tuning column provides phenyl selectivity in addition to true boiling point separation. This helps resolve certain aromatics from alkanes and alkenes

Figure 1.
Analysis of ASTM D6730 Components by GC-FID on ZB-DHA-PONA & ZB-DHA-PONA-TUNE GC Column




Sample #	Analyte	Time min.
1.	Ethanol	4.039
2.	C5 (n-pentane)	4.696
3.	2-Methylbutane	5.064
4.	Tert-Butanol	6.174
5.	2,3-Dimethylbutane	6.276
6.	Methyl tert-butyl ether (MTBE)	7.446
7.	C6 (n-hexane)	9.708
8.	1-Methylcyclopentene	9.823
9.	Benzene	10.106
10.	Cyclohexane	11.493
11.	3-Ethylpentane	11.568
12.	trans-1,2-Dimethylcyclopentane	12.278
13.	C7 (n-heptane)	16.167
14.	2,3,3-Trimethylpentane	16.355
15.	Toluene	21.131
16.	C8 (n-octane)	28.485
17.	Ethylbenzene	30.269
18.	2,3-Dimethylheptane	39.526
19.	p-Xylene	48.182
20.	C9 (n-nonane)	51.611
21.	5-Methylnonane	54.958
22.	1-Methyl-2-ethylbenzene	55.011
23.	C10 (n-decane)	59.164
24.	C11 (undecane)	59.628
25.	1,2,3,5-Tetramethylbenzene	63.09
26.	Naphthalene	63.568
27.	C12 (dodecane)	64.984
28.	1-Methylnaphthalene	70.146
29.	C13 (Tridecane)	70.331

GC Conditions for Analysis:

- Column 1 (Tuning):** Zebtron™ ZB-DHA-PONA-TUNE
- Phase:** 5% Phenyl 95% Dimethylpolysiloxane
- Dimensions:** 5 meter x 0.25 mm x 1.0 µm
- Part No.:** [7AG-G042-22](#)
- Column 2:** Zebtron™ ZB-DHA-PONA
- Phase:** 100% Dimethylpolysiloxane
- Dimensions:** 100 meter x 0.25 mm x 0.50 µm
- Part No.:** [7MG-G042-17](#)
- Recommended Column Union:** [AG0-4716](#)
- Injection:** Split 150:1 @ 200 °C, 0.2 µL
- Recommended Liner:** ZebtronPLUS Straight Z-Liner™
- Part No.:** [AG2-0A03-05](#)
- Carrier Gas:** Hydrogen @ 2 mL/min (Constant Flow Mode)
- Oven Program:** 30 °C for 8.5 min, to 48 °C at 22 °C/min, for 27 min, to 141 °C at 3 °C/min, for 1 min, to 275 °C at 1 °C/min, for 2 min
- Detection:** Flame Ionization (FID) @ 275 °C

Figure 2.
Benefits of ZB-DHA-PONA GC Column



Flexible Polyimide

Provides stability & flexibility at different temperatures

Specially Deactivated Fused Silica Surface

Gives excellent peak shape for polar & nonpolar compounds

Highly Selective Stationary Phase

- High efficiency dimension and consistent film thickness delivers excellent separation of paraffins, isoparaffins, olefins, naphthenes, aromatics, and polar compounds
- Special designed Engineered Self Cross-linking™ (ESC) stationary phase for low bleed with exceptional peak shape, resolution, and separation

Results and Discussion

The spark ignition fuel matrix contains a variety of compounds including paraffins, isoparaffins, olefins, naphthenes, and related compounds which makes the DHA separation complex. This technical work describes the use of a Zebron™ ZB-DHA-PONA GC column in combination with a 5 meter ZB-DHA-PONA-TUNE column which are specially designed to provide excellent peak symmetry and resolution of hydrocarbon components and oxygenates. The primary ZB-DHA-PONA column's low polar 100 % Dimethylpolysiloxane phase provides very good separation of low and high boiling impurities of complex mixtures in the spark ignition fuels while the 5m ZB-DHA-PONA-TUNE GC tuning column provides a phenyl selectivity that helps the effective separation of aromatics such as benzene and toluene from hydrocarbons.

The DHA of a 29-component mixture was performed for this technical piece (**Figure 1**). The column efficiency of the Zebron ZB-DHA-PONA far exceeded the ASTM-6730-01 required minimum 400,000 plates for the pentane peak. Zebron ZB-DHA-PONA has a special deactivated fused silica surface that provides excellent peak shape for polar and nonpolar compounds, plus it's station-

ary phase has an Engineered Cross-linking™ stationary phase designed for lower bleed, which together are wrapped in a flexible polyimide outer coating that can withstand up to 340/360 °C maximum temperature (**Figure 2**). The smaller 0.25mm inner diameter along with the 100-meter column length provides an efficient mass transfer of the complex analytes for great separation and resolution of critical pairs.

Conclusion

The Zebron ZB-DHA-PONA GC column provides the optimal separation of DHA critical pairs with symmetric peaks and it is suitable for ASTM-D6730-01. The Zebron ZB DHA-PONA GC column with the 100 % Dimethylpolysiloxane phase in combination with ZB-DHA-PONA-TUNE is a great selectivity to provide separation of paraffins, isoparaffins, aromatics, naphthalene, olefins and Oxygenates.

Australia

t: +61 (0)2-9428-6444
auiinfo@phenomenex.com

Austria

t: +43 (0)1-319-1301
anfrage@phenomenex.com

Belgium

t: +32 (0)2 503 4015 (French)
t: +32 (0)2 511 8666 (Dutch)
beinfo@phenomenex.com

Canada

t: +1 (800) 543-3681
info@phenomenex.com

China

t: +86 400-606-8099
cninfo@phenomenex.com

Denmark

t: +45 4824 8048
nordicinfo@phenomenex.com

Finland

t: +358 (0)9 4789 0063
nordicinfo@phenomenex.com

France

t: +33 (0)1 30 09 21 10
franceinfo@phenomenex.com

Germany

t: +49 (0)6021-58830-0
anfrage@phenomenex.com

India

t: +91 (0)40-3012 2400
indiainfo@phenomenex.com

Ireland

t: +353 (0)1 247 5405
eireinfo@phenomenex.com

Italy

t: +39 051 6327511
italiainfo@phenomenex.com

Luxembourg

t: +31 (0)30-2418700
nlinfo@phenomenex.com

Mexico

t: 01-800-844-5226
tecnicomx@phenomenex.com

The Netherlands

t: +31 (0)30-2418700
nlinfo@phenomenex.com

New Zealand

t: +64 (0)9-4780951
nzinfo@phenomenex.com

Norway

t: +47 810 02 005
nordicinfo@phenomenex.com

Portugal

t: +351 221 450 488
ptinfo@phenomenex.com

Singapore

t: +65 800-852-3944
sginfo@phenomenex.com

Spain

t: +34 91-413-8613
espinfo@phenomenex.com

Sweden

t: +46 (0)8 611 6950
nordicinfo@phenomenex.com

Switzerland

t: +41 (0)61 692 20 20
swissinfo@phenomenex.com

Taiwan

t: +886 (0) 0801-49-1246
twinfo@phenomenex.com

United Kingdom

t: +44 (0)1625-501367
ukinfo@phenomenex.com

USA

t: +1 (310) 212-0555
info@phenomenex.com

**All other countries
Corporate Office USA** 

t: +1 (310) 212-0555
info@phenomenex.com

Ordering Information**Zebtron ZB-DHA-PONA GC Columns**

ID(mm)	df(μm)	Temp. Limits °C	Part No.
5-Meter			
0.25	1.0	-60 to 340/360	7AG-G042-22
50-Meter			
0.20	0.50	-60 to 360/370	7GE-G042-17
100-Meter			
0.25	0.50	-60 to 360/370	7MG-G042-17
150-Meter			
0.25	1.0	-60 to 340/360	7QG-G042-22

**Find Your Liner
Online!**

Easily search by part numbers, applications, injection mode, or system manufacturer for quick selection in under 1 minute!

www.phenomenex.com/FindLiner

**Inlet Consumables
Are Available Online!**

Need inlet seals, septa, or syringes? Explore hundreds of available parts online at:

www.phenomenex.com/InletGC



www.phenomenex.com

Phenomenex products are available worldwide. For the distributor in your country, contact Phenomenex USA, International Department at international@phenomenex.com

Our Technical Nerds are available nearly 24/7.

Live Chat today!

► www.phenomenex.com/chat ◀

BE-HAPPY™
guarantee

Your happiness is our mission. Take 45 days to try our products. If you are not happy, we'll make it right.

www.phenomenex.com/behappy

Terms and Conditions

Subject to Phenomenex Standard Terms and Conditions which may be viewed at www.phenomenex.com/TermsAndConditions.

Trademarks

Zebtron, Z-Liner, Engineered Self Cross-linking (ESC), and Be-Happy are trademarks of Phenomenex.

FOR RESEARCH USE ONLY. Not for use in clinical diagnostic procedures.

© 2019 Phenomenex, Inc. All rights reserved.