

# APPLICATIONS

## A Means of Improving the Resolution Between Acetaldehyde and Methanol in the Limit Test for 96% Ethanol (Ph.Eur.1317)

Zachary Woodward and Heiko Behr

Phenomenex, Inc., 411 Madrid Avenue., Torrance, CA 90501 USA

Data courtesy of WZZ Herbapol S, Katarzyna Piechota, Senior Assistant of Laboratory Quality Control

### Introduction

The European Pharmacopoeia has published a gas chromatography method for the purity analysis of 96% ethanol (Ph. Eur. Monograph 1317; Supplement 10), which stipulates that the peaks for residual amounts of acetaldehyde and methanol be resolved from one another with a minimum resolution factor of 1.5. Acetaldehyde and methanol can be challenging to resolve from one another within the confines of this method, however baseline resolution is imperative for the purity analysis of ethanol. The Ph. Eur. method specifies the use of a poly[(cyanopropylphenyl)-(dimethyl)siloxane capillary GC column, for which the Zebtron™ ZB-624 has been demonstrated to suitably resolve acetaldehyde from methanol.

Acetaldehyde is a precursor to ethanol during fermentation and is also classified as an irritant and carcinogen in humans. The ingestion of relatively low amounts of methanol (as little as 3 – 20 mL) can result in acute central nervous system toxicity in humans. This European Pharmacopoeia (Ph. Eur.) monograph was developed as a limit test for the determination of low levels of acetaldehyde and methanol in 96% ethanol solutions.


Operating conditions for GC methods that are published by the European Pharmacopoeia may be adjusted within certain allowable limits to further improve the existing method, or to overcome limitations of the instrument used during analysis. These allowances include adjusting the oven temperature to within  $\pm 10\%$  and the flow rate to within  $\pm 50\%$  of their respective published values. The ZB-624 column was demonstrated to suitably resolve acetaldehyde and methanol upon making adjustments within these allowable limits.

### Materials and Methods

#### Sample Preparation:

“Reference Solution b” was prepared as per Ph. Eur. Monograph 1317; Supplement 10, for which 50  $\mu\text{L}$  of acetaldehyde and 50  $\mu\text{L}$  of anhydrous methanol were diluted to 50 mL with the ethanol sample. A 100  $\mu\text{L}$  aliquot of this spiking solution was transferred to a 10 mL volumetric flask and diluted to volume with ethanol for analysis.

**Figure 1.**  
Chromatogram for Ethanol Analysis on ZB-624 using Modified Method Parameters



#### GC-FID Method Parameters - Modification of Ph. Eur. Monograph 1317; Supplement 10

**Column:** Zebtron ZB-624  
**Dimensions:** 30 meter x 0.32 mm x 1.80  $\mu\text{m}$   
**Part No.:** 7HM-G005-31  
**Injection:** Split 20:1 @ 200 °C, 1  $\mu\text{L}$

**Recommended Liner:** Zebtron PLUS Straight Z-Liner

**Liner Part No.:** AG2-0A03-05 (for Agilent® & Thermo Scientific® systems)  
AG2-4B03-05 (for Shimadzu® systems)


**Carrier Gas:** Helium @ 25 cm/sec (linear velocity)

**Oven Program:** 36 °C for 12 min, to 260 °C at 10 °C/min for 15 min

**Detector:** FID @ 280 °C

**Sample:** 10 ppm additional spike in 96% ethanol sample  
1) Acetaldehyde  
2) Methanol

**Figure 2.**  
Original Chromatogram for Ethanol Analysis on ZB-624



#### GC-FID Method Parameters – Original Ph. Eur. Monograph 1317; Supplement 10

**Column:** Zebtron ZB-624  
**Dimensions:** 30 meter x 0.32 mm x 1.80  $\mu\text{m}$   
**Part No.:** 7HM-G005-31  
**Injection:** Split 20:1 @ 200 °C, 1  $\mu\text{L}$

**Recommended Liner:** Zebtron PLUS Straight Z-Liner

**Liner Part No.:** AG2-0A03-05 (for Agilent® & Thermo Scientific® systems)  
AG2-4B03-05 (for Shimadzu® systems)

**Carrier Gas:** Helium @ 35 cm/sec (linear velocity)

**Oven Program:** 40 °C for 12 min, to 260 °C at 10 °C/min for 15 min

**Detector:** FID @ 280 °C

**Sample:** 10 ppm additional spike in 96% ethanol sample  
1) Acetaldehyde  
2) Methanol

# APPLICATIONS

## Results and Discussion

The ZB-624 afforded a resolution factor of 1.82 between acetaldehyde and methanol upon making suitable adjustments to Ph. Eur. Monograph 1317; Supplement 10 within the range of allowed adjustments, as shown in **Figure 1**. The original method was only able to generate resolution factors up to 1.23 across a selection of equivalent phases that shared the same column dimensions and film thickness, as shown in **Figure 2**.

“Reference Solution b” is spiked to approximate 10 ppm each of acetaldehyde and methanol with respect to ethanol. The original Ph. Eur. method calls for an initial isothermal hold of 40 °C for 12 minutes, while a lower isothermal hold of 36 °C allows for increased retention of the targeted impurities, and therefore increased separation. The original Ph. Eur. method also calls for a linear velocity of 35 cm/sec, whereas a lower flow-rate of 25 cm/sec affords more opportunities for mass transfer of the target impurities into and out of the stationary phase. Furthermore, the lower flow-rate preserves the integrity of the narrow analytical bands for each analyte as they transfer between the carrier gas and stationary phase at a lower temperature.

The lower isothermal hold of 36 °C is within the allowable  $\pm 10\%$  adjustment to the original temperature of 40 °C. Likewise, the lower flow-rate is well within the allowable  $\pm 50\%$  adjustment to the original 35 cm/sec flow of the published Ph. Eur. method.

## Conclusion

The selectivity of the Zebron ZB-624 column is particularly adept at resolving acetaldehyde from methanol when making allowable adjustments to the Ph. Eur. Monograph 1317; Supplement 10 analytical method. The resulting resolution is well above the minimum value that is stipulated within the monograph, affording the analyst leeway for fluctuations in GC instrument performance, and for routine cutting of each column end during repeated installations of the column. The combined levels of acetaldehyde and other acetal impurities within the ethanol sample may not exceed 10 ppm, and the ZB-624 column affords the resolution required to separate acetaldehyde from methanol for quantitative determination.

## References

1. European Pharmacopeia 10, Section 2.2.46, (2019)

**BE-HAPPY™**  
guarantee

Your happiness is our mission. Take 45 days to try our products. If you are not happy, we'll make it right.  
[www.phenomenex.com/behappy](http://www.phenomenex.com/behappy)

### Terms and Conditions

Subject to Phenomenex Standard Terms and Conditions which may be viewed at [www.phenomenex.com/TermsAndConditions](http://www.phenomenex.com/TermsAndConditions)

### Trademarks

Zebron is a trademark of Phenomenex. Agilent is a registered trademark of Agilent Technologies, Inc. Thermo-Scientific is a registered trademark of Thermo Fisher Scientific, Inc. Shimadzu is a registered trademark of Shimadzu Corporation.

### Disclaimers

Comparative separations may not be representative of all applications. Phenomenex is not affiliated with Agilent, Thermo-Scientific, or Shimadzu.

FOR RESEARCH USE ONLY. Not for use in clinical diagnostic procedures

© 2019 Phenomenex, Inc. All rights reserved.

## Ordering Information

Zebron ZB-624 GC Columns			
ID(mm)	df( $\mu$ m)	Temp. Limits °C	Part No.
<b>20-Meter</b>			
0.18	1.00	-20 to 260	7FD-G005-22
<b>30-Meter</b>			
0.25	1.40	-20 to 260	7HG-G005-27
0.32	1.80	-20 to 260	7HM-G005-31
0.53	3.00	-20 to 260	7HK-G005-36
<b>60-Meter</b>			
0.25	1.40	-20 to 260	7KG-G005-27
0.32	1.80	-20 to 260	7KM-G005-31
0.53	3.00	-20 to 260	7KK-G005-36
<b>75-Meter</b>			
0.53	3.00	-20 to 260	7LK-G005-36
<b>105-Meter</b>			
0.53	3.00	-20 to 260	7NK-G005-36

Note: If you need a 5 in. cage, simply add a (-B) after the part number, e.g., 7HG-G005-27-B. Some exceptions may apply. Agilent 6850 and some SRI and process GC systems use only 5 in. cages.

### Australia

t: +61 (0)2-9428-6444  
auinfo@phenomenex.com

### Austria

t: +43 (0)1-319-1301  
anfrage@phenomenex.com

### Belgium

t: +32 (0)2 503 4015 (French)  
t: +32 (0)2 511 8666 (Dutch)  
beinfo@phenomenex.com

### Canada

t: +1 (800) 543-3681  
info@phenomenex.com

### China

t: +86 400-606-8099  
cninfo@phenomenex.com

### Denmark

t: +45 4824 8048  
nordicinfo@phenomenex.com

### Finland

t: +358 (0)9 4789 0063  
nordicinfo@phenomenex.com

### France

t: +33 (0)1 30 09 21 10  
franceinfo@phenomenex.com

### Germany

t: +49 (0)6021-58830-0  
anfrage@phenomenex.com

### India

t: +91 (0)40-3012 2400  
indiainfo@phenomenex.com

### Ireland

t: +353 (0)1 247 5405  
eireinfo@phenomenex.com

### Italy

t: +39 051 6327511  
italiainfo@phenomenex.com

### Luxembourg

t: +31 (0)30-2418700  
nlinfo@phenomenex.com

### Mexico

t: 01-800-844-5226  
tecnicomx@phenomenex.com

### The Netherlands

t: +31 (0)30-2418700  
nlinfo@phenomenex.com

### New Zealand

t: +64 (0)9-4780951  
nzinfo@phenomenex.com

### Norway

t: +47 810 02 005  
nordicinfo@phenomenex.com

### Portugal

t: +351 221 450 488  
ptinfo@phenomenex.com

### Singapore

t: +65 800-852-3944  
sginfo@phenomenex.com

### Spain

t: +34 91-413-8613  
espinfo@phenomenex.com

### Sweden

t: +46 (0)8 611 6950  
nordicinfo@phenomenex.com

### Switzerland

t: +41 61 692 20 20  
swissinfo@phenomenex.com

### United Kingdom

t: +44 (0)1625-501367  
ukinfo@phenomenex.com

### USA

t: +1 (310) 212-0555  
info@phenomenex.com

### All other countries Corporate Office USA

t: +1 (310) 212-0555  
info@phenomenex.com

## www.phenomenex.com

Phenomenex products are available worldwide. For the distributor in your country, contact Phenomenex USA, International Department at [international@phenomenex.com](mailto:international@phenomenex.com)