

Application Report 403


Author: Katherine Stenerson
Acquisition System: 9499
Notebook Reference: 1569-040

Aroclors on the 20m x 0.18 mm I.D., 0.36 μ m SLB-5ms

Aroclors are commercial mixtures of PCB congeners. Due to their stability, they have been used in many different industrial and commercial applications. Because of their toxicity and ability to bioaccumulate, production of these materials ceased in 1977. This application demonstrates the separation of two common Aroclor mixtures, Aroclor 1016 and Aroclor 1260 on a 20 m x 0.18 mm I.D., 0.36 μ m SLB-5ms. The higher efficiency of the 0.18 mm I.D. allowed for good resolution of the PCB congeners, and subsequent pattern recognition of the Aroclor mixtures, while maintaining an analysis time of < 20 minutes.

Key Words

Aroclor, PCB, SLB-5ms, US EPA Method 608, US EPA Method 8082, 28576-U, 46846-U



G003582

Conditions

column: SLB-5ms, 20 m x 0.18 mm I.D., 0.36 μ m (28576-U)
oven: 100 °C (2 min.), 15 °C/min. to 325 °C (3 min.)
inj.: 250 °C
det.: micro-ECD, 325 °C
carrier gas: helium, 0.5 mL/min., constant flow
injection: 1.0 μ L, splitless (0.75 min.)
liner: 4 mm I.D., single taper
sample: Aroclor standard mix 1 (46846-U) diluted to 500 ppb/50 ppb (Aroclors/surrogates) in n-hexane

Peak IDs

1. Tetrachloro-m-xylene (surr.)
2. Aroclor 1016
3. Aroclor 1260
4. Decachlorobiphenyl (surr.)