

Application Report 478


Author: Katherine Stenerson
Acquisition System: 5973 GC-MS
Notebook Reference: 1569-084

US EPA Method OLM04.2 SVOA (CLP Semivolatiles) by GC-MS on the Supelco SLB™-5ms

The Contract Laboratory Program (CLP) is a nationwide network consisting of testing laboratories, contractors, and EPA employees. Its primary function is to provide analytical data in support of EPA's Superfund program. To ensure the data is of consistent and documented quality, labs participating in the program must analyze samples using methods described in the CLP Statement of Work (SOW). This application illustrates the separation of semivolatile compounds, on a 30 m x 0.25 mm I.D., 0.25 µm SLB-5ms, including those commonly analyzed as part of the CLP.

Key Words

Semivolatiles, CLP, SOW, US EPA Method OLM04.2 SVOA, SLB-5ms, 28471-U, 506508, 47508, 47509-U, 47514, 40060, 40026, 46955-U, 47960-U



G003745

Conditions

column: SLB-5ms, 30 m x 0.25 mm I.D., 0.25 µm (28471-U)
oven: 40 °C (2 min.), 15 °C/min. to 240 °C, 10 °C/min. to 330 °C (2 min.)
inj.: 250 °C
MSD interface: 330 °C
scan range: m/z 40-450
carrier gas: helium, 1.0 mL/min. (11 min.), ramp to 1.5 mL/min. (constant for remainder of run)
injection: 1 µL, pulsed (15 psi until 0.10 min.) splitless (0.75 min.)
liner: 4 mm I.D., single taper
sample: 50 ppm of an 80 component semivolatile standard, plus 6 internal standards (at 40 ppm) in methylene chloride

27. 2,4-Dichlorophenol	59. Azobenzene
28. 1,2,4-Trichlorobenzene	60. 2,4,6-Tribromophenol (surr.)
29. Naphthalene-d ₈ (I.S.)	61. 4-Bromophenyl phenyl ether
30. Naphthalene	62. Hexachlorobenzene
31. 5-Chloroaniline	63. Atrazine
32. Hexachlorobutadiene	64. Pentachlorophenol
33. Caprolactam	65. Phenanthrene-d ₁₀ (I.S.)
34. 4-Chloro-3-methylphenol	66. Phenanthrene
35. 2-Methylnaphthalene	67. Anthracene
36. Hexachlorocyclopentadiene	68. Carbazole
37. 2,4,6-Trichlorophenol	69. Di-n-butyl phthalate
38. 2,4,5-Trichlorophenol	70. Fluoranthene
39. 2-Fluorobiphenyl (surr.)	71. Pyrene
40. 1,1'-Biphenyl	72. Terphenyl-d ₁₄ (surr.)
41. 2-Chloronaphthalene	73. Butyl benzyl phthalate
42. 2-Nitroaniline	74. 3,3'-Dichlorobenzidine
43. Dimethyl phthalate	75. Benzo(a)anthracene
44. 2,6-Dinitrotoluene	76. Chrysene-d ₁₂ (I.S.)
45. Acenaphthylene	77. Chrysene
46. 3-Nitroaniline	78. Bis(2-ethylhexyl)phthalate
47. Acenaphthene-d ₁₀ (I.S.)	79. Di-n-octyl phthalate
48. Acenaphthene	80. Benzo(b)fluoranthene
49. 2,4-Dinitrophenol	81. Benzo(k)fluoranthene
50. 4-Nitrophenol	82. Benzo(a)pyrene
51. 2,4-Dinitrotoluene	83. Perylene-d ₁₂ (I.S.)
52. Dibenzofuran	84. Indeno(1,2,3-cd)pyrene
53. Diethyl phthalate	85. Dibenzo(a,h)anthracene
54. 4-Chlorophenyl phenyl ether	86. Benzo(g,h,i)perylene
55. Fluorene	
56. 4-Nitroaniline	
57. 2-Methyl-4,6-dinitrophenol	
58. n-Nitrosodiphenylamine	

Peak IDs

1. 2-Fluorophenol (surr.)	14. 2-Methylphenol
2. Benzaldehyde	15. Bis(2-chloroisopropyl)ether
3. Phenol-d ₆ (surr.)	16. 4-Methylphenol
4. Phenol	17. n-Nitroso-di-n-propylamine
5. Bis(2-chloroethyl)ether	18. Acetophenone
6. 2-Chlorophenol-d ₄ (surr.)	19. Hexachloroethane
7. 2-Chlorophenol	20. Nitrobenzene-d ₅ (surr.)
8. 1,3-Dichlorobenzene	21. Nitrobenzene
9. 1,4-Dichlorobenzene-d ₄ (I.S.)	22. Isophorone
10. 1,4-Dichlorobenzene	23. 2-Nitrophenol
11. Benzyl alcohol	24. 2,4-Dimethylphenol
12. 1,2-Dichlorobenzene (surr.)	25. Bis(2-chloroethoxy)methane
13. 1,2-Dichlorobenzene	26. Benzoic acid