

"<u>Overview</u> of KnowItAll Software for Unknown Identifications in EI GC-MS Analyses"

James Little Mass Spec Interpretation Services

Handouts for Videos: Website: "Little Mass Spec and Sailing" <u>https://littlemsandsailing.wordpress.com</u>


Note: In depth training videos/handouts on this new mass spec software on my website!


5 Primary Software Functions


ProcessIt:Manually process nominal and accurate mass EI GC-MS filesSearchIt:Manually search EI spectra, structures, names, etc.Minelt:Display EI manual search resultsMS Expert:Automatic deconvolution of nominal and accurate mass EI GC-MS filesChemWindow:Structure drawing program

Also, the new novel Adaptive MS search greatly extends the identification capabilities for compounds *not* present in EI libraries!
 Discussed in detail in another video


Left Side Bar for Selecting 5 Functions








ProcessIt for Manually Processing El GC-MS Data




Sending Selected El Spectrum to SearchIt Window



Left click on "Start a new search"

SearchIt Window to Define Search Parameters Such as Libraries, Types of Search, etc.



Search Results Shown in Minelt >1.2 Million Spectra Searched



MS Expert for Automatic Deconvolution and Library Searching of Complex Mixtures in GC-MS

- Chemical background removed
- Coeluting species resolved
- Both nominal and accurate EI GC-MS data
- ➢ 50 components characterized in ~6 secs

KnowItAll Informatics System 2023, Analytical Edition File Edit View Analysis License Help 😓 📲 🚽 🚽 Transfer to: 🚸 Reportit 🎕 Searchit 🚯 Minelt Database 💉 QC Expert 🗞 Assignit NMR 📀 Processit 🐼 Analyzeit Spectral Processing - # MS Expert 👺 🏡 🗁 🖳 🐚 🥑 📢 Analysis Profiles: <no profile> · 🎛 💕 Basics Components Data TIC 1e+06-Spectral Processing Comp. RT [m... # Match Sco... HQI R.H. 40 0.7241 99.22 99.17 99.68 I Argon 44 0.9833 1 Carbamic acid, monoammonium salt 99.95 99.95 99.95 800000. ProcessIt 69 2 0 2 9 7 Methane, chloro-98.75 98.62 99.92 H 1 In the 43 Methane, oxybis-3.1198 H 1 97.68 97.68 97.68 119 99.61 99.61 99.61 600000-3.2225 I 1 Methyl bromide **MS Expert** 131 4.2028 ⊞ 1 Ethene, chloro-97.51 97.51 97.54 91 5.3466 ⊞ 1 Ethyl chloride 99.82 99.82 99.82 400000-57 7.6648 H 1 Methylene chloride 99.79 99.79 99.80 8.5025 H 1 Acetone 95.43 95.23 97.25 141 9,1442 ⊞ 1 Dithioxomethane 98.82 98.82 98.82 51 200000. 9.7463
H 1 Trichloromonofluoromethane 99.85 99.85 99.87 169 10.4328 II 1 Ethene 1.1-dichloro-99.27 99.26 99.30 181 10.7473 1 METHANE, DIMETHOXY 96.45 96.39 97.03 39 11.0188
H 1 Methane, bromochloro-99.38 99.38 99.39 83 11.9073 1 Ethane, 1,1-dichloro-98.99 98.99 99.00 63 mir 99.64 99.63 99.69 12.6332 II Ethylene 1.2-dichloro- (E) 85 13.3568 I 1 METHANE, TRICHLORO 99.53 99.52 99.53 70 14 0243 II 1 2-Dichloroethane-d4 93.36 93.36 93.39 100-- Extracted Spectrum Raw Spectrum Reference lons(s): 78 - 83 74 76 77 61 63 Bh Parameters: Contributing lons: 78 Energy: 0.0054 Base Peak Vertex: 55755.5 75 00 35 50 60 65 ROR m/z 78 - Extracted Spectrum Automatic 100-MTX #40910 BENZENE 74 76 77 50 52 53 61 63 55 70 75 90 60 65 m/z V50 - James Little.cdf × Spectral Analysis

<u>More Detailed</u> Training of Newly Released Wiley KnowItAll Mass Spec Software on My Website

James Little Website: "Little Mass Spec and Sailing" <u>https://littlemsandsailing.wordpress.com</u>