


Analysis of Soil Samples using the HTD

The analysis of complex environmental materials such as soils and sludges for semivolatile pollutants may be simplified using rapid thermal desorption techniques, in which the sample material is desorbed directly to the gas chromatograph. This process eliminates solvent extraction by volatilizing the organics of interest rapidly from the sample matrix and transferring them immediately to the injection port of the GC without trapping or dilution.

To achieve such rapid analysis, it is important that the instrument be capable of heating samples which weigh from 5 to 200 mg in a rapid, yet controlled manner. If the samples are heated too slowly, an additional collection or trapping step will be required; if heated too rapidly, or to an excessive temperature, unwanted byproducts may be generated. The CDS High Temperature Desorber (HTD) is controlled by the CDS Analytical Pyroprobe temperature controller for exacting supervision of heating rate and temperature. Using the Pyroprobe 2000, multiple step controls may be programmed for even more versatile temperature management.

Figure 1.



Samples for the HTD are placed into a quartz tube 25 mm in length and 6 mm inside diameter. This tube is then inserted into a large platinum heating coil, which is controlled by the Pyroprobe. The sample coil is placed into the GC carrier gas stream in a heated interface, then rapidly heated to desorb the volatile and semivolatile organic contaminants.

Figure 1 compares soil contaminated with diesel fuel to a run of blank soil. In each case, a 5 mg sample was placed into the quartz tube and held in place with glass wool. The HTD coil was programmed to heat the sample to 200°C for 50 seconds, with an initial 3 second burst to 1000° to shorten the heating time. The gas chromatograph run was initiated immediately since no collection step is required, so the entire analysis takes only as long as the GC run itself.

Equipment:

DESORPTION

CDS HTD platinum coil, programmed to 1000°C for 3 seconds, then 200° for 50 seconds. Interface temperature: 100°C isothermal.

CHROMATOGRAPHY

30 m x 0.53 mm SE-54 megabore capillary column operated at 7 ml/min. helium flow. Program: Initial temperature 40°C for 2 minutes, then 6°/min. to 225°C. Detection: Flame Ionization Detector.

FOR MORE INFORMATION
CONCERNING THIS AND RELATED
APPLICATIONS, WE RECOMMEND
THE FOLLOWING READING:

Systems approach to automated cryofocusing in purge and trap, headspace and pyrolytic analysis, T.P. Wampler, W. Bowe, J. Higgins, American Laboratory, 17, 8 (1985) 82-87.

A dedicated purge and trap system for environmental analysis, J.W. Washall, T. Wampler, W. Bowe, K. Kristunas, CDS Application paper #155.

© courtesy of CDS Analytical Inc.

Headquarters

JSB International
Tramstraat 15
5611 CM Eindhoven
T +31 (0) 40 251 47 53
F +31 (0) 40 251 47 58

Zoex Europe
Tramstraat 15
5611 CM Eindhoven
T +31 (0) 40 257 39 72
F +31 (0) 40 251 47 58

Sales and Service

Netherlands
Apolloweg 2B
8239 DA Lelystad
T +31 (0) 320 87 00 18
F +31 (0) 320 87 00 19

Belgium
Grensstraat 7
Box 3 1831 Diegem
T +32 (0) 2 721 92 11
F +32 (0) 2 720 76 22

Germany
Max-Planck-Strasse 4
D-47475 Kamp-Lintfort
T +49 (0) 28 42 9280 799
F +49 (0) 28 42 9732 638

UK & Ireland
Cedar Court,
Grove Park Business Est.
White Waltham, Maidenhead
Berks, SL6 3LW
T +44 (0) 16 288 220 48
F +44 (0) 70 394 006 78

info@go-jsb.com
www.go-jsb.com

