

Qualitative analysis of 4 impurities in enalapril using liquid chromatography-ion trap time of flight hybrid mass spectrometry

IMSC 2012 PMo-004

Feng Ji, Yueqi Li, Jing Dong
Analytical Applications Center, Shimadzu
International Trading (Shanghai) Co., Ltd., 6F No. 16
Chao Yang Men Wai Street, China Life Tower,
Beijing 100020

Qualitative analysis of 4 impurities in enalapril using liquid chromatography-ion trap time of flight hybrid mass spectrometry

Introduction

Enalapril Maleate is the maleate salt of enalapril, the ethyl ester of a long acting angiotensin converting enzyme inhibitor. Enalapril is a pro-drug; it is bioactivated by hydrolysis of the ethyl ester to enalaprilat, which is the active angiotensin converting enzyme inhibitor. Enalapril has side effects: fainting; urinating more; chest pain, etc.

However, by-products, unknown impurities and their structures are a big concern. Currently, high-performance liquid chromatography (HPLC) has been the essential method for analysis of enalapril. Atmospheric pressure ionization sources, such as electrospray ionization (ESI) are routinely used in drug impurities analysis for sensitivity.




Fig. 1 Molecular structure of Enalapril and 4 impurities.

A liquid chromatography/electrospray ion trap time-of-flight mass spectrometry (LCMS-IT-TOF) was used in this study. It has the ability to provide multistage tandem spectra (MS^n) with accurate masses (error < 5 ppm) in both MS and MS^n modes. ToF analyzer allows fast acquisition of full spectra with high sensitivity and elevated mass resolution (12000 FWHM). With this technique, we can acquire reliable structural information about the ionized

molecule and the product ions, based on the exact mass and the multistage tandem mass spectrometric analyses (MS^n). In addition, as the mass/charge ratio gets smaller, the number of candidates gets smaller, too. Since mass/charge ratios of the product ions in the highest order spectra are usually less than 200 Da and the neutral losses are usually less than 100 Da, their formulae could be confirmed with significantly higher confidence.




Fig. 2 chromatogram of 4 impurities.

Qualitative analysis of 4 impurities in enalapril using liquid chromatography-ion trap time of flight hybrid mass spectrometry




Fig. 3 LCMS-IT-TOF basic configuration diagram.

Experimental

Enalapril samples were separated on a Shimadzu Shim-pack XR-ODS column (2.0 mm I.D. x 75 mm L., 2.2 μm) with isocratic elution using 10 mM ammonium acetate and acetonitrile as mobile phase. High concentration of enalapril (300 mg/L) was injected to concentrate impurities. Flow path switching valve was used in the instrument, cutting the enalapril main constituent to waste line in order

not to contaminate MS detector. Chromatogram is shown in Fig. 2. Four impurities were separated and detected by LCMS-IT-TOF. MSⁿ spectra (n ≥ 4) of each impurity were obtained in both positive and negative ion mode. Formulae, chemical structures of impurities in an Enalapril sample were suggested with supporting results on the probable fragmentation pathways are shown in Fig. 4-5.




Fig. 4 Probable fragmentation pathways of Enalapril.

Qualitative analysis of 4 impurities in enalapril using liquid chromatography-ion trap time of flight hybrid mass spectrometry

Conclusions

To our knowledge, this is the first complete characterization of the fragmentation of enalapril impurities, mainly are its by-products, using LCMS-IT-TOF tandem mass spectrometry with accurate mass measurements.