

Application of a sensitive liquid chromatography-tandem mass spectrometric method to determination of octreotide in human plasma

ASMS 2015 MP 291

Yuling Song, Jinting Yao, Hongyuan Hao, Qiang Li,
Taohong Huang, Shin-ichi Kawano, Yuki Hashi
Shimadzu Global COE, Shimadzu (China) Co., Ltd.

Application of a sensitive liquid chromatography-tandem mass spectrometric method to determination of octreotide in human plasma

Overview

A fast and sensitive method for octreotide in human plasma by ultra high performance liquid chromatography-tandem mass spectrometry was developed. The lower limit of quantitation was 5.0 pg/mL and the results of other parameters for method validation were good.

Introduction

Analysis of therapeutic peptides in bio-matrices is very important in the pharmaceutical industry especially in the area of oncology therapy. Octreotide (Fig. 1) is a synthetic linc-acting cyclic octapeptide. It is a somatostatin analog that has a longer half life and more selectivity. Radio-immunoassay methods have been used for determination of octreotide in human blood, urine and saliva after dosing octreotide. A LC/MS method has been developed for determination of octreotide in human

plasma with a low limit of quantitation (LLOQ) of 500 pg/mL. A LC-MS/MS method has been used for determination of octreotide with LLOQ of 50 pg/mL. However, this method required 400 μ L plasma and laborious sample pretreatment. In this study, a rapid, sensitive method using a simple μ Elution solid phase extraction and ultra high performance liquid chromatography-tandem mass spectrometry was developed and validated.




Fig.1 Structure of octreotide

Methods

Sample Preparation

- (1) Add 200 μ L of plasma into the polypropylene tube, add 50 μ L of internal standard (IS) working solution (1.0 ng/mL leuprolide), add 200 μ L of 4% phosphoric acid, seal the tube and vortex for approximately 1 minute.
- (2) Condition the μ Elution WCX SPE plate with 200 μ L methanol then equilibrate with 200 μ L of water.
- (3) Load the sample onto the SPE plate, wash with 400 μ L of 5% ammonium hydroxide, then with 400 μ L of 20% acetonitrile, and finally elute the sample with 50 μ L of 1% TFA in 75:25 acetonitrile/water (v/v) into a 96-position, 2.0 mL, polypropylene plate.
- (4) Dilute the sample with 150 μ L of water and vortex for 3 minutes.
- (5) Inject the sample into the UHPLC-MS/MS system.

Application of a sensitive liquid chromatography-tandem mass spectrometric method to determination of octreotide in human plasma

LC-MS/MS Analysis

The analysis was performed on a Shimadzu Nexera UHPLC system (Kyoto, Japan). A triple quadrupole mass spectrometer (Shimadzu LCMS-8050, Kyoto, Japan) was connected to the UHPLC instrument via an ESI interface.

HPLC (Nexera UHPLC system)	
Column	: InertSustain Swift C18 (2.1 mm.i.d. x 100 mmL, 1.9 µm, GL Sciences)
Mobile phase A	: water with 0.1% formic acid
Mobile phase B	: acetonitrile
Flow rate	: 0.4 mL/min
Column temperature	: 40 °C
Mode	: Gradient Elution
MS (LCMS-8050 triple quadrupole mass spectrometer)	
Ionization	: ESI
Polarity	: Positive
Ionization voltage	: +4.0 kV
Nebulizing gas flow	: 3.0 L/min
Heating gas flow	: 15.0 L/min
Drying gas flow	: 5.0 L/min
Interface temperature	: 400 °C
Heat block temperature	: 400 °C
DL temperature	: 200 °C
Mode	: MRM

Table 1 MRM parameters

Compound	Precursor <i>m/z</i>	Product <i>m/z</i>	Dwell time (ms)	Q1 Pre bias (V)	CE (V)	Q3 Pre bias (V)
Octreotide	510.30	120.00	100	-36.0	-31.0	-21.0
Leuprolide	605.40	249.05	100	-22.0	-29.0	-26.0

Results

Human plasma samples containing octreotide ranging from 5 to 2000 pg/mL were prepared, extracted by µSPE and the final extracts were analyzed by LC-MS/MS. The typical MRM chromatogram of octreotide is presented below (Fig.2). The linear regression for octreotide was found to be >0.999. Excellent precision and accuracy were maintained, demonstrating a linear dynamic range suitable for real-world applications. The LOQ for octreotide was 5 pg/mL, which meet the criteria for

precision and accuracy within ±20% both within run and between run. The intra- and inter-day precision and accuracy of the assay were investigated by analyzing QC samples. Intra-day precision values at three concentration levels 15 pg/mL, 150 pg/mL and 1500 pg/mL were between 2.32~7.43% and inter-day precision was between 2.44~8.99%. accuracy was within ±15.0%.

Application of a sensitive liquid chromatography-tandem mass spectrometric method to determination of octreotide in human plasma

The recoveries of octreotide were $70.6 \pm 7.8\%$, $71.5 \pm 1.8\%$ and $70.1 \pm 1.9\%$ at three concentration levels individually. The matrix effect factors for three concentration levels were all larger than 80% and

internal standard normalized factors for each levels were larger than 90%. There was no significant carryover after analysis of highest concentration of calibration (Fig. 4).




Figure 1 Representative MRM chromatograms of spiked (5 pg/mL) human plasma (left: octreotide, right: IS)




Figure 2 Representative MRM chromatograms of blank plasma (left: octreotide, right: IS)

Application of a sensitive liquid chromatography-tandem mass spectrometric method to determination of octreotide in human plasma




Figure 3 Calibration curve of octreotide in human plasma

Table 2 Calibration curve information

Compound	Calibration curve	Linear range (ng/mL)	Accuracy (%)	r
Octreotide	$Y = (2.89 \times 10^{-3})X + (-6.34 \times 10^{-3})$	5~2000	89.6~112.3%	0.9993

Table 3 Accuracy and precision for the analysis of octreotide in human plasma (in pre-study validation, n=3 days, six replicates per day)

Added conc. (ng/mL)	Intra-day precision (%RSD)	Inter-day precision (%RSD)	Accuracy (%)
15	7.43	8.99	87.5~109.4
150	2.41	4.37	95.4~112.1
1500	2.32	2.44	95.4~103.1

Table 4 Recovery for QC samples (n=6)

QC Level	Concentration (ng/mL)	Recovery (%)
LQC	15	70.6
MQC	150	71.5
HQC	1500	70.1

Application of a sensitive liquid chromatography-tandem mass spectrometric method to determination of octreotide in human plasma

Table 5 Matrix effect for QC samples (n=6)

QC Level	Added conc. (ng/mL)	Matrix factor	IS-Normalized matrix factor
LQC	15	82.1%	90.5%
MQC	150	87.6%	96.6%
HQC	1500	89.7%	98.9%




Figure 4 MRM chromatograms of blank plasma after ULOQ sample (Carryover test, left: octreotide, right: IS)

Conclusions

Results of parameters for method validation such as dynamic range, linearity, LLOQ, intra-day precision, inter-day precision, and matrix effect factors were excellent. The sensitive LC-MS/MS technique provides a powerful tool for the high-throughput and highly selective analysis of octreotide in DMPK study.