

Analysis of residual Ethylene Oxide in Intraocular lenses using HS-GCMS

ASMS 2013 ThP11-200

Dheeraj Handique, Ankush Bhone, Durvesh Sawant,
Prashant Hase, Sanket Chiplunkar, Ajit Datar,
Jitendra Kelkar, Pratap Rasam
Shimadzu Analytical (India) Pvt. Ltd., 1 A/B Rushabh
Chambers, Makwana Road, Marol, Andheri (E),
Mumbai-400059, Maharashtra, India.

1. Introduction

It is imperative to employ effective agents to disinfect and sterilize instruments and equipments, used for patient care and healthcare staff. There are many methods of disinfection and it is very important to confirm the residues of disinfectant in finished product. In this poster we discuss the quantitation of Ethylene Oxide (EtO) as disinfectant residue in intraocular lenses. Ethylene Oxide sterilization is mainly used to sterilize medical and pharmaceutical products that cannot support conventional high temperature steam sterilization, such as devices that incorporate electronic components, plastic containers and intraocular lenses^[1]. Boiling point of pure EtO is 10.73°C at atmospheric pressure (The molecular structure of EtO is given in Fig. 1). It is highly volatile and can be easily




Fig. 1 Structure of Ethylene Oxide

removed from product. As it is highly volatile there are less chances of finding it in product as an impurity. EtO is a potential carcinogen. Since it is still used in sterilization of medical products, it is necessary to quantify EtO precisely at very low level. Occupational Safety and Health Administration (OSHA) has set the permissible exposure limit for EtO as 1.0 ppm^[2]. The objective of this study is to quantify EtO at very low concentration as the intraocular lenses are permanently implant in human eye. (Picture of intraocular lens is given in Fig. 2). The analysis is carried out by using Shimadzu Headspace-Gas chromatograph - Mass Spectrometer (GCMS-QP2010 Ultra coupled with HS-20).



Fig. 2 Intraocular lens

2. Method of Analysis

2-1: Extraction of Ethylene Oxide from Intraocular Lenses

Intraocular lenses were procured from local medical store. Standard stock solution of EtO (500 ppm in Dimethyl sulfoxide) was procured from Sigma Aldrich. Dimethyl sulfoxide was used as diluent for further sample and

standard solution preparations. HS-GCMS technique was used for quantitation of EtO at very low concentration. Solutions were prepared as follows,

- 1) Blank Solution – 5 mL diluent was added to headspace vial and was crimped tightly with automated crimper.
- 2) Standard Solutions – 5 mL of each Std. EtO in the range of 2.5 to 3000 ppb was added to headspace vial and was crimped tightly with automated crimper.
- 3) Sample Solution – A piece of intraocular lens was transferred to headspace vial, 5 mL diluent was added and vial was crimped tightly with automated crimper

Partial method validation was carried out by performing Reproducibility, Linearity, LOD-LOQ determination and Recovery study. For validation, solutions of different concentrations were prepared using standard stock solution of EtO (500 ppm) as mentioned in Table 1.

Table 1 Method validation parameters

| Parameter | Concentration |
|------------------------|-------------------------------|
| Reproducibility | 1000 ppb |
| Linearity | 5 levels – 100 ppb – 3000 ppb |
| Accuracy/Recovery | 3 levels – 100 ppb – 500 ppb |
| Precision at LOQ level | 2.5 ppb |

Analysis of residual Ethylene Oxide in Intraocular lenses using HS-GCMS

2-2. HS-GCMS Analytical Conditions

Samples were analyzed using HS-20 coupled with GCMS-QP2010 Ultra (Fig. 3) as per the conditions given in Table 2.



Fig. 3 HS-20 coupled with GCMS-QP2010 Ultra by Shimadzu

Table 2 HS-GCMS analytical parameters

Headspace parameters

| | |
|------------------------|------------|
| Mode | : Loop |
| Oven Temp | : 100°C |
| Sample Line Temp | : 110°C |
| Transfer Line Temp | : 120°C |
| Equilibrating Time | : 30.0 min |
| Pressurizing Time | : 1.0 min |
| Pressure Equilib. Time | : 0.10 min |
| Load Time | : 0.50 min |
| Load Equilib. Time | : 0.10 min |
| Injection Time | : 1.0 min |
| Needle Flush Time | : 5.0 min |
| GC Cycle Time | : 26.0 min |

Chromatographic parameters

| | | | |
|-------------------|--------------------------------------|----------------|-----------------|
| Column | : Rtx-624 (60 m × 0.53 mm × 3.00 μm) | | |
| Injection Mode | : Split | | |
| Split Ratio | : 5.0 | | |
| Carrier Gas | : Helium | | |
| Flow Control Mode | : Linear Velocity | | |
| Linear Velocity | : 59.8 cm/sec | | |
| Pressure | : 12.1 kPa | | |
| Column Flow | : 5.50 mL/min | | |
| Total Flow | : 33.0 mL/min | | |
| Total Run Time | : 16.0 min | | |
| Column Oven Temp | : Rate °C /min | Temperature °C | Hold time (min) |
| | 30.0 | 40.0 | 5.0 |
| | | 220.0 | 5.0 |

Mass Spectrometry parameters

| | |
|-----------------|-----------------|
| Ion Source Temp | : 200°C |
| Interface Temp | : 220°C |
| Ionization Mode | : EI |
| Event Time | : 0.30 sec |
| Mode | : SIM |
| m/z | : 29, 43 and 44 |
| Start Time | : 2.0 min |
| End Time | : 4.0 min |

3. Results

3-1. Fragmentation of Ethylene Oxide

According to fragmentation of Ethylene Oxide, molecular ion peak was m/z -44 with base peak at m/z -29, which was used for quantitation where as fragment ion with m/z -43 and m/z -44 were used as reference ions. Mass chromatograms of 1000 ppb EtO standard solution with

m/z 29,43 and 44 are shown in Fig. 4, Mass spectrum of EtO is shown in Fig. 5. Validation data is summarized in Table 3. Fig. 6 and 7 shows, overlay mass chromatograms for m/z -29 at different concentrations and calibration curve for linearity levels, respectively.




Fig. 4 Mass chromatograms of 1000 ppb Ethylene Oxide standard solution




Fig. 5 Mass spectrum of Ethylene Oxide

3-2. Summary of validation results

Table 3 Summary of results for validation parameter

| ID | Compound Name | Parameter | Concentration | Result |
|----|----------------|------------------|--------------------|--|
| 1 | Ethylene Oxide | Reproducibility | 1000 ppb | %RSD is 0.8 for area (n=6) |
| 2 | | Linearity | 100 ppb - 3000 ppb | Correlation Coefficient is 0.9999* |
| 3 | | LOD | 2.5 ppb - 10 ppb | 0.8 ppb** |
| 4 | | LOQ | | 2.5 ppb** |
| 5 | | Precision at LOQ | 2.5 ppb | Average of S/N ratio is 16 (n=6) %RSD is 9.6 for area (n=6) |

*Linearity level– 100 ppb, 250 ppb, 500 ppb,1000 ppb and 3000 ppb. For linearity, refer Fig. 6 and Fig. 7.

**As per software calculations.

Analysis of residual Ethylene Oxide in Intraocular lenses using HS-GCMS




Fig. 6 Linearity levels overlay mass chromatogram




Fig. 7 Linearity levels calibration curve

3-3. Quantification of Ethylene Oxide in intraocular lens sample

Analysis of intraocular lens samples was done as per the given method. Recovery studies were carried out by spiking 100 ppb, 250 ppb and 500 ppb of standard solutions in

sample of intraocular lenses. Fig. 8 gives overlay mass chromatogram of spiked and unspiked samples. Table 4, gives the summary of results.




Fig. 8 Spiked and unspiked sample solution overlay mass chromatograms

Table 4 Summary of results for sample analysis

| ID | Sample Name | Parameter | Concentration | Result |
|----|---|-----------|----------------|-----------------|
| 1 | Unspiked Sample | Precision | NA | Below LOQ level |
| 2 | Intraocular lens samples spiked with different linearity levels standards | Recovery | 100 ppb spiked | Recovery - 87% |
| | | | 250 ppb spiked | Recovery - 92% |
| | | | 500 ppb spiked | Recovery - 97% |

Analysis of residual Ethylene Oxide in Intraocular lenses using HS-GCMS

4. Conclusion

- HS-GCMS method was developed for quantitation of residual EtO present in intraocular lenses sample. Part method validation was performed successfully. Results obtained for Reproducibility, Linearity, LOQ and Recovery studies were well within limit, as per ICH guidelines^[3].
- With “Low Carryover” – The characteristics feature of HS-20 headspace, reproducibility even at very low concentration level could be achieved easily.
- High speed scan rate 20,000 u/sec is the characteristic feature of GCMS-QP2010 Ultra mass spectrometer, useful for quantitation of residual EtO at very low level (ppb level) with high sensitivity.

5. References

- [1] “Biological Evaluation of Medical Devices-Part 7: Ethylene Oxide Sterilization Residuals,” ANSI/AAMI/ISO 10993-7, Arlington, VA, Association for the Advancement of Medical Instrumentation, 1995.
- [2] Potter, Wayne “OSHA Method No. 30, Ethylene Oxide ”, OSHA Analytical Laboratory, Salt Lake City, Utah 84115, August 1981.
- [3] ICH guidelines Q2(R1), Validation of Analytical Procedure Text And Methodology, 2005