

It Isn't Always the Column: Troubleshooting Your HPLC Separation

Jean Lane
Applications Engineer
LC Columns and Consumables Technical Support
July 25, 2023



You've Recognized that There Is a Problem

Ask Questions:

When did the system and chromatography last function properly?

Has anything been changed?

For the method, was the procedure followed correctly?

Are the instrument settings correct?

What is the problem that is being seen?

Where might the problem be?


- Pump
- Injector/autosampler
- Column
- Detector
- Data system
- Mobile phase
- Sample
- Tubing/fittings
- User

The problem could be one, some, many, or all of these independently or together.

How's that for a challenge?

Common Symptoms and Problems



- Pressure
- Peak shape
- Retention
- Baseline



- Increased pressure
- Low pressure
- Leaks
- Pressure fluctuations
- Tailing
- Peak splitting and doubling
- Fronting
- Broadening
- Changing retention time
- Loss of resolution
- Noisy baseline
- Drifting baseline
- Reduced intensity or sensitivity

Understand Your HPLC System

Know your flow path



Changes in System Pressure

Causes of Increase in Back Pressure


- Particles leading to blockage can come from sources located both *outside* and *inside* the LC system:
 - Solvent, buffer
 - Microbial growth in solvent reservoirs
 - The sample
 - Wear of LC components – piston seals, autosampler valve
- Debris will either be captured on the a filter, frit, or inline filter (inexpensive replaceable frit) or a column frit (column = expensive)

Reduce LC problems by eliminating most common sources of flow blockage:
Preventing this is the key


Filter, filter, filter

Blockages and Clogging

Characteristics	
Parts affected	<p>Blockages:</p> <ul style="list-style-type: none"> • Capillaries, needle and needle seat • Detector flow cells <p>Clogging:</p> <ul style="list-style-type: none"> • Filter frits (inline filter, column filter)
Characteristic	●
Identification	<ul style="list-style-type: none"> • Start by disconnecting the capillary at the column inlet • Install a test setup with restriction capillary • Continue disconnecting capillaries, one-by-one, moving back toward the pump
Possible root cause	<ul style="list-style-type: none"> • Debris from mechanically worn parts (needle seat material, rotor seal at injection valve) • Coring of vial septa material
Instant action/first aid	<ul style="list-style-type: none"> • Backflush affected part • Replace part
Preventive measures	<ul style="list-style-type: none"> • Replace wear parts in time; apply proper preventive maintenance schedules • Use high-quality septa • Install inline filters



Blockages: instant pressure increase step



Clogging: constant pressure increase over time

Microbial Growth

- Potential problems
 - **Increased system pressure or pressure fluctuations**
 - **Increased column pressure, premature column failure**
 - Can mimic application problems
 - Gradient inaccuracies
 - Ghost peaks
- Prevent and/or Reduce Microbial Growth
 - Use freshly prepared mobile phase
 - Filter
 - Do not leave mobile phase in instrument for days without flow
 - Always discard “old” mobile phase
 - Do not add fresh mobile phase to old. No “topping off”
 - Use an amber solvent bottle for aqueous mobile phase
 - If possible, can add:
 - 5% organic added to water can be used to reduce bacterial growth
 - Use a few mg/L sodium azide



PN 3150-0577
Solvent filter/degasser assembly



Glass solvent inlet filter (20 mm), PN 5041-2168
Stainless Steel solvent inlet filter. PN 01018-60028
Amber solvent bottle 1L, PN 9301-6526

Why Filter the LC Sample


Capillaries, frits, and the column inlet are less likely to end up with blockages

Less wear and tear of injection and switching valves

Less downtime

Agilent Syringe Filter Selector tool

[Captiva Syringe Filter Selector | Agilent](#)



What is the Particle Size of Your LC Column?




Applications

Type of Filtration	Recommended	Alternatives
HPLC • UHPLC • LC/MS • GC	RC	PTFE or Nylon
ICP-MS	PTFE	Glass Fiber/PTFE (High Particle Samples)
CE	RC	Nylon
Undiluted Organic Solvents	PTFE	Nylon
Protein Analysis • Samples with Biomolecules – Buffers	PES	RC or CA
Tissue Culture Media	PES	RC or CA
High Particle-Load Samples – Organic Solvents	Glass Fiber/PTFE	-
High Particle-Load Samples – Aqueous Solutions	Glass Fiber/Nylon	-


Filtration

Captiva premium syringe filters

Column lifetime test



Impact of filtering a 0.3 µm latex-bead suspension on lifetime of a sub-2 µm column.



Filter cleanliness comparison of the Agilent Captiva Premium PES syringe filter with non-Agilent PVDF syringe filters using LC/MS under positive mode.

Captiva syringe filters guide: [5991-1230EN](#)

Changes in System Pressure


Low pressure

Potential Cause	Recommended Action
● Leak in high-pressure flow path	<ul style="list-style-type: none">• Visual inspection of the flow path• Instrument diagnostic tests ^{LA}
● Wrong mobile phase	<ul style="list-style-type: none">• Check for correct mobile phase• Check solvent reservoir and tube connections

Helpful Troubleshooting Tool





LA: With its advanced diagnostic and maintenance capabilities, **Agilent Lab Advisor SW** helps you to keep your Agilent analytical instruments in top condition. Agilent Lab Advisor is independent of the chromatography software you are using.



Characteristics

Parts affected	<ul style="list-style-type: none"> Potentially all parts in the flow path High potential at frequently operated fitting connections (for example, column inlet) and parts with high mechanical stress (rotor seal, needle, and needle seat)
Characteristic	<ul style="list-style-type: none"> Lower pressure Potentially impacting retention times and peak shape
Identification	<ul style="list-style-type: none"> Drops of solvent or residues of salt System diagnostic tests ^{LA}
Possible root cause	<ul style="list-style-type: none"> Loose or bad fitting connections Cracked capillaries Worn needle and needle seat
Instant action/first aid	<ul style="list-style-type: none"> Replace affected parts Renew or redo fitting connections
Preventive measures	<ul style="list-style-type: none"> Use proper fitting connections Replace fittings and wear parts in time




Changes in System Pressure

Pressure fluctuations

Potential Cause	Recommended Action
● Air in the system	<ul style="list-style-type: none">• Prime and flush instrument• Check for sufficient solvent supply• Check for correct plumbing (SSV/MCGV)• Check for correct degassing
● Malfunctions at pump head	<ul style="list-style-type: none">• Perform pump head diagnostic tests LA• Replace defective parts• Implement proper maintenance schedule
● Cavitation effects	<ul style="list-style-type: none">• Check for flow restrictions (solvent bottle to pump head inlet)• Clean or replace parts• Verify that solvent supply is positioned above the pump inlet

In addition

Pressure fluctuations will typically also impact the UV-signal due to refractive index effects.



Peak Shape Changes

Changes in Peak Shape

What is typically seen

Peak tailing

Peak splitting/doublets

Peak fronting



Peak broadening

Changes in Peak Shape

Peak tailing

If applicable to some peaks	Recommended Action
Secondary interactions	<ul style="list-style-type: none"> • Change pH • Change stationary phase
Small peak eluting on tail of larger peak	<ul style="list-style-type: none"> • Change selectivity (column, mobile phase) • Switch to methods with higher resolution (UHPLC, 2D-LC)

If applicable to all peaks	Recommended Action
Silica based – column degradation	<ul style="list-style-type: none"> • Use specialty, polymeric or sterically protected column
Silica based – basic interactions with stationary phase	<ul style="list-style-type: none"> • Use stronger mobile phase or add appropriate base (e.g. TEA)
Poor tubing connections; high dispersion volume	<ul style="list-style-type: none"> • Minimize number of connections • check connections / fitting condition and proper seat of fittings • use fittings with spring-load function



InfinityLab Quick Connect and Quick Turn Fittings


- Spring-loaded design
- Easy-to-use
- Works for all column types
- Reusable
- Consistent ZDV connection

Quick Connect Fitting


- Finger tight up to 1300 bar
- Hand tighten the nut, then depress the lever

Quick Turn fitting

- Finger tight up to 400 bar
- Up to 1300 bar with a wrench
- Compact design




Capillary for Quick Connect fitting



InfinityLab Quick Turn fitting

Capillary for Quick Turn fitting



Brochure: 5991-5164EN


Peak Tailing

Injector seal failure

Column: Bonus-RP, 4.6 x 75 mm, 3.5 mm
Temperature: R.T. Detection: UV 254 nm


Mobile phase: 30% H₂O : 70% MeOH Flow rate: 1.0 mL/min
Sample: 1. Uracil 2. Phenol 3. N,N-Dimethylaniline

<u>Plates</u>	<u>USP TF (5%)</u>
1. 2235	1.72
2. 3491	1.48
3. 5432	1.15



Before

<u>Plates</u>	<u>USP TF (5%)</u>
1. 3670	1.45
2. 10457	1.09
3. 10085	1.00




After replacing rotor seal
and isolation seal

Changes in Peak Shape

Peak splitting/doubling


If Applicable to Some Peaks	Recommended Action
Partially plugged column frit	<ul style="list-style-type: none">• Backflush the column (if applicable)• Use an inline filter• Use a guard column
Column void	<ul style="list-style-type: none">• Use a guard column• Use less aggressive mobile phase conditions• Replace the column
Sample volume overload	<ul style="list-style-type: none">• Use a smaller injection volume
Sample solvent incompatibility with mobile phase	<ul style="list-style-type: none">• Use mobile phase or weaker miscible solvent as injection solvent
Issues with injection valve	<ul style="list-style-type: none">• Check injector valve parts• Replace worn parts




Split Peaks from Injection Solvent Effects

Column: Agilent StableBond SB-C8, 4.6 x 150 mm, 5 mm Mobile phase: 82% H₂O:18% ACN
Injection volume: 30 mL Sample: 1. Caffeine 2. Salicylamide

A. Injection solvent
100% acetonitrile



B. Injection solvent
mobile phase




Injecting in a solvent stronger than the mobile phase can cause peak shape problems, such as peak splitting or broadening.

Keep organic concentration in sample solvent \leq mobile phase


Strong Sample Solvent Can Compromise Peak Shape

Isocratic method



Strong Sample Solvent Can Compromise Peak Shape




Gradient analysis



Changes in Peak Shape

Fronting


Potential Cause	Recommended Action
Channeling in column	<ul style="list-style-type: none">• Replace the column• Use guard columns
Column overload	<ul style="list-style-type: none">• Use a higher capacity column (increase length, diameter or change to high-capacity material)• Decrease sample amount



Changes in Peak Shape


Peak broadening

Potential Cause	Recommended Action
Injection volume is too large	<ul style="list-style-type: none">Decrease the injection volume
Long retention times	<ul style="list-style-type: none">Use gradient elution or stronger mobile phase
System settings	<ul style="list-style-type: none">Check data collection rate:Adjust the detector setting or time constant to the fastest possible value without compromising signal-to-noise.
Viscosity of the mobile phase is too high	<ul style="list-style-type: none">Increase the column temperature
Detector cell volume is too large	<ul style="list-style-type: none">Use the smallest possible cell volume
Improper fittings/connections	<ul style="list-style-type: none">Ensure that your fitting connections are correct
Extra tubing volume on system	<ul style="list-style-type: none">Ensure that the tubing is narrow and as short as possible to avoid extra volume.
Sample diluent too strong	<ul style="list-style-type: none">Reduce diluent strength




Changes in Peak Shape

Influence of data rate on appearance



Peak Shape

Extracolumn dispersion (Volume)




Increasing extracolumn volume

- Use short, small internal diameter tubing between the injector and the column and between the column and the detector.
- Make certain all tubing connections are made with matched fittings.
- Use a low-volume detector cell.
- Inject small sample volumes.



Length	10mm	50mm	100mm	150mm
Tubing ID	Volume	Volume	Volume	Volume
0.17mm (green)	0.227 uL	1.1uL	2.27 uL	3.3 uL
0.12mm (red)	0.113 uL	0.55uL	1.13 uL	1.65 uL

Peak Shape

Negative peaks



Normal



Negative

Causes

- Absorbance of the sample is less than the mobile phase
- Equilibrium disturbance when the sample solvent passes through the column
- Normal with Refractive Index Detectors
- Indirect UV detection

Retention

What is the Specific Issue?

- Retention times of all peaks shift
- Retention time of only one peak shifts


Increase the detail to be more specific:

- Retention time of all peaks shift
- Retention time of all peaks is shifting earlier
- Retention time of all peaks is shifting to earlier times and the extent of the shift appears to be the same


Changes in Separation

Retention time changing

Potential Cause	Recommended Action
Flow rate changing	Check "Pressure fluctuation", pump flow rate
Inconsistent online mobile phase mixing	Ensure gradient system is delivering constant composition check vs. manual preparation of mobile phase
Column temperature varying	Thermostat column and ensure constant lab temperature
Equilibration time insufficient with gradient run or change in isocratic mobile phase	Flush with at least 10 column volumes after solvent change or gradient conclusion
Selective evaporation of mobile phase component	Keep solvent reservoirs covered Prepare fresh mobile phase
Buffer capacity insufficient	Use >20 mM concentration of buffer
Contamination buildup	Occasionally flush the column with strong solvent to remove contaminants
First few injections – adsorption on active sites	Condition the column by initial injection of a concentrated sample
Column overloaded with sample	Decrease injection volume or concentration
Active sites on silica packing	Add competing base to mobile phase
Mobile phase composition changing	Follow 'best practices'



Change in Volatile Buffer Concentration and Shift in Retention Time and Peak Shape




Tip: The definition of “volatile” is “evaporating rapidly” or “passing off rapidly in the form of vapor”

Changes in Separation

Ghost peaks, carry over

	Potential Cause	Recommended Action
●	Peaks from previous injections	<ul style="list-style-type: none"> Flush the column to remove contaminants Check with blank injection
●	Specific interaction with metal surfaces	<ul style="list-style-type: none"> Passivate instrument Use InfinityLab Deactivator Additive Use bio-inert LC equipment
●	Contamination or unknown interferences in samples	<ul style="list-style-type: none"> Proper sample cleanup
●	Ion pair – disequilibrium	<ul style="list-style-type: none"> Prepare sample in actual mobile phase to minimize disturbance


**BIO
INERT**



[PN 5191-4506 | Deactivator Additive 50 ml](#)
[PN 5191-3940 | Deactivator Additive 25 ml](#)

Solvent Contamination

Injections on Agilent 1100



Solvent: acetonitrile

- Solvent source 1
- Solvent source 2
- Solvent source 3
- Solvent source 4, lot 1
- Solvent source 4, lot 2

Column Cleaning

Do what's recommended for your column



Tips for cleaning columns:

- Flush with stronger solvents than your mobile phase
- Make sure the detector is taken out of the flow path
- Do not add your organic solvent directly to the buffer, as this may cause the buffer salts to precipitate out and lead to more backpressure

For reversed phase:

Use at least 10 column volumes of each solvent for analytical columns

1. Start with your mobile phase without buffer salts (water/organic)
2. 100% organic (MeOH or ACN)
3. Check the pressure to see if it has returned to normal; if not, then
4. Discard the column or consider more drastic conditions: 75% acetonitrile/25% isopropanol
5. 100% isopropanol
6. 100% methylene chloride, solvent wash for very nonpolar compounds
7. Hexane

*Always see your specific column user guide for instructions



[LC Column User Guides | Agilent](#)

Baselines and Detection

Changes in Detection

Noisy baseline



	Potential Cause	Recommended Action
●	Gas bubbles in the mobile phase	<ul style="list-style-type: none">• Apply degassing• Check degasser performance
●	Low difference between sample and mobile phase absorbance	<ul style="list-style-type: none">• Check absorbance values of sample vs. mobile phase
●	Contamination	<ul style="list-style-type: none">• Use degassed HPLC-grade solvents• Flush system• Clean up sample
●	Detector optics	<ul style="list-style-type: none">• Perform intensity test• Check signal with flow cell removed if possible• Replace lamp
	Pressure instability	<ul style="list-style-type: none">• Check “Pressure fluctuation”



Changes in Detection


Drifting baseline

Potential Cause	Recommended Action
● Contamination in the mobile phase	<ul style="list-style-type: none">• Make up new mobile phase• If running a gradient, might need to adjust modifier
● Low difference between sample and mobile phase absorbance	<ul style="list-style-type: none">• Check absorbance values of sample vs. mobile phase
● Contamination	<ul style="list-style-type: none">• Use degassed HPLC-grade solvents• Flush system• Clean up sample
● Detector	<ul style="list-style-type: none">• Check temperature stability• Check for leaks• Replace lamp
Pressure instability	<ul style="list-style-type: none">• Check "Pressure fluctuation"




Reference Wavelength

How settings can affect the baseline



Gradient: 0.02 m KH₂PO₄/ACN, from 12% ACN to 45% ACN in 12 min

Influence of Data Collection Rate on Noise



Column Considerations

LC Columns Are Not Indestructible

- Columns are packed using hydraulic pressure and can be damaged by it.
- Silica dissolves (slowly) at higher pH
- Acid hydrolysis of bonded phase can occur at low pH
- Column failure
 - Void
 - Contamination


Columns must be stored properly

- Check your user guide

Important to Do:

- Know the technical specifications for your column
- Choose a mobile phase that is right for your column
- Keep record/history of your column

[LC Column User Guides | Agilent](#)



Columns:	4.6 x 150 mm, 5 µm
Purge:	50% ACN / 50% 0.02 M K ₂ HPO ₄ , pH 11
Flow Rate:	1.5 mL / min
Temperature:	25°C
Detection:	Silicate concentration by silicomolybdate color reaction

Choice of Your Column

Low and high pH can cause column failure

The InfinityLab Poroshell 120 portfolio offers choices for low and high pH

Best all around	Best for low pH mobile phases	Best for high pH mobile phases	Best for alternative selectivity	Best for more polar analytes	Chiral
EC-C18 1.9 µm, 2.7 µm, 4 µm	SB-C18 1.9 µm, 2.7 µm, 4 µm	HPH-C18 1.9 µm, 2.7 µm, 4 µm	Bonus-RP 2.7 µm	SB-Aq 1.9 µm, 2.7 µm, 4 µm AQ-C18 2.7 µm	Chiral-V 2.7 µm
EC-C8 1.9 µm, 2.7 µm, 4 µm	SB-C8 2.7 µm	HPH-C8 2.7 µm, 4 µm	PFP 1.9 µm, 2.7 µm, 4 µm	EC-CN 2.7 µm	Chiral-T 2.7 µm
Phenyl-Hexyl 1.9 µm, 2.7 µm, 4 µm		CS-C18 2.7 µm ← →		HILIC 1.9 µm, 2.7 µm, 4 µm	Chiral- CD 2.7 µm
				HILIC-Z 1.9 µm, 2.7 µm, 4 µm	Chiral-CF 2.7 µm
				HILIC- OH5 2.7 µm	


Every new column should be tested on your instrument

Performance verification based on Agilent checkout

- Run Agilent checkout before use
 - Record difference between your instrument and performance report (use as base value)
- Perform again if column seems to lose performance
 - Compare with results from first run

Performance verification based on in-house checkout

- Run in house checkout before use
 - Record key specifications such as tailing factor, plates, and backpressure
- Perform again if column seems to lose performance
 - Compare with results from first run



LC Troubleshooting Poster Available

LC Troubleshooting Guide

Your guide to solving common problems and staying productive

Agilent
InfinityLab

Places to Start

Solvents

- Use brown borosilicate bottles to avoid algae growth
- Prepare solvent volume to be used up within 1 to 2 days
- Use only HPLC-grade solvents filtered through 0.2 µm filters

Preparing and powering up the pump

- Inspect solvent bottles and inlet filters for damage or coloring
- Always use seal wash when installed and purge the pump
- Use the appropriate system conditioning method

Daily tasks

- Replace aqueous and organic mobile phases every second day
- Check seal wash solvent
- Flush the system with the composition of your application

Weekly tasks


- Change seal wash solvent and bottle and inspect solvent filters
- Check system backpressure and change filters if necessary

Pump shutdown

- Flush all channels to remove salt deposits and particulate matter
- Flush the system with appropriate storage solvent and power down the system

Handling of acetonitrile

- If possible, use 5 to 10% of water in your mobile phase
- Be sure to avoid ACN evaporation
- Don't leave ACN on the system for more than 2 to 3 days
- Perform a periodic warm water wash (60 to 70 °C) if you face problems



Maintenance

Agilent Lab Advisor software helps you manage your Agilent LC instruments to achieve high-quality chromatographic results in the most efficient way by ensuring high instrument performance, productivity, and reliability. It is available free-of-charge.

- Diagnostic tests to evaluate performance
- Easier maintenance of all Agilent LC modules
- Comprehensive reports generated to ease communication with Agilent service

Retention Time Drift



Possible Cause	Solution
Inconsistent online mobile phase mixing	Ensure gradient system delivers constant composition; compare with manual preparation of mobile phase

Possible Cause	Solution
Variation in column temperature	Thermostat or insulate column; ensure constant lab temperature

Possible Cause	Solution
Inefficient equilibration time with gradient run or change in isocratic mobile phase	Make sure at least 10 column volumes pass through column after sample run

Possible Cause	Solution
Selective evaporation of mobile phase component	Less vigorous helium sparging; keep solvent reservoirs covered; prepare fresh mobile phase

Possible Cause	Solution
Contamination buildup	Occasionally flush column with strong solvent

Possible Cause	Solution
Column overloaded with sample	Decrease injection volume or concentration

Pressure Fluctuation



Possible Cause	Solution
Leak in the system	Identify the channel and clean or replace check valve; replace pump seals

Possible Cause	Solution
Buildup of particulates	Filter sample and mobile phase

Possible Cause	Solution
Bubble in pump	Perform solvent degassing; sparge solvent with helium

Pressure Increase



Possible Cause	Solution
System blockage	Check flowpath (needle seat, capillaries, filter and frits)

Possible Cause	Solution
Water/organic systems: buffer precipitation	Test buffer/organic mixtures to ensure compatibility

High Column Backpressure



Possible Cause	Solution
Column blockage	Better sample cleanup; use guard column

Possible Cause	Solution
Mobile phase viscosity too high	Use lower viscosity solvents or higher temperature

Possible Cause	Solution
Particle size too small	Use larger d_p packing

Possible Cause	Solution
Plugged inlet frit	Replace column

Drifting Baseline




Possible Cause	Solution
Positive/negative direction: contaminant buildup/dilution	Flush column; clean up sample; use pure solvents

Possible Cause	Solution
Positive/negative: difference in refractive index of injection solvent	Use mobile phase for sample solvent

Possible Cause	Solution
Temperature changes	Insulate and thermostat column and tubing

Noisy Baseline



Possible Cause	Solution
Contamination	Use degassed HPLC-grade solvents; flush system; clean up sample

Possible Cause	Solution
Detector problems	Check number of hours of UV lamp; replace UV lamp or flow cell

Ghost Peaks



Possible Cause	Solution
Peaks from previous injection	Flush column to remove contaminants; check with blank injection

Possible Cause	Solution
Contamination; unknown interference in samples	Proper sample cleanup; Prepare sample in actual mobile phase to minimize disturbance

Possible Cause	Solution
Contaminated mobile phase	Check your mobile phase

Possible Cause	Solution
Bubbles in solvent	Check and degas your solvents

Peak Tailing



Possible Cause	Solution
Unwiped dead volumes	Minimize number of connections; ensure injector seal is tight; ensure fittings are properly seated

Possible Cause	Solution
Column performance	Change mobile phase; replace column

Possible Cause	Solution
Silica-based column degradation	Use specialty, polymeric, or sterically protected column

Possible Cause	Solution
Silica-based: basic interactions with stationary phase	Use stronger mobile phase or add appropriate base (e.g., TEA)

Peak Broadening



Possible Cause	Solution
Injection volume too large	Decrease injection volume or solvent strength of injection solvent; use gradient method

Possible Cause	Solution
Low sampling rate of data system	Increase data rate

Possible Cause	Solution
Detector cell volume too large	Use smallest possible cell volume

Possible Cause	Solution
Injection volume too large	Decrease injection volume

Sensitivity Problems



Possible Cause	Solution
Peaks are outside of sensitivity range of detector	Dilute/concentrate sample to bring into linear region

Possible Cause	Solution
Sample-related losses during preparation	Use internal standard during sample preparation; optimize sample preparation method

Leaks



Possible Cause	Solution
White powder at fitting/ loose fitting	Tighten fittings; replace capillaries

Possible Cause	Solution
System leak	Identify location checking leak sensors/sensors; check flow cell

Discover more best practices for using an Agilent LC system:
<https://www.agilent.com/chem/lc-best-practices>



Training courses are available at:
<https://www.agilent.com/crosslab/university>



Get answers. Share insights. Join the Agilent Community at:
<https://community.agilent.com>



For Lab Advisor software, please visit:
<https://www.agilent.com/chem/lab-advisor>



Request yours today at

www.agilent.com/chem/troubleshootLC

Contact Agilent Chemistries and Supplies Technical Support

Available in the U.S. and Canada, 8-5 all time zones

1-800-227-9770 option 3, option 3:

Option 1 for GC and GC/MS columns and supplies

Option 2 for LC and LC/MS columns and supplies

Option 3 for sample preparation, filtration, and QuEChERS

Option 4 for spectroscopy supplies

Option 5 for chemical standards

Option 6 for Prozyme products



gc-column-support@agilent.com

lc-column-support@agilent.com

spp-support@agilent.com

spectro-supplies-support@agilent.com

chem-standards-support@agilent.com

pzi.info@agilent.com



Agilent

Trusted Answers