

Analysis of Odors in Surface Water with SPME-GC-MS

Application Note

Authors

Johan Kuipers
Agilent Technologies, Inc.

Annegret Tischendorf
Zentrallabor der Südsachsen Wasser
GmbH

Introduction

The presence of algae, especially green-blue algae, together with some bacteria in surface water, can cause problems for drinking water supplies. This is due to the release of certain organic chemicals from algal or bacteria cells occurring either during normal metabolic processes or following cell rupture after the organism dies. These chemicals may cause off-odors and affect taste with implications for the drinking water treatment. Their development is dependent on temperature, sunlight and the age of the algae. The occurrence of trichlorophenol in the water is a result of the use of chlorine for disinfection purposes.

The analysis of these relatively volatile compounds in surface water can be performed by Headspace-SPME/GC/MS using a polyethylene glycol liquid phase column, such as the Agilent J&W FactorFour VF-WAXms, for optimal polar selectivity.

Components	RT (min)	Mode	Detection Limit (ng/L)
Camphen	8.17	SIM	5
Limonen	11.30	FS	100
Cineol	11.60	FS	100
Isopropyl-methoxypyrazin	17.61	SIM	2
Isobutyl-methoxypyrazin	17.73	SIM	2
Campher	17.74	FS	100
2-Methyl-isoborneol	18.72	SIM	6
Beta-cyclocitral	19.16	SIM	10
Geranylacetone	21.48	FS	100
2,4,6-Trichloranisole	21.57	SIM	1
Geosmin	21.84	SIM	3
Beta-ionone	22.88	SIM	5




Agilent Technologies

Conditions

Column: VF-WAXms, 60 m x 0.25 mm x 0.25 μ m (part number CP9207) Detection: MS-Quadrupol, full scan + SIM mode, TL=250 $^{\circ}$ C
Sample Size: 10 mL + 3 g NaCl in 22 mL Extraction: 60 min rapid stirring, 60 $^{\circ}$ C
 Headspace-Vial Headspace
Carrier Gas: Helium, 1.5 mL/min, constant flow Fiber: 2 cm StableFlex 50/30 μ m DVB/
Injector: SPME-Faser, splitless 3 min, 250 $^{\circ}$ C Carboxen/PDMS
Temperature: 60 $^{\circ}$ C/2 min - 4 $^{\circ}$ C/min-90 $^{\circ}$ C - 10 $^{\circ}$ C/min - 250 $^{\circ}$ C

RT: 7.50 - 24.00



www.agilent.com/chem

This information is subject to change without notice.
© Agilent Technologies, Inc. 2010
Published in UK, October 08, 2010
SI-01548

Analysis of odors in surface water using SPME-GC-MS
(data courtesy of *Zentrallabor der Südsachsen Wasser GmbH, Chemnitz*)



Agilent Technologies