A new instrument design and methodology for trace-level headspace analysis of polycyclic aromatic hydrocarbons is described.

 Using a 1 mL sample of water, quantitation from 5 ppt – 10 ppb is demonstrated for a wide range of polycyclic aromatic hydrocarbons from Naphthalene (BP 218°C) to Benzo(g,h,i) pyrene (BP 550°C).

Recovery of a 10 ppb PAH mixture from 10 mL water

Introduction

- Since its conception in 1989 (Belardi & Pawliszyn, 1989), solid phase headspace (HS) sampling coupled to thermal desorption (TD) gas chromatography-mass spectrometry (GC-MS) has become a powerful technique for both targeted and untargeted chemical analysis of volatile and semi volatile organic compounds (VOCs & SVOCs).
- We recently developed a HS sampling approach, termed vacuum assisted sorbent extraction (VASE), in which the sample i positioned in a disposable sample vial and placed in a reduced pressure environment in the presence of a packed headspace sorbent pen




- Following a controlled, diffusive headspace extraction with optional heat and agitation, the HSP containing adsorbed analytes is thermally desorbed directly onto the GC column via an inline sorbent pen thermal desorption unit (SPDU).
- Here we present a configuration and operating mode of the SPDU specifically designed for tracelevel analysis of SVOCs, and we apply this operating mode in conjunction with VASE to analyze polycyclic aromatic hydrocarbons (PAHs) in water.
- In comparison to traditional operating modes of the SPDU, this dual-column setup enables simultaneous split-mode analysis of VOCs (which are unretained on column 1) and spitless analysis
- Using 1 mL of sample volume, quantitation from 5 ppt 10 ppb is demonstrated for a wide range of polycyclic aromatic hydrocarbons from Naphthalene (BP 218°C) to Benzo(g,h,i)pyrene (BP 550°C).

Method Development: Adjusting the Split/Splitless Cutoff


Variation in Desorption Time & Benzo(g,h,i)Perylene Anthracene (m/z 178)

Variation in Desorption Flow



Headspace Sorbent Pen™

was used.



sample size is less than 1 mL. Schematic of the Headspace Sorbent Pen (HSP), showing the Micro-QT vacuum seal and the internal sorbent bed. The HSP can be packed with a variety of sorbents, including multi-component beds of varying physical properties. For the purposes of this work, ~80 mg of Tenax TA (35/60 mesh)

3 Water Management bed and HSP.

are available however the typical

Step3:Thesample vials are placed onto a chilled block to remove water from the sorbent

The Vacuum Assisted Sorbent Extraction Process

2 Diffusive Headspace Extraction

Step 2: The samples are heated and agitated while under vacuum. Diffusive headspace extraction conditions promote adherence of analyte to the front of the sorbent bed, enabling flash desorption and narrow chromatographic peaks. Typical extractions require 1–24 h.

4 Storage




Step 4: Extracted Sorbent Pens are placed in isolation sleeves, where they can typically be stored for up to 1 week prior to TD-GC-MS.

Thermal Desorption GC-MS

Methods

Column 1

Column 2



Sorbent Pen Thermal Desorption (SPDU) and Sample Preparation Rail (SPR) positioned on top of a Thermo GC-MS. The SPR transfers Sorbent Pens between the SPDU and air-tight isolation sleeves. The SPDU is connected directly to a capillary GC column held within the GC oven, which helps to minimize the flow path and mitigate analyte loss and/or carryover. Important optimization parameters include flow rate, preheat temperature & duration, desorption temperature & duration, and bakeout temperature

Split-Mode TD for VOCs & SVOCs

For analysis of the widest range of compounds (i.e., VOCs & **SVOCs**), the SPDU is arranged using an inert Silonite coated metal tubing (column 1) as an expansion loop and a capillary column (column 2) for analytical separation. The system is operated in split mode for adequate flow of desorption gas through the sorbent bed.

> During preheat, the bypass and split 2 valves (V2 & V4) are open while the desorb and split 1 valves (V1 & V3) are closed to allow heating without flow.

During desorption, the bypass valve is closed

and the desorb valve is opened so that the

desorb gas can flow through the HSP and onto column 1 and column 2. During bakeout, the split 1 valve is opened,

enable effective flushing.


GC-MS: Thermo TRACE 1310 ISQ (50-300 amu), 15:1 Split, Entech inert expansion loop

(Column 1), Agilent DB-5MS (30m×0.25mm×0.5µm) (Analytical Column)

For trace level analysis of **SVOCs**, the SPDU is arranged using a

short megabore capillary precolumn (column 1) to retain analytes

Dual-Column Spitless TD for trace SVOCs



 During bakeout, the split 1 valve is opened and the **split 2** valve is closed to allow SVOCs to pass splitlessly onto the analytical column.

column 1. SVOCs are retained on column

and lighter compounds are split at the oven

Step 3: Bakeout


ByPass GC EPC Inj1 or Inj2 ByPass Desorb Split Control GC EPC Split Vent Valve V2 V3 Split 1/V4 Split 2

GC-MS: Thermo TRACE 1310 ISQ (50-300 amu), 15:1 Split, Quadrex UAC-1MS/HT (7m×0.53mm×0.25μm) (Column 1), Agilent DB-5MS (30m×0.25mm×0.5μm) (Analytical Column)

Results: Quantitative Recovery of Parts-Per-Trillion Level 6-Ring PAHs from 1 mL Water

TIC (10 ppb) **EIC** (10 ppb)


Quantitative Recovery of PAHs from Water Using VASE & Splitless-Mode TD-GC-MS




2.6 TIC (10 ppb) ਰੌ 0.65 $R^2 = 0.991$ $R^2 = 0.995$ $R^2 = 0.990$ 7.5E+08 5.0E+08 2.0E+08 2.5E+08 Concentration (ppb) 1.2E+09 Pyrene R² = 0.997 2.0E+09 Benzo(a)Pyrene R² = 0.998 Fluranthene 0.5 5 Concentration (ppb)

Quantitative Recovery of PAHs from Water Using VASE & Split-Mode TD-GC-MS

Reproducibility (5 ppb replicates, n = 5)



Carryover: HSP Reruns in Split and Splitless Modes



Conclusions

- A dual column split/splitless operation of the sorbent pen thermal desorption unit (SPDU) is
- The split/splitless operating mode enables lower limits of detection for SVOCs, with a cutoff point dependent on desorption time, desorption flow, and the physical properties of the pre-column.
- In combination with VASE, limits of detection as low as 5 ppt are demonstrate for PAHs extracted from 1 mL of water.
- With 5 technical replicates, RSDs of less than 10% are demonstrated for the full range of standards examined (BP 218°C to 550°C).

References

R.G. Belardi & J. Pawliszyn, "The application of chemically modified fused silica fibers in the extraction of organics from water matrix samples and their rapid transfer to capillary columns", Water Pollut. Res. J. Can. 24 (1989) 179-191

Acknowledgements

The authors express their gratitude to the entire Entech staff, including the excellent engineers, machinists, software developers, assemblers, and graphic designers.