

Application Note

No.2

Multi-Residue Analysis of Pesticides in Green Tea Using Caffeine Removal Pretreatment

Y. Okamura



Food

1. Introduction

Green tea is becoming a popular beverage worldwide. Table 1-1 through Table 1-3 show the survey results for worldwide green tea production, and import and export quantities. With about 4 million tons produced worldwide, which is about half that of coffee bean production, China boasts the greatest rate of green tea production, followed by India, Kenya and Sri Lanka.

Sri Lanka is the greatest exporter of green tea, followed by Kenya, India and China. Domestic consumption is very high in India and China, with about 80 % of the production consumed in those countries. Japan also is a high-producing country, but due to even higher consumption, Japan also imports an amount which is equivalent to 50 % of its own production level.

The 27 European Union (EU) countries are the greatest importers, followed by the Russian Federation and the United Kingdom, with very high consumption clearly occurring in Europe.

Due to recent concern regarding food safety among consumers, advances in analytical methods for detecting and quantifying pesticide residues now permit the inspection of many crops for the presence of residual pesticides. With an increasing number of pesticides becoming subject to inspection every year, mass spectrometers are the instrument of choice for conducting simultaneous analyses targeting multiple pesticide residues.

The multi-residue analysis of pesticides in teas has become common worldwide. Caffeine, which is typically present in large quantities, can interfere with detection and quantitation of pesticides and other tea constituents, and is also a source of contamination in analytical instrumentation.

The development of an analytical method for multi-residue analysis of pesticides in green tea by gas chromatography with mass spectrometry (GCMS) is reported in this Application Note. A novel technique was employed to easily and efficiently eliminate caffeine to avoid any adverse effect on pesticide recoveries.

Table 1-1 Tea Production (2009)

Rank	Area	Production (tonnes)
1	China	1375780
2	India	972700
3	Kenya	314100
4	Sri Lanka	290000
5	Turkey	198601
6	Viet Nam	185700
7	Indonesia	146440
8	Japan	86000
9	Argentina	71715
10	Thailand	63707
11	Bangladesh	59500
12	Malawi	52559
13	Uganda	48663
14	Iran (Islamic Republic of)	165717
15	United Republic of Tanzania	32000
16	Myanmar	30500
17	Zimbabwe	20862
18	Rwanda	20000

Table 1-2 Tea Export (2009)

Rank	Area	Production (tonnes)
1	Sri Lanka	288528
2	Kenya	331594
3	China	305352
4	India	203863
5	EU(27)ex.int	29882
6	United Kingdom	27741
7	Germany	25301
8	Indonesia	92304
9	United Arab Emirates	23681
10	Viet Nam	82416
11	Malawi	47356
12	Belgium	7859
13	Argentina	69816
14	United Republic of Tanzania	30438
15	Russian Federation	9713
16	Netherlands	18158
17	Poland	8609
18	Uganda	44446

Table 1-3 Tea Import (2009)

Rank	Area	Production (tonnes)
1	EU(27)ex.int	249930
2	Russian Federation	182149
3	United Kingdom	145960
4	United States of America	110861
5	United Arab Emirates	75255
6	Egypt	80304
7	Pakistan	96932
8	Iran (Islamic Republic of)	51733
9	Japan	43301
10	Saudi Arabia	20331
11	Syrian Arab Republic	30651
12	Germany	44267
13	Canada	17353
14	France	17695
15	Poland	41784
16	Morocco	54400
17	Ukraine	26915
18	Netherlands	29982

Reference : Food and Agriculture Organization of The United Nations , FAOSTAT
<http://faostat.fao.org/default.aspx>

2. Maximum Pesticide Residue Levels in Tea and Analytical Method

The levels of pesticide residues in food are established in various countries around the world using Maximum Residue Levels (MRL) and Tolerance values. Methods of regulation vary depending on the country, but the applicable pesticides are recorded, their appropriate usage conditions are specified, and their MRL values are set for foods. Although these have been determined based on impact assessments on the human body, reference values are set in consideration of the various types of food products, as intake varies depending on the type of food. Reference values are set by the EU and Japan for many pesticides used in tea production. In the EU, MRL are specified for each agricultural product under Regulation (EC) No. 396/2005 Annexes. In Japan, the Ministry of Health, Labour and Welfare establishes MRL values for each product, and these can be viewed on the Ministry's home page. The analytical method follows the Official Method of the AOAC

INTERNATIONAL (AOAC). The United States Food and Drug Administration (FDA) publishes the Pesticide Analytical Manual (PAM), which specifies multi-residue methods as well as methods for individual pesticide compounds. The PAM multi-residue simultaneous analysis methods include methods for Non-fatty Foods and for Fatty Foods. Japan's Ministry of Health, Labour and Welfare categorizes the multi-residue analytical method based on whether the product is a cereal or a fruit, and a test method for tea is also indicated. Fig. 1 illustrates the analysis flow specified in the test method indicated by Japan's Ministry of Health, Labour and Welfare. In Japan, the system by which pesticides are regulated in foods is referred to as a "positive list system," which is referred to below as Japan's Positive List Test Method. In this investigation, 250 target pesticide compounds were analyzed following Japan's Positive List Test Method, and using the Shimadzu GCMS QP-2010 Plus shown below. Fig. 2 shows GCMS chromatogram of a 1 mg/L (ppm) standard mixture of the 250 pesticides. Analytical conditions are shown in Table 2.

Regulation (EC) No 396/2005

http://ec.europa.eu/food/plant/protection/pesticides/community_legislation_en.htm

Pesticide Analytical Manual (PAM)

<http://www.fda.gov/Food/ScienceResearch/LaboratoryMethods/PesticideAnalysisManualPAM/default.htm>

Japan's Ministry of Health, Labour and Welfare

<http://www.mhlw.go.jp/english/topics/foodsafety/positivelist060228/index.html>




Fig. 1 Flow Diagram of Japan's Positive List Test Method




Fig. 2 Total Ion Chromatogram (TIC) of 250 Pesticides Analyzed by GCMS

Table 2 Analytical Conditions

Instrument	: Shimadzu GCMS-QP2010 Plus	Oven program	: 50 °C (1 minute)
Inlet	: 1- μ L injection volume		: 25 °C/minute to 125 °C (0 minute)
	: High-pressure splitless mode (250 kPa, 1.5 minute)		: 10 °C/minute to 300 °C (10 minutes)
	: 250 °C	Interface	: 250 °C
Column	: Rtx-5MS, 30 mL. \times 0.25 mm I.D, df 0.25 μ m	MS operation	: Electron Impact (EI) ionization
	: Helium carrier gas, constant linear velocity (47.0 cm/second)		: Full scan mode, m/z 45-550

3. Analysis of Tea

Analysis of pesticide residues in commercially available green tea was conducted using the method conditions described above. After pretreatment according to Japan's Positive List Test Method, analysis by GCMS indicated that none of the 250 target pesticides were detected in the real-world sample. The Total Ion Chromatogram (TIC) obtained from analysis of the green tea extract is shown in Fig. 3. Ideally, method verification should be performed using a sample known to contain one or more of the target pesticides within the calibration range of the method. Since an actual tea sample contaminated with the target compounds was not available, a spike and recovery test was used to verify detection of pesticides and validate the method. A standard solution of pesticides was added to the green tea at a known concentration during the homogenization step. The extract was analyzed by GCMS and individual pesticide peaks were quantified against a calibration curve to verify recovery of the spiked pesticides.

The standard mixture of pesticides was added to the sample at a concentration of 100 µg/L (ppb). The recovery rate (%) obtained in this test was used to assess the test method. Table 3 shows the pesticides for which good recovery was obtained, at 70 % - 120 %.

Analytical interferences such as pigments, proteins, waxes, and other high molecular weight materials are co-extracted from the analytical sample along with the pesticides. Despite the absence of pesticide peaks in the chromatogram of Fig. 3, many peaks were detected. Because caffeine is present in large quantities in tea, its peak, which is seen to elute in the retention time range of 8 – 10 minutes, interferes with the detection of several pesticides in this portion of the chromatogram. Caffeine, which is present at high concentrations in coffee and various teas, including green tea and black tea, behaves much like the targeted pesticide compounds

during extraction and cleanup, and is therefore difficult to eliminate using the pretreatment process that was used here.

Depending on the type of solid phase cartridge used as an extract cleanup step, caffeine can be retained, thereby allowing its elimination from the sample solution. Fig. 4 shows an example of caffeine reduction through the use of a Florisil® column. However, due to the similar characteristics of caffeine and pesticides, they are both likely to remain in the cartridge with this processing. Thus, a separate step would be required to elute the pesticides from the cartridge, which would increase the pretreatment time. It would also require two separate GCMS analyses for each sample of green tea. Thus, finding a simplified procedure for caffeine removal was a major priority with respect to the analysis of pesticide residues in tea.



Shimadzu GCMS-QP2010 Plus Used in This Study




Fig. 3 GCMS Chromatogram of a Green Tea Extract




Fig. 4 Caffeine Reduction Due to Florisil® Column

Table 3 Green Tea Sample Spike and Recovery Test Results

Pesticide Name	Recovery (%)	Pesticide Name	Recovery (%)	Pesticide Name	Recovery (%)	Pesticide Name	Recovery (%)
EPTC	104	Vinclozolin	93	Fenamiphos	97	Thenylchlor	111
Mevinphos	86	Parathion-methyl	101	Flutolanil	114	Diflufenican	99
Etridiazole	102	Chlorpyrifos-methyl	98	Hexaconazole	113	Propargite	101
Chloroneb	97	Tolclofos-methyl	107	Imazalil	92	Piperonyl butoxide	102
XMC	106	Carbaryl	106	Isoprothiolane	105	Zoxamide	70
Fenobucarb	88	Alachlor	116	Profenofos	110	Mefenpyl-diethyl	101
Tecnazene	92	Heptachlor	75	Tribufos	114	Iprodione	114
Propoxur	93	Prometryn	100	Pretilachlor	109	Pyridaphenthion	110
Propachlor	100	Metalaxyl	86	Uniconazole P	102	Bifenthrin	105
Diphenylamine	103	Spiroxamin-2	89	p,p'-DDE	104	Bromopropylate	119
Ethoprophos	103	Terbutryn	114	Oxadiazon	113	Phosmet	93
Chloroprotham	103	Malathion	108	Dieldrin	86	EPN	92
Ethalfuralin	103	Thiobencarb	111	Oxyfluorfen	112	Tebufenpyrad	112
Trifluralin	113	Chlorpyrifos	109	Flamprop-methyl	82	BifenoX	112
Bendiocarb	92	Diethofencarb	110	Myclobutanil	98	Anilofos	91
Benfluralin	111	Aldrin	92	Buprofezin	84	Phenothrin-2	102
Cadusafos	100	Metolachlor	95	Imibenconazole-debenzyl	70	Tetradifon	102
alpha-BHC	99	Fenpropimorph	98	Flusilazole	95	Phosalone	93
Hexachlorobenzene	96	Fenthion	104	Thifluzamide	106	Pyriproxyfen	96
Dicloran	103	(Z)-Dimethylvinphos	100	Bupirimate	88	Cyhalothrin-1	104
Dimethoate	90	Parathion	120	Kresoxim-methyl	101	Cyhalofop-butyl	100
Carbofuran	78	Triadimefon	110	Isoxathion	117	Mefenacet	90
Atrazine	72	Isofenphos oxon	111	Cyproconazole	92	Cyhalothrin-2	120
Propazine	79	Chlorthal-dimethyl	109	Chlorfenapyr	115	Fenarimol	95
beta-BHC	96	Nitrothal-isopropyl	106	Fenoxanil	115	Pyrazophos	105
gamma-BHC	96	Bromophos	98	Chlorobenzilate	106	Pyraclofos	90
Propetamphos	101	Fthalide	105	beta-Endosulfan	91	Fenoxaprop-ethyl	77
Terbufos	101	Diphenamid	91	Fensulfothion	94	Bitertanol-1	89
Cyanophos	111	Fosthiazate-2	99	(Z)-Pyriminobac-methyl	120	trans-Permethrin	109
Quintozene	110	E-Chlorfenvinphos	106	p,p'-DDD	107	Bitertanol-2	101
Pyroquilon	84	Dimethametryn	103	o,p'-DDT	104	cis-Permethrin	106
Pyrimethanil	93	Penconazole	100	Mepronil	106	Pyridaben	97
Diazinone	100	Heptachlor epoxide (A)	93	Fluacrypyrim	115	Cyfluthrin-1	93
Phosphamidon-1	73	Oxy-Chlordane	96	Triazophos	91	Cafenstrole	83
Prohydrojasmon-1	103	(Z)-PyrifenoX	102	Benalaxyl	99	Fenbuconazole	88
Tefluthrin	102	Heptachlor epoxide (B)	84	Edifenphos	87	Halfenprox	112
delta-BHC	103	alpha-Chlorfenvinphos	117	Quinoxifen	74	Flucythrinate-1	104
Triallate	109	Diclocymet-1	110	Propiconazole-1	104	Flucythrinate-2	99
Iprobenfos	111	Quinalphos	100	Trifloxystrobin	115	Fenvalerate-1	91
Pirimicarb	90	Phenthoate	102	Norflurazon	88	Fluvalinate-1	114
Benoxacor	107	Zoxamide deg.	106	Lenacil	90	Fenvalerate-2	108
Benfuresate	111	Procyimidone	100	Endsulfan sulfate	91	Fluvalinate-2	98
Dichlofenthion	96	trans-Chlordane	105	p,p'-DDT	109	Difenoconazole-1	87
Propanil	89	Methidathion	70	Propiconazole-2	95	Difenoconazole-2	88
Bromobutide	108	Diclocymet-2	105	(E)-Pyriminobac-methyl	103	Flumiclorac-pentyl	96
Spiroxamin-1	82	(E)-PyrifenoX	92	Tebuconazole	97	Tolfenpyrad	93
Acetochlor	102	Tetrachlorvinphos	105	Diclofop-methyl	115	Imibenconazole	91

4. Investigation of Caffeine Removal Prior to Detection of Pesticide Residues in Tea by GCMS

The presence of a large amount of caffeine in tea not only interferes with the detection of pesticides by GCMS, it is a source of contamination of the injection port liner and the GC column. Furthermore, it can sometimes affect the analysis results, shifting the retention times of pesticides in the same chromatographic region as caffeine. The method which employs a solid phase cartridge, as shown in Fig. 4, is time-consuming, and increases the number of analyses. For this study caffeine was removed using a simple procedure that exploits the physical properties of caffeine.

4.1 Caffeine Removal Study

Fig. 5 shows the structural formula of caffeine. The high solubility of caffeine in polar solvents permits large amounts of caffeine to be dissolved in the extraction solvent, e.g. acetone. Moreover, as the temperature of the solution rises, even greater amounts of caffeine are dissolved.

Japan's Positive List Test Method for GCMS analysis of pesticides specifies an extraction solvent having a 1:1 ratio of acetone and hexane. Here we considered the physical property of polarity, with acetone as a polar solvent, and hexane a non-polar solvent. Because the polar solvent accounts for half of the ratio of the solution, this solvent mixture is thought to permit caffeine to dissolve easily. Therefore, it was decided to use hexane alone as the solvent, eliminating the use of acetone.

In addition, a lower solution temperature was considered to provide the benefit of impeding the dissolution of caffeine, and therefore the sample extracts were stored in a freezer. This freezing of the solution is referred to below as "freeze processing."

When the sample solution was freeze processed, deposits became suspended in solution (Fig. 6). By applying centrifugation, these deposits were precipitated, and analysis of the supernatant was equivalent to analysis without most of the caffeine. This process was applied to a green tea sample for confirmation of the effect.

To confirm the effect on pesticide recovery using hexane as the extraction solvent along with the application of freeze processing, pesticides were added to a green tea sample solution that was previously subjected to the described pretreatment, and then we applied the change in solvent (using only hexane) in addition to freeze processing. No adverse effect on recovery was observed.

C₈H₁₀N₄O₂
Mol. Wt.:194.19




Fig. 5 Structure of Caffeine




Fig. 6 Status at Vial Tip Before and After Centrifugation

4.2 Caffeine Removal Study Results

After preparing a caffeine-saturated hexane solution, the effect of freeze processing on caffeine removal was evaluated. The chromatograms of caffeine generated before and after processing are shown in Fig. 7, and the caffeine peak area ratio comparison is shown in Table 4. The caffeine content was reduced by 63.5 % as a result of freeze processing.




Fig. 7 Effect of Freeze Processing

Table 4 Comparison of Caffeine Area Ratios Before and After Freeze Processing

Before Freeze Processing	After Freeze Processing	Reduction Rate (%)
5209488	1906330	63.5

When freeze processing is conducted, deposits become suspended in the solution. Since centrifugal separation is known to effectively remove these suspended particles, centrifugation time was investigated to determine its effect on the separation. However, any rise in temperature that would occur during longer centrifugation would lessen the effect of freeze processing. Fig. 8 shows the results of the study of centrifugal separation time, allowing for precipitation of deposits following freeze processing. It took at least 1 minute to attain centrifugal separation, but as the centrifugal separation time increased beyond 1 minute, the caffeine area values increase accordingly, until they become constant after 3 minutes. Thus, at 3 minutes, it is thought that the effect of freeze processing would become counterproductive. From this result, we determined that the optimum time required for centrifugal separation is 1 minute.




Fig. 8 Result of Investigation of Centrifugal Separation Time

4.3 Applying Caffeine Removal Operation to an Actual Sample

The caffeine removal operations described in the above study were employed to remove caffeine from a real-world tea sample. After processing commercially available tea according to Japan's Positive List Test Method, the solvent that had been specified for use in the obtained final solution was replaced with hexane, followed by freeze processing (-20 °C) and centrifugal separation (1 minute). The obtained supernatant was then analyzed by GCMS. Fig. 9 shows a flow chart of Japan's Positive List Test Method + caffeine removal operations.

Fig. 10 shows the TIC chromatograms obtained before and after caffeine removal by freeze processing. Prior to the removal processing, the caffeine peak is detected as a broad peak that exceeds the column load capacity, but after caffeine removal, a sharp caffeine peak within

the column load range is seen, indicating that most of the caffeine was removed. As for the stability of the caffeine removal operation, the repeatability of caffeine area values following caffeine removal operations is shown in Table 5. Caffeine removal by freeze processing can thus be considered to be a stable process that provides good repeatability.

Caffeine removal by solvent replacement and freeze processing was thus confirmed, however actual target pesticides might be removed along with the caffeine. Therefore, spike and recovery testing was performed with respect to the above caffeine removal operation. A test solution prepared using Japan's Positive List Test Method was spiked with a standard mixture of pesticides, and after subjecting this solution to the above described caffeine removal processing, the pesticide recovery rates were obtained.




Fig. 9 Flow Diagram of Japan's Positive List Test Method + Caffeine Removal




Fig. 10 Effect of Caffeine Removal by Freeze Processing

Table 5 Repeatability of Caffeine Removal Effect

	1	2	3	Average Value	CV%
Caffeine Area Values	4963206	5010227	4996268	4989900	0.48

Fig. 11 shows the pesticide spike and recovery test results obtained after caffeine removal processing. Recovery for most compounds fell between 80 and 120 %, and it was confirmed that switching to a hexane solvent and use of freeze processing did not result in significant loss of pesticides.

Not only was the caffeine peak drastically reduced as a result of the caffeine removal processing, other contaminant substances were eliminated. Removal of these contaminants lessened the adverse effects

on the pesticide peak shapes. Fig. 12 shows examples of improved peak detection.

In the case of Mevinphos and Fosthiazate, the interfering peaks before these compounds were removed to allow clear detection of these pesticides. As for Propoxur, the m/z 152 peak shape was very different from the m/z 110 peak, but this improved as a result of freeze processing. Carbofuran, which co-eluted with another peak, became a single Carbofuran peak.




Fig. 11 Pesticide Spike and Recovery Test Results Due to Caffeine Removal Processing




Fig. 12 Examples of Effect of Caffeine Removal Processing (green tea 0.1 ppm)

5. Conclusion

Pesticide residue analysis in tea using Japan's Positive List Test Method for GCMS simultaneous analysis yields good results, but during pesticide analysis by GCMS, the large amount of caffeine that remains in these tea samples can contaminate the GC injection port and column. Similarly, the presence of caffeine can co-elute with and mask the presence of pesticide residues in the same region of the chromatogram. To eliminate the interfering effects of caffeine, the solvent was changed from acetone : hexane (1:1) to 100 % hexane followed by freeze processing, exploiting the physical properties of caffeine, thereby efficiently decreasing the presence of caffeine in the sample. Further, utilizing this removal process, a method was developed for removing caffeine while retaining good recovery and minimal loss of the target pesticides.

*This document is based on information valid at the time of publication. It may be changed without notice.

First Edition: November, 2012

