

Peptide Mapping

Hardware and Column Optimization



Peptide Mapping Hardware Challenges

Very complex samples require long, shallow gradients

Mobile phase components require mixing, however delay volume within the pump needs to be low to deliver accurate gradients


Overall system dispersion needs to be kept to a minimum to keep peaks separated

Compositional and flow accuracy and precision are needed for reproducibility of complex runs

Why is delay/dwell volume important

1. Different dwell volumes result in a RT time shift
(the time for the mobile phase to reach the column head)
2. Different dwell volume could effect resolution
(peaks spends different time under isocratic/gradient conditions)
3. additionally, the dwell volume effects the gradient shape
*(dispersion effects, flush out behavior
=> the programmed gradient becomes deteriorated)*
4. Therefore even with the same “geometrical” delay volume the chromatograms could look different on different systems
5. The dwell volume has an big impact for narrow bore applications, especially combined with fast gradient


Dwell volume determination



W. Dolan LCGC Vol 24, No 5, 458-466

The system setup is simple. Use water for the A-solvent and water spiked with 0.1% acetone for the B-solvent. Replace the column with ≈ 1 m of 0.005-in. i.d. tubing, set the detector to 265 nm, and set the flow rate so that there is sufficient backpressure for reliable check valve operation (for example, 2 mL/min). Run a 0–100% B gradient in 20 min. The data system output should be a curve similar to the one below. You can measure the dwell time (t_D) by drawing a tangent to the main part of the gradient curve (dashed line in Figure 4) and extend the baseline to intersect this tangent. The time it takes from the start of the program to this intersection is the dwell time. Multiply by the flow rate to get the dwell volume.

Delay Volume Profiles



Delay Volume Variability within Agilent Systems

1200 Infinity Series

Configuration	Delay V* (μ l)
1290 Pump	10
1290 Pump + Fixed Loop ¹ (for MS)	20
1290 Pump + Jet Weaver + Fixed ¹ Loop	55
1290 Pump + ALS (for MS)	75
1290 Pump + Jet Weaver ¹ + ALS	110
1200 RRLC (low delay volume)	260
1200 RRLC (standard delay volume)	740-940



Agilent 1290 Infinity Quaternary Pump

Specifications & Benefits

Power Range

- For any kind of analysis

Composition Accuracy and Precision


- < 0.15 % RSD or 0.02 min SD
- ± 0.4 % (1-99 % Composition B)

- High RT precision in gradient runs

Flow Accuracy and Precision

- < 0.07 % RSD or 0.01 min SD
- ± 1.0 % or 10 µL

- High RT precision in isocratic runs



Composition Range

1-99 %

- Wide analytical range


Delay Volume

< 350 µL

- For fast quaternary gradients


Composition Accuracy and Precision

Comparison with Quaternary UHPLC pump from other vendor




Composition Range

Shallow gradient at low organic concentration



Composition Range


Shallow gradient at low organic concentration




Precision & Accuracy

1st Quaternary pump with Bin-like performance

Different mixtures (H₂O, MeOH, ACN) at different pressures and different flow rates




1290 Quaternary Pump Flow Diagram




Multipurpose Valve

Optional Extra Mixing Volume for lowest baseline ripple (TFA applications)



TFA Mixing Noise: Effect of JetWeaver Mixer: Lower Short Term Noise



Innovation: 1290 Infinity Quaternary Pump

Quaternary Pump with Binary Pump-like Performance




1290 Infinity Quaternary Pump

- **Boost Performance**
Highest accuracy and precision for composition and flow, with exceptionally low delay volume
- **Save Time**
Accelerates transfer of existing HPLC methods to UHPLC
- **Reduce Costs**
Outstanding UHPLC performance at HPLC-like costs of ownership

1290 Infinity Quaternary Pump - BlendAssist

Simple tool for online-dilution of modifiers and gradient set-up

You need different concentrations of modifiers in your analysis, would like to have just one stock-solution and do online dilution to profite from the quaternary mixing capability of your pump? Here is a simple tool – **BlendAssist!**



Desired method conditions - example:


1. 5 to 95% gradient of ACN with 0.1% TFA in Water and 0.08% TFA in ACN
2. 20 – 80% gradient of ACN with 0.5% TFA in Water and 0.4% TFA in ACN

Without BlendAssist you need to either pre-mix the required solvents or by using stock-solutions of TFA in Water and ACN to program complex gradients (%A, B, C, D).

With BlendAssist: just program your binary organic/aqueous gradient and define the dilution factor!

1290 Infinity Quaternary Pump - BlendAssist

Simple tool for online-dilution of modifiers and gradient set-up

Quat. Pump (G4204A) 

Flow: 0.000 mL/min

Solvents

Enable Blend Assist

Solvent	Used	%	Name
Water/TFA (0.1%)	<input checked="" type="checkbox"/>	57.0	
ACN/TFA (0.01mM)	<input checked="" type="checkbox"/>	43.0	

Advanced

Timetable (empty)

Blend Assist

Channel	Type	Calibration	Name	Stock concentration	Final concentration	Conc. unit
A	Solvent 1	100.0 % Water V.02	Water	1.00	1.00	%
B	Solvent 1 Additive	100.0 % Water V.02	TFA	1.00	0.10	%
C	Solvent 2	100.0 % Acetonitrile V.02	ACN	1.00	1.00	%
D	Solvent 2 Additive	100.0 % Acetonitrile V.02	TFA	0.10	0.01	mM

Pressure Limits

Min: 0.00 bar Max: 1,200.00 bar

Stoptime

As Injector/No Limit
 1.00 min

Posttime

Off
 1.00 min

Peak Dispersion in HPLC

- Dispersion is the sample peak broadening or dilution which occurs in connecting tubing, sample valves, flow cells and in column end-fittings. It begins with the injector and ends at the **last detector** in the system
- As column internal diameter and length decrease the potential peak broadening in a non-optimized LC system increases.
- Higher efficiency in the column can only be realized if the system dispersion does not substantially degrade the column performance.
- As particle size decreases, resolution increases as a result of narrower peak widths.
- Narrower peaks are more susceptible to extra-column dispersion.
- The smaller the column dimension, too, the smaller the expected peak volume. Thus the small particle size columns used in low volume configurations require the greatest attention to plumbing in the LC system.

How to optimize your LC-system

Non-column sources of peak broadening:

- **General**

- Connecting tubing (internal diameter too big, tubing too long)
- Connectors (unions, tees, bulkhead fittings)
- Switching valves for automated SPE or alternating column regeneration

- **Sampler**

- Sample aspirating needle and loading/transfer port
- Sampler switching valve(s) contacting sample


- **Detection**

- Inlet heat exchangers, flow cell volume and geometry
- Incorrect data rate selection and data filtering effects in high speed applications


Tips for minimizing dispersion in LC systems

Minimize interconnection volume from the injector to detector, with minimal junctions and smallest i.d. and length tubing


Make fittings carefully, using appropriate connectors, and do not re-use seated stainless steel ferrules in different locations, including different ports on switching valves



LC System -- Variable Configurations for Dispersion Volume and Delay/Dwell Volume



Time-aligned overlay of HPLC dispersion tests



Dispersion Volume within Flow Cell

Simple example: dispersion in a straight tube

Dispersion:

$$\sigma^2 = \frac{\pi \cdot r^4 \cdot F \cdot L}{24 \cdot D_m}$$

Volumetric peak variance


F: Flow rate

D_m : Molecular diffusion coefficient

L: Capillary length

Derived from Aris, Taylor and Golay equation

Same geometrical volume ($V_1 = V_2$), but totally different dispersion volumes



R. Tijssen, "Mechanism and Importance of Zone-Spreading, in Handbook of HPLC, Vol 78, 1998, Marcel Dekker, NY


Why do we need a new cell technology?

Short 2.1 mm ID column




+

10 mm pathlength
13 μ l volume





=



How to achieve smaller cell volume?

short path-length (3 mm, 2 μ l)





Low S/N

High light transmission =>

low noise

long path-length (10 mm, 0.5 μ l)




Low S/N

Low light transmission => high noise

Solution:

Optofluidic Waveguides



- Long path length
- Small cell volume
- High light transmission



Highest S/N

Max-Light Cartridge Cell

- *Optofluidic waveguides*



High Light Transmission due to Total-Internal Reflection (TIR) principle
(~ 100 % Light efficiency)

Benefits:

- *Highest Sensitivity (S/N) with small cell volumes (dispersion effects)*
- *More reliable and robust peak integration (automated) due to nearly no Refractive Index and thermal effects (solvent temperature)*
- *Coating free fused silica (no special care instructions or smiling baseline effects)*
- *Easy cell selection (one cell for all major applications)*
- *Cartridge design for ease of use*

Materials Needed to Optimize a 1290 Infinity LC for Ultra-Low Dispersion

Ultra-Low Dispersion Kit for 1290 Infinity LC System

Includes:

0.075 x 220 mm capillary

0.075 x 340 mm capillary


1.0 uL heat exchanger

0.075 mm id needle seat


Ultra-Low Dispersion Max-Light Cartridge Flow Cell

1290 Infinity LC System Set-Up


Default 1290



Optimized 1290




Isocratic Analyses of Alkylphenones




A: H₂O; B: CH₃CN; 0.4 mL/min
 Isocratic, 60% B
 1 μ L injection of RRLC Checkout
 Sample (PN 5188-6529) spiked w/ 50
 μ L 2 mg/mL thiourea in
 water/acetonitrile
 TCC: 26 °C
 DAD: Sig = 254, 4 nm; Ref = Off
 Agilent ZORBAX RRHD Eclipse Plus
 C18, 2.1 mm x 50 mm, 1.8 μ m
 Sample:

1. Thiourea (v_0 marker)
2. Acetanilide
3. Acetophenone
4. Propiophenone
5. Butyrophenone
6. Benzophenone
7. Valerophenone
8. Hexanophenone
9. Heptanophenone
10. Octanophenone

% Improvement from Default to Optimized Agilent 1290 Infinity LC System with an Isocratic Analysis




Premixed versus Blend Assist; 0.05%TFA-Water/0.04%TFA ACN BSA Tryptic Digest Peptide Map



Duplicate Runs Premixed


0.05%TFA-Water/0.04%TFA ACN

BSA Tryptic Digest Peptide Map



Triplicate Runs Blend Assist: Stock 1%TFA; 1200 Infinity Series


0.05%TFA Water/0.04%TFA ACN BSA Tryptic Digest Peptide Map



Overlay of Blend Assist TFA percent Differences: 1200 Infinity Series

Blue: 0.05%TFA-Water/0.04%TFA-ACN
Red: 0.1%TFA-Water/0.08%TFA-ACN

BSA Tryptic Digest Peptide Map



OpenLAB CDS MatchCompare

Comparison tool for complex
Chromatograms

Problem / Solution

Problem

I need to compare a chromatogram against a reference and obtain **objective** data


Solution

OpenLAB CDS Match Compare Add-on Software

Matching Chromatograms

Reference defined

- Area percent tolerance defined by peak
- RT (RI) tolerance defined in time (Index)
- Initial shift allows for injection delays



Matching Chromatograms

Reference defined

- Area percent tolerance defined by peak
- RT (RI) tolerance defined in time (Index)
- Initial shift allows for injection delays

Parameters:

Temporal tolerance: [min]

Initial shift: [min]

Allow to change the shift sign

Filter small peaks Minimum area: [%]

Hide identical peaks

Summary - Area Example 40 ppm.merg

#	Components	Ref Time	Ref Area	40 ppm.cdf Time	Area
9		5.44	0.00	5.60	0.00
10		5.60	0.00		
11		5.72	0.01	5.73	0.02
12		6.20	0.01	6.19	0.01
13		8.75	0.03	8.75	0.03
14				8.88	0.01
15		11.30	0.00	11.30	0.00
16	analyte 4ds	12.66	0.02	12.66	0.02
17		12.97	0.00	12.98	0.01
18		13.31	0.01	13.33	0.01
19		14.23	0.01	14.24	0.01
20		15.72	0.01	15.74	0.03
21		15.77	0.00	15.78	0.01
22	analyte 32al	15.96	0.02	15.97	0.02
23		16.05	0.00	16.05	0.00
24		19.24	0.01	19.25	0.02
25	analyte 34ai	22.03	0.02	22.03	0.02
26		22.17	0.01	22.18	0.02
27	analyte 231	24.34	0.02	24.36	0.03
28		24.38	0.01	24.41	0.02
29		25.06	0.02	25.06	0.02
30		26.10	0.00	26.10	0.00
31		26.62	0.00	26.62	0.00

Agilent OpenLAB MatchCompare - Reference


Reference name: Area Example 40 ppm.ref Nb of peaks:30

A comparison will be successful if and only if:

- The sum of the missing or additional peaks area is under %
- and all errors are under %

Success criterion Tolerances

Time shifts



Peak results

Comparison Summary

The screenshot shows the 'Agilent OpenLAB MatchCompare - Comparison' window. The 'Results' section is highlighted with a red box and contains the following summary:

- 79.17 % Identical
- 12.50 % Out of tolerance
- 4.17 % Ref. only
- 4.17 % Samp. only

Additional text in the summary: 'stand for 0.02 % of total area. The 4.17 % of unknown peaks in sample stand for 0.01 % of total area.' The similarity score is 0.9711.

Parameters:

- Temporal tolerance: 0.100 [min]
- Initial shift: 0.170 [min]
- Allow to change the shift sign
- Filter small peaks Minimum area: 0.05 [%]
- Hide identical peaks


Buttons: Calculate, Stop

Name	Rt Samp	Rt Ref [min]	DT	% Samp	% Ref	% Error	Tol [%]	Info
	5.44	5.29	0.15	0.0000	0.0000	13.18	100.00	Id.
	5.60	5.44	0.15	0.0001	0.0000	140.73	100.00	Id. OT
	-----	5.60	-----	-----	0.0001	0.00	100.00	Ref.
analyte 32a	5.73	5.72	0.01	0.0211	0.0107	97.03	100.00	Id., Tail peak
	6.19	6.20	0.00	0.0120	0.0114	5.34	100.00	Id.
analyte 32	8.75	8.75	0.00	0.0304	0.0283	7.29	100.00	Id., Tail peak
	8.88	-----	-----	0.0075	-----	-----	-----	Samp.
	11.30	11.30	0.00	0.0001	0.0001	10.37	100.00	Id.
32a	12.66	12.66	0.00	0.0212	0.0201	5.57	100.00	Id., Tail peak
	12.98	12.97	0.01	0.0083	0.0041	103.85	100.00	Id. OT

Legend:

- Identical (Green)
- Identical out of tolerance (Cyan)
- Reference only (Grey)
- Sample only (White)

Area comparison



OpenLAB CDS Match Compare Reports

Comparison Report

Comparison between Area Example 40 ppm.ref and SIGNAL01.cdf Processing date: 4/25/2013 2:47:54 PM

Similarity: 0.9781
 87.50% Identical
 4.17% Out of tolerance
 4.17% In reference
 4.17% In sample

4.17 % of peaks out of tolerance represents 0.00 % of total area
 4.17 % of unknown peaks represents 0.00 % of total area

Validated by: _____ Signature: _____

Results table

Code	Name	Rt Samp (min)	Rt Ref (min)	DT	% Samp	% Ref	% Error	Tol (%)	Info
0.00	5.29	-----	-----	0.0000	-----	-----	-----	-----	Samp
0.00	5.44	5.29	5.29	0.15	0.0000	0.0000	25.72	100.00	id
0.00	5.60	5.44	5.44	0.16	0.0001	0.0000	163.96	100.00	id, OT
0.00	5.80	-----	-----	-----	-----	0.0001	0.00	100.00	Ref
0.00	5.72	5.72	0.00	0.0107	0.0107	0.00	100.00	id	Tail peak
0.00	6.20	6.20	0.00	0.0114	0.0114	0.00	100.00	id	Tail peak
0.00	8.75	8.75	0.00	0.0283	0.0283	0.00	100.00	id	Tail peak
0.00	11.30	11.30	0.00	0.0021	0.0021	0.00	100.00	id	
0.00	12.66	12.66	0.00	0.0021	0.0021	0.00	100.00	id	Tail peak
0.00	12.97	12.97	0.00	0.0041	0.0041	0.00	100.00	id	
0.00	13.31	13.31	0.00	0.0065	0.0065	0.00	100.00	id	
0.00	14.23	14.23	0.00	0.0062	0.0062	0.00	100.00	id	
0.00	15.72	15.72	0.00	0.0030	0.0030	0.00	100.00	id	Tail peak
0.00	15.77	15.77	0.00	0.0035	0.0035	0.00	100.00	id	
0.00	15.96	15.96	0.00	0.0024	0.0024	0.00	100.00	id	Tail peak
0.00	16.05	16.05	0.00	0.0001	0.0001	0.00	100.00	id	
0.00	19.24	19.24	0.00	0.0100	0.0100	0.00	100.00	id	Tail peak
0.00	22.03	22.03	0.00	0.0028	0.0028	0.00	100.00	id	Tail peak
0.00	22.17	22.17	0.00	0.0082	0.0082	0.00	100.00	id	Tail peak
0.00	24.34	24.34	0.00	0.0163	0.0163	0.00	100.00	id	Tail peak
0.00	24.38	24.38	0.00	0.0126	0.0126	0.00	100.00	id	Tail peak
0.00	25.06	25.06	0.00	0.0022	0.0022	0.00	100.00	id	Tail peak
0.00	26.10	26.10	0.00	0.0001	0.0001	0.00	100.00	id	
0.00	26.62	26.62	0.00	0.0003	0.0003	0.00	100.00	id	

Summary - Area Example 40 ppm.merg

File Edit Help

#	Area Example 40 ppm.ref			40 ppm.cdf	
	Components	Ref Time	Ref Area	Time	Area
9		5.44	0.00	5.60	0.00
10		5.60	0.00		
11		5.72	0.01	5.73	0.02
12		6.20	0.01	6.19	0.01
13		8.75	0.03	8.75	0.03
14				8.88	0.01
15		11.30	0.00	11.30	0.00
16	analyte 4ds	12.66	0.02	12.66	0.02
17		12.97	0.00	12.96	0.01
18		13.31	0.01	13.33	0.01
19		14.23	0.01	14.24	0.01
20		15.72	0.01	15.74	0.03
21		15.77	0.00	15.78	0.01
22	analyte 32al	15.96	0.02	15.97	0.02
23		16.05	0.00	16.05	0.00
24		19.24	0.01	19.25	0.02
25	analyte 34ai	22.03	0.02	22.03	0.02
26		22.17	0.01	22.18	0.02
27	analyte 231	24.34	0.02	24.36	0.03
28		24.38	0.01	24.41	0.02
29		25.06	0.02	25.06	0.02
30		26.10	0.00	26.10	0.00
31		26.62	0.00	26.62	0.00

Overview of the principals of Match Compare

- Comparison of two chromatograms
 - Pattern matching, not pattern recognition
 - Matching based on retention time (or retention index) and area percent
- Comparison parameters can be individually tailored
- Report can be on all peaks or just those falling outside the limits



Live Demonstration of MatchCompare Using Peptide Maps



Appendix

1200 Infinity
Series



Analyzing protein drug by RP-HPLC peptide map - pyroglutamate formation

Reverse-Phase Chromatography/ Mass Spectrometry Analysis of Reduced Monoclonal Antibodies in Pharmaceuticals

Douglas Rehder, Thomas Dillon, Gary Pipes, and Pavel Bondarenko, Journal of ChromatographyA, 1102 (2006), p.164-175

Analyzing protein drug by RP-HPLC Peptide map 1- amino acid substitution

Identification of a Glu > Lys substitution in the activation segment of human pepsinogen A-3 and -5 isozymogens by peptide mapping using endoproteinase Lys-C

Ruud A. Bank, Bart C. Crusius, Toon Zwieters, Stephan G.M. Meuwissen*, Fre Arwert and Jan C. Pronk Institute of Human Genetics and *Department of Gastroenterology, Free University, Amsterdam. The Netherlands, Volume 238, number 1, 105-108 FEB 06339 September 1988

Analyzing protein drug by RP-HPLC Peptide map-mirror image

LC/ESI-MS/MS analysis of recombinant IgG2 mAb after Lys-C digest

Wypych J et al. J. Biol. Chem. 2008;283:16194-16205

Analyzing protein drug by RP-HPLC Peptide map-deamidation

Quantification and characterization of antibody deamidation by peptide mapping with mass spectrometry

Weijie Wang, Andrea R. Meeler, Luke T. Bergerud, Mark Hesselberg, Michael Byrne, Zhuchun Wu, Analytical Sciences Department, Human Genome Sciences, Inc., 14200 Shady Grove Road, Rockville, MD 20850, United States

Analyzing protein drug by RP-HPLC Peptide map - disulfide bonds

Disulfide Linkage Analysis of IgG1 using an Agilent 1260 infinity Bio-inert LC System with an Agilent Zorbax RRHD Diphenyl sub-2um Column

M. Sundaram Palaniswamy, Agilent Technologies, Inc., Bangalore, India

Analyzing protein drug by RP-HPLC Peptide map - PEGylation

Toward Top-Down Determination of PEGylation Site Using MALDI In-Source Decay MS Analysis

Chul Yoo, Detlev Suckau, Volker Sauerland, Michael Ronk, and Minhui Maa, Analytical R&D, Amgen Inc., Thousand Oaks, California, USA

Questions

