

Understanding Synthetic Biology using the Q Exactive GC Orbitrap GC-MS and a High Resolution Accurate Mass Metabolomics Library for Untargeted Metabolomics

Cristian Cojocariu¹, Maria Vinaxia², Mark Dunstan², Adrian J. Jervis², Paul Silcock¹ and Nicholas J. W. Rattray²

¹Thermo Fisher Scientific, Runcorn, UK; ²SYNBIOCHEM – Manchester Synthetic Biology Research Centre for Fine and Specialty Chemicals, Manchester Institute of Biotechnology, University of Manchester, Manchester, UK

ABSTRACT

In this study, the application of untargeted metabolomics was used to understand the metabolic effects of inducer type and concentration on the metabolic fingerprint of engineered bacteria under different growth conditions. To track the metabolic shifts caused by promotor induction in the bacteria, a large number of metabolites from the cellular exo-metabolome were detected and identified using the Thermo Scientific™ Q Exactive™ GC Orbitrap™ mass spectrometer. The untargeted workflow used in this study involved data acquisition of randomized biological samples and quality controls. Compound identification was made using both NIST2017 nominal mass library and the Thermo Scientific™ Orbitrap™ GC-MS HRAM Metabolomics Library, the first commercially available high resolution accurate mass metabolomics library for electron ionization GC-MS.

INTRODUCTION

Meeting the demand for specialty chemical compounds for the pharmaceutical, agricultural, and manufacturing industries is one of the grand challenges of the modern chemical industry. This demand must be met under increasing regulatory scrutiny using environmentally friendly methodologies. Biotechnological approaches, powered by the techniques and concepts of synthetic biology, have the potential to deliver the necessary sustainable solutions. At its core, synthetic biology applies a design-build-test framework¹ to the redesigning of natural biological systems for beneficial purposes. By inserting and fine-tuning genetic information within microbial bio-factories (such as *Escherichia coli*² and *Streptomyces* spp.³), it is possible to assemble complex enzymatic pathways for rapid and diverse chemical production. Inducible bacterial promoters, such as the isopropyl β-D-1-thiogalactopyranoside (IPTG)-inducible lac promoters or tetracycline-inducible *tet*-promoters, are commonly used to strongly activate gene transcription, switching on engineered biosynthetic pathways. Yet, it is often not fully understood how these inducible systems interact with the global bacterial metabolome, potentially with toxic side effects. The aim of this study was to investigate how the application of untargeted metabolomics can be used to understand the metabolic effects of inducer type and concentration on the metabolic fingerprint of engineered bacteria (*E. coli* DH5α harboring an IPTG-inducible, red fluorescent protein expression plasmid pBbA1a-RFP) under different growth conditions. This phenotypic information has the potential to inform upstream genetic strategies while at the same time better defining the most efficient use of this promotor for biochemical pathway expression.

MATERIALS AND METHODS

Growth conditions: *Escherichia coli* DH5α from glycerol stocks were inoculated onto Lysogeny Broth (LB) agar plates followed by inoculation into either Terrific Broth (TB) or LB with 0.4% glucose and incubated overnight at 37 °C with shaking. Cultures were inoculated (1/100) into 1 mL of fresh broth in a 24-well plate and grown to mid-logarithmic phase, whereupon they were induced using a Hamilton® Multistar robotic system with variable levels of IPTG (25, 50, and 100 μM final concentration) and incubated for a further 24 h at 37 °C with shaking.

Sample preparation and derivatization: Following incubation, samples were quenched with 1 mL of cold methanol (−48 °C) to halt any enzymatic action within the bacteria and centrifuged for 15 min at 12,225 RCF to remove cellular debris from the media. Then, 100 μL of supernatant was filtered using a 0.45 μm syringe filter, combined with 100 μL of a 100 μg/mL internal standard solution of D-glucose and L-alanine-d7, and dried down under vacuum. Lyophilized pellets were then subjected to a common two-step sample derivatization method carried out by the initial addition of 50 μL of a 20 mg/mL methoxyamine/pyridine solution to enable the methoximation of any potentially labile ketone groups. Incubation at 65 °C for 40 min was followed by silylation in which 50 μL of MSTFA +1% TMCS (N-methyl-N-(trimethylsilyl)trifluoroacetamide + 1% trimethylchlorosilane) was added. Subsequent heating at 65 °C for 40 min afforded volatility to any labile hydroxyl and amine groups and the addition of the TMS (trimethylsilyl) moiety. The TMCS acted as a catalyst to ensure optimal TMS addition.

GC-MS analysis: In all experiments, a Q Exactive GC-MS/MS Orbitrap mass spectrometer was used. Sample injection into a hot split/splitless injector (280 °C) was performed using a Thermo Scientific™ TriPlus RSH™ autosampler, and chromatographic separation was obtained with a Thermo Scientific™ TRACE™ 1310 GC system and a Thermo Scientific™ TraceGOLD™ TG-5SIIIMS 30 m × 0.25 mm I.D. × 0.25 μm film capillary column (P/N 26096-1425). A total GC run time of 33 min per sample was used. Additional details of instrument parameters are shown in Table 1 and Table 2.

Data processing workflow for unknown metabolite detection and identification: Full-scan, lock mass corrected data were imported into Thermo Scientific™ Compound Discoverer™ software and subjected to a qualitative untargeted workflow that involved retention time alignment, normalization, and statistical analysis (principal component analysis and differential analysis).⁴ Compound identification was achieved using Thermo Scientific™ TraceFinder™ software following spectral deconvolution and using the Orbitrap GC-MS HRAM metabolomics library. In addition to this, the NIST 2017 nominal mass library was used to further extend the number of annotations assigned to putatively detected metabolites.

Table 1. GC and injector conditions.

TRACE 1310 GC System Parameters	
Injection Volume:	1.0 mL
Liner:	Single taper (P/N: 453A1345)
Inlet:	280 °C
Inlet Module and Mode:	SSL/SL, split 40:1
Carrier Gas:	He, 1.2 mL/min.
Oven Temperature Program:	
Temperature 1:	70 °C
Hold Time:	2 min.
Temperature 2:	325 °C
Rate:	10 °C/min.
Hold Time:	6 min.

Table 2. Mass spectrometer parameters.

Q Exactive GC Mass Spectrometer Parameters	
Transfer Line (°C):	280 °C
Ionization Type:	Electron Ionization (EI)
Ion Source:	250 °C
Electron Energy (eV):	70 eV
Acquisition Mode:	Full-scan
Mass Range:	50-550 Da
Mass Resolution:	60,000 FWHM
Lockmass:	207.03235 m/z

RESULTS

Relative metabolite levels were determined from *E. coli* culture media (LB and TB) following incubation in the presence and absence of IPTG at various concentration levels. In addition to these samples, pooled quality control (QC) samples were also analyzed. Raw data files were imported into Compound Discoverer software and grouped according to the treatment (IPTG) and media (LB and TB) (Figure 1). Data processing in Compound Discoverer involved a retention time alignment step to compensate for small differences in the retention times of the components in the sequence (Figure 1).




Figure 1. Example of retention time alignment in Compound Discoverer software for several peaks detected in *E. coli* DH5α cultures induced with IPTG and grown in LB media.

The component extraction (unknown detection using a ±5 ppm extraction window and a signal intensity threshold of 500,000 peak area counts) step was followed by data normalization to correct for potential batch effects. To identify class differences, data was subjected to principal component analysis (PCA) (Figure 2). In this case, the first two principal components explain 54% of the variance within the dataset, with PC1 (30%) dominated by differences between the two media types, and PC2 (24%) by differences between induced and uninduced cultures. By using such an approach, comparison of PCA loadings against blank media allows the identification of bacterial metabolites that differ between sample classes.

As complete group separation was noted within the PCA, a wholly unsupervised approach was adopted. The next three data processing steps were designed to select significant features that contributed to group differences, in this case LB control vs. LB IPTG treated. An analysis of variance test (ANOVA) was performed alongside a subsequent multiple comparison Tukey Honest Significance Difference (T-HSD) test. This supplied an adjusted p-value of compound significance that was subsequently used as input, alongside associated compound fold change, into the volcano plot tool available in Compound Discoverer software. This tool plotted log2 fold change vs. -log10 p-value and identified compounds that were important in group discrimination and also had a suitable large fold change (Figure 3). Significant compounds (2045 ions corresponding to 212 compounds selected based on p-values < 0.05 and log2 fold change values > 1) were then selected and sent to TraceFinder software for attempted identification using spectral matches against libraries/databases.




Figure 2. Centered and log₂-scaled Principal Component Analysis (PCA) scores plot (top) and loading plot (bottom). Data points within specific ellipses represent *E. coli* DH5α growing in either TB or LB media (Blue - TB media control, LB media control) and after the inclusion of various levels of the pBbA1a-RFP (IPTG) plasmid, under different media conditions (Orange - TB treated, LB treated). QCs (n=8, pooled samples) were also analyzed to test instrument and method performance.

Figure 3. Discriminatory analysis (Volcano plot) generated for LB IPTG 100 (green) and LB control (red) samples showing the compounds that significantly contributed to group difference to the left and right sides. The x-axis represents the log₂ of the fold change between the two sample groups, and the y-axis represents the log₁₀ of the adjusted ANOVA p-value. The top-ranking ions in each group are highlighted in blue.

Compound Identification using Orbitrap-GC HRAM Metabolomics Library

The overall goal of untargeted GC-MS metabolomics studies is to detect and annotate (identify) the metabolites responsible for group differences. This is usually accomplished by comparing the measured spectra against in-house standard databases or unit mass spectral libraries such as NIST or Wiley. Statistically significant features were sent to TraceFinder software and identified using both NIST 2017 and the Thermo Scientific Orbitrap GC-MS HRAM metabolomics library and retention time index derived from a C10–C19 alkane mix. This HRAM metabolomics library was created using pure metabolite standards analyzed on the Orbitrap-GC, and it contains ~850 unique metabolite spectra (each with retention time index, CAS numbers, and PubChem identifiers) acquired in EI using 70 eV and 60,000 resolution. An example of spectral matching is shown in Figure 4 for glycine trimethylsilyl ester (glycine 3TMS).




Figure 4. Glycine 3TMS identification using the Orbitrap GC-MS HRAM metabolomics library. Forward and reverse search indices in addition to accurate mass information add to the confidence in compound identification.

Compound annotation was achieved using a search index dot product value of >750, a total score 5 of >80, and a maximum retention time index difference (ΔRI) of 100 (measured versus expected). An example of the TraceFinder software deconvolution browser showing GABA 3TMS identified based on the criteria stated above is shown in Figure 5.




Figure 5. Example of metabolite identification in the TraceFinder software deconvolution browser showing a list of compounds (a), identified based on a total (average) score and retention index information (b), across the retention time aligned media samples (c), spectral match against the HRAM Orbitrap metabolomics library (c) as well as the deconvoluted spectrum (d) for GABA 3TMS are shown.

Following this process, 39 significant metabolites were confidently identified from the ANOVA and volcano plot analysis, and their corresponding peak area fold-changes in each sample were calculated (Table 3).

Table 3. Table of fold change of detected metabolites that significantly contributed to group differences in the LB IPTG 25, LB IPTG 50, and LB IPTG 100 groups. Green to red color gradients indicate the fold change of associated metabolite upon comparison to blank media. All comparisons are made between LB media control and each of the IPTG-treated samples (e.g., putrescine is upregulated 5-fold in LB IPTG 25 sample as compared to LB media control). The color and intensity of the boxes is used to represent changes of fold change. In the example below, red represents up-regulated metabolites, and green represents down-regulated metabolites.

RT	putative ID	m/z	Base peak	Average Score	LB IPTG 25LB	IPTG 50LB	IPTG 100LB
5.15	propylene glycol 2TMS	73.04691	97.3	-2.3	-2.0	-2.3	
6.06	(R)-3-hydroxybutyric acid 2TMS	147.0657	84.1	3.2	3.6	2.5	
7.89	D-isoleucine, N-acetyl TMS	86.09647	84.6	-0.8	-0.7	-0.6	
9.22	2-[Bis(trimethylsilyl)amino]ethyl bis(trimethylsilyl) phosphate	299.0715	98.9	1.5	-1.4	-1.7	
9.25	Leucine 2TMS	158.1358	98.9	3.5	3.1	3.4	
9.55	L-allo-isoleucine 2TMS	158.1358	96.5	6.9	7.0	7.1	
9.62	Proline 2TMS	142.10467	98.7	-2.1	-1.9	-1.7	
9.75	glycine 3TMS	174.11307	99.5	-0.7	0.3	0.3	
9.82	succinic acid 2TMS	147.0657	99.4	-1.7	-1.8	-2.3	
10	glyceric acid 3TMS	147.0657	99.1	-0.4	-0.3	-0.3	
10.1	uracil 2TMS	241.08238	95.5	1.4	1.6	1.6	
10.3	fumaric acid 2TMS	245.06578	96	5.2	5.3	5.2	
10.5	serine 3TMS	204.12357	98.4	-6.2	-6.0	-6.4	
11.2	L-methionine 3TMS	56.0496	94.6	-1.3	-1.0	-1.2	
11.8	cadaverine 4TMS	174.11308	92.6	-2.8	-3.1	-3.1	
11.9	malonic acid 2TMS	147.06569	94.6	-1.1	-0.7	-0.8	
12.5	putrescine 4TMS	174.11307	81	5.0	5.2	5.0	
12.6	pyroglutamic acid 2TMS	156.08385	98.6	-2.1	-1.3	-1.8	
12.7	GABA 3TMS	174.11292	93.4	-4.7	-3.8	-3.3	
12.9	Phenylalanine TMS	120.08079	91.5	-0.8	-0.6	-0.6	
13.7	glutamic acid 3TMS	246.13425	98.7	-1.6	-1.2	-1.0	
16	L-arginine 3TMS	157.11508	96.7	-4.1	-3.5	-3.7	
16.5	sorbose methoxyamine isomer 1	217.1075	97.6	-2.4	-2.3	-2.2	
16.6	sorbose methoxyamine isomer 2	217.1075	96.5	-1.4	-1.1	-1.2	
16.8	d-Mannose, 2,3,4,5,6-pentakis-O-(trimethylsilyl)-, o-methylxyme, (1Z)	147.06561	96.9	-2.0	-1.8	-1.7	
17	glucose 5TMS isomer 2	73.04685	97.5	-2.8	-2.5	-1.4	
17	L-lysine 4TMS	156.12053	98.4	-1.4	-1.2	-1.1	
17.2	tyrosine 3TMS	218.1029	98.7	-4.3	-3.8	-3.9	
17.3	ribose 4TMS	217.10722	84.8	-1.4	-1.2	-1.4	
17.6	mannose 5TMS isomer 1	204.09943	99	-2.3	-2.0	-1.9	
18.6	myoinositol 6 TMS	217.10748	99	-0.5	-0.2	-0.3	
18.9	guanine 3TMS	352.1441	95	1.8	2.0	1.8	
19.8	L-tryptophan 3TMS	202.10474	97.3	-1.7	-2.0	-1.9	
20.1	sorbose isomer 1	73.04684	90.6	-3.9	-3.8	-4.2	
21.7	uridine 2TMS	217.10728	84	-3.7	-3.5	-4.2	
23.2	adenosine 4TMS	230.11531	92.8	-4.8	-5.3	-5.8	
24	trehalose 8TMS	361.16836	97.4	-9.0	-8.8	-8.2	
24.1	guanoside 3TMS	324.13055	91.2	-1.8	-2.1	-2.1	
24.5	melibiose 3TMS isomer 1	204.09943	92.8	3.9	2.8	3.3	

CONCLUSIONS

The data obtained in this study clearly show significant changes in metabolism arising from IPTG induction and subsequent protein production on the metabolic fingerprints of *E. coli* cells. However, it is important to note that due to the proof-of-principle nature of this research, it is difficult to determine if induction or protein production and subsequent extracellular transportation caused the associated metabolic changes observed. A larger secondary experiment is required to fully validate and confirm the biological findings.

- The untargeted analysis leads to the identification of 39 significant metabolites that showed statistically significant differences between the exo-metabolome of *E. coli* grown in LB media and exposed to various levels of IPTG and the corresponding control.
- Metabolites identification was simplified by using a dedicated HRAM metabolomics library retention time index information, significantly increasing the confidence in compound identification, one of the most critical steps in metabolomics.
- Compound Discoverer software and TraceFinder software streamline data interrogation, proving qualitative and quantitative information and increase the confidence in the results.

Taken together, using the Q Exactive GC Orbitrap GC/MS/MS system operated in full-scan mode at high resolving power allows confident metabolite detection and identification, and eases the ability to discriminate between various biological sample groups.

REFERENCES

- Carbonell, P., et al. SYNBIOCHEM – a SynBio foundry for the biosynthesis and sustainable production of fine and specialty chemicals. *Biochem. Soc. Trans.* **2016**, *44* (3), 675–677.
- Toogood, H.S., et al. Enzymatic Menthol Production: One-Pot Approach Using Engineered *Escherichia coli*. *ACS Synth. Biol.* **2015**, *4* (10), 1112–1123.
- Muhamadali, H., et al. Metabolomics investigation of recombinant mTNFα production in *Streptomyces lividans*. *Microb. Cell Fact.* **2015**, *14* (1), 157.
- Sumner, L. W., et al. Proposed minimum reporting standards for chemical analysis Chemical Analysis Working Group (CAWG) Metabolomics Standards Initiative (MSI). *Metabolomics* **2007**, *3* (3), 211–221.
- Application of GC Orbitrap mass spectrometry for untargeted metabolomics of pathogenic microorganisms. Thermo Fisher Scientific Application Note 10532, 2016.
- Carneiro, S., et al. Metabolic footprint analysis of recombinant *Escherichia coli* strains during fed-batch fermentations. *Molec. BioSyst.* **2011**, *7* (3), 899–910.
- Muhamadali, H., et al. Metabolomic analysis of riboswitch containing *E. coli* recombinant expression system. *Molec. BioSyst.* **2016**, *12* (2), 350–361.

TRADEMARKS/LICENSING

© 2018 Thermo Fisher Scientific Inc. All rights reserved. All trademarks are the property of Thermo Fisher Scientific and its subsidiaries. Hamilton is a registered trademark of Hamilton Company. This information is presented as an example of the capabilities of Thermo Fisher Scientific products. It is not intended to encourage use of these products in any manner that might infringe the intellectual property rights of others.

PO10619