

Rapid and Comprehensive Proteome Analysis Using LC-FAIMS-MS/MS

Alexander S Hebert¹, Satendra Prasad², Michael W Belford², Derek J Bailey², Susan E Abbatiello², Romain Huguet², Graeme C McAlister², Eloy R Wouters²; Jean-Jacques Dunyach², Dain R Brademan⁴, Michael S Westphall¹, and Joshua J Coon^{3, 4, 5} ¹Genome Center of Wisconsin, University of Wisconsin-Madison, WI; ²Thermo Fisher Scientific, San Jose, CA; ³Department of Biomolecular Chemistry, University of Wisconsin-Madison, WI; ⁵Morgridge Institute for Research, Madison, WI

Conclusion

- 1) Single shot experiments with FAIMS produce >8,000 protein and >100,000 peptide identifications
- 2) 1,000 more proteins guantified with addition of FAIMS.
- 3) FAIMS gas phase fractionation works similarly or better than LC fractionation.



, or from LC fractions without FAIMS (blue), in 4 hours (�) or 6 hour (•) of instrument time. The top point consists of 2x3 hour runs, each with two CV internal steppir 45V|-75V and -60V|-90V). The top protein result for each method type is annotated



FAIMS sources transmit ions through inner and outer electrodes based on their difference in mobility when in a high or low electric field is generated from an asymmetric waveform. Since this field causes the ion to disperse, the maximum peak amplitude of the asymmetric waveform is referred to as the dispersion voltage (DV). Ions with a large enough difference in mobility between the high and low field migrate towards the electrodes, while ions with no or limited difference in mobility are transmitted. The trajectory of an ion may be altered by the addition of a DC voltage. Termed as the compensation voltage (CV), the selection of an appropriate DC level will compensate for the drift of a specific ion or group of ions, allowing them to pass through the device. By changing the CV (a), alternate groups of ions will pass. Thus, the CV provides a handle by which one can control which population of ions are traversing the FAIMS device. The mobility, and thus the appropriate (to use cannot be easily predicted and must be determined empirically. (b) Schematic of device used in these





All results are from 60 minute analyses of human cell line K562 derived tryptic peptides, analyzed on an Orbitrap Fusion Lumos except where noted.





