Reduction of Emissions and Operating Costs of GCMS Usage with a New Ecology Mode

BHIMADZU

Nicole M. Lock, Clifford M. Taylor, Richard R. Whitney, Ph.D. and Zhuangzhi 'Max' Wang, Ph.D. Shimadzu Scientific Instruments, Inc.

Introduction

For the first time in US history, the government requires greenhouse gas emissions to be reported to the EPA in an effort to reduce the United States total carbon footprint. The EPA 40 CFR 98 states that facilities that emit 25,000 metric tons or more per year of greenhouse gas (GHG) emissions will be required to submit an annual report. Not only is it the duty of large industrial leaders to reduce their GHG emissions, but reduction should also take place at all levels of industry.

One way to assist in the reduction of operating costs and emissions previously was by manually turning down temperatures, turning off unnecessary power to selected parts of analytical instrumentation, and reducing gas consumption in standby mode. In most cases, these could only be done either by downloading additional complex methods or by manually adjusting individual parameters.

A new feature in GCMS software called "Eco-mode" automatically performs the power- and gas-saving functions at the end of an analytical batch. Eco-mode is used to reduce gas flow and turn off unnecessary heated zones, projecting an estimated 35% in energy savings as well as reduction of carbon emissions by up to 1.1 tons per year.

In this work, daily electrical power consumption was measured and compared to that of typical GCMS operation. Measured power correlated to predicted energy savings and emission reductions very closely.

Ecology Simulation

Ecology Mode

Ecology mode = Eco mode

Rather than running a traditional 'Shut-down


method', we press a single button to send

maximum savings depending on the type of standby parameters they choose.

the instrument into a energy-saving standby

Shimadzu Corporation has developed a Comparison calculator that will allow you to compare power consumption of previous GCMS systems to a new GCMS-QP2010 Ultra using Eco Mode.

http://www.shimadzu.com/eco_sim/eco.htm



The Eco mode can be customized depending on the users' needs; however, GCMSsolution software

comes with a default Eco mode method. Each of the method parameters can be adjusted to give users the

Downloads the following parameters to the instrument when performing shutdown. ntrols the following units when performing OK Cancel Help In addition to the method setup, the parameters and instrument control allow more customization. Once these parameters are set, Eco mode can be manually started or can be set up to start at the end of an analytical batch. DK Cancel Help DK Cancel Help m 🗖 🥦


Power Consumption Results: Start-up

A GCMS system was powered on and started.

 Power was turned on to both the GC and MS

Parameters


- Auto start sequence was initiated
- System was allowed to pump down
- Heat was turned on and method was downloaded



Instrument Control

Example Analysis Conditions

Basic Analysis Setup



GCMS parameters were set up to represent a typical analysis.


- Split injection port temperature was set
- High-pressure injection was turned on, and the split ratio was set at 10:1
- Column oven ramp was programmed from 40°C to 330°C and back to start

GC Program Time: 22.40 min

Power Consumption Results: Post Analysis


System was tuned and a method was set up to run.

- System was auto tuned
- A method was started
- Seventeen samples were run over the course of two days



Eco Mode Conditions


Default Eco Mode Set-up



Power Consumption Results: Post Eco Mode

Eco mode was set up to be in effect all weekend.

- Default 'Shutdown' method
- parameters were used Eco mode was started
- Friday afternoon Eco mode was stopped Monday afternoon
- (roughly 3 days) System was returned to
- standard settings and allowed to run



Conclusions

The new GCMS-QP2010 Ultra's Eco mode is a valuable tool in today's economy. The system's ability to automatically enter a standby mode upon the completion of a batch allows for maximum savings. The flexibility of customized settings allows for turning on Eco mode at any time; therefore, the QP2010 Ultra saves energy more efficiently.

It was estimated that the overall consumption of electricity is 6,370 kWh/yr, and CO2 emission is

After checking the power usage in start-up, stand by, sampling and Eco mode, it was determined you will save on the energy consumption and overall operating costs

| CO2 emissions | 3,554 | kg/year |
|-----------------|-------|-----------------|
| Electricity | 6,370 | kWh / year |
| Including a PC | 573 | US\$/year |
| Sas consumption | 1.56 | Cylinder / year |
| | 468 | US\$/year |

The GC used three times less power while in Eco mode verses normal operation, while the MS used one and a half times less power during a weekend Eco mode state.

From the Laboratory

QP-2010 Ultra:

GC on average draws: 360 W in stand-by

1000 W while heating up 90 W in Eco Mode


MS on average draws: 450 W in stand-by

500 W while Tuning 420 W in Eco Mode




Eco mode operation TOTAL TIME: 90 HR

to the Home



all electronics consume power and add to our carbon foot print



Making popcorn on High: draws 979 W TOTAL TIME PLUGGED IN: 10 MIN

One cup coffee maker: Off (clock on): draws 4.7W On and heating up: 1376W Brewing (dispense hot water): 461W TOTAL TIME PLUGGED IN: 2 HR KWH = 0.12

On with Low (hot/warm): 413 W On with High (hot/warm): 914W TOTAL TIME PLUGGED IN: 15 MIN

Normal operation TOTAL TIME: 163 HR KWH (for GC): 60 KWH (for MS): 77

KWH (for MS): 38

at 225°C Interface temperature was set at MS to start data collection after 2min. with a scan speed of 3333 OK Cancel Help

00000

Cancel Help