

Determination of inborn errors of metabolism

Quantification of amino acids and acylcarnitines in dried blood spots by FIA-MS/MS for clinical research

Author: Mariana Barcenas,
Thermo Fisher Scientific, Les Ulis, France

Keywords: Amino acids, acylcarnitines,
dried blood spots, non-derivatized, mass
spectrometry, TSQ Fortis MS

Goal

Implementation of a flow injection analysis-tandem mass spectrometry (FIA-MS/MS) method for simultaneous quantification of 13 amino acids and 13 acylcarnitines in dried blood spots

Introduction

Newborn screening (NBS) is a public health program provided by most of the countries around the world aimed at screening newborns for a list of serious genetic and metabolic disorders.¹⁻⁴ Determination of the amino acids and acylcarnitines profile is an integral part of the newborn screening programs. In the last decade, significant improvements in mass spectrometry (MS) technology enabled widespread use of MS/MS for newborn screening. MS can offer differentiation between true and false positive samples and also facilitate screening for untreatable disorders—all with a single injection.



In this report, an MS-based analytical method for determination of amino acids and acylcarnitines in dried blood spots (DBS) for inborn errors of metabolism (IEM) research is reported. FIA was performed on a Thermo Scientific™ Vanquish™ Flex Binary UHPLC system. Detection was performed on a Thermo Scientific™ TSQ Fortis™ triple-stage quadrupole mass spectrometer with heated electrospray ionization (HESI) by selected reaction monitoring (SRM). Reagents and controls were obtained from the MassChrom™ LC-MS/MS Complete Kit for Amino Acids and Acylcarnitines in Dried Blood Spots (Ref 57000) from Chromsystems Instruments & Chemicals GmbH (Gräfelfing, Germany).

Table 1. List of analytes

| Amino acids | Amino acid internal standards | Acylcarnitines | Acylcarnitine internal standards |
|--------------------|--|------------------------------|----------------------------------|
| Alanine | d ₄ -Alanine | Carnitine (C0) | d ₉ -C0 |
| Arginine | d ₇ -Arginine | Acetylcarnitine (C2) | d ₃ -C2 |
| Aspartic acid | d ₃ -Aspartic acid | Propionylcarnitine (C3) | d ₃ -C3 |
| Citrulline | d ₂ -Citrulline | Butyrylcarnitine (C4) | d ₃ -C4 |
| Glutamic acid | d ₅ -Glutamic acid | Isovalerylcarnitine(C5) | d ₉ -C5 |
| Glycine | ¹³ C ₂ ¹⁵ N-Glycine | Glutarylacarnitine (C5DC) | d ₆ -C5DC |
| Leucine/Isoleucine | d ₃ -Leucine | Hexanoylcarnitine (C6) | d ₃ -C6 |
| Methionine | d ₃ -Methionine | Octanoylcarnitine (C8) | d ₃ -C8 |
| Ornithine | d ₆ -Ornithine | Decanoylcarnitine (C10) | d ₃ -C10 |
| Phenylalanine | d ₅ -Phenylalanine | Dodecanoylcarnitine (C12) | d ₃ -C12 |
| Proline | d ₇ -Proline | Tetradecanoylcarnitine (C14) | d ₃ -C14 |
| Tyrosine | d ₄ -Tyrosine | Hexadecanoylcarnitine (C16) | d ₃ -C16 |
| Valine | d ₈ -Valine | Octadecanoylcarnitine (C18) | d ₃ -C18 |

Experimental

Target analytes

Target analytes and their respective internal standards are summarized in Table 1.

Sample preparation

Reagents from Chromsystems included an extraction reagent and two controls, as well as the internal standards.

Samples were prepared as described by Chromsystems using the following procedure:

1. The internal standards were reconstituted using the extraction reagent.
2. 3.2 mm punches were placed in a well plate, and 100 µL of extraction reagent (containing the internal standards) were added to each well.
3. The plate was shaken at 600 rpm for 20 minutes.
4. Extracted samples were transferred to a clean plate.

Liquid chromatography

Flow injection analysis was performed on a Vanquish Flex Binary UHPLC system, using an injection volume of 10 µL of extracted sample. The mobile phase was provided by Chromsystems. Details of the analytical method are reported in Table 2. Total runtime was 1.5 minutes.

Table 2. LC method description

| Time (min) | Flow rate (mL/min) | A (%) |
|------------|--------------------|-------|
| 0.00 | 0.09 | 100 |
| 1.23 | 0.09 | 100 |
| 1.25 | 0.30 | 100 |
| 1.50 | 0.09 | 100 |

Mass spectrometry

Analytes and internal standards were detected by SRM on a TSQ Fortis triple-stage quadrupole mass spectrometer with heated electrospray ionization operated in positive ionization mode. A summary of the MS source conditions is reported in Table 3. SRM transitions with optimized collision energy and tube lens values are summarized in Table 4.

Method evaluation

The robustness, reliability, and reproducibility of the method were evaluated in terms of intra- and inter-assay precision and accuracy for all analytes.

Table 3. MS settings

| | |
|-------------------------------|---------------------------------------|
| Source type | Heated electrospray ionization (HESI) |
| Vaporizer temperature | 150 °C |
| Capillary temperature | 300 °C |
| Spray voltage (positive mode) | 3500 V |
| Sheath gas | 15 AU |
| Sweep gas | 0 AU |
| Auxiliary gas | 5 AU |
| Data acquisition mode | Selected-reaction monitoring (SRM) |
| Source fragmentation | 10 V |
| Collision gas pressure | 1.5 mTorr |
| Cycle time | 0.800 s |
| Q1 mass resolution (FWMH) | 0.7 |
| Q3 mass resolution (FWMH) | 0.7 |

Table 4. SRM transitions, collision energies, and tube lens values

| Analyte | Precursor ion | Product ion | Internal standard | Precursor ion | Product ion | Collision energy (V) | Tube lens (V) |
|---------------|---------------|-------------|--|---------------|-------------|----------------------|---------------|
| Alanine | 90.038 | 44.071 | d ₄ -Alanine | 94.075 | 48.125 | 12 | 84 |
| Arginine | 175.088 | 70.054 | d ₇ -Arginine | 182.162 | 77.179 | 24 | 85 |
| Aspartic acid | 134.038 | 116.125 | d ₃ -Aspartic acid | 137.175 | 119.196 | 8 | 83 |
| Citrulline | 176.162 | 113.125 | d ₂ -Citrulline | 178.162 | 115.125 | 17 | 60 |
| Glutamic acid | 148.088 | 130.054 | d ₅ -Glutamic acid | 153.138 | 135.125 | 10 | 75 |
| Glycine | 76.088 | 30.179 | ¹³ C ₂ ¹⁵ N-Glycine | 79.088 | 32.196 | 13 | 130 |
| Leucine | 132.125 | 86.125 | d ₃ -Leucine | 135.138 | 89.125 | 10 | 83 |
| Methionine | 150.050 | 132.982 | d ₃ -Methionine | 153.088 | 136.125 | 10 | 84 |
| Ornithine | 133.038 | 70.054 | d ₆ -Ornithine | 139.162 | 76.196 | 18 | 68 |
| Phenylalanine | 166.112 | 120.125 | d ₅ -Phenylalanine | 171.162 | 125.125 | 14 | 82 |
| Proline | 116.112 | 70.125 | d ₇ -Proline | 123.175 | 77.196 | 16 | 89 |
| Tyrosine | 182.162 | 136.125 | d ₄ -Tyrosine | 186.162 | 140.125 | 14 | 88 |
| Valine | 118.125 | 72.125 | d ₈ -Valine | 126.162 | 80.250 | 11 | 84 |
| C0 | 162.125 | 85.054 | d ₉ -C0 | 171.162 | 85.054 | 21 | 77 |
| C2 | 204.112 | 85.071 | d ₃ -C2 | 207.125 | 85.071 | 20 | 64 |
| C3 | 218.125 | 85.054 | d ₃ -C3 | 221.162 | 85.054 | 20 | 71 |
| C4 | 232.162 | 85.071 | d ₃ -C4 | 235.212 | 85.054 | 20 | 77 |
| C5 | 246.162 | 85.054 | d ₉ -C5 | 255.162 | 85.054 | 21 | 78 |
| C5DC | 276.175 | 85.071 | d ₆ -C5DC | 282.175 | 85.054 | 24 | 90 |
| C6 | 260.212 | 85.054 | d ₃ -C6 | 263.175 | 85.125 | 22 | 94 |
| C8 | 288.212 | 85.125 | d ₃ -C8 | 291.212 | 85.125 | 23 | 99 |
| C10 | 316.212 | 85.125 | d ₃ -C10 | 319.212 | 85.125 | 24 | 104 |
| C12 | 344.212 | 85.125 | d ₃ -C12 | 347.250 | 85.125 | 25 | 99 |
| C14 | 372.300 | 85.054 | d ₃ -C14 | 375.300 | 85.125 | 25 | 110 |
| C16 | 400.350 | 85.125 | d ₃ -C16 | 403.350 | 85.125 | 27 | 116 |
| C18 | 428.350 | 85.125 | d ₃ -C18 | 431.400 | 85.054 | 28 | 118 |

Intra-assay precision for each day was evaluated in terms of percentage coefficient of variation (%CV) using the controls at two different concentration levels in replicates of five (n=5). Inter-assay precision was evaluated as the %CV on the full set of samples (control samples at two levels in replicates of five prepared and analyzed on three different days, n=15). Analytical accuracy was evaluated in terms of percentage bias between nominal and average calculated concentrations using quality control samples at two different levels provided by Chromsystems (0192 and 0193 batches #1519 and #2518).

Data analysis

Data were acquired and processed using Thermo Scientific™ TraceFinder™ 4.1 software. Quantification of the analytes is done by comparison with the corresponding isotopically labeled internal standard, using the formula $Conc = A Area/IS Area \times IS conc$.

Results and discussion

Representative chromatograms for phenylalanine, acetylcarnitine, and the corresponding internal standards are reported in Figure 1.




Figure 1. FIA-MS/MS profiles for (a) phenylalanine, (b) d₅-phenylalanine, (c) acetylcarnitine, and (d) d₃-acetylcarnitine

The reported method showed good reproducibility, with the maximum intra- and inter-assay precision below 9% and 12.5%, respectively, for all analytes. C5DC showed slightly higher values (14.9% and 17.7% for intra- and inter-assay precision, respectively). Results for intra- and inter-assay precision are reported in Tables 5 and Table 6, respectively.

Analytical accuracy for all analytes was always within the acceptance range provided by the supplier, with values between -12.8% and 11.8% (Table 7).

Table 5. Analytical intra-assay precision for controls 0192 and 0193 batches #1519 and #2518

| Analyte | Control I | | | Control II | | |
|---------------|-----------|--------|--------|------------|--------|--------|
| | Day 1 | Day 2 | Day 3 | Day 1 | Day 2 | Day 3 |
| | CV (%) | CV (%) | CV (%) | CV (%) | CV (%) | CV (%) |
| Alanine | 5.4 | 3.5 | 5.3 | 4.7 | 5.3 | 5.5 |
| Arginine | 6.2 | 4.3 | 3.4 | 2.5 | 3.0 | 4.4 |
| Aspartic acid | 7.8 | 6.0 | 5.4 | 3.7 | 4.7 | 5.3 |
| Citrulline | 7.8 | 4.5 | 6.3 | 5.3 | 5.4 | 5.5 |
| Glutamic acid | 5.5 | 4.0 | 4.2 | 4.4 | 5.0 | 5.3 |
| Glycine | 8.3 | 7.8 | 8.5 | 6.7 | 6.3 | 6.8 |
| Leucine | 5.5 | 2.2 | 5.7 | 2.2 | 5.0 | 5.1 |
| Methionine | 8.0 | 6.5 | 7.9 | 6.1 | 4.8 | 6.1 |
| Ornithine | 5.4 | 5.9 | 5.5 | 2.4 | 3.9 | 3.9 |
| Phenylalanine | 5.8 | 2.6 | 5.9 | 2.6 | 4.3 | 4.6 |
| Proline | 5.2 | 2.5 | 4.9 | 3.5 | 4.3 | 4.7 |
| Tyrosine | 5.0 | 3.6 | 4.6 | 2.7 | 3.0 | 4.3 |
| Valine | 6.0 | 2.5 | 6.0 | 2.8 | 5.3 | 5.1 |
| C0 | 5.6 | 3.0 | 4.1 | 3.6 | 4.9 | 5.5 |
| C2 | 6.1 | 2.7 | 6.7 | 3.0 | 4.4 | 5.2 |
| C3 | 6.6 | 4.5 | 5.9 | 3.1 | 5.5 | 5.4 |
| C4 | 7.6 | 5.5 | 5.4 | 4.6 | 5.0 | 6.6 |
| C5 | 7.0 | 6.8 | 7.3 | 3.9 | 6.1 | 5.6 |
| C5DC | 14.9 | 12.8 | 12.3 | 9.4 | 10.6 | 11.9 |
| C6 | 8.7 | 5.6 | 8.5 | 5.1 | 6.3 | 6.0 |
| C8 | 7.4 | 5.8 | 8.7 | 3.6 | 4.6 | 6.2 |
| C10 | 8.2 | 5.0 | 8.1 | 6.2 | 5.3 | 5.8 |
| C12 | 5.2 | 4.1 | 7.8 | 3.3 | 5.6 | 6.9 |
| C14 | 6.3 | 3.3 | 7.5 | 3.0 | 5.1 | 6.2 |
| C16 | 6.3 | 2.8 | 6.1 | 2.9 | 5.9 | 6.1 |
| C18 | 6.2 | 3.6 | 6.5 | 3.5 | 4.2 | 5.4 |

Table 6. Analytical inter-assay precision results for controls 0192 and 0193 batches #1519 and #2518

| Analyte | Control I CV (%) | Control II CV (%) |
|----------------|-------------------------|--------------------------|
| Alanine | 8.1 | 5.7 |
| Arginine | 7.6 | 11.2 |
| Aspartic acid | 7.1 | 6.6 |
| Citrulline | 7.5 | 6.0 |
| Glutamic acid | 5.6 | 6.0 |
| Glycine | 8.7 | 7.7 |
| Leucine | 8.7 | 4.7 |
| Methionine | 12.1 | 5.9 |
| Ornithine | 5.9 | 8.3 |
| Phenylalanine | 8.3 | 4.6 |
| Proline | 7.5 | 5.4 |
| Tyrosine | 6.4 | 4.5 |
| Valine | 10.5 | 5.2 |
| C0 | 8.5 | 4.8 |
| C2 | 8.0 | 4.8 |
| C3 | 8.9 | 4.9 |
| C4 | 7.9 | 5.5 |
| C5 | 9.6 | 5.7 |
| C5DC | 17.7 | 13.2 |
| C6 | 10.6 | 5.6 |
| C8 | 9.9 | 5.1 |
| C10 | 12.3 | 5.4 |
| C12 | 11.3 | 5.6 |
| C14 | 11.0 | 5.5 |
| C16 | 9.9 | 5.8 |
| C18 | 7.6 | 6.1 |

Table 7. Accuracy for controls 0192 and 0193 batches #1519 and #2518

| Analyte | Control I (μmol/L) | Calculated conc. (μmol/L) | Bias (%) | Control I range (μmol/L) | Outcome | Control II (μmol/L) | Calculated conc. (μmol/L) | Bias (%) | Control II range (μmol/L) | Outcome |
|---------------|--------------------|---------------------------|----------|--------------------------|---------|---------------------|---------------------------|----------|---------------------------|---------|
| Alanine | 358 | 373 | 4.2 | 164–552 | Normal | 736 | 645 | -12.3 | 323–1149 | Normal |
| Arginine | 98.0 | 90.4 | -7.7 | 36.0–160 | Normal | 225 | 196 | -12.8 | 115–335 | Normal |
| Aspartic acid | 144 | 134 | -6.7 | 93.0–195 | Normal | 261 | 286 | 9.7 | 173–350 | Normal |
| Citrulline | 67.0 | 74.3 | 11 | 47.0–87.0 | Normal | 300 | 273 | -8.9 | 221–379 | Normal |
| Glutamic acid | 444 | 402 | -9.5 | 286–603 | Normal | 730 | 794 | 8.7 | 502–958 | Normal |
| Glycine | 218 | 221 | 1.2 | 158–278 | Normal | 649 | 611 | -5.8 | 456–842 | Normal |
| Leucine | 302 | 275 | -8.8 | 166–438 | Normal | 504 | 544 | 8.0 | 335–673 | Normal |
| Methionine | 60.0 | 55.4 | -7.7 | 19.0–101 | Normal | 191 | 206 | 7.8 | 76.0–306 | Normal |
| Ornithine | 199 | 179 | -10 | 117–281 | Normal | 547 | 488 | -10.8 | 343–751 | Normal |
| Phenylalanine | 149 | 158 | 5.8 | 95.0–203 | Normal | 436 | 456 | 4.6 | 269–603 | Normal |
| Proline | 311 | 304 | -2.3 | 229–393 | Normal | 774 | 846 | 9.3 | 475–1074 | Normal |
| Tyrosine | 160 | 153 | -4.7 | 108–212 | Normal | 556 | 515 | -7.4 | 381–731 | Normal |
| Valine | 239 | 225 | -5.7 | 150–328 | Normal | 424 | 435 | 2.7 | 278–570 | Normal |
| C0 | 51.8 | 48.5 | -6.4 | 28.9–74.7 | Normal | 101 | 112 | 10.9 | 63.0–139 | Normal |
| C2 | 23.2 | 25.9 | 11.8 | 15.9–30.5 | Normal | 60.1 | 63.4 | 5.4 | 37.2–83.0 | Normal |
| C3 | 5.03 | 5.00 | -0.7 | 2.91–7.15 | Normal | 13.0 | 12.4 | -5.0 | 8.44–17.6 | Normal |
| C4 | 0.940 | 0.860 | -8.0 | 0.440–1.44 | Normal | 3.29 | 3.60 | 9.1 | 1.92–4.66 | Normal |
| C5 | 0.540 | 0.560 | 4.6 | 0.280–0.800 | Normal | 2.17 | 2.23 | 2.8 | 1.22–3.12 | Normal |
| C5DC | 0.550 | 0.570 | 3.2 | 0.130–0.970 | Normal | 2.60 | 2.48 | -4.6 | 1.20–4.00 | Normal |
| C6 | 0.460 | 0.500 | 9.0 | 0.270–0.650 | Normal | 2.12 | 2.26 | 6.8 | 1.35–2.89 | Normal |
| C8 | 0.510 | 0.480 | -5.9 | 0.270–0.750 | Normal | 2.17 | 2.22 | 2.4 | 1.22–3.12 | Normal |
| C10 | 0.450 | 0.460 | 2.5 | 0.260–0.630 | Normal | 1.96 | 2.06 | 5.1 | 1.08–2.84 | Normal |
| C12 | 0.450 | 0.440 | -3.3 | 0.200–0.700 | Normal | 2.09 | 2.06 | -1.7 | 1.37–2.81 | Normal |
| C14 | 0.460 | 0.450 | -2.9 | 0.240–0.680 | Normal | 2.09 | 2.02 | -3.5 | 1.24–2.94 | Normal |
| C16 | 4.60 | 4.42 | -4.0 | 2.79–6.41 | Normal | 13.2 | 12.3 | -6.6 | 8.08–18.3 | Normal |
| C18 | 2.61 | 2.43 | -7.0 | 1.45–3.77 | Normal | 8.28 | 8.09 | -2.3 | 4.47–12.1 | Normal |

Conclusion

The MassChrom LC-MS/MS Complete Kit for Amino Acids and Acylcarnitines in DBS from Chromsystems was developed by an FIA-MS/MS workflow and implemented on a Vanquish Flex Binary UHPLC system coupled to a TSQ Fortis triple-stage quadrupole mass spectrometer. SRM data acquisition was used to provide high selectivity and sensitivity. The method proved to be robust and reliable, and met the required sensitivity requirements typically demanded by clinical research laboratories.

References

1. Lemonde, H. *Pediatr. Child Health*, **2014**, *25*, 103–107.
2. Marca, G.J. *Pharm. Biomed. Anal.*, **2014**, *101*, 174–182.
3. Lund, A.M. *Mol. Genet. Metab.*, **2012**, *107*, 281–293.
4. Sandlers, Y. *Transl. Res.*, **2017**, *189*, 65–75.

Find out more at thermofisher.com/clinicalresearch