

Explore what's new in Empower 3.8.0

Empower Software offers new features to enhance compliance and quality, simplify work for analysts and administrators, and increase productivity and efficiency, while providing new technology stack support for IT.

Reduce risk and variability in your chromatographic workflow, while simplifying work for your administrators with Empower Software.

Elevate your quality



Comply with USP Chapter <621> and Ph. Eur. 11th Edition formula and terminology requirements



New assignable privileges and policy updates



Record review of result Audit Trails



Quickly identify permission levels for users



Streamline and optimize your method validation workflow







Unique features that work specifically with the Alliance™ iS HPLC System and Waters eConnect™ HPLC Columns solve common challenges you face every day:





Prevent Errors with Method Matching: Empower Software intelligently determines which methods are applicable to your LC system.



Increase Traceability

with Column Tracking and History: **Empower Software** is your single electronic source of column-use data.



Prevent and Respond to Failures with **Pre-run Checks:** Empower Software and

the Alliance iS HPLC System work together to identify and alert you of common pre-analysis errors.



Increase Traceability with Empower System **Audit Trail:**

Log installation of and modifications to your Alliance iS HPLC System hardware and software.



Maximize Lab Productivity Anytime, Anywhere: waters_connect™ System Monitoring provides real-time information about your lab's Empower Software-controlled systems.

Simplify work for analysts and increase efficiency with Empower 3.8.0

