



# EXTEND THE LIFE OF YOUR NEW UHPLC COLUMN



Particulate matter in biological, food, environmental, and other “dirty” samples can cause premature wear to your valuable UHPLC columns.

You can extend the life of your analytical columns – **while maintaining fast LC performance** – by using an inexpensive guard column specifically designed for ZORBAX RRHD 1.8 µm (1200 bar) columns and ZORBAX RRHT 1.8 µm (600 bar) columns. Other advantages include:

- ▶ Less *time* and *money* spent replacing analytical columns
- ▶ Easy installation with no special tools or hardware required



# Guard columns can prevent untimely analytical column failure

## ORDERING INFORMATION

These Fast Guards for UHPLC work with:

- ZORBAX RRHD columns, 1.8  $\mu\text{m}$  (1200 bar)
- ZORBAX RRHT columns, 1.8  $\mu\text{m}$  (600 bar)



Size (mm)	Eclipse Plus C18	Eclipse XDB-C18	SB-C18	SB-C8
2.1	821725-901	821725-903	821725-902	821725-904
3.0	823750-901	823750-903	823750-902	823750-904
4.6	820750-901	820750-903	820750-902	820750-904

Get proof, and more details about how easy Fast Guards are to use, at [agilent.com/chem/fastguardsuhplc](http://agilent.com/chem/fastguardsuhplc) (look under “literature”)

This information is subject to change without notice.

© Agilent Technologies, Inc., 2013

Printed in the USA, July 30, 2013

820120-005



**Agilent Technologies**