

TraceFinder 4.1 SP4 System Requirements

The following table lists the hardware and software configurations for Thermo TraceFinder™ 4.1 SP4 operation.

Minimum System Requirements

Table 1. Minimum system requirements

System	Requirements
Hardware	<ul style="list-style-type: none">• 3.6 GHz dual-core processor enabled• 16GB RAM with system-managed memory enabled• DVD drive• Resolution display 1280 × 1024 (SXGA)• 5 GB available on drive C• NTFS format
Software	<ul style="list-style-type: none">• Microsoft™ Windows™ 7 SP1 Professional or Windows 10 (64-bit) English only• Microsoft .NET Framework 4.6.2• Adobe™ Reader™ 10.1• Thermo Scientific™ software:<ul style="list-style-type: none">– Xcalibur™ data system 4.1– Foundation™ platform 3.1 SP4– mzVault™ 1.0 SP1– FreeStyle™ 1.3– LC Devices 3.1 or later– Exactive™ Series 2.7 or later– TSQ Endura™ and TSQ Quantiva™ 2.0 or later– TSQ™ 2.6 or later– Fusion™ 1.2 SP1 or later– Aria™ MX 2.2 or later– GC Devices 3.2 or later– ISQ™ 3.2 SP1 or later– SII 1.3 for Xcalibur or later– TSQ 8000 3.2 or later– LTQ™ Series 2.7 SP4 or later

Trademarks

Aria, Exactive, Foundation, FreeStyle, Fusion, ISQ, mzVault, TraceFinder, and TSQ are trademarks; and LTQ, Thermo Scientific, TSQ Endura, TSQ Quantiva, and Xcalibur are registered trademarks of Thermo Fisher Scientific Inc. in the United States.

The following are registered trademarks in the United States and other countries: Microsoft and Windows are registered trademarks of Microsoft Corporation. Adobe and Reader are registered trademarks of Adobe Systems Inc.