

595 North Harrison Road
Bellefonte, PA 16823-0048 USA
Telephone 800-247-6628 ● 814-359-3441
Fax 800-447-3044 ● 814-359-3044

email: supelco@sial.com sigma-aldrich.com/supelco

Product Information

TO-14 Gas Standards

These high pressure, transportable gas calibration blends for US EPA Method TO-14 ensure that laboratories can perform Method TO-14 and Clean Air Act analyses at a fraction of their former costs. These blends are available only through Supelco.

Scott TO-14 gas blends are packaged in SCOTTY® HP high pressure, transportable, laboratory-size cylinders. A SCOTTY HP cylinder measures only 12¼" high and 3¼" in diameter (31 cm x 8.3 cm) , but contains 104 liters of gas at 1800 psi – enough gas to deliver adequate calibration for 80 or more uses. These products provide the reliability and performance of high pressure gas mixes without the maintenance and safety problems associated with large high pressure cylinders. SCOTTY HP cylinders are equipped with CGA 110/180 fittings, to plumb directly to your gas sampling system or through an appropriate regulator.

Scott TO-14 gas calibration blends are traceable to NIST reference mixtures (NIST SRM 1804), and are certified by Scott Gases for superior stability and accuracy.

Four blends are available. All products contain 100 ppb each analyte in nitrogen.

Ordering Information:

Gas Calibration Blends

TO-14 Calibration Mix 1

100 ppb each component in nitrogen.

Benzene Halocarbon 12
Bromomethane Halocarbon 113
Carbon tetrachloride Halocarbon 114

Chlorobenzene Hexachloro-1,3-butadiene Chloroform Methylene chloride

Chloromethane Styrene

1,2-Dibromoethane 1,1,2,2-Tetrachloroethane

1,2-Dichlorobenzene Tetrachloroethene

1,3-Dichlorobenzene Toluene

1,4-Dichlorobenzene 1,2,4-Trichlorobenzene
1,1-Dichloroethane 1,1,1-Trichloroethane
1,2-Dichloroethane 1,1,2-Trichloroethane
1,1-Dichloroethene Trichloroethene

cis-1,2-Dichloroethene 1,2,4-Trimethylbenzene 1,2-Dichloropropane 1,3,5-Trimethylbenzene

p-Xylene

cis-1,3-Dichloropropene Vinyl chloride trans-1,3-Dichloropropene m-Xylene Ethyl chloride o-Xylene

Ethylbenzene Halocarbon 11

SCOTTY HP cylinder 41900-U

Ordering information continued on reverse side of this sheet.







TO-14 Chlorinated Hydrocarbons Subset Mix

100 ppb each component in nitrogen.

Carbon tetrachloride Ethyl chloride
Chloroform Hexachloro-1,3-butadiene
Chloromethane Methylene chloride
1,1-Dichloroethane Tetrachloroethene

1,2-Dichloroethane1,1,2,2-Tetrachloroethane1,1-Dichloroethene1,1,1-Trichloroethanecis-1,2-Dichloroethene1,1,2-Trichloroethane1,2-DichloropropaneTrichloroethenecis-1,3-DichloropropeneVinyl chloride

trans-1,3-Dichloropropene

SCOTTY HP cylinder 41902

TO-14 CFC/HCFC Subset Mix

100 ppb each component in nitrogen.

Halocarbon 11 Halocarbon 113 Halocarbon 12 Halocarbon 114

SCOTTY HP cylinder 41903

TO-14 Aromatics Subset Mix

100 ppb each component in nitrogen.

Benzene Toluene
Chlorobenzene 1,2,4-Trichlorobenzene
1,2-Dichlorobenzene 1,3-Dichlorobenzene 1,3,5-Trimethylbenzene

1,4-Dichlorobenzenem-XyleneEthylbenzeneo-XyleneStyrenep-XyleneSCOTTY HP cylinder41901

TO-14 Reactive Subset Mix

100 ppb each component in nitrogen.

1,3-Butadiene 4-Ethyltoluene

3-Chloropropene

SCOTTY HP cylinder 41911

TO-14 GC/MS Tuning Standard

2 ppm bromofluorobenzene in nitrogen.SCOTTY HP cylinder 41913

Accessories For Use with SCOTTY HP Cylinders



994-0243

Model 6 Series Regulator with Gauge

Accepts inlet pressures up to 4000 psig. **501115**

Model 26 Pressure Regulator with Gauge

For use with Scotty HP cylinders, plus high pressure gauge. 41910-U

High Pressure Gauge

For use with Model 26 regulator. 41908

Stanc

Ensures your SCOTTY HP cylinder will be stable on a benchtop or other surface. **41909**

SCOTTY is a registered trademark of Scott Specialty Gases, Inc.

For expert answers to your questions contact our Technical Service Department:

Phone **800-359-3041**, **814-359-3041**Fax **800-359-3044**, **814-359-5468**E-mail **techservice@sial.com**

To download Supelco's free technical literature visit us at

sigma-aldrich.com/supelco-literature

Order/Customer Service 800-247-6628, 800-325-3010 Fax 800-325-5052 ● E-mail supelco@sial.com

Technical Service 800-359-3041, 814-359-3041

Fax 800-359-3044, 814-359-5468 • E-mail techservice@sial.com

SUPELCO ● 595 North Harrison Road, Bellefonte, PA 16823-0048 ● 814-359-3441

We are committed to the success of our Customers, Employees and Shareholders through leadership in Life Science, High Technology and Service.

The SIGMA-ALDRICH Family

