



# Analysis of Chlorinated Pesticides by GCxGC- $\mu$ ECD

Life Science and Chemical Analysis Centre  
LECO Corporation, Saint Joseph, Michigan USA



## Columns:

10 m x 0.18 mm x 0.2  $\mu$ m Rtx-5  
1.1 m x 0.10 mm x 0.1  $\mu$ m DB-17

## Detector:

ECD at 50 Hz

For questions on this analysis  
e-mail us at: [life\\_science@leco.com](mailto:life_science@leco.com)

For a complete listing of snapshots  
and application notes visit us on the  
web at [www.leco.com](http://www.leco.com)

