

For TOC Analyzers Shimadzu 40 mL Vial LE



The reliability of the vials in direct contact with the sample has a direct impact on the reliability of daily measurement results. Use the Shimadzu 40 mL vial LE (Low Elution) for TOC analyzers, which aims for the manufacturing quality standards of the manufacturers of TOC analyzers.

Optimal for Shimadzu TOC Analyzers for Laboratories

This vial can be used in autosamplers (ASI-L or ASI-V) for laboratory TOC analyzers (the TOC-L or TOC-V series).



TOC-L



TOC-Vw

Strict Quality Management

Our proprietary glass manufacturing technology and strict examinations have achieved a level of quality matching the requirements of pharmaceutical manufacturers.



Specifications

Product Name	Shimadzu 40 mL Vial LE Set
P/N	227-34200-01
Content	Shimadzu 40 mL Vial LE, 72PCs. Product Specification Sheet, and Vial identification Label provided. Cap 72PCs.; Septum 72PCs.

Product Name	Shimadzu 40 mL Vial LE, 72 pcs
P/N	227-34200-02
Content	Shimadzu 40 mL Vial LE, 72PCs. Product Specification Sheet, and Vial identification Label provided. Note: Caps and Septums not included.

Product Name	Cap for Shimadzu 40 mL Vial LE
P/N	227-34200-03
Content	Cap 72PCs.

Product Name	Septum for Shimadzu 40 mL Vial LE
P/N	227-34200-04
Content	Septum 72PCs.

Note: The vials in the Shimadzu 40 mL Vial LE set are not pre-cleaned.

Before using the vials, they must be cleaned.

If pre-cleaned vials with a TOC guarantee (<10 µgC/L) are required, contact Shimadzu separately.



Shimadzu GLC Ltd.

<https://solutions.shimadzu.co.jp/glc>
gssupport@glc.shimadzu.co.jp



Shimadzu Corporation

www.shimadzu.com/an/

For Research Use Only. Not for use in diagnostic procedures.

This publication may contain references to products that are not available in your country. Please contact us to check the availability of these products in your country.

Company names, products/service names and logos used in this publication are trademarks and trade names of Shimadzu Corporation, its subsidiaries or its affiliates, whether or not they are used with trademark symbol "TM" or "®".

Third-party trademarks and trade names may be used in this publication to refer to either the entities or their products/services, whether or not they are used with trademark symbol "TM" or "®".

Shimadzu disclaims any proprietary interest in trademarks and trade names other than its own.

Shimadzu does not assume any responsibility or liability for any damage, whether direct or indirect, relating to the use of this publication.