


Analysis of Gasoline by GCxGC-FID

Life Science and Chemical Analysis Centre
LECO Corporation, Saint Joseph, Michigan USA



GCxGC-FID separation of gasoline showing clear separation of compound classes.

Columns:

50 m x 0.2 mm x 0.5 μm J&W PONA
0.5 m x 0.1 mm x 0.1 μm BP-20

Detector:

FID

For questions on this analysis
e-mail us at: life_science@leco.com

For a complete listing of snapshots
and application notes visit us on the
web at www.leco.com

