



Detailed Analysis of Wine

Application Note

Food and Agriculture

Author

Laurent Maignial
Bureau National Interprofessionnel
du Cognac (BNIC)
Station Viticole
Cognac
France

Introduction

Wine can be injected without sample preparation on a very polar polyethylene glycol-type of stationary phase, Agilent J&W CP-Wax 57 CB. An almost complete separation of a wide range of alcohols, aldehydes, and esters is possible with this column. A reasonable separation of isomers 2-methyl-1-butanol and 3-methyl-1-butanol is also achievable, which is not possible on standard polyethylene glycol columns.

Conditions


Technique:	GC-capillary
Columns:	Agilent J&W CP-Wax 57 CB, 0.25 mm × 50 m, df = 1.2 µm, fused silica WCOT (Part no. CP97723)
Temperature:	35 °C (5 min) → 220 °C, 4 °C/min, 220 °C (10 min)
Carrier gas:	H ₂ , 140 kPa (1.4 bar, 20 psi)
Injector:	Split, T = 220 °C
Detector:	FID, T = 220 °C
Sample size:	1.0 µL
Concentration range:	10 to 1,000 ppm
Solvent sample:	ethanol/water



Agilent Technologies

Peak identification

1. acetaldehyde (ethanal)
2. ethyl acetate
3. methanol
4. 1-propanol
5. isobutanol
6. 1-butanol
7. 4-methyl-2-pentanol
8. 2-methyl-1-butanol
9. 3-methyl-1-butanol
10. acetoin
11. ethyl lactate
12. 1-hexanol
13. cos-3-hexen-1-ol (not detected)
14. S 2,3-butanediol
15. meso-2,3-butanediol
16. 2-phenylethanol
17. 1,5-pentanediol
18. glycerol



For More Information

These data represent typical results. For more information on our products and services, visit our Web site at www.agilent.com/chem.

www.agilent.com/chem

Agilent shall not be liable for errors contained herein or for incidental or consequential damages in connection with the furnishing, performance, or use of this material.

Information, descriptions, and specifications in this publication are subject to change without notice.

© Agilent Technologies, Inc., 2014
Printed in the USA
February 7, 2014
5991-4009EN