


Analysis of Pump Diesel by GCxGC-FID Using Consumable-Free Modulator

Life Science and Chemical Analysis Centre
LECO Corporation, Saint Joseph, Michigan USA



A Highly Structured GCxGC-FID Contour Plot Showing a Pump Diesel Sample. The inset shows an expanded 3D surface plot which displays the excellent peak shape attainable using the consumable-free modulator (CFM).

Columns:

50 m x 0.2 mm x 0.5 μ m DB-PONA
1.7 m x 0.10 mm x 0.10 μ m DB-Wax

Detector:

FID

For questions on this analysis
e-mail us at: life_science@leco.com

For a complete listing of snapshots
and application notes visit us on the
web at www.leco.com



Application Snapshot

LECO