


## Analysis of herbicides

BPX5

Column part number	054101		
Phase	BPX5	Rate 2	5°C/min to 260°C
Column	30 m x 0.25 mm x 0.25 µm	Final temperature	260°C, 10 min
Initial temperature	90°C, 1 min	Detector	NPD
Rate 1	30°C/min to 180°C	Carrier gas	He, 10 psi


 For more information visit [www.trajanscimed.com](http://www.trajanscimed.com) or contact [techsupport@trajanscimed.com](mailto:techsupport@trajanscimed.com)