


Analysis of polychlorinated dibenzofurans


BPX5

| Column Part Number | 054114 | | |
|---------------------|--------------------------|-------------------|------------------|
| Phase | BPX5 | Rate 3 | 5°C/min to 330°C |
| Column | 50 m x 0.22 mm x 0.25 µm | Final temperature | 330°C |
| Initial temperature | 180°C, 2 min | Carrier gas | He, 30 psi |
| Rate 1 | 4°C/min to 220°C, 12 min | Injection mode | Splitless, 270°C |
| Rate 2 | 5°C/min to 235°C, 7 min | Detector | MS |


Hexachlorodibenzofurans




Tetrachlorodibenzofurans



Heptachlorodibenzofurans




Pentachlorodibenzofurans



Analysis of polychlorinated dibenzofurans

BPX5

Octachlorodibenzofurans



For more information visit www.trajanscimed.com or contact techsupport@trajanscimed.com