


Analysis of triglyceride standards

HT5

| | | | |
|---------------------------|------------------------|----------------|----------------------------|
| Column part number | 054655 | | |
| Phase | HT5 | Carrier gas | H ₂ , 2 psi |
| Column | 6 m x 0.53 mm x 0.1 μm | Detector | FID |
| Initial temperature | 60°C, 0 min | Sensitivity | 32 x 10 ⁻¹² AFS |
| Rate | 10°C/min | Injection mode | On-column |
| Final temperature | 370°C, 5 min | | |



For more information visit www.trajanscimed.com or contact techsupport@trajanscimed.com